

www.hoggi.de

Made in Germany
by HOGGI 

BINGO EVOLUTION

Innovative and modern rehab pushchair.



Designed and hand built in Germany.

BINGO EVOLUTION *Mini*

Rehab pushchair in a proven concept.

BINGO Evolution Mini: The Rehab pushchair for early care in a proven concept.

- A feeling of security from the start
- Adaptable in proven manner
- Large growth capacity
- Modular mobility base

Finally it was logical, to offer **BINGO Evolution** as a version for our smallest ones. To provide security was the top priority without losing the need for a good positioning. This issue has been taken account by a strong lateral side support and a comfortable and soft upholstery, according to the lower user weight.

Adaptability and **growth capacity** are equivalent to the proven features of **BINGO Evolution**. A modular mobility base allows maximum economy. By the replacement of only a few small items, the next bigger seat unit - size 1 - can be applied which makes the mobility base a faithful companion over a long period.



Seat width	13 - 25 cm
Seat depth	12 - 27 cm
Tilt in space seat unit	
In driving direction	45°
Opposite to the driving direction	45°
Backrest height	37 - 51 cm
Backrest angle	80° - 180° in 10° steps
Shoulder height	25,5 - 39,5 cm
Lower leg length	16 - 30 cm
Knee angle	85° - 180°, stepless adjustable
Weight limit	25 kg