

BINGO Evolution Mini

Release date: February 2026 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address see invoice address

Company: _____

Street: _____

City, Zip code: _____



- Basic configuration**
- Individual Upholstery System (IUS)
 - Thigh supports
 - Footrest hanger angle adjustable
 - Angle adjustable Footplate
 - Side covers

BINGO Evolution Mini Trendline



- Trendline package**
- Upholstery in reed green
 - Canopy in reed green
 - Frame colour black
 - Hand loop
 - Whitewall rims (optional)



Frame colour black

Upholstery colour reed green

Frame colour



Black dull matt

Upholstery colour



Black



Berry soft



Berry soft // Black




Grey soft



Grey soft // Black



Turquoise soft



Turquoise soft // Black

Presettings

Seat width, Mini

enter seat width: (13 - 25 cm)

Back height, Mini

enter back height: (37 - 51 cm)

Seat depth, Mini

enter seat depth: (12 - 27 cm)

Lower leg length, Mini


enter lower leg length: (16 - 30 cm)



Notice:

If dimensions are missing,
we adjust to medium
dimensions.

Seat unit



BINGO Evolution Mini seat unit
incl. footrest hanger
angle adjustable

Headrest pads



Headrest pads,
low version



Headrest pads



Anatomical neck-headrest with
length compensator

Supports



Trunk support,
stepless height
adjustable (8 x 10 x 2 cm)



Thigh supports, wide- and
depth adjustable, removable



Hip supports, wide- and
depth adjustable, removable



Calf pad support plate
and foam insert

Notice:
Minimum seat width is 16 cm.

Calf pad support plate, Footplates



Footplate with heel stop

- 1=27 cm**
- 2=29 cm**
- 3=31 cm**


for short lower leg length

- 1=27 cm**
- 2=29 cm**
- 3=31 cm**



Footplate in boat form, width approx. 30 cm

for short lower leg lengths



Alu-footplate „footprint“ with heel stop

width **approx. 30 cm**

Footrest flip up lock, Steel footplate-carrier



Foot rest bracket angle adjustable, flip up type

with flip up lock

with flip up lock, for short lower leg lengths (reduction by 4 cm)



Reinforced steel footplate-carrier (about 0,3 kg heavier) angle adjustable, flip up type

with flip up lock

Belting



Lap belt



4-point lap belt



Five-point harness

size XS sw 18 - 24 cm



Groin strap



Foot straps



Ankle hugger

size 1 ø 4 - 6,5 cm



Fixation jacket

size 0

size 1




Chest-shoulder harness

size 0

size 1

size 2




Chest-shoulder harness
with 4 buckels

size 1

size 2

Abduction block



Abduction block, swing away and removable



Abduction block, fixed

Retainer for accessories, Grab rail



Retainer for accessories, grab rail or therapy tray



Grab rail with imitation leather, pivotable, black

Notice:
Only possible with retainer for accessories.


Therapy tray




Therapy tray, transparent



Therapy tray, transparent, height adjustable



Therapy tray, anthracite



Therapy tray, anthracite, height adjustable

Notice:
Only possible with retainer for accessories / with seat in the direction of travel.
Hip supports only possible in combination with the height adjustable therapy tray.

Winter warmer, Rain cape, Canopy



Winter warmer

size Mini
(up to sd 20)
size 1
(up to sd 18)



Sheep skin insert

!
Notice:
Possible from a seat
depth of 18 cm.




Rain cape, black




Canopy incl. adapter
and rain cover

!
Notice:
If the Trendline package
is selected, the canopy is
supplied in reed green.


Mobility bases



Rigid



Swivelling



Swivelling for additional equipment



Cross

Side covers




Side covers matching the upholstery color



Side covers matching the frame color

Frame suspension system

Without suspension system



Frame suspension system, fourfold adjustable



Please select the desired spring set:

Spring set white

Spring set black

Spring set matching the upholstery color

for load up to 15 kg

for load up to 25 kg

Front forks, Swivel lock (selectable for the Swivelling and the Cross mobility base)



Plastic front forks



Swivel lock for Plastic front forks




Aluminum front forks



Swivel lock for Aluminum front forks

Notice:
Cross mobility base only possible with aluminum front forks. These are included with the Cross mobility base.


Tires (selectable for the swivelling, the rigid mobility base)



Standard wheels



Drum brake wheels with tyres and with locking function



Polyurethan wheels with 7-spoke rims



Polyurethan wheels with whitewall rims

Pneumatic wheels

Polyurethan wheels

Pneumatic wheels

Polyurethan wheels

Notice:
Wheel removable via quick release.

Notice:
Only possible with the swivelling mobility base.
Wheel can be removed via quick-release fastener.

Notice:
Only possible with the **Trendline** version of the swivelling mobility base.
Wheel can be removed via quick-release fastener.

Tires (only selectable for the Cross mobility base)



Drum brake wheels with tyres on black 5-spoke rims, with locking function

Pneumatic wheels

Polyurethan wheels



Notice:

Wheel removable via quick release.

Brake variant



Attendant brake, two hand use



Notice:

Not for drum brake wheels.

Anti tip, Tie down kit



Anti tip left and right



Tie down kit according to ISO 7176-19



Notice:

Can only be selected for the Cross mobility base. Not possible with the mobility base for additional equipment.



Notice:

Only possible in conjunction with a Chest-shoulder harness or Fixation jacket as a positioning aid.

Accessory bags



Soft accessory bag (max. 3 kg)



Accessory bag (max. 8 kg)

Notice:
Not possible with tie down kit & only possible with seat facing the direction of travel.

Tray for additional equipment, Oxygen bottle holder



Tray for additional equipment



Oxygen bottle holder (d=100 mm)



Oxygen bottle holder for contoured oxygen devices

Notice:
Only possible with the swivelling mobility base to accommodate the tray for additional equipment.

for bottles up to 75 cm

for 2 liter bottles (short version)

Attachment left

Attachment right

for bottles up to 75 cm

for 2 liter bottles (short version)


Attachment left

Attachment right

Infusions bottle holder, Reinforced rear axle, Seat adapter



Infusion pole with 2-armed holder, articulated joint



Reinforced rear axle (10 kg more load capacity)



Seat adapter 10 cm heightened

Connecting parts



Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base

**Notice:**

Only possible with 2 mobility bases of one size.
Not possible in combination with aluminium front forks & mobility base Cross.



Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution swivelling or swivelling mobility base for additional equipment rear

**Notice:**

Only possible with the swivelling mobility base for additional equipment.
Not possible in combination with aluminium front forks.

Sun shade, Changing backpack, Cosies



Sun shade, black



Changing backpack, black



„Cosies“ Thermofleece hand warmer for companion

Buggy-Board



Buggy-Board

BASIC home underframe



BASIC home underframe

**with 75 mm
twin castors**

**with 125 mm
twin castors**



Side covers

in standard color white

**matching the upholstery
colour**

in black



Push handle for BASIC
home under frame

COBRA indoor base



COBRA indoor base

**with 75 mm
twin castors**

**with 125 mm
twin castors**

Presetting of the gas spring height adjustment:

up to 20 kg

up to 35 kg

up to 50 kg

Remark / Customer request / Individual special custom-made request

All prices are subject to change and non-binding. Errors, misprints, and price changes are reserved.
Binding prices are only those stated in our written offer or order confirmation.