



BINGO Evolution Mini rehab pushchair

Size Mini for children at the age of approx. 1 - 3 years

Order form incl. Exportprices (€) valid until revoked, release date 02/2025

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable, rotatable
- Maximum load 60 kg

Seat unit

- IUS - Individual Upholstery System
- Seat depth 12 - 27 cm
- Seat width 13 - 25 cm
- Lower leg length 16 - 30 cm
- Back height 37 - 51 cm
- Should height 25,5 - 39,5 cm
- Recline 80° - 180°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 25 kg

Minimum information (Without this information your order can not be handled.)

Frame colour:	<input type="checkbox"/> Black dull matt 032
Upholstery single-coloured:	<input type="checkbox"/> Black 346 <input type="checkbox"/> Berry soft 347 <input type="checkbox"/> Grey soft 348
with accessories colour black	<input type="checkbox"/> Turquoise soft 349
Upholstery two-coloured:	<input type="checkbox"/> Berry soft / Black 350 <input type="checkbox"/> Grey soft / Black 351 <input type="checkbox"/> Turquoise soft / Black 352
with accessories colour black	

Seat depth: (12 - 27 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height: (37 - 51 cm)

Lower leg length: (16 - 30 cm)

Basic variations - seat unit

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3230-2000	032	BINGO Evolution Mini seat unit		
			incl. footrest hanger angle adjustable	1 pc.	
with	<input type="checkbox"/> 3230-2900		High supports, size 1, wide- and depth-adjustable	01 1 pair	incl.
			* minimum seat width 16 cm		
or	<input type="checkbox"/> 3230-2600		Hip supports, size 1, wide- and depth-adjustable, removable	(upcharge) 02 1 pair	
with	<input type="checkbox"/> 3230-5200		Upholstery, size 1, adaptable (IUS)	03 1 set	incl.
	<input type="checkbox"/> 3230-5190		Calf pad support plate and foam insert	04 1 set	
with	<input type="checkbox"/> 3231-249x		Footplate with heel stop	05 1 pc.	incl.
			width: <input type="checkbox"/> 1=27 cm, <input type="checkbox"/> 2=29 cm, <input type="checkbox"/> 3=31 cm		
or	<input type="checkbox"/> 3231-2551		Footplate in boat form, width: approx. 30 cm	(upcharge) 05 1 pc.	
or	<input type="checkbox"/> 3231-2463	015	Alu-footplate "footprint"with heel stop, width 30cm	(upcharge) 06 1 pc.	
with	<input type="checkbox"/> 3240-2950		Footrest flip up lock	(upcharge) 07 1 pc.	
with	<input type="checkbox"/> SA-0430		Reinforced steel footplate-carrier (about 0,3 kg heavier)	08 1 pc.	incl.

Accessories - seat unit

	Orderno.	Colour	Description	Qty.	
	<input type="checkbox"/> 3230-7601		Headrest pads low version size 1	09 1 set	
	<input type="checkbox"/> 3230-7602		Headrest pads, size 1	10 1 set	
	<input type="checkbox"/> 3230-6900		Anatomical neck-headrest with length compens, size 1	11 1 set	
	<input type="checkbox"/> 3230-7700		Trunk supports, size 1, stepless height adjustable (8 x 10 x 2 cm)	12 1 pair	
	<input type="checkbox"/> 3230-6100		Lap belt Mini	13 1 pc.	
	<input type="checkbox"/> 3231-6701		4-point lap belt, size XS (for seat width 18 - 24cm)	14 1 pc.	
	<input type="checkbox"/> 3230-6200		Five-point harness Mini	15 1 pc.	
	<input type="checkbox"/> 3201-6300		Foot straps	16 1 pair	
	<input type="checkbox"/> 3201-6400		Ankle huger, size 1	17 1 pair	
	<input type="checkbox"/> 3230-5300		Groin strap, Mini	18 1 pc.	
	<input type="checkbox"/> 3201-5800		Fixation jacket, size 0	19 1 pc.	
	<input type="checkbox"/> 3201-5400		Fixation jacket, size 1	1 pc.	
	<input type="checkbox"/> 3230-5600		Chest-shoulder harness, size 0	20 1 pc.	
	<input type="checkbox"/> 3201-5600		Chest-shoulder harness, size 1	1 pc.	
	<input type="checkbox"/> 3201-5601		Chest-shoulder harness, size 1 with 4 buckels	21 1 pc.	
	<input type="checkbox"/> 3201-5700		Chest-shoulder harness, size 2	1 pc.	
	<input type="checkbox"/> 3201-5701		Chest-shoulder harness, size 2 with 4 buckels	1 pc.	
	<input type="checkbox"/> 3231-7800		Abduction block, size 1, swing away and removable	22 1 pc.	
	<input type="checkbox"/> 3231-7801		Abduction block, size 1, fixed	23 1 pc.	
	<input type="checkbox"/> 3230-2205	032	Retainer for accessories, grab rail or therapy tray	24 1 pair	
	<input type="checkbox"/> 3230-8102	306	Grab rail with imitation leather, pivotable, size Mini, black ¹⁾	25 1 pc.	
	<input type="checkbox"/> 3230-8400		Therapy tray, size 1, transparent ^{1) 2) 3)}	26 1 pc.	
	<input type="checkbox"/> 3230-8402		Therapy tray, size 1, transparent, height adjustable ^{1) 2)}	27 1 pc.	
	<input type="checkbox"/> 3230-8401		Therapy tray, size 1, grey ^{1) 2) 3)}	28 1 pc.	
	<input type="checkbox"/> 3230-8403		Therapy tray, size 1, grey, height adjustable ^{1) 2)}	1 pc.	
			¹⁾ only possible with retainer for accessories		
			²⁾ only possible with seat in driving direction		
			³⁾ not possible when using hip supports		
	<input type="checkbox"/> 3230-8300	306	Winter warmer Mini, black (up to SD 20)	29 1 pc.	
	<input type="checkbox"/> 3201-8300	306	Winter warmer, size 1, black (from SD 18)	1 pc.	
	<input type="checkbox"/> 7000-0018		Sheep skin insert, size 1	30 1 pc.	
	<input type="checkbox"/> 3201-9200	306	Rain cape, size 1, black	31 1 pc.	
	<input type="checkbox"/> 3230-9300	306	Canopy incl. adapter and rain cover, black	32 1 pc.	
	<input type="checkbox"/> 3002-0012	032	BASIC home underframe size Mini with 75 mm twin castors	33 1 pc.	
	<input type="checkbox"/> 3002-0013	032	BASIC home underframe size Mini with 125 mm castors	34 1 pc.	
	<input type="checkbox"/> 3002-9500	022	Side covers in standard color white	35	
	<input type="checkbox"/> 3002-9500		Side covers matching the upholstery colour	36	
	<input type="checkbox"/> 3002-9500	032	Side covers in black		
	<input type="checkbox"/> 3002-5400		Push handle for BASIC home under frame	1 pc.	
	<input type="checkbox"/> 3004-0002	032	COBRA indoor base, size Mini, with 75 mm twin castors	37 1 pc.	
	<input type="checkbox"/> 3004-0003	032	COBRA indoor base, size Mini, with 125 mm castors	34 1 pc.	
			Presetting of the gas spring height adjustment: <input type="checkbox"/> up to 20kg, <input type="checkbox"/> up to 35kg, <input type="checkbox"/> up to 50kg		

Basic variations - mobility base swivelling front wheels

	Orderno.	Colour	Description	Qty.	
	<input type="checkbox"/> 3230-1000	032	BINGO Evolution Mini mobility base, swivelling front wheels	1 pc.	
	<input type="checkbox"/> 3238-1011	032	BINGO Evolution Mini mobility base, size 1, swivelling front wheels,for respiratory equipment	1 pc.	
with	<input type="checkbox"/> 1281-1027		Plastic front forks	38	1 set incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	39	1 set
or	<input type="checkbox"/> 1281-1028		Aluminum front forks	(upcharge) 40	1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks		1 set
with	<input type="checkbox"/> 3231-1301		Without suspension system		1 set incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable	(upcharge) 41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	42	1 pair incl.
or	<input type="checkbox"/> 3231-9500	032	Side covers matching the frame color	43	1 pair incl.
with	<input type="checkbox"/> 1219-0004S		Pneumatic wheels with 6-spoke rims	44	1 set incl.
or	<input type="checkbox"/> 1219-0030		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge) 45	1 set
or	<input type="checkbox"/> 1229-0004S		Polyurethan wheels with 6-spoke rims	46	1 set incl.
or	<input type="checkbox"/> 1229-0009		Polyurethan wheels with 7-spoke rims	47	1 set incl.
or	<input type="checkbox"/> 1229-0030		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge) 48	1 set
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use	49	1 set
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	50	1 pc. incl.
or	<input type="checkbox"/> 3230-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge) 51	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3230-1501		Without tie down kit		1 set incl.
or	<input type="checkbox"/> 3230-7300		Tie down kit according to ISO 7176-19 (from 12 kg body weight) *		1 set
			* only possible with chest-shoulder harness or fixation jacket as positioning aid		

Basic variations - mobility base CROSS with Aluminium front forks

	Orderno.	Colour	Description	Qty.	
	<input type="checkbox"/> 3230-1030	032	BINGO Evolution CROSS Mini mobility base	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system		1 set incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable	(upcharge) 41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	42	1 pair incl.
or	<input type="checkbox"/> 3231-9500	032	Side covers matching the frame color	43	1 pair incl.
with	<input type="checkbox"/> 1219-0050		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	44	1 set incl.
or	<input type="checkbox"/> 1229-0050		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	48	1 set incl.
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	50	1 pc. incl.
or	<input type="checkbox"/> 3230-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge) 51	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3230-1501		Without tie down kit		1 set incl.
or	<input type="checkbox"/> 3230-7300	032	Tie down kit according to ISO 7176-19 (from 12 kg body weight) *		1 set
			* only possible with chest-shoulder harness or fixation jacket as positioning aid		
with	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks		1 set
with	<input type="checkbox"/> 3234-7100		Anti tip left and right *		1 set
			* not possible in combination with frame for respiratory equipment		

! BINGO Evolution Cross: Custom-made only on request !

Basic variations - mobility base rigid front wheels

	Orderno.	Colour	Description	Qty.	
	<input type="checkbox"/> 3231-1011	032	BINGO Evolution Mini mobility base, rigid front wheels	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	42	1 pair incl.
or	<input type="checkbox"/> 3231-9500	032	Side covers matching the frame color	43	1 pair incl.
with	<input type="checkbox"/> 1219-0025S		Pneumatic wheels with 6-spoke rims	44	1 set incl.
or	<input type="checkbox"/> 1219-0029		Drum brake wheels, pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	45	1 set
or	<input type="checkbox"/> 1229-0021S		Polyurethan wheels with 6-spoke rims	46	1 set incl.
or	<input type="checkbox"/> 1229-0029		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	48	1 set
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use	49	1 set
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	50	1 pc. incl.
or	<input type="checkbox"/> 3230-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge)	51	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3230-1501		Without tie down kit		1 set incl.
or	<input type="checkbox"/> 3230-7300	032	Tie down kit according to ISO 7176-19 (from 12 kg body weight) *		1 set
			* only possible with chest-shoulder harness or fixation jacket as positioning aid		

Basic variations - mobility base TRIKE (custom-made product)

	Orderno.	Colour	Description	Qty.	
	<input type="checkbox"/> 3236-1011	032	BINGO Evolution Mini mobility base TRIKE	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	42	1 pair incl.
or	<input type="checkbox"/> 3231-9500	032	Side covers matching the frame color	43	1 pair incl.
with	<input type="checkbox"/> 1219-0011S		Pneumatic wheels with 6-spoke rims	44	1 set incl.
or	<input type="checkbox"/> 1219-0031		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	45	1 set
or	<input type="checkbox"/> 1229-0011S		Polyurethan wheels with 6-spoke rims	46	1 set incl.
or	<input type="checkbox"/> 1229-0031		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	48	1 set
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use	49	1 set
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	50	1 pc. incl.
or	<input type="checkbox"/> 3230-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge)	51	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input checked="" type="checkbox"/> 3230-1501		Without tie down kit		1 set incl.

Accessories - mobility base

	Orderno.	Colour	Description	Qty.	
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	39	1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks		1 set
	<input type="checkbox"/> 3231-7400E		Frame suspension system, fourfold adjustable, retrofit (upcharge)	41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg		
	<input type="checkbox"/> 3230-9402E		Accessory bag, for "mini" mobility bases (max. 8 kg) *	51	1 pc.
	<input type="checkbox"/> 3230-7300E	032	Tie down kit according to ISO 7176-19 (from 12 kg body weight) *		1 set
			* only possible with chest-shoulder harness or fixation jacket as positioning aid		
	<input type="checkbox"/> 3231-9500E		Side covers matching the upholstery colour or the frame color		1 pair

Custom-made				
	Orderno.	Colour	Description	Qty.
<input type="checkbox"/>	3202-4100		Tray for additional equipment, size 1 (mobility base 3238-1000 only)	52 1 pc.
<input type="checkbox"/>	3231-4200		Oxygen bottle holder (d=100 mm) for bottles up to 75 cm *	53 1 pc.
			* assembly on the right side only	
<input type="checkbox"/>	SA-0225-MP		Adjustable Angle brackets for contoured oxygen units (upcharge)	54 1 set
<input type="checkbox"/>	3231-4201		Oxygen bottle holder (d=100 mm) for short bottles up to 40 cm *	1 pc.
			* assembly on the right side only	
<input type="checkbox"/>	SA-0031-KP		Infusions bottle holder	55 1 pc.
<input type="checkbox"/>	SA-0605-MP		Footrest carrier for short LLL	1 pc.
<input type="checkbox"/>	SA-0016/3-MP		Reinforced rear axle for mobility base max. load 70 kg (upcharge)	56 1 pc.
<input type="checkbox"/>	SA-0298-MP		Seat adapter 10 cm heightened (upcharge)	1 set
<input type="checkbox"/>	SA-0321-MP		Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base * (upcharge)	57 1 set
			* only possible with 2 mobility bases of one size, not for base 3230-1000	
			* mobility base BINGO Evolution Mini only possible at the front	
			* not possible in combination with aluminium front forks & mobility base Cross	
<input type="checkbox"/>	SA-0441-MP		Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution mobility base 3238-1011 rear * (upcharge)	1 set
			* not possible in combination with aluminium front forks & mobility base Cross	

Extra equipment				
	Orderno.	Colour	Description	Qty.
<input type="checkbox"/>	7000-0010		Buggy-Board	58 1 pc.
<input type="checkbox"/>	7000-0012	306	Sun shade, black	59 1 pc.
<input type="checkbox"/>	7000-0006		Headrest speakers, for headrest pads 3201-7601	1 set
<input type="checkbox"/>	7000-0202	306	Changing backpack, black	60 1 pc.
<input type="checkbox"/>	7000-0199		"Cosies" Thermofleece hand warmer for companion	61 1 pair

Pre-assembly of pads, upholstery parts & belts				
	Orderno.	Colour	Description	Qty.
<input type="checkbox"/>	6999-1012		BINGO Evolution Pre-assembly service (without settings)	1 pc.


CONFIGURATIONS // ACCESSORIES // BINGO Evolution Mini




Thigh supports



Hip supports




IUS - Upholstery



Calf pad support plate and foam insert



Footplate with heel stop or boat form



Alu-footplate „foot-print“ with heel stop



Footrest flip up lock



Reinforced steel footplate-carrier



Headrest pads low version



Headrest pads high version



Anatomical neck-headrest



Trunk supports



Lap belt



4-point lap belt



Five-point harness



Foot straps



Ankle huger



Groin strap



Fixation jacket



Chest-shoulder harness



Chest-shoulder harness with 4 buckles



Abduction block, swing away and removable




Abduction block, fixed




Retainer for accessories



Grap rail



Therapy tray



Therapy tray, height adjustable



Therapy tray, grey



Winter warmer



Sheep skin insert

CONFIGURATIONS // ACCESSORIES // BINGO Evolution Mini



Rain cape



Canopy



BASIC home underframe,
75 mm twin castors



125 mm castors



Side covers in standard
color white



Side covers matching the
upholstery colour



COBRA indoor base,
75 mm twin castors



Plastic front forks



Swivel lock



Aluminum front forks



Frame suspension
system



Side covers matching the
upholstery color



Side covers matching the
frame color



Pneumatic wheel,
front wheel



Pneumatic wheel,
rear wheel



Pneumatic tyres front
wheel, Drum brake wheel



Pneumatic tyres rear
wheel, Drum brake wheel



Polyurethan wheel,
front wheel



Polyurethan wheel,
rear wheel



Polyurethan wheel front
wheel, black rim



Polyurethan wheel rear
wheel, black rim



Polyurethan tyres front
wheel, Drum brake wheel



Polyurethan tyres rear
wheel, Drum brake wheel




Attendant brake,
two hand use



Soft accessory bag



Accessory bag



Tray for additional
equipment



Oxygen bottle holder




Adjustable Angle
brackets



Infusions bottle holder

CONFIGURATIONS // ACCESSORIES // BINGO Evolution Mini

56



Reinforced rear axle

57



Connecting parts for mobility bases

58



Buggy-Board

59



Sun shade, grey

60



Changing backpack

61



„Cosies“ Thermofleece hand warmer