



BINGO Evolution Mini rehab pushchair

Size Mini for children at the age of approx. 1 - 3 years

Order form incl. Exportprices (€) valid until revoked, release date 04/2024

| | |
|--------------------------------|--------------------------------|
| <input type="checkbox"/> Offer | <input type="checkbox"/> Order |
| Order number / reference | Date |
| Invoice address | |
| Company | Customer number |
| Street | City, Zip code |
| E-Mail | Telephone |
| Fax | |
| Delivery address | |
| Company | |
| Street | City, Zip code |



Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable, rotatable
- Maximum load 60 kg

Seat unit

- IUS - Individual Upholstery System
- Seat depth 12 - 27 cm
- Seat width 13 - 25 cm
- Lower leg length 16 - 30 cm
- Back height 37 - 51 cm
- Should height 25,5 - 39,5 cm
- Recline 80° - 180°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 25 kg

Minimum information (Without this information your order can not be handled.)

| | |
|------------------------------------|---|
| Frame colour: | <input type="checkbox"/> Black dull matt 032 |
| Upholstery single-coloured: | <input type="checkbox"/> Black 346 <input type="checkbox"/> Berry soft 347 <input type="checkbox"/> Grey soft 348 |
| with accessories colour black | <input type="checkbox"/> Turquoise soft 349 |
| Upholstery two-coloured: | <input type="checkbox"/> Berry soft / Black 350 <input type="checkbox"/> Grey soft / Black 351 <input type="checkbox"/> Turquoise soft / Black 352 |
| with accessories colour black | |

Seat depth: (12 -27 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height: (37 - 51 cm)

Lower leg length: (16 - 30 cm)

Basic variations - seat unit

| | Orderno. | Colour | Description | Qty. | Price |
|------|------------------------------------|--------|---|-----------------------------|-------|
| | <input type="checkbox"/> 3230-2000 | 032 | BINGO Evolution Mini seat unit, size 1 | 1 pc. | |
| with | <input type="checkbox"/> 3230-2900 | | Thigh supports, size 1, wide- and depth-adjustable * minimum seat width 16 cm | 01 1 pair | incl. |
| or | <input type="checkbox"/> 3230-2600 | | Hip supports, size 1, wide- and depth-adjustable, removable | (upcharge) 02 1 pair | |
| with | <input type="checkbox"/> 3230-5200 | | Upholstery, size 1, adaptable (IUS) | 03 1 set | incl. |
| | <input type="checkbox"/> 3230-5190 | | Calf pad support plate and foam insert | 04 1 set | |
| with | <input type="checkbox"/> 3231-249x | | Footplate with heel stop width: <input type="checkbox"/> 1=27 cm, <input type="checkbox"/> 2=29 cm, <input type="checkbox"/> 3=31 cm | 05 1 pc. | incl. |
| or | <input type="checkbox"/> 3231-2551 | | Footplate in boat form, width: approx. 30 cm | (upcharge) 05 1 pc. | |
| or | <input type="checkbox"/> 3231-2463 | 015 | Alu-footplate "footprint" with heel stop, width 30cm | (upcharge) 06 1 pc. | |
| with | <input type="checkbox"/> 3240-2950 | | Footrest flip up lock | (upcharge) 07 1 pc. | |
| with | <input type="checkbox"/> SA-0430 | | Reinforced steel footplate-carrier (about 0,3 kg heavier) | 08 1 pc. | incl. |

Accessories - seat unit

| | Orderno. | Colour | Description | Qty. | |
|--|------------------------------------|--------|---|------------------|--|
| | <input type="checkbox"/> 3230-7601 | | Headrest pads low version size 1 | 09 1 set | |
| | <input type="checkbox"/> 3230-7602 | | Headrest pads, size 1 | 10 1 set | |
| | <input type="checkbox"/> 3230-6900 | | Anatomical neck-headrest with length compens, size 1 | 11 1 set | |
| | <input type="checkbox"/> 3230-7700 | | Trunk supports, size 1, stepless height adjustable (8 x 10 x 2 cm) | 12 1 pair | |
| | <input type="checkbox"/> 3230-6100 | | Lap belt Mini | 13 1 pc. | |
| | <input type="checkbox"/> 3231-6701 | | 4-point lap belt, size XS (for seat width 18 - 24cm) | 14 1 pc. | |
| | <input type="checkbox"/> 3230-6200 | | Five-point harness Mini | 15 1 pc. | |
| | <input type="checkbox"/> 3201-6300 | | Foot straps | 16 1 pair | |
| | <input type="checkbox"/> 3201-6400 | | Ankle huger, size 1 | 17 1 pair | |
| | <input type="checkbox"/> 3230-5300 | | Groin strap, Mini | 18 1 pc. | |
| | <input type="checkbox"/> 3201-5800 | | Fixation jacket, size 0 | 19 1 pc. | |
| | <input type="checkbox"/> 3201-5400 | | Fixation jacket, size 1 | 1 pc. | |
| | <input type="checkbox"/> 3230-5600 | | Chest-shoulder harness, size 0 | 20 1 pc. | |
| | <input type="checkbox"/> 3201-5600 | | Chest-shoulder harness, size 1 | 1 pc. | |
| | <input type="checkbox"/> 3201-5601 | | Chest-shoulder harness, size 1 with 4 buckels | 21 1 pc. | |
| | <input type="checkbox"/> 3201-5700 | | Chest-shoulder harness, size 2 | 1 pc. | |
| | <input type="checkbox"/> 3201-5701 | | Chest-shoulder harness, size 2 with 4 buckels | 1 pc. | |
| | <input type="checkbox"/> 3231-7800 | | Abduction block, size 1, swing away and removable | 22 1 pc. | |
| | <input type="checkbox"/> 3231-7801 | | Abduction block, size 1, fixed | 23 1 pc. | |
| | <input type="checkbox"/> 3230-2201 | 032 | Retainer for accessories, grab rail or therapy tray | 24 1 pair | |
| | <input type="checkbox"/> 3230-8100 | | Grab rail with upholstery, size 1 ¹⁾ | 25 1 pc. | |
| | <input type="checkbox"/> 3230-8200 | | Therapy tray, size 1, transparent ^{1) 2) 3)} | 26 1 pc. | |
| | <input type="checkbox"/> 3230-8202 | | Therapy tray, size 1, transparent, height adjustable ^{1) 2)} | 27 1 pc. | |
| | <input type="checkbox"/> 3230-8201 | | Therapy tray, size 1, grey ^{1) 2) 3)} | 1 pc. | |
| | <input type="checkbox"/> 3230-8203 | | Therapy tray, size 1, grey, height adjustable ^{1) 2)} | 1 pc. | |
| | | | ¹⁾ only possible with retainer for accessories | | |
| | | | ²⁾ only possible with seat in driving direction | | |
| | | | ³⁾ not possible when using hip supports | | |
| | <input type="checkbox"/> 3230-8300 | 306 | Winter warmer Mini, black (up to SD 20) | 28 1 pc. | |
| | <input type="checkbox"/> 3201-8300 | 306 | Winter warmer, size 1, black (from SD 18) | 1 pc. | |
| | <input type="checkbox"/> 7000-0018 | | Sheep skin insert, size 1 | 29 1 pc. | |
| | <input type="checkbox"/> 3201-9200 | 306 | Rain cape, size 1, black | 30 1 pc. | |
| | <input type="checkbox"/> 3230-9300 | 306 | Canopy incl. adapter and rain cover, black | 31 1 pc. | |
| | <input type="checkbox"/> 3002-0012 | 032 | BASIC home underframe size Mini with 75 mm twin castors | 32 1 pc. | |
| | <input type="checkbox"/> 3002-0013 | 032 | BASIC home underframe size Mini with 125 mm castors | 33 1 pc. | |
| | <input type="checkbox"/> 3002-9500 | 022 | Side covers in standard color white | 34 | |
| | <input type="checkbox"/> 3002-9500 | | Side covers matching the upholstery colour | 35 | |
| | <input type="checkbox"/> 3002-9500 | 032 | Side covers in black | | |
| | <input type="checkbox"/> 3002-5400 | | Push handle for BASIC home under frame | 1 pc. | |
| | <input type="checkbox"/> 3004-0002 | 032 | COBRA indoor base, size Mini, with 75 mm twin castors | 36 1 pc. | |
| | <input type="checkbox"/> 3004-0003 | 032 | COBRA indoor base, size Mini, with 125 mm castors | 33 1 pc. | |
| | | | Presetting of the gas spring height adjustment: <input type="checkbox"/> up to 20kg, <input type="checkbox"/> up to 35kg, <input type="checkbox"/> up to 50kg | | |

Basic variations - mobility base swivelling front wheels

| | Orderno. | Colour | Description | Qty. |
|------|-------------------------------------|--------|--|-----------------|
| | <input type="checkbox"/> 3230-1000 | 032 | BINGO Evolution Mini mobility base, swivelling front wheels | 1 pc. |
| | <input type="checkbox"/> 3238-1011 | 032 | BINGO Evolution Mini mobility base, size 1, swivelling front wheels,for respiratory equipment | 1 pc. |
| with | <input type="checkbox"/> 1281-1027 | | Plastic front forks | 37 1 set incl. |
| | <input type="checkbox"/> 3202-7100 | | Swivel lock for Plastic front forks | 38 1 set |
| or | <input type="checkbox"/> 1281-1028 | | Aluminum front forks (upcharge) | 39 1 set |
| | <input type="checkbox"/> 3232-7100 | | Swivel lock for Aluminum front forks | 1 set |
| with | <input type="checkbox"/> 3231-1301 | | Without suspension system | 1 set incl. |
| or | <input type="checkbox"/> 3231-7400 | | Frame suspension system, fourfold adjustable (upcharge) | 40 1 set |
| | <input type="checkbox"/> 1131-9017 | 022 | spring set (white) for load up to 15 kg | |
| | <input type="checkbox"/> 1131-9017 | | spring set (in upholstery color) foar load up to 15 kg | |
| | <input type="checkbox"/> 1131-9018 | 022 | spring set (white) for load up to 25 kg | |
| | <input type="checkbox"/> 1131-9018 | | spring set (in upholstery color) foar load up to 25 kg | |
| with | <input type="checkbox"/> 3231-9500 | | Side covers matching the upholstery color (Standard) | 41 1 pair incl. |
| or | <input type="checkbox"/> 3231-9500 | 032 | Side covers matching the frame color | 42 1 pair incl. |
| with | <input type="checkbox"/> 1219-0004S | | Pneumatic wheels with 6-spoke rims | 43 1 set incl. |
| or | <input type="checkbox"/> 1219-0030 | | Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge) | 44 1 set |
| or | <input type="checkbox"/> 1229-0004S | | Polyurethan wheels with 6-spoke rims | 45 1 set incl. |
| or | <input type="checkbox"/> 1229-0009 | | Polyurethan wheels with 7-spoke rims | 46 1 set incl. |
| or | <input type="checkbox"/> 1229-0030 | | Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge) | 47 1 set |
| | <input type="checkbox"/> 3231-7201 | | Attendant brake, two hand use | 48 1 set |
| with | <input type="checkbox"/> 3232-9401 | | Soft accessory bag, size 1 (max. 3 kg) | 49 1 pc. incl. |
| or | <input type="checkbox"/> 3230-9402 | | Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge) | 50 1 pc. |
| | | | * do not combine with tie down kit, use with seat in driving direction only | |
| with | <input type="checkbox"/> 3230-1501 | | Without tie down kit | 1 set incl. |
| or | <input type="checkbox"/> 3230-7300 | | Tie down kit according to ISO 7176-19 (from 12 kg body weight) * | 1 set |
| | | | * only possible with chest-shoulder harness or fixation jacket as positioning aid | |

Basic variations - mobility base CROSS with Aluminium front forks

| | Orderno. | Colour | Description | Qty. |
|------|------------------------------------|--------|--|-----------------|
| | <input type="checkbox"/> 3230-1030 | 032 | BINGO Evolution CROSS Mini mobility base | 1 pc. |
| with | <input type="checkbox"/> 3231-1301 | | Without suspension system | 1 set incl. |
| or | <input type="checkbox"/> 3231-7400 | | Frame suspension system, fourfold adjustable (upcharge) | 40 1 set |
| | <input type="checkbox"/> 1131-9017 | 022 | spring set (white) for load up to 15 kg | |
| | <input type="checkbox"/> 1131-9017 | | spring set (in upholstery color) foar load up to 15 kg | |
| | <input type="checkbox"/> 1131-9018 | 022 | spring set (white) for load up to 25 kg | |
| | <input type="checkbox"/> 1131-9018 | | spring set (in upholstery color) foar load up to 25 kg | |
| with | <input type="checkbox"/> 3231-9500 | | Side covers matching the upholstery color (Standard) | 41 1 pair incl. |
| or | <input type="checkbox"/> 3231-9500 | 032 | Side covers matching the frame color | 42 1 pair incl. |
| with | <input type="checkbox"/> 1219-0050 | | Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function | 43 1 set incl. |
| or | <input type="checkbox"/> 1229-0050 | | Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function | 47 1 set incl. |
| with | <input type="checkbox"/> 3232-9401 | | Soft accessory bag, size 1 (max. 3 kg) | 49 1 pc. incl. |
| or | <input type="checkbox"/> 3230-9402 | | Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge) | 50 1 pc. |
| | | | * do not combine with tie down kit, use with seat in driving direction only | |
| with | <input type="checkbox"/> 3230-1501 | | Without tie down kit | 1 set incl. |
| or | <input type="checkbox"/> 3230-7300 | 032 | Tie down kit according to ISO 7176-19 (from 12 kg body weight) * | 1 set |
| | | | * only possible with chest-shoulder harness or fixation jacket as positioning aid | |
| with | <input type="checkbox"/> 3232-7100 | | Swivel lock for Aluminum front forks | 1 set |
| with | <input type="checkbox"/> 3234-7100 | | Anti tip left and right * | 1 set |
| | | | * not possible in combination with frame for respiratory equipment | |

! BINGO Evolution Cross: Custom-made only on request !

Basic variations - mobility base rigid front wheels

| | Orderno. | Colour | Description | Qty. | |
|------|-------------------------------------|--------|--|--------------|--------------|
| | <input type="checkbox"/> 3231-1011 | 032 | BINGO Evolution Mini mobility base, rigid front wheels | 1 pc. | |
| with | <input type="checkbox"/> 3231-1301 | | Without suspension system | 1 set | incl. |
| or | <input type="checkbox"/> 3231-7400 | | Frame suspension system, fourfold adjustable (upcharge) | 40 | 1 set |
| | <input type="checkbox"/> 1131-9017 | 022 | spring set (white) for load up to 15 kg | | |
| | <input type="checkbox"/> 1131-9017 | | spring set (in upholstery color) foar load up to 15 kg | | |
| | <input type="checkbox"/> 1131-9018 | 022 | spring set (white) for load up to 25 kg | | |
| | <input type="checkbox"/> 1131-9018 | | spring set (in upholstery color) foar load up to 25 kg | | |
| with | <input type="checkbox"/> 3231-9500 | | Side covers matching the upholstery color (Standard) | 41 | 1 pair incl. |
| or | <input type="checkbox"/> 3231-9500 | 032 | Side covers matching the frame color | 42 | 1 pair incl. |
| with | <input type="checkbox"/> 1219-0025S | | Pneumatic wheels with 6-spoke rims | 43 | 1 set incl. |
| or | <input type="checkbox"/> 1219-0029 | | Drum brake wheels, pneumatic tyres on black 5-spoke rims, with locking function (upcharge) | 44 | 1 set |
| or | <input type="checkbox"/> 1229-0021S | | Polyurethan wheels with 6-spoke rims | 45 | 1 set incl. |
| or | <input type="checkbox"/> 1229-0029 | | Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge) | 47 | 1 set |
| | <input type="checkbox"/> 3231-7201 | | Attendant brake, two hand use | 48 | 1 set |
| with | <input type="checkbox"/> 3232-9401 | | Soft accessory bag, size 1 (max. 3 kg) | 49 | 1 pc. incl. |
| or | <input type="checkbox"/> 3230-9402 | | Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge) | 50 | 1 pc. |
| | | | * do not combine with tie down kit, use with seat in driving direction only | | |
| with | <input type="checkbox"/> 3230-1501 | | Without tie down kit | | 1 set incl. |
| or | <input type="checkbox"/> 3230-7300 | 032 | Tie down kit according to ISO 7176-19 (from 12 kg body weight) * | | 1 set |
| | | | * only possible with chest-shoulder harness or fixation jacket as positioning aid | | |

Basic variations - mobility base TRIKE (custom-made product)

| | Orderno. | Colour | Description | Qty. | |
|------|---|--------|--|--------------|--------------|
| | <input type="checkbox"/> 3236-1011 | 032 | BINGO Evolution Mini mobility base TRIKE | 1 pc. | |
| with | <input type="checkbox"/> 3231-1301 | | Without suspension system | 1 set | incl. |
| or | <input type="checkbox"/> 3231-7400 | | Frame suspension system, fourfold adjustable (upcharge) | 40 | 1 set |
| | <input type="checkbox"/> 1131-9017 | 022 | spring set (white) for load up to 15 kg | | |
| | <input type="checkbox"/> 1131-9017 | | spring set (in upholstery color) foar load up to 15 kg | | |
| | <input type="checkbox"/> 1131-9018 | 022 | spring set (white) for load up to 25 kg | | |
| | <input type="checkbox"/> 1131-9018 | | spring set (in upholstery color) foar load up to 25 kg | | |
| with | <input type="checkbox"/> 3231-9500 | | Side covers matching the upholstery color (Standard) | 41 | 1 pair incl. |
| or | <input type="checkbox"/> 3231-9500 | 032 | Side covers matching the frame color | 42 | 1 pair incl. |
| with | <input type="checkbox"/> 1219-0011S | | Pneumatic wheels with 6-spoke rims | 43 | 1 set incl. |
| or | <input type="checkbox"/> 1219-0031 | | Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge) | 44 | 1 set |
| or | <input type="checkbox"/> 1229-0011S | | Polyurethan wheels with 6-spoke rims | 45 | 1 set incl. |
| or | <input type="checkbox"/> 1229-0031 | | Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge) | 47 | 1 set |
| | <input type="checkbox"/> 3231-7201 | | Attendant brake, two hand use | 48 | 1 set |
| with | <input type="checkbox"/> 3232-9401 | | Soft accessory bag, size 1 (max. 3 kg) | 49 | 1 pc. incl. |
| or | <input type="checkbox"/> 3230-9402 | | Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge) | 50 | 1 pc. |
| | | | * do not combine with tie down kit, use with seat in driving direction only | | |
| with | <input checked="" type="checkbox"/> 3230-1501 | | Without tie down kit | | 1 set incl. |

Accessories - mobility base

| | Orderno. | Colour | Description | Qty. | |
|--|-------------------------------------|--------|--|------|--------|
| | <input type="checkbox"/> 3202-7100 | | Swivel lock for Plastic front forks | 38 | 1 set |
| | <input type="checkbox"/> 3232-7100 | | Swivel lock for Aluminum front forks | | 1 set |
| | <input type="checkbox"/> 3231-7400E | | Frame suspension system, fourfold adjustable, retrofit (upcharge) | 40 | 1 set |
| | <input type="checkbox"/> 1131-9017 | 022 | spring set (white) for load up to 15 kg | | |
| | <input type="checkbox"/> 1131-9017 | | spring set (in upholstery color) foar load up to 15 kg | | |
| | <input type="checkbox"/> 1131-9018 | 022 | spring set (white) for load up to 25 kg | | |
| | <input type="checkbox"/> 1131-9018 | | spring set (in upholstery color) foar load up to 25 kg | | |
| | <input type="checkbox"/> 3230-9402E | | Accessory bag, for "mini" mobility bases (max. 8 kg) * | 50 | 1 pc. |
| | <input type="checkbox"/> 3230-7300E | 032 | Tie down kit according to ISO 7176-19 (from 12 kg body weight) * | | 1 set |
| | | | * only possible with chest-shoulder harness or fixation jacket as positioning aid | | |
| | <input type="checkbox"/> 3231-9500E | | Side covers matching the upholstery colour or the frame color | | 1 pair |

| Custom-made | | | | |
|--------------------------|-----------------|---------------|---|-----------------|
| | Orderno. | Colour | Description | Qty. |
| <input type="checkbox"/> | 3202-4100 | | Tray for additional equipment, size 1 (mobility base 3238-1000 only) | 51 1 pc. |
| <input type="checkbox"/> | 3231-4200 | | Oxygen bottle holder (d=100 mm) for bottles up to 75 cm * | 52 1 pc. |
| | | | * assembly on the right side only | |
| <input type="checkbox"/> | SA-0225-MP | | Adjustable Angle brackets for contoured oxygen units (upcharge) | 53 1 set |
| <input type="checkbox"/> | 3231-4201 | | Oxygen bottle holder (d=100 mm) for short bottles up to 40 cm * | 1 pc. |
| | | | * assembly on the right side only | |
| <input type="checkbox"/> | SA-0031-KP | | Infusions bottle holder | 54 1 pc. |
| <input type="checkbox"/> | SA-0605-MP | | Footrest carrier for short LLL | 1 pc. |
| <input type="checkbox"/> | SA-0016/3-MP | | Reinforced rear axle for mobility base max. load 70 kg (upcharge) | 55 1 pc. |
| <input type="checkbox"/> | SA-0298-MP | | Seat adapter 10 cm heightened (upcharge) | 1 set |
| <input type="checkbox"/> | SA-0321-MP | | Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base * (upcharge) | 56 1 set |
| | | | * only possible with 2 mobility bases of one size, not for base 3230-1000 | |
| | | | * mobility base BINGO Evolution Mini only possible at the front | |
| | | | * not possible in combination with aluminium front forks & mobility base Cross | |
| <input type="checkbox"/> | SA-0441-MP | | Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution mobility base 3238-1011 rear * (upcharge) | 1 set |
| | | | * not possible in combination with aluminium front forks & mobility base Cross | |

| Extra equipment | | | | |
|--------------------------|-----------------|---------------|---|------------------|
| | Orderno. | Colour | Description | Qty. |
| <input type="checkbox"/> | 7000-0010 | | Buggy-Board | 57 1 pc. |
| <input type="checkbox"/> | 7000-0012 | 306 | Sun shade, black | 58 1 pc. |
| <input type="checkbox"/> | 7000-0006 | | Headrest speakers, for headrest pads 3201-7601 | 1 set |
| <input type="checkbox"/> | 7000-0202 | 306 | Changing backpack, black | 59 1 pc. |
| <input type="checkbox"/> | 7000-0199 | | "Cosies" Thermofleece hand warmer for companion | 60 1 pair |

| Pre-assembly of pads, upholstery parts & belts | | | | |
|---|-----------------|---------------|--|-------------|
| | Orderno. | Colour | Description | Qty. |
| <input type="checkbox"/> | 6999-1012 | | BINGO Evolution Pre-assembly service (without settings) | 1 pc. |

CONFIGURATIONS // ACCESSORIES // BINGO Evolution Mini




01 Thigh supports



02 Hip supports




03 IUS - Upholstery



04 Calf pad support plate and foam insert



05 Footplate with heel stop or boat form



06 Alu-footplate „foot-print“ with heel stop



07 Footrest flip up lock



08 Reinforced steel footplate-carrier



09 Headrest pads low version



10 Headrest pads high version



11 Anatomical neck-headrest



12 Trunk supports



13 Lap belt



14 4-point lap belt



15 Five-point harness



16 Foot straps



17 Ankle huger



18 Groin strap



19 Fixation jacket



20 Chest-shoulder harness




21 Chest-shoulder harness with 4 buckles



22 Abduction block, swing away and removable




23 Abduction block, fixed




24 Retainer for accessories



25 Grap rail with upholstery



26 Therapy tray



27 Therapy tray, height adjustable



28 Winter warmer



29 Sheep skin insert



30 Rain cape

CONFIGURATIONS // ACCESSORIES // BINGO Evolution Mini



Canopy



BASIC home underframe, 75 mm twin castors



125 mm castors



Side covers in standard color white



Side covers matching the upholstery colour



COBRA indoor base, 75 mm twin castors



Plastic front forks



Swivel lock



Aluminum front forks



Frame suspension system



Side covers matching the upholstery color



Side covers matching the frame color



Pneumatic wheel, front wheel



Pneumatic wheel, rear wheel



Pneumatic tyres front wheel, Drum brake wheel



Pneumatic tyres rear wheel, Drum brake wheel



Polyurethan wheel, front wheel



Polyurethan wheel, rear wheel



Polyurethan wheel front wheel, black rim



Polyurethan wheel rear wheel, black rim



Polyurethan tyres front wheel, Drum brake wheel



Polyurethan tyres rear wheel, Drum brake wheel




Attendant brake, two hand use



Soft accessory bag



Accessory bag



Tray for additional equipment




Oxygen bottle holder



Adjustable Angle brackets



Infusions bottle holder



Reinforced rear axle

CONFIGURATIONS // ACCESSORIES // BINGO Evolution Mini

56



Connecting parts for mobility bases

57



Buggy-Board

58



Sun shade, grey

59



Changing backpack

60



„Cosies“ Thermofleece hand warmer