

BINGO Evolution

Release date: April 2025 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address

see invoice address

Company: _____

Street: _____

City, Zip code: _____



Basic configuration

- Individual Upholstery System (IUS)
- Thigh supports
- Footrest hanger 90°
- Angle adjustable Footplate
- Side covers

Frame colour



Black dull matt



Silver


Upholstery colour




Black




Berry soft




Berry soft // Black




Grey soft



Grey soft // Black



Turquoise soft



Turquoise soft // Black



Presettings // size 1

Seat width, size 1

enter seat width: (18 - 30 cm)

Back height, size 1

enter back height: (45 - 63 cm)

Seat depth, size 1

enter seat depth: (16 - 34 cm)

Lower leg length, size 1

enter lower leg length: (16 - 37 cm)

**Notice:**

If dimensions are missing,
we adjust to medium
dimensions.

Presettings // size 2

Seat width, size 2

enter seat width: (23 - 35 cm)

Back height, size 2

enter back height: (54 - 75 cm)*

Seat depth, size 2

enter seat depth: (22 - 40 cm)

Lower leg length, size 2

enter lower leg length: (16 - 40 cm)

**Notice:**

If dimensions are missing,
we adjust to medium
dimensions.

* From a back height of
65 cm, the long push
handle is recommended.

Presettings // size 2 XL ! Considered as custom-made with corresponding accessories !

Seat width, size 2 XL

enter seat width: (23 - 35 cm)

Back height, size 2 XL

enter back height: (62 - 83 cm)**

Seat depth, size 2 XL

enter seat depth: (30 - 45 cm)

Lower leg length, size 2 XL

enter lower leg length: (16 - 40 cm)

**Notice:**


If dimensions are missing,
we adjust to medium
dimensions.

** From a back height
of 75 cm, the long push
handle is recommended.

BINGO Evolution Pre-assembly

Pre-assembly service (without settings)

Seat unit



BINGO Evolution seat unit

su size 1**su size 2****su size 2 XL**

Headrest pads

Headrest pads,
low version

Headrest pads

Headrest pads with spacer
fabric, for speakersAnatomical neck-headrest
with length compens

Supports



Trunk supports

8 x 10 x 2 cm**10 x 12 x 2 cm****12 x 13 x 2 cm**Thigh supports, wide- and depth-
adjustable, removableHip supports, wide- and depth-
adjustable, removable**Notice:**Size 1: Minimum seat width
is 20 cm.Size 2 / Size 2 XL: Minimum
seat width is 25 cm.

Calf pad support plate, Footrest hanger

Calf pad support plate
and foam insert


Footrest hanger 90°

Footrest hanger angle
adjustable 85 - 170°**Notice:**Minimum lower leg length of
16 cm possible.**Notice:**Minimum lower leg length of
19 cm possible.

**Footrest hanger short LLL
with bracket and footrest
flip up lock,**
Lower leg length from
15-36 cm possible

**Footrest hanger short LLL
with bracket and footrest
flip up lock,**
Lower leg length from
15-36 cm possible

Footplates



Footplate with heel stop

1=27 cm

2=29 cm

3=31 cm

4=33 cm


5=35 cm



Footplate in boat form

su width approx. 30 cm
size 1su width approx. 35 cm
size 2
size 2 XL

Alu foot plate in boat form

su width approx. 30 cm
size 1su width approx. 35 cm
size 2
size 2 XLAlu-footplate „footprint“
with heel stopsu width approx. 30 cm
size 1su width approx. 34 cm
size 2
size 2 XL**Notice:**The width of 5=35 cm is
only possible for seat unit
size 2 and size 2 XL.

Footrest flip up lock, Steel footplate-carrier



Footrest flip up lock



Reinforced steel footplate-carrier
(about 0,3 kg heavier)

Belting



Lap belt



4-point lap belt



Five-point harness

size XS
sw 18 - 24 cm

size S
sw 22 - 30 cm

size M
sw 28 - 36 cm



Groin strap



Foot straps



Ankle hugger

su **size 1**
size 1 ø 4 - 6,5 cm

su **size 2**
size 2 ø 6 - 8,5 cm
size 2 XL

Belting



Fixation jacket

size 0

size 1

size 2

size 3



Chest-shoulder harness

size 1

size 2

size 3



Chest-shoulder harness
with 4 buckels

size 1

size 2

size 3

Abduction block



Abduction block,
swing away and removable

size 1

size 2

size 3



Abduction block, fixed

size 1

size 2

size 3



Notice:

The abduction block in size 3 is only for the seat unit in size 2 and size 2 XL possible.

Retainer for accessories, Grab rail



Retainer for accessories,
grab rail or therapy tray



Grab rail with imitation leather,
pivotable, black




Notice:

Only possible with retainer for accessories.


Therapy tray



Therapy tray, transparent

Therapy tray, transparent,
height adjustable

Therapy tray, anthracite

Therapy tray, anthracite,
height adjustable**Notice:**

Only possible with retainer for accessories / with seat in the direction of travel.
Hip supports only possible in combination with the height adjustable therapy tray.

Winter warmer, Rain cape, Canopy



Winter warmer

size 1**size 2****size 3**


Sheep skin insert

size 1**size 2****size 3**


Rain cape, black

Canopy incl. adapter
and rain cover


Mobility bases



Rigid



Swivelling




Swivelling for additional equipment



Cross



Cross for additional equipment



Trike

**Notice:**

The mobility bases **Rigid** and **Trike** are only available in **size 1** and **2**.

Side covers



Side covers matching
the upholstery color



Side covers matching
the frame color

Frame suspension system

Without suspension system



Frame suspension system, fourfold adjustable



Please select the desired spring set:

Spring set white

for load up to 15 kg

Spring set black

for load up to 25 kg

Spring set matching the upholstery color

for load above to 25 kg



Notice:

Only available for the mobility base in **size 1**.

Front forks, Swivel lock (selectable for the Swivelling and the Cross mobility base)



Plastic front forks



Swivel lock for Plastic front forks



Aluminum front forks



Swivel lock for Aluminum front forks



Notice:

Cross mobility base only possible with aluminum front forks. These are included with the Cross mobility base.

Tires (selectable for the swivelling, the rigid and the Trike mobility base)

Standard wheels

Pneumatic wheels**Polyurethan wheels****Notice:**

Wheel removable via quick release.

Drum brake wheels with
tyres and with locking function**Pneumatic wheels****Polyurethan wheels**Polyurethan wheels with
7-spoke rims**Notice:**Only possible with the swivel-
ling mobility base.Wheel can be removed via
quick-release fastener.**Tires** (only selectable for the Cross mobility base)Drum brake wheels with tyres on
black 5-spoke rims, with locking
function**Pneumatic wheels****Polyurethan wheels****Notice:**Wheel removable via quick
release.**Brake variant**

Attendant brake, two hand use

**Notice:**

Not for drum brake wheels.

Anti tip, Push handle, Tie down kit



Anti tip left and right



Push handle extended by 7 cm



Tie down kit according to ISO 7176-19

**Notice:**

Can only be selected for the Cross mobility base. Not possible with the mobility base for additional equipment

**Notice:**

Not possible with the Trike mobility base.

Accessory bags



Soft accessory bag (max. 3 kg)




Accessory bag (max. 8 kg)

**Notice:**

Not possible with tie down kit & only possible with seat facing the direction of travel.

Tray for additional equipment, Oxygen bottle holder



Tray for additional equipment



Oxygen bottle holder (d=100 mm)



Oxygen bottle holder for
contoured oxygen devices

for bottles up to 75 cm

for 2 liter bottles
(short version)

for bottles up to 75 cm

for 2 liter bottles
(short version)

Attachment left

Attachment right

Attachment left

Attachment right




Notice:

Only possible with the swivelling
& CROSS mobility base to accom-
modate the tray for additional
equipment.

Infusions bottle holder, Reinforced rear axle, Seat adapter



Infusion pole with 2-armed holder,
articulated joint



Reinforced rear axle
(10 kg more load capacity)



Seat adapter 10 cm heightened

Mobility base and seat wider (only possible for the 2 XL mobility base)



Mobility base and seat 5 cm wider



Notice:

Only possible in black frame
color.
Additional costs apply for all
widened accessories.

Connecting parts



Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base



Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution swivelling or swivelling mobility base for additional equipment rear



Connecting parts for coupling a BINGO Evolution mobility base and a HOGGI wheelchair

**Notice:**

Only possible with 2 mobility bases of one size.
not possible in combination with aluminium front forks & mobility base Cross.

**Notice:**

Not possible in combination with aluminium front forks & mobility base Cross.

SWINGBO VTi**CLEO**

Sun shade, Changing backpack, Headrest speakers, Cosies



Sun shade, black



Changing backpack, black



Headrest speakers, for headrest pads




„Cosies“ Thermofleece hand warmer for companion

Buggy-Board, Buggypod io




Buggy-Board



temporarily
not available

Buggypod io



Preparation for mounting a “Buggypod”

**Notice:**

Preparation is required to use the Buggypod io.

BASIC home underframe



BASIC home underframe

**with 75 mm
twin castors**

**with 125 mm
twin castors**



Side covers

in standard color white

**matching the upholstery
colour**

in black

Push handle for BASIC
home under frame

COBRA indoor base



COBRA indoor base

**with 75 mm
twin castors**

**with 125 mm
twin castors**

Presetting of the gas spring height adjustment:

up to 20 kg

up to 35 kg

up to 50 kg

Remark / Customer request / Individual special custom-made request