



## DURO / DURO Cross rehab buggy

Size 1 for children at the age of approx. 3 - 8 years

Crash tested according to ISO 7176-19

Order form incl. Exportprices (€) valid until revoked, release date 01/2025

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<b>Invoice address</b>	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
<b>Delivery address</b>	
Company	
Street	City, Zip code



### Basic configuration

#### Stable and compact rehab buggy

- very stable aluminum construction
- Compact folding
- knee angle adjustment from 80°-150° without A-Block
- knee angle adjustment from 80°-120° with A-Block
- angle adjustable footrest 80°-100°
- ergonomic lap belt position
- IUS - Individual Upholstery System

- Seat depth 26 - 40 cm
- Seat width 23 - 34 cm
- Seat tilt 16 or 20°
- Lower leg length 20 - 44 cm
- Back height 45 - 63 cm
- back tilt 90° - 120°
- Maximum load 75 kg

#### Minimum information (Without this information your order can not be handled.)

<b>Frame colour:</b>	<input type="checkbox"/> Black <b>052</b>
<b>Upholstery single-coloured:</b>	<input type="checkbox"/> Black <b>346</b> <input type="checkbox"/> Berry soft <b>347</b> <input type="checkbox"/> Grey soft <b>348</b>
with accessories colour black	<input type="checkbox"/> Turquoise soft <b>349</b>
<b>Seat depth:</b>	<input type="checkbox"/> (26 -40 cm)

#### Presettings (At missing specifications we will adjust to medium sizes.)

**Back height:**  (45 - 63 cm) **Lower leg length:**  (20 - 44 cm) **seat tilt:**  16°,  20° (recommend)

**Basic variations - seat unit**

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.
	<input type="checkbox"/> 3211-0001	052	<b>DURO size 1</b>	<b>1 pc.</b>
with	<input type="checkbox"/> 1281-1007		Plastic front forks	<b>01</b> 1 set incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	<b>02</b> 1 set
or	<input type="checkbox"/> 1281-1008		Aluminum front forks	<b>(upcharge) 03</b> 1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set
with	<input type="checkbox"/> 1219-0004S		Pneumatic tyres with 6-spoke rims	<b>04</b> 1 set incl.
or	<input type="checkbox"/> 1229-0004S		Polyurethan tyres with 6-spoke rims	<b>05</b> 1 set incl.
or	<input type="checkbox"/> 1229-0009		Polyurethan wheels with 7-spoke rims	<b>06</b> 1 set incl.
	<input type="checkbox"/> 3211-7201		Attendant brake, two hand use	<b>07</b> 1 set
or	<input type="checkbox"/> 1229-0045		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	<b>(upcharge) 08</b> 1 set
or	<input type="checkbox"/> 1219-0045		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	<b>(upcharge) 09</b> 1 set
	<input type="checkbox"/> 3214-0001	052	<b>DURO Cross size 1</b>	<b>1 pc.</b>
with	<input type="checkbox"/> 1281-1007		Plastic front forks	<b>01</b> 1 set incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	<b>02</b> 1 set
or	<input type="checkbox"/> 1281-1008		Aluminum front forks	<b>(upcharge) 03</b> 1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set
with	<input type="checkbox"/> 1229-0044		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	<b>08</b> 1 pc. incl.
or	<input type="checkbox"/> 1219-0044		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	<b>09</b> 1 pc. incl.
	<input checked="" type="checkbox"/> 3211-5100		Upholstery, size 1, adaptable (IUS)	<b>10</b> 1 set incl.
	<input type="checkbox"/> 3232-5190		Calf pad support plate and foam insert	<b>11</b> 1 set
with	<input type="checkbox"/> 3231-2900		Thigh supports, size 1, wide- and depth-adjustable <b>* minimum seat width 25 cm</b>	<b>12</b> 1 pair incl.
or	<input type="checkbox"/> 3231-2600		Hip supports, size 1, wide- and depth-adjustable, removable	<b>(upcharge) 13</b> 1 pair
with	<input type="checkbox"/> 3212-1500	022	Supporting brackets in standard color signal white	<b>14</b> 1 pair incl.
or	<input type="checkbox"/> 3212-1500		Supporting brackets in upholstery color	<b>15</b> 1 pair incl.
with	<input type="checkbox"/> 3245-4302		Foot rest bracket angle adjustable, flip up type	1 pair incl.
or	<input type="checkbox"/> 3245-4303		Foot rest bracket angle adjustable, with interlocking	<b>(upcharge)</b> 1 pair
with	<input type="checkbox"/> 3211-249x		Footplate with heel stop, (width: <input type="checkbox"/> 2=30cm, <input type="checkbox"/> 3=34 (cm))	<b>16</b> 1 pair incl.
or	<input type="checkbox"/> 3211-9401		Soft accessory bag, size 1 (max. 3 kg)	<b>17</b> 1 pc.
with	<input type="checkbox"/> 3211-9401_SA-0279		Soft accessory bag, size 1 (max. 3 kg) for use with Outdoor front frame	1 pc.
	<input type="checkbox"/> 3211-7300	052	Tie down kit according to ISO 7176-19	<b>18</b> 1 pair

**Accessories**

	Orderno.	Colour	Description	Qty.
	<input type="checkbox"/> 3231-7600		Headrest pads low version size 1	<b>19</b> 1 set
	<input type="checkbox"/> 3201-7600		Headrest pads size 1	<b>20</b> 1 set
	<input type="checkbox"/> 3201-7601		Headrest pads with spacer fabric for speakers, size 1	<b>21</b> 1 set
	<input type="checkbox"/> 3231-6900		Anatomical neck-headrest with length compens, size 1	<b>22</b> 1 set
	<input type="checkbox"/> 3231-7700		Trunk supports, size 1, stepless height adjustable (8 x 10 x 2 cm)	<b>23</b> 1 pair
	<input type="checkbox"/> 3233-7700		Trunk supports, size 2, stepless height adjustable (10 x 12 x 2 cm)	1 pair
	<input type="checkbox"/> 3201-6100		Lap belt	<b>24</b> 1 pc.
	<input type="checkbox"/> 3231-6701		4-point lap belt, size XS (for seat width 18 - 24)	<b>25</b> 1 pc.
	<input type="checkbox"/> 3231-6702		4-point lap belt, size S (for seat width 22 - 30)	1 pc.
	<input type="checkbox"/> 3201-6200		Five-point harness	<b>26</b> 1 pc.
	<input type="checkbox"/> 3201-6300		Foot straps	<b>27</b> 1 pair
	<input type="checkbox"/> 3201-6400		Ankle huger, size 1	<b>28</b> 1 pair
	<input type="checkbox"/> 3231-5300		Groin strap, size 1	<b>29</b> 1 pc.
	<input type="checkbox"/> 3201-5800		Fixation jacket, size 0	<b>30</b> 1 pc.
	<input type="checkbox"/> 3201-5400		Fixation jacket, size 1	1 pc.
	<input type="checkbox"/> 3201-5500		Fixation jacket, size 2	1 pc.
	<input type="checkbox"/> 3201-5600		Chest-shoulder harness, size 1	<b>31</b> 1 pc.
	<input type="checkbox"/> 3201-5700		Chest-shoulder harness, size 2	1 pc.
	<input type="checkbox"/> 3231-7800		Abduction block, size 1, swing away and removable	<b>32</b> 1 pc.
	<input type="checkbox"/> 3231-7801		Abduction block, size 1, fixed	<b>33</b> 1 pc.
	<input type="checkbox"/> 3231-7900		Abduction block, size 2, swing away and removable	1 pc.
	<input type="checkbox"/> 3231-7901		Abduction block, size 2, fixed	1 pc.

**Accessories**

	<b>Orderno.</b>	<b>Colour</b>	<b>Description</b>	<b>Qty.</b>	
<input type="checkbox"/>	3211-2201	052	Retainer for accessories, grab rail or therapy tray	<b>34</b>	1 pair
<input type="checkbox"/>	3231-8100		Grap rail with upholstery, size 1 *	<b>35</b>	1 pc.
<input type="checkbox"/>	3211-8200		Therapy tray, size 1, transparent <sup>1) 2)</sup>	<b>36</b>	1 pc.
<input type="checkbox"/>	3211-8202		Therapy tray, size 1, transparent, height adjustable <sup>1)</sup>	<b>37</b>	1 pc.
<input type="checkbox"/>	3211-8201A		Therapy tray, size 1, anthracite <sup>1) 2)</sup>		1 pc.
<input type="checkbox"/>	3211-8203A		Therapy tray, size 1, anthracite, height adjustable <sup>1)</sup>		1 pc.
	<sup>1)</sup> only in combination with retainer for accessories				
	<sup>2)</sup> not recommended with hip supports, please use height adjustable trays				
<input type="checkbox"/>	3201-8300	306	Winter warmer, size 1, black	<b>38</b>	1 pc.
<input type="checkbox"/>	3201-8400	306	Winter warmer, size 2, black		1 pc.
<input type="checkbox"/>	7000-0018		Sheep skin insert, size 1	<b>39</b>	1 pc.
<input type="checkbox"/>	7000-0019		Sheep skin insert, size 2		1 pc.
<input type="checkbox"/>	3201-9200	306	Rain cape, size 1, black	<b>40</b>	1 pc.
<input type="checkbox"/>	3211-9300	306	Canopy incl. adapter and rain cover, black, size 1	<b>41</b>	1 pc.































**Custom-made**

	<b>Orderno.</b>	<b>Colour</b>	<b>Description</b>	<b>Qty.</b>	<b>Price</b>
<input type="checkbox"/>	SA-0516-KP		Tray for additional equipment made of alu	<b>42</b>	1 pc.
<input type="checkbox"/>	SA-0279/8.1-KP		Outdoor front end for <b>DURO / DURO CROSS</b> size 1	<b>43</b>	1 pc.

**Extra equipment**

	<b>Orderno.</b>	<b>Colour</b>	<b>Description</b>	<b>Qty.</b>	
<input type="checkbox"/>	7000-0010		Buggy-Board	<b>44</b>	1 pc.
<input type="checkbox"/>	7000-0012	306	Sun shade, black	<b>45</b>	1 pc.
<input type="checkbox"/>	7000-0006		Headrest speakers, for headrest pads 3201-7601	<b>46</b>	1 set
<input type="checkbox"/>	7000-0202	306	Changing backpack, black	<b>47</b>	1 pc.
<input type="checkbox"/>	7000-0199		"Cosies" Thermofleece hand warmer for companion	<b>48</b>	1 pair

# CONFIGURATIONS // ACCESSORIES // DURO

				
01	02	03	04	
Plastic front forks	Swivel lock	Aluminum front forks	Pneumatic front wheel	Pneumatic rear wheel
				
05	05	06	06	07
Polyurethan front wheel	Polyurethan rear wheel	Polyurethan wheels with black rim front wheel	Polyurethan wheels with black rim rear wheel	Attendant brake, two hand use
				
08	08	09	09	10
Polyurethan front wheel for drum brake wheels	Polyurethan drum brake rear wheels, 5-spoke rims	Pneumatic front wheel for drum brake wheels	Polyurethan drum brake rear wheels, 5-spoke rims	IUS - Individual upholstery system
				
11	12	13	14	15
Calf pad support plate	High supports	Hip supports	Supporting brackets in color signal white	Supporting brackets in upholstery color
				
16	17	18	19	20
Footplate with heel stop	Accessory bag	Tie down kit according to ISO 7176-19	Headrest pads, low version	Headrest pads, high version
				
21	22	23	24	25
Headrest pads with spacer fabric for speakers	Anatomical neck-headrest	Trunk supports	Lap belt	4-point lap belt

# CONFIGURATIONS // ACCESSORIES // DURO



Five-point harness



Foot straps



Ankle huger




Groin strap



Fixation jacket




Chest-shoulder harness



Abduction block, swing away and removable



Abduction block, fixed



Retainer for accessories



Grasp rail



Therapy tray



Therapy tray, height adjustable



Winter warmer with sheep skin insert



Sheep skin



Rain cape



Canopy




Buggy-Board



Sun shade



Headrest speakers for headrest pads



Tray for additional equipment



Outdoor front end



Changing backpack



„Cosies“ Thermofleece hand warmer