

HOGGI®

Partnering together to mobilise kids



DURO / DURO Cross rehab buggy

Size 2 for children at the age of approx. 4 - 10 years

Crash tested according to ISO 7176-19

Order form incl. Exportprices (€) valid until revoked, release date 08/2024

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<input type="text"/>	<input type="text"/>
Invoice address	
Company	Customer number
<input type="text"/>	<input type="text"/>
Street	City, Zip code
<input type="text"/>	<input type="text"/>
E-Mail	Telephone
<input type="text"/>	<input type="text"/>
Fax	
<input type="text"/>	
Delivery address	
Company	
<input type="text"/>	
Street	City, Zip code
<input type="text"/>	<input type="text"/>



Basic configuration

Stable and compact rehab buggy

- very stable aluminum construction
- Compact folding
- knee angle adjustment from 80°-150° without A-Block
- knee angle adjustment from 80°-120° with A-Block
- angle adjustable footrest 80°-100°
- ergonomic lap belt position
- IUS - Individual Upholstery System

- Seat depth 31 - 45 cm
- Seat width 28 - 39 cm
- Seat tilt 16 or 20°
- Lower leg length 20 - 44 cm
- Back height 62 - 83 cm
- back tilt 90° - 120°
- Maximum load 75 kg

Minimum information (Without this information your order can not be handled.)

- Frame colour:** Black **052**
-
- Upholstery single-coloured:** Black **346** Berry soft **347** Grey soft **348**
- with accessories colour black Turquoise soft **349**
- Seat depth:** (31 - 45 cm)

Presettings (At missing specifications we will adjust to medium sizes.)

Back height: (62 - 83 cm) **Lower leg length:** (20 - 44 cm) seat tilt: 16°, 20° (recommend)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.
	<input type="checkbox"/> 3212-0001	052	DURO size 2	1 pc.
with	<input type="checkbox"/> 1281-1007		Plastic front forks	01 1 set incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	02 1 set
or	<input type="checkbox"/> 1281-1008		Aluminum front forks	(upcharge) 03 1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set
with	<input type="checkbox"/> 1219-0004S		Pneumatic tyres with 6-spoke rims	04 1 set incl.
or	<input type="checkbox"/> 1229-0004S		Polyurethan tyres with 6-spoke rims	05 1 set incl.
or	<input type="checkbox"/> 1229-0009		Polyurethan wheels with 7-spoke rims	06 1 set incl.
	<input type="checkbox"/> 3211-7201		Attendant brake, two hand use	07 1 set
or	<input type="checkbox"/> 1229-0045		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge) 08 1 set
or	<input type="checkbox"/> 1219-0045		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge) 09 1 set
	<input type="checkbox"/> 3215-0001	052	DURO Cross size 2	1 pc.
with	<input type="checkbox"/> 1281-1007		Plastic front forks	01 1 set incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	02 1 set
or	<input type="checkbox"/> 1281-1008		Aluminum front forks	(upcharge) 03 1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set
with	<input type="checkbox"/> 1229-0044		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	08 1 pc. incl.
or	<input type="checkbox"/> 1219-0044		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	09 1 pc. incl.
	<input checked="" type="checkbox"/> 3212-5100		Upholstery, size 2, adaptable (IUS)	10 1 set incl.
	<input type="checkbox"/> 3234-5190		Calf pad support plate and foam insert	11 1 set
with	<input type="checkbox"/> 3233-2900		Thigh supports, size 2, wide- and depth-adjustable * minimum seat width 25 cm	12 1 pair incl.
or	<input type="checkbox"/> 3233-2600		Hip supports, size 2, wide- and depth-adjustable, removable	(upcharge) 13 1 pair
with	<input type="checkbox"/> 3212-1500	022	Supporting brackets in standard color signal white	14 1 pair incl.
or	<input type="checkbox"/> 3212-1500		Supporting brackets in upholstery color	15 1 pair incl.
with	<input type="checkbox"/> 3245-4302		Foot rest bracket angle adjustable, flip up type	1 pair incl.
with	<input type="checkbox"/> 3245-4303		Foot rest bracket angle adjustable, with interlocking	(upcharge) 1 pair
with	<input type="checkbox"/> 3211-249x		Footplate with heel stop, (width: <input type="checkbox"/> 3=34 cm <input type="checkbox"/> 4=38 cm)	16 1 pair incl.
or	<input type="checkbox"/> 3212-9401		Soft accessory bag, size 2 (max. 3 kg)	17 1 pc.
with	<input type="checkbox"/> 3212-9401_SA-0279		Soft accessory bag, size 2 (max. 3 kg) for use with Outdoor front frame	1 pc.
	<input type="checkbox"/> 3211-7300		Tie down kit according to ISO 7176-19	18 1 pair

Accessories

	Orderno.	Colour	Description	Qty.
	<input type="checkbox"/> 3233-7600		Headrest pads low version size 2	19 1 set
	<input type="checkbox"/> 3203-7600		Headrest pads size 2	20 1 set
	<input type="checkbox"/> 3203-7601		Headrest pads with spacer fabric for speakers, size 2	21 1 set
	<input type="checkbox"/> 3233-6900		Anatomical neck-headrest with length compens, size 2	22 1 set
	<input type="checkbox"/> 3233-7700		Trunk supports, size 2, stepless height adjustable (10 x 12 x 2 cm)	23 1 pair
	<input type="checkbox"/> 3234-7700		Trunk supports, size 3, stepless height adjustable (12 x 13 x 2 cm)	1 pair
	<input type="checkbox"/> 3201-6100		Lap belt	24 1 pc.
	<input type="checkbox"/> 3231-6702		4-point lap belt, size S (for seat width 22 - 30)	25 1 pc.
	<input type="checkbox"/> 3231-6703		4-point lap belt, size M (for seat width 28 - 36)	1 pc.
	<input type="checkbox"/> 3201-6200		Five-point harness	26 1 pc.
	<input type="checkbox"/> 3201-6300		Foot straps	27 1 pair
	<input type="checkbox"/> 3203-6400		Ankle huger, size 2	28 1 pair
	<input type="checkbox"/> 3233-5300		Groin strap, size 2	29 1 pc.
	<input type="checkbox"/> 3201-5500		Fixation jacket, size 2	30 1 pc.
	<input type="checkbox"/> 3203-5400		Fixation jacket, size 3	1 pc.
	<input type="checkbox"/> 3201-5700		Chest-shoulder harness, size 2	31 1 pc.
	<input type="checkbox"/> 3203-5600		Chest-shoulder harness, size 3	1 pc.
	<input type="checkbox"/> 3231-7900		Abduction block, size 2, swing away and removable	32 1 pc.
	<input type="checkbox"/> 3231-7901		Abduction block, size 2, fixed	33 1 pc.
	<input type="checkbox"/> 3233-7800		Abduction block, size 3, swing away and removable	1 pc.
	<input type="checkbox"/> 3233-7801		Abduction block, size 3, fixed	1 pc.

Accessories

	Orderno.	Colour	Description	Qty.	
<input type="checkbox"/>	3211-2201	052	Retainer for accessories, grab rail or therapy tray	34	1 pair
<input type="checkbox"/>	3233-8100		Grap rail with upholstery, size 2 *	35	1 pc.
<input type="checkbox"/>	3212-8200		Therapy tray, size 2, transparent ^{1) 2)}	36	1 pc.
<input type="checkbox"/>	3212-8202		Therapy tray, size 2, transparent, height adjustable ¹⁾	37	1 pc.
<input type="checkbox"/>	3212-8201A		Therapy tray, size 2, anthracite ¹⁾²⁾		1 pc.
<input type="checkbox"/>	3212-8203A		Therapy tray, size 2, anthracite, height adjustable ¹⁾		1 pc.
			¹⁾ only in combination with retainer for accessories ²⁾ not recommended with hip supports, please use height adjustable trays		
<input type="checkbox"/>	3201-8400	306	Winter warmer, size 2, black	38	1 pc.
<input type="checkbox"/>	3203-8300	306	Winter warmer, size 3, black		1 pc.
<input type="checkbox"/>	7000-0019		Sheep skin insert, size 2	39	1 pc.
<input type="checkbox"/>	7000-0020		Sheep skin insert, size 3		1 pc.
<input type="checkbox"/>	3203-9200	306	Rain cape, size 2, black	40	1 pc.
<input type="checkbox"/>	3212-9300	306	Canopy incl. adapter and rain cover, black, size 2	41	1 pc.































Custom-made

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	SA-0516-KP		Tray for additional equipment made of alu	42	1 pc.
<input type="checkbox"/>	SA-0279/8.2-KP		Outdoor front end for DURO / DURO CROSS size 2	43	1 pc.

Extra equipment

	Orderno.	Colour	Description	Qty.	
<input type="checkbox"/>	7000-0010		Buggy-Board	44	1 pc.
<input type="checkbox"/>	7000-0012	306	Sun shade, black	45	1 pc.
<input type="checkbox"/>	7000-0006		Headrest speakers, for headrest pads 3203-7601	46	1 set
<input type="checkbox"/>	7000-0202	306	Changing backpack, black	47	1 pc.
<input type="checkbox"/>	7000-0199		"Cosies" Thermofleece hand warmer for companion	48	1 pair

CONFIGURATIONS // ACCESSORIES // DURO

<p>01</p> 	<p>02</p> 	<p>03</p> 	<p>04</p> 	
<p>Plastic front forks</p>	<p>Swivel lock</p>	<p>Aluminum front forks</p>	<p>Pneumatic front wheel</p>	<p>Pneumatic rear wheel</p>
<p>05</p> 		<p>06</p> 		<p>07</p> 
<p>Polyurethan front wheel</p>	<p>Polyurethan rear wheel</p>	<p>Polyurethan wheels with black rim front wheel</p>	<p>Polyurethan wheels with black rim rear wheel</p>	<p>Attendant brake, two hand use</p>
<p>08</p> 		<p>09</p> 		<p>10</p> 
<p>Polyurethan front wheel for drum brake wheels</p>	<p>Polyurethan drum brake rear wheels, 5-spoke rims</p>	<p>Pneumatic front wheel for drum brake wheels</p>	<p>Polyurethan drum brake rear wheels, 5-spoke rims</p>	<p>IUS - Individual upholstery system</p>
<p>11</p> 	<p>12</p> 	<p>13</p> 	<p>14</p> 	<p>15</p> 
<p>Calf pad support plate</p>	<p>Thigh supports</p>	<p>Hip supports</p>	<p>Supporting brackets in color signal white</p>	<p>Supporting brackets in upholstery color</p>
<p>16</p> 	<p>17</p> 	<p>18</p> 	<p>19</p> 	<p>20</p> 
<p>Footplate with heel stop</p>	<p>Accessory bag</p>	<p>Tie down kit according to ISO 7176-19</p>	<p>Headrest pads, low version</p>	<p>Headrest pads, high version</p>
<p>21</p> 	<p>22</p> 	<p>23</p> 	<p>24</p> 	<p>25</p> 
<p>Headrest pads with spacer fabric for speakers</p>	<p>Anatomical neck-headrest</p>	<p>Trunk supports</p>	<p>Lap belt</p>	<p>4-point lap belt</p>

CONFIGURATIONS // ACCESSORIES // DURO



Five-point harness



Foot straps



Ankle hugger




Groin strap



Fixation jacket




Chest-shoulder harness



Abduction block, swing away and removable



Abduction block, fixed



Retainer for accessories



Grip rail



Therapy tray



Therapy tray, height adjustable



Winter warmer with sheep skin insert




Sheep skin



Rain cape



Canopy



Tray for additional equipment



Outdoor front end



Buggy-Board



Sun shade



Headrest speakers for headrest pads



Changing backpack



„Cosies“ Thermofleece hand warmer