



DURO / DURO Cross rehab buggy

Size 2 for children at the age of approx. 4 - 10 years

Crash tested according to ISO 7176-19

Order form incl. Exportprices (€) valid until revoked, release date 04/2024

☐ Offer	☐ Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



Basic configuration

Stable and compact rehab buggy

- very stable aluminum construction
- Compact folding
- knee angle adjustment from 80°-150° without A-Block
- knee angle adjustment from 80°-120° with A-Block
- angle adjustable footrest 80°-100°
- ergonomic lap belt position
- IUS - Individual Upholstery System

- Seat depth 31 - 45 cm
- Seat width 28 - 39 cm
- Seat tilt 16 or 20°
- Lower leg length 20 - 44 cm
- Back height 62 - 83 cm
- back tilt 90° - 120°
- Maximum load 75 kg

Minimum information (Without this information your order can not be handled.)

Frame colour:	☐ Black 052
Upholstery single-coloured:	☐ Black 346 ☐ Berry soft 347 ☐ Grey soft 348
with accessories colour black	☐ Turquoise soft 349
Seat depth:	(31 - 45 cm)

Presettings (At missing specifications we will adjust to medium sizes.)

Back height: (62 - 83 cm) Lower leg length: (20 - 44 cm) seat tilt: ☐ 16°, ☐ 20° (recommend)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	
	☐ 3212-0001	052	DURO size 2		1 pc.	
with	☐ 1281-1007		Plastic front forks		01 1 set	incl.
	☐ 3202-7100		Swivel lock for Plastic front forks		02 1 set	
or	☐ 1281-1008		Aluminum front forks	(upcharge)	03 1 set	
	☐ 3232-7100		Swivel lock for Aluminum front forks		1 set	
with	☐ 1219-0004S		Pneumatic tyres with 6-spoke rims		04 1 set	incl.
or	☐ 1229-0004S		Polyurethan tyres with 6-spoke rims		05 1 set	incl.
or	☐ 1229-0009		Polyurethan wheels with 7-spoke rims		06 1 set	incl.
	☐ 3211-7201		Attendant brake, two hand use		07 1 set	
or	☐ 1229-0045		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge)	08 1 set	
or	☐ 1219-0045		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge)	09 1 set	
	☐ 3215-0001	052	DURO Cross size 2		1 pc.	
with	☐ 1281-1007		Plastic front forks		01 1 set	incl.
	☐ 3202-7100		Swivel lock for Plastic front forks		02 1 set	
or	☐ 1281-1008		Aluminum front forks	(upcharge)	03 1 set	
	☐ 3232-7100		Swivel lock for Aluminum front forks		1 set	
with	☐ 1229-0044		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function		08 1 pc.	incl.
or	☐ 1219-0044		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function		09 1 pc.	incl.
	☒ 3212-5100		Upholstery, size 2, adaptable (IUS)		10 1 set	incl.
	☐ 3234-5190		Calf pad support plate and foam insert		11 1 set	
with	☐ 3233-2900		Thigh supports, size 2, wide- and depth-adjustable		12 1 pair	incl.
			* minimum seat width 25 cm			
or	☐ 3233-2600		Hip supports, size 2, wide- and depth-adjustable, removable	(upcharge)	13 1 pair	
with	☐ 3212-1500	022	Supporting brackets in standard color signal white		14 1 pair	incl.
or	☐ 3212-1500		Supporting brackets in upholstery color		15 1 pair	incl.
with	☐ 3245-4302		Foot rest bracket angle adjustable, flip up type		1 pair	incl.
with	☐ 3245-4303		Foot rest bracket angle adjustable, with interlocking	(upcharge)	1 pair	
with	☐ 3211-249x		Footplate with heel stop, (width: ☐ 3=34 cm ☐ 4=38 cm)		16 1 pair	incl.
or	☐ 3212-9401		Soft accessory bag, size 2 (max. 3 kg)		17 1 pc.	
with	☐ 3212-9401_SA-0279		Soft accessory bag, size 2 (max. 3 kg) for use with Outdoor front frame		1 pc.	
	☐ 3211-7300		Tie down kit according to ISO 7176-19		18 1 pair	

Accessories

	Orderno.	Colour	Description		Qty.	
	☐ 3233-7600		Headrest pads low version size 2		19 1 set	
	☐ 3203-7600		Headrest pads size 2		20 1 set	
	☐ 3203-7601		Headrest pads with spacer fabric for speakers, size 2		21 1 set	
	☐ 3233-6900		Anatomical neck-headrest with length compens, size 2		22 1 set	
	☐ 3233-7700		Trunk supports, size 2, stepless height adjustable (10 x 12 x 2 cm)		23 1 pair	
	☐ 3234-7700		Trunk supports, size 3, stepless height adjustable (12 x 13 x 2 cm)		1 pair	
	☐ 3201-6100		Lap belt		24 1 pc.	
	☐ 3231-6702		4-point lap belt, size S (for seat width 22 - 30)		25 1 pc.	
	☐ 3231-6703		4-point lap belt, size M (for seat width 28 - 36)		1 pc.	
	☐ 3201-6200		Five-point harness		26 1 pc.	
	☐ 3201-6300		Foot straps		27 1 pair	
	☐ 3203-6400		Ankle huger, size 2		28 1 pair	
	☐ 3233-5300		Groin strap, size 2		29 1 pc.	
	☐ 3201-5500		Fixation jacket, size 2		30 1 pc.	
	☐ 3203-5400		Fixation jacket, size 3		1 pc.	
	☐ 3201-5700		Chest-shoulder harness, size 2		31 1 pc.	
	☐ 3203-5600		Chest-shoulder harness, size 3		1 pc.	
	☐ 3231-7900		Abduction block, size 2, swing away and removable		32 1 pc.	
	☐ 3231-7901		Abduction block, size 2, fixed		33 1 pc.	
	☐ 3233-7800		Abduction block, size 3, swing away and removable		1 pc.	
	☐ 3233-7801		Abduction block, size 3, fixed		1 pc.	

Accessories

	Orderno.	Colour	Description	Qty.	
☐	3211-2201	052	Retainer for accessories, grab rail or therapy tray	34	1 pair
☐	3233-8100		Grap rail with upholstery, size 2 *	35	1 pc.
☐	3212-8200		Therapy tray, size 2, transparent ^{1) 2)}	36	1 pc.
☐	3212-8202		Therapy tray, size 2, transparent, height adjustable ¹⁾	37	1 pc.
☐	3212-8201		Therapy tray, size 2, grey ^{1) 2)}		1 pc.
☐	3212-8203		Therapy tray, size 2, grey, height adjustable ¹⁾		1 pc.
¹⁾ only in combination with retainer for accessories					
⁴⁾ not recommended with hip supports, please use height adjustable trays					
☐	3201-8400	306	Winter warmer, size 2, black	38	1 pc.
☐	3203-8300	306	Winter warmer, size 3, black		1 pc.
☐	7000-0019		Sheep skin insert, size 2	39	1 pc.
☐	7000-0020		Sheep skin insert, size 3		1 pc.
☐	3203-9200	306	Rain cape, size 2, black	40	1 pc.
☐	3212-9300	306	Canopy incl. adapter and rain cover, black, size 2	41	1 pc.


Custom-made

	Orderno.	Colour	Description	Qty.	Price
☐	SA-0516-KP		Tray for additional equipment made of alu	42	1 pc.
☐	SA-0279/8.2-KP		Outdoor front end for DURO / DURO CROSS size 2	43	1 pc.

Extra equipment

	Orderno.	Colour	Description	Qty.	
☐	7000-0010		Buggy-Board	44	1 pc.
☐	7000-0012	306	Sun shade, black	45	1 pc.
☐	7000-0006		Headrest speakers, for headrest pads 3203-7601	46	1 set
☐	7000-0202	306	Changing backpack, black	47	1 pc.
☐	7000-0199		"Cosies" Thermofleece hand warmer for companion	48	1 pair

CONFIGURATIONS // ACCESSORIES // DURO



CONFIGURATIONS // ACCESSORIES // DURO



Five-point harness



Foot straps



Ankle hugger



Groin strap



Fixation jacket




Chest-shoulder harness



Abduction block, swing away and removable



Abduction block, fixed



Retainer for accessories



Grasp rail



Therapy tray



Therapy tray, height adjustable



Winter warmer with sheep skin insert




Sheep skin



Rain cape



Canopy



Tray for additional equipment



Outdoor front end



Buggy-Board



Sun shade



Headrest speakers for headrest pads



Changing backpack



„Cosies“ Thermofleece hand warmer