

www.hoggi.de

Made in Germany
by HOGGI 

DIJRO

Solid rehab buggy for every challenge.



Designed and hand built in Germany.

DURO

Product description.

DURO - Solid rehab buggy for everyday challenges.

HOGGI extends its range of rehab buggies for kids and adolescents with the addition of **DURO** - the new solid rehab buggy meeting the **challenging needs of daily applications**. Furthermore the version **DURO Cross** convinces with its extreme stable back frame and large rear wheels for all outdoor situations.

The drum brake (basic configuration of **DURO Cross**) ensures reliable control even on inclined terrain. Thanks to the **compact construction** of the mobility base, **DURO** provides extreme agility - an advantage in crowded pedestrian zones as well as in narrows of the terrain. The **perfect individual seating position** is achieved by choosing a seat angle between 16° or 20°. The approved accessories of **BINGO Evolution** are compatible with **DURO** and provide **excellent support and positioning possibilities**. Footrest, knee and backrest angle can be adjusted as well.

Enormous growth capacity:

- **Seat depth** growth of up to 10 cm (without any additive elements):
seat depth range size 1: 30 - 40 cm, seat depth range size 2: 35 - 45 cm
- By **moving the backrest**, **shorter seat depths** can be achieved:
seat depth range size 1: 26 - 36 cm, seat depth range size 2: 31 - 41 cm
- **Lower leg length** from 20 - 44 cm, infinitely adjustable
- **Backrest height** from 45 - 63 cm (size 1), respectively 62 - 83 cm (size 2)



DURO.



DURO Cross.


DURO

Specifications.

| Specifications | DURO (Size 1) | DURO (Size 2) | DURO Cross (Size 1) | DURO Cross (Size 2) |
|------------------|---------------|---------------|---------------------|---------------------|
| Seat width | 23 - 34 cm | 28 - 39 cm | 23 - 34 cm | 28 - 39 cm |
| Seat depth | 26 - 40 cm | 31 - 45 cm | 26 - 40 cm | 31 - 45 cm |
| Seat angle | 16° or 20° | 16° or 20° | 16° or 20° | 16° or 20° |
| Backrest height | 45 - 63 cm | 62 - 83 cm | 45 - 63 cm | 62 - 83 cm |
| Backrest angle | 78° - 120° | 78° - 120° | 78° - 120° | 78° - 120° |
| Knee angle | 85° - 160° | 85° - 160° | 85° - 160° | 85° - 160° |
| Lower leg length | 20 - 44 cm | 20 - 44 cm | 20 - 44 cm | 20 - 44 cm |
| Wheel diameter | | | | |
| front | 20 cm | 20 cm | 20 cm | 20 cm |
| rear | 30 cm | 30 cm | 40 cm | 40 cm |
| Weight limit | 75 kg | 75 kg | 75 kg | 75 kg |



Recline mechanism.



Seat angle 16° or 20°.



Tie down kit.



Compact folding size.






Solid rehab buggy for every challenge.

HOGGI

Partnering together to mobilise kids.

HOGGI GmbH

 Eulerstraße 27
D-56235 Ransbach-Baumbach

 Fon: (+49) 2623 / 92 499-0
 Fax: (+49) 2623 / 92 499-99

 info@hoggi.de

www.hoggi.de