



ZIP® rehab buggy

Size 1 for children at the age of approx. 3 - 8 years

Crash tested according to ISO 7176-19 and ANSI/RESNA WC/Vol.1. - Section 19

Order form incl. Exportprices (€) valid until revoked, release date 05/2024

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<input type="text"/>	<input type="text"/>
Invoice address	
Company	Customer number
<input type="text"/>	<input type="text"/>
Street	City, Zip code
<input type="text"/>	<input type="text"/>
E-Mail	Telephone
<input type="text"/>	<input type="text"/>
Fax	
<input type="text"/>	
Delivery address	
Company	
<input type="text"/>	
Street	City, Zip code
<input type="text"/>	<input type="text"/>



Basic configuration

Foldable compact rehab-buggy in a light aluminium construction

- Reinforced seat- and back cover
- Height adjustable footrest, angle adjustable and foldable
- Ball bearing rear and front wheels
- Wheel lock

Measurements:

Seat width	35 cm
Seat depth	25 / 30 / 35 cm
Back height	60 / 66 cm
Back recline	90° - 100°
Lower leg length	20 - 39 cm
Load capacity	50 kg
Weight	12,5 kg

Minimum information (Without this information your order can not be handled.)

Frame colour:	<input type="checkbox"/> Black 052
Nylon cover and Canopy:	<input type="checkbox"/> Grey 304 <input type="checkbox"/> Black 320
Upholstery: with accessory colour black	<input type="checkbox"/> Black 346 <input type="checkbox"/> Berry soft 347 <input type="checkbox"/> Grey soft 348 <input type="checkbox"/> Turquoise soft 349

Presettings (At missing specifications we will adjust to medium sizes.)

Seat depth:	<input type="checkbox"/> 25 cm	<input type="checkbox"/> 30 cm	<input type="checkbox"/> 35 cm
Back height:	<input type="checkbox"/> 60 cm	<input type="checkbox"/> 66 cm	
Lower leg length:	<input type="checkbox"/> (20 - 39 cm)		

Basic variations

with/or: Please choose one version, only!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3210-0001		ZIP® size 1 with back height 60 cm	1 pc.	
and	<input checked="" type="checkbox"/> 3210-2101		Reinforced nylon cover, back height 60 cm	01	1 set incl.
	<input type="checkbox"/> 3210-0002		ZIP® size 1 with back height 66 cm	1 pc.	
and	<input checked="" type="checkbox"/> 3210-2102		Reinforced nylon cover, back height 66 cm		1 set incl.
	<input type="checkbox"/> 3210-8501		Recline angle extension 120°	02	1 pc.
and	<input type="checkbox"/> 3210-2103		Reinforced nylon cover (120°) back height 60 cm		1 set incl.
and	<input type="checkbox"/> 3210-2104		Reinforced nylon cover (120°) back height 66 cm		1 set incl.
with	<input type="checkbox"/> 1281-1007		Plastic front forks	03	1 set incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	04	1 set
or	<input type="checkbox"/> 1281-1008		Aluminum front forks	(upcharge) 05	1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks		1 set
with	<input type="checkbox"/> 1229-0008S		Polyurethan tires with 6-spoke rims	06	1 set incl.
or	<input type="checkbox"/> 1229-0010		Polyurethan wheels with 7-spoke rims	07	1 set incl.
or	<input type="checkbox"/> 1219-0004S		Pneumatic tires with 6-spoke rims	08	1 set incl.
with	<input type="checkbox"/> 3210-2410	015	Foot board assembly height adjustable, angle adjustable, foldable	09	1 set incl.
or	<input type="checkbox"/> 3210-2510	015	Foot board assembly straight angle adjustable, foldable, swing away	(upcharge) 10	1 set
with	<input type="checkbox"/> 3210-7300		Tie down kit ISO 7176-19	11	1 set

Accessories

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3202-7100		Retrofit kit Swivel lock for Plastic front forks	04	1 set
	<input type="checkbox"/> 3232-7100		Retrofit kit Swivel lock for Aluminum front forks		1 set
	<input type="checkbox"/> 3210-7201		Attendant brake, two hand use	12	1 set
	<input type="checkbox"/> 3210-7300E		Tie down kit, retrofit		1 set
	<input type="checkbox"/> 3210-5000		Foldable back rest upholstery, height 60 cm	13	1 pc.
	<input type="checkbox"/> 3210-5100		Foldable back rest upholstery, height 66 cm		1 pc.
	<input type="checkbox"/> 3210-5200		Foldable seat upholstery	14	1 pc.
	<input type="checkbox"/> 3210-7600		Headrest	15	1 pc.
	<input type="checkbox"/> 3210-6900		Anatomical neck headrest, size 1	16	1 set
	<input type="checkbox"/> 3210-7800		Abduction block, size 2, swing away and detachable	17	1 pc.
	<input type="checkbox"/> 3213-7800		Abduction block, size 3, swing away and detachable		1 pc.
	<input type="checkbox"/> 3210-5300		Groin strap, size 1	18	1 pc.
	<input type="checkbox"/> 3210-5400		Groin strap, size 2		1 pc.
	<input type="checkbox"/> 3201-6100		Lap belt	19	1 pc.
	<input type="checkbox"/> 3213-6200		Five-point harness	20	1 pc.
	<input type="checkbox"/> 3201-6400		Ankle hugger, size 1	21	1 pair
	<input type="checkbox"/> 3203-6400		Ankle hugger, size 2		1 pair
	<input type="checkbox"/> 3210-8100		Grap rail with upholstery, size 1	22	1 pc.
	<input type="checkbox"/> 3210-5500	306	Calf strap with upholstery, black	23	1 pc.
	<input type="checkbox"/> 3210-8203		Therapy tray, size 1, grey, height adjustable	24	1 pc.
	<input type="checkbox"/> 3210-8202		Therapy tray, size 1, transparent, height adjustable	25	1 pc.
	<input type="checkbox"/> 3201-8400	306	Winter warmer, size 2, black	26	1 pc.
	<input type="checkbox"/> 3203-8300	306	Winter warmer, size 3, black		1 pc.
	<input type="checkbox"/> 7000-0019		Sheep skin insert, size 2	27	1 pc.
	<input type="checkbox"/> 7000-0020		Sheep skin insert, size 3		1 pc.
	<input type="checkbox"/> 3203-9200	306	Rain cape, size 2, black	28	1 pc.
	<input type="checkbox"/> 3210-9300		Canopy incl. adapter and rain cover *	29	1 pc.
	<input type="checkbox"/> 3210-9302		Canopy incl. adapter and rain cover, for back rest angle 120° *		1 pc.

* color of the canopy matches the color of the nylon cover

Custom-made

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> SA-0260-KP		Tip assist	30	1 pc.

Extra equipment

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 7000-0010		Buggy-Board	31	1 pc.
	<input type="checkbox"/> 7000-0012	306	Sun shade with fasteners, black	32	1 pc.
	<input type="checkbox"/> 3210-9500		Nylon bag for transport	33	1 pc.
	<input type="checkbox"/> 7000-0202	306	Changing backpack, black	34	1 pc.
	<input type="checkbox"/> 7000-0199		"Cosies" Thermofleece hand warmer for companion	35	1 pair

CONFIGURATIONS // ACCESSORIES // ZIP



01 Reinforced nylon cover



02 Recline angle extension 120°



03 Plastic front forks



04 Swivel lock



05 Aluminum front forks



06 Polyurethan tires with 6-spoke rims



07 Polyurethan tires with 6-spoke rims



08 Polyurethan wheels with 7-spoke rims



09 Polyurethan wheels with 7-spoke rims



10 Pneumatic tires with 6-spoke rims



11 Pneumatic tires with 6-spoke rims



12 Foot board assembly



13 Foot board assembly straight



14 Tie down kit



15 Attendant brake, two hand use



16 Back rest upholstery



17 Seat upholstery



18 Headrest



19 Anatomical neck headrest



20 Abduction block



21 Groin strap



22 Lap belt




23 Five-point harness




24 Ankle huger



25 Crap rail with upholstery



26 Calf strap with upholstery



27 Therapy tray, grey



28 Therapy tray, transparent



29 Winter warmer



30 Sheep skin

CONFIGURATIONS // ACCESSORIES // ZIP



Rain cape



Canopy incl. transparent rain cover



Tip assist



Buggy-Board



Sun shade



Nylon bag for transport



Changing backpack



„Cosies“ Thermofleece hand warmer