



CESA Active-wheelchair for kids & adolescents

Basic configuration inclusive:

- Footrest
- Adjustable seat depth
- Seat- and back panel
- Wheelbase adjustment
- HOGGI light wheels + brake

Order form incl. Exportprices (€) valid until revoked, release date 04/2024

Price:

☐ Offer	☐ Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



Basic configuration

- Boxframe-Active wheelchair. Aluminum construction powder coated with growth seat depth
- **3 frame sizes:**
 - Size 1: Seat width (SW) 24 - 30 cm in 2 cm steps
 - Size 2: Seat width (SW) 26 - 36 cm in 2 cm steps
 - Size 3: Seat width (SW) 30 - 44 cm in 2 cm steps
 - Convertible to other seat width by part replacement
 - Seat height 36 – 49 cm
 - Seat angle 0°-12°
- seat frame angle 100°
- Active degree change 1 – 8 cm
- back angle 80 - 105°
- Clothes cover
- **HOGGI** light rear wheels including push rims
- **HOGGI** light wheel locks
- Camber 3°, 6°, 9° or 12°
- Load capacity: 60 kg (size1), 80 kg (size 2), 100 kg (size 3)

Minimum information (Without this information your order can not be handled.)

Frame colour		(Please see table „frame colours“)
Seat width	☐ 24 ☐ 26 ☐ 28 ☐ 30 ☐ 32 ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44	
Camber	☐ 3° (3274-1291) ☐ 6° (3274-1292) ☐ 9° (3274-1293) ☐ 12° (3274-1294)	
Seat height	(36-46)	Seat angle* (0° - 12°)

*Please see attachment of “seat measurements”

Presettings

Seat depth (SD) ¹	(24 – 48 cm)	Lower leg length (LLL) ⁴	(15 - 46 cm)
Backrest angle ²	(80 - 105°)		
Active degree ³	(1 - 8 cm)		

¹ Please see table seat depth

² Angle compared to the vertical line

At missing specification we will adjust to middle sized measurements

³ Please see drawing "active degree"

⁴ Please see table Lower leg length

Basic variation

	Order no.	Colour	Description	Qty.
	☐ 3220-00SW		CESA size 1 (20")	1 pc.
	☐ 3221-00SW		CESA size 2 (22")	1 pc.
	☐ 3222-00SW		CESA size 3 (24")	1 pc.
with	☐ 3220-21SW	052	Seat panel	01 1 set incl.
or	☐ 3220-2199		Without Seat panel	1 set
	☐ 3220-2110		Mounting set for custom made seat panel	1 set
with	☐ 3220-31SW		Back base "rigid" (back height 20)	02 1 set incl.
or	☐ 3220-32SW		Back base "rigid" (back height 25)	1 set incl.
or	☐ 3220-33SW		Back base "rigid" (back height 30)	1 set incl.
or	☐ 3220-34SW		Back base "rigid" (back height 35)	1 set incl.
or	☐ 3220-35SW		Back base "foldable & angle adjustable" (back height 20) (Upcharge)	03 1 set
or	☐ 3220-35SW		Back base "foldable & angle adjustable" (back height 25) (Upcharge)	1 set
or	☐ 3220-37SW		Back base "foldable & angle adjustable" (back height 30) (Upcharge)	1 set
or	☐ 3220-38SW		Back base "foldable & angle adjustable" (back height 35) (Upcharge)	1 set
(back angle 80-105 ° in 5 cm steps)				
with	☐ 3220-51SW		Back panel in frame color (back height 20) SW24-36 only	04 1 set incl.
or	☐ 3220-52SW		Back panel in frame color (back height 25)	1 set incl.
or	☐ 3220-53SW		Back panel in frame color (back height 30)	1 set incl.
or	☐ 3220-54SW		Back panel in frame color (back height 35)	1 set incl.
or	☐ 3220-56SW		Nylon back cover (back height 25-35*)	05 1 set incl.
or	☐ 3220-57SW		Nylon back cover (back height 30-40*)	1 set incl.
or	☐ 3220-58SW		Nylon back cover (back height 35-40*)	1 set incl.
* height adjustment only with 3251-2175 or 3261-5102				
or	☐ 3220-5199		Without back panel or back cover	1 set
	☐ 3220-5110		Mounting set for custom made back panel	1 set
with	☐ 3221-48SW		Pipe bracket footrest standard *	06 1 set incl.
or	☐ 3221-41SW		Footrest system standard, one piece, angle adjustable, fixed	07 1 set incl.
or	☐ 3221-43SW		Footrest system, one piece, angle adjustable, foldable** (Upcharge)	08 1 set
or	☐ 3221-45SW		Single foot plate footrest system, angle adjustable, foldable** (Upcharge)	09 1 pair
or	☐ 3221-46SW		Footrest system, one piece, angle adjustable, for short LLL*** (Upcharge)	10 1 set
or	☐ 3221-47SW		Footrest system, one piece, angle adjustable, for short LLL*** (Upcharge)	1 set
* from SW 30, ** from SW 28, *** only up to SW 36				
or	☐ 3221-40SW		Footrest system, one piece foot plate, fold backwards	11 1 set
	☐ 3220-880x	346	Calf strap with upholstery, length according to seat width	12 1 set
	☐ 3220-7300		Tie down kit (minimum back rest height 30) *	13 1 set
	☐ 3274-7301		Back rest extension (with back panel only)	1 set
	☐ 3274-7370	346	Cushion for back rest extension (Upcharge)	1 pc
* Use with back rest extension if there is no support for shoulder and head mounted in the transport vehicle.				
with	☐ 3275-52SW		Push-bar	14 1 pc
or	☐ 3275-54SW		Push bar "wide"	15 1 pc
or	☐ 3275-53SW		Push-handles	16 1 pc
or	☐ 3281-5200		Push handles displaced to the back, height adjustable	17 1 pc
or	☐ 3281-5203		Single push handles, foldable	18 1 pc
	☐ 3251-2175		Teleskopik tubes, plugged into in back tubes for back height adjustment (from BH 25) *	19 1 set
* only in connection with push-bar, push-handles or one hand push-handle				
or	☐ 3261-5102		Single push handle, plugged into in back tubes (from BH 25)	20 1 set
or	☐ 3281-5201		Standard single push handles 90°, size 1	21 1 set
or	☐ 3281-5202		Standard single push handles 90°, size 2	1 set
or	☐ 3275-5502		One hand push-handle delta form (max. SW 36)	22 1 pc

Brakes, fender

	Order no.	Colour Description	Qty.	
with ☐	3275-1750	HOGGI light sport brake	23	1 set incl.
or ☐	3221-1700	Standard brakes CESA		1 set incl.
with ☐		Clothes cover		1 set incl.
or ☐	3275-160x	Fender (Upcharge)	24	1 set
or ☐	3275-169x	Fender with integrated wheel lock levers *	(Upcharge) 25	1 set
with ☐	3276-160x	Fender with height adjustable and detachable arm rests (Upcharge)	26	1 set
or ☐	3276-169x	Fender with integrated brake system * and with height adjustable and detachable arm rests (Upcharge)	27	1 set
with ☐	7000-0051	PU-Armrest pads, 235 mm long		1 pair incl.
or ☐	7000-0052	PU-Armrest pads, 270 mm long		1 pair incl.
☐	3240-8202	Therapy tray size 2, SW 24 and 26 cm **	28	1 pc
☐	3240-8203	Therapy tray size 3, SW 28, 30 and 32 cm **		1 pc
☐	3240-8204	Therapy tray size 4, SW 34, 36 and 38 cm **		1 pc
☐	3240-8205	Therapy tray size 5, SW 40, 42 and 44 cm **		1 pc
* for 22 ", SA-0512 must be selected for SH 43 + 45				
* for 24", SA-0512 must be selected for SH 47 + 49				
** only in conjunction with arm rests				

Aluminium fender

	Order no.	Colour Description	Qty.	
with ☐	3277-1604	Wheel guards aluminum in frame color, 20-24"	29	1 set
or ☐	3278-1604	Wheel guards aluminum in frame color, 20-24" with height adj. and detachable arm rests	30	1 set
or ☐	3277-1605	Wheel guards aluminum in frame color, 20-24", wide for 3° camber		1 set
or ☐	3278-1605	Wheel guards aluminum in frame color 20-24" with height adj. a. detach. arm rests, wide for 3° camber		1 set
or ☐	3277-1694	Wheel guards aluminum in frame color 20-24" with integrated wheel lock levers *	31	1 set
or ☐	3278-1694	Wheel guards alu. frame color 20-24"with integrated brake system, height adjust. detachable Armrests *	32	1 set
or ☐	3277-1695	Wheel guards aluminum in frame color 20-24" with integrated wheel lock levers, wide for 3° camber *		1 set
or ☐	3278-1695	Wheel guards ALU fr-color 20-24"with integ. brake sys., height adjust. detach. Armrests, wide for 3° camber *		1 set
* for 22 " SA-0512 must be selected for SH 41, 43 + 45				
* for 24" SA-0512 must be selected for SH 47 + 49				

Rear wheels

	Order no.	Colour Description	Qty.	
with ☐	1296-0xxx	Rear wheel HOGGI light	33	1 set incl.
or ☐	1297-1xxx	Rear wheel HOGGI light with drum brakes* (Upcharge)	34	1 set
*in combination with push-bar or push-handles				
with ☐	129x-x2xx	20" wheel diameter		1 set incl.
or ☐	129x-x3xx	22" wheel diameter		1 set incl.
or ☐	129x-x4xx	24" wheel diameter		1 set incl.
with ☐	129x-xx1x	Pneumatic tires 1"		1 set incl.
or ☐	129x-xx2x	Flat free tires (20x1 3/8", 22x1" and 24x1")		1 set incl.
with ☐	129x-xxx0	Without push-rims		
or ☐	129x-xxx1	Push rims aluminium anodised, "standard" mounting	35	1 set incl.
or ☐	129x-xxx2	Push rims aluminium anodised, "high" mounting	36	1 set incl.
or ☐	129x-xxx3	Push rims stainless steel, "standard" mounting (Upcharge)	37	1 set
or ☐	129x-xxx4	Push rims stainless steel, "high" mounting (Upcharge)	38	1 set
with ☐		Push-rim mounting "narrow"		
or ☐		Push-rim mounting "wide"		

Rear wheels

	Order no.	Colour	Description		Qty.
or	☐	1204-1155	Closed ergonomic hand rims 22 ", black anodized, for Hoggi-light wheels	(Upcharge)	67 1 Set
or	☐	1204-1156	Closed ergonomic hand rims 24 ", black anodized, for Hoggi-light wheels	(Upcharge)	1 Set
with	☐	1204-106X	"HOGGI Grip Plus" hand rim inserts black, for closed ergonomic hand rims	(Upcharge)	68 1 Set
	☐	1204-114X	Aluminum hand rim with plastic anti-slip coating, black, standard mounting *	(Upcharge)	65 1 pair
	☐	1204-114X	Aluminum hand rim with plastic anti-slip coating, black, high mounting *	(Upcharge)	1 pair
* only in combination with HOGGI light rear wheels					
with	☐	1296-4211	Pneumatic wheel 20x1" HOGGI Extra light *	(Upcharge)	39 1 set
with	☐	1296-4311	Pneumatic wheel 22x1" HOGGI Extra light *	(Upcharge)	1 set
or	☐	1296-4411	Pneumatic wheel 24x1" HOGGI Extra light *	(Upcharge)	1 set
* incl. black alu push-rim					

Spoke guards // Old motifs // Sale

	Order no.	Colour	Description		Qty.
	☐	1204-101x	Spoke guard, transparent *		40 1 set
	☐	1204-101x	Spoke guard, transparent, for drum brakes wheels *		1 set
	☐	1204-101x	Spoke guard, printed **	(Sale)	41 1 set
	☐	1204-101x	Spoke guard, printed, for drum brakes **	(Sale)	1 set
* Motifs 2021 only possible in combination with 1204-2001-XXX					
** please enter code of your print design see attachment or brochure "spoke guard" // Only while stocks last!					
*** printed spoke guards are not available for 20" wheels					

Spoke guards // New motifs 2021 (with 10% dealer discount)

	Order no.	Colour	Description	Qty.	Price €
	☐	1204-2001-000	Spoke guard Special interior foiling, single-coloured according to customer requirements ^{1 2}	1 set	
	☐	1204-2001-2XX	Spoke guard Special interior foiling, printed ^{1 3}	1 set	
	☐	1204-2001-SA	Spoke guard Special interior foiling, printed according to customer requirements ¹		Price on request
¹ only in combination with transparent spoke guard					
² please specify RAL tone					
³ Spoke guard motif selectable in "Order form spoke-guards" (available on the website)					

Front castors

with/or: Please select only one version!

	Order no.	Colour	Description		Qty.
with	☐	1250-1002S	Wheel with rubber tire, black 100 x 32 mm *	42	1 set incl.
or	☐	1250-1003S	Wheel with rubber tire, black 125 x 32 mm		1 set incl.
or	☐	1250-1005S	Wheel with rubber tire, black 140 x 40 mm		1 set incl.
or	☐	1249-1005	Skater wheel with LED 100 x 24 mm *	(Upcharge) 43	1 set
or	☐	1249-1006	Skater wheel with LED 125 x 24 mm	(Upcharge)	1 set
or	☐	1221-1003	PU-wheel with plastic design rim 100 x 34 mm*	(Upcharge) 44	1 set
or	☐	1221-1004	PU-wheel with plastic design rim 125 x 34 mm	(Upcharge)	1 set
* 100 mm wheels not with camber 3°					

Accessories

	Order no.	Colour	Description	Qty.
☐	3220-1801		Anti tip, swingaway, assembly left	45 1 pc.
☐	3220-1802		Anti tip, swingaway, assembly right	1 pc.
☐	3275-1803		Tip assist left	46 1 pc.
☐	3275-1804		Tip assist right	1 pc.
☐	3275-1910		Bumper wheels *	47 1 pc.
* not in combination with outdoor front end				
☐	1204-101x		Push-rim cover, silicone black fluted	48 1 set
☐	3222-1xSW		Seat width extension kit of 2 cm, "rigid" back base, back height	1 set
			☐ 6 = BH 20, ☐ 7 = BH 25, ☐ 8 = BH 30, ☐ 9 = BH 35	
☐	3229-1xSW		Seat width extension kit of 2 cm, "foldable" back base, back height	1 set
			☐ 6 = BH 20, ☐ 7 = BH 25, ☐ 8 = BH 30, ☐ 9 = BH 35	
☐	3201-6400		Ankle hugger, size 1	49 1 pair
☐	3203-6400		Ankle hugger, size 2	1 pair
☐	3220-6100		Lap belt, black for rigid back base only	50 1 pc.
☐	3220-6101		Lap belt, black for angle adj. back base only	1 pc.
☐	3220-6702		4-Point lap belt, size S (SW 24 - 30)	51 1 pc.
☐	3220-6703		4-Point lap belt, size M (SW 28 - 36)	1 pc.
☐	3220-6704		4-Point lap belt, size L (SW 34 - 42)	1 pc.
☐	3251-7600		Edge protector for back rest tubes with nylon back cover	52 1 Set

Accessoires for Seat- and Back panel

	Order no.	Colour	Description	Qty.
☐	3231-7800	346	Abduction block size 1, swing away, detachable, (9 x 6 x 9 cm)*	53 1 set
☐	3231-7900	346	Abduction block size 2 swing away, detachable, (11 x 8 x 10 cm)*	1 set
☐	3233-7800	346	Abduction block size 3 swing away, detachable, (12 x 9 x 10 cm)*	1 set
* only with seat panel 3220-21SW				
☐	3220-71SW	346	Seat cushion; 3 cm thick	54 1 pc.
☐	3276-7100	346	Back cushion for back panel, 2,5 cm thick	55 1 pc.
with ☐	2240-1505		MDP pad holder swing away and detachable for attachment to round tube D25 mm, left side	56 1 pc.
with ☐	2240-1506		MDP pad holder swing away and detachable for attachment to round tube D25 mm, right side	1 pc.
or ☐	2240-1405		MDP pad holder fixed for attachment to round tube D25 mm, left/right side	57 1 pc.
with ☐	2240-1461	346	MDP trunk support size 1 (8 x 10 x 2 cm), cover black	58 1 pc.
or ☐	2240-1462	346	MDP trunk support size 2 (10 x 12 x 2 cm), cover black	1 pc.
or ☐	2240-1463	346	MDP trunk support size 3 (12 x 13 x 2 cm), cover black	1 pc.
or ☐	2240-1471		MDP trunk support curved size 1 (6 x 11,5 x 2,5 cm), PU w/o cover	59 1 pc.
or ☐	2240-1472		MDP trunk support curved size 2 (8 x 14 x 2,5 cm), PU w/o cover	1 pc.
☐	3201-3830		Universal headrest bracket (3245-7800 or third-party product) *	60 1 pc.
☐	3245-7800	306	Headrest with occiput upholstery	61 1 pc.
* starting from backrest height 30 cm				
☐	3200-0401		Belt routing as a T-strap above the shoulders	62 1 pc.

Custom-made

	Order no.	Colour	Description	Qty.
☐	3220-80SW	346	Contoured seat cushion with cover, black, SW 24 - 32	63 1 pc.
☐	3220-80SW	346	Contoured seat cushion with cover, black, SW 34 - 44	1 pc.
☐	3275-7200	346	"HOGGI align" Contoured back cushion, black	64 1 pc.
☐	SA-0286/2-KP		Abduction block 30 mm height extension kit *	1 set
* only in combination with abduction block and contoured seat cushion				
☐	SA-0210-MP		Special Frame colour, RAL code only (Upcharge)	1 pc.
☐	6999-0014	060	Special Frame colour, Copper Gold (Upcharge)	1 pc.
☐	6999-0017		Castor forks and camber adapter powder coated according to customer requirements * (no glitter) (Upcharge)	1 set
☐	SA-0279/5-KP		CESA outdoor front end *	1 pc.
* additional costs for retrofitting				
☐	SA-0390		Upholstery in special colour *	
* please enter colour code or upholstery code respectively				
☐	SA-0572-F-MP		Plastic anti-slip coating in custom colour: Pure white (RAL 9010), Silver grey (RAL 7001), Sky blue (RAL 5015), Traffic red (RAL 3020), Pure orange (RAL 2004)* (Upcharge)	66 1 pair
* only in combination with HOGGI light rear wheels				

Custom-made


	Order no.	Colour Description		Qty.
☐	SA-0093-MP	Maxgrepp Original cell rubber coated push rims	(Upcharge) 69	1 set
☐	SA-0094-MP	Maxgrepp Ergo Para (closed push rims)	(Upcharge) 70	1 set
☐	SA-0095-MP	Maxgrepp Ergo Para CP (advanced closed push rims)	(Upcharge) 71	1 set
☐	SA-0528-MP	Preparation for electric drive Alber "e-fix"	(Upcharge)	1 set
☐	SA-0528/1-MP	Tie down kit in combination with auxiliary drives	(Upcharge)	1 set
☐	SA-0528/2-MP	Preparation for electric drive Alber "Twion, e-motion"	(Upcharge) 72	1 set
☐	SA-0315/7-MP	Reconstruction CESA for ProActiv One hand drive	(Upcharge) 73	1 pc.
☐	SA-0505-MP	Locking of the integrated brake	74	1 pc.
☐	SA-0594	Carbon fender	(Upcharge) 75	Price on request

Frame colours

001 Black (RAL 9005)	020 Yellow green (RAL 6018)	028 Pearl violet (RAL 4011)	056 Turquoise blue (RAL 5018)
002 Silver	021 Grey aluminium (RAL 9007)	029 Pearl beige (RAL 1035)	057 Iron sparkling dark light
003 Pearl ruby red (RAL 3032)	022 White (RAL 9003)	032 Black dull matt (RAL 9005M)	058 Salmon orange (RAL 2012)
004 Pearl night blue (RAL 5026)	023 Pure orange (RAL 2004)	039 Signal violet (RAL 4008)	059 Pastel blue (RAL 5024)
005 Pearl gentiant blue (RAL 5025)	024 Sparkling iron effect	040 Night blue (RAL 5022)	061 Light pink (RAL 3015)
017 Traffic yellow (RAL 1023)	025 Signal blue (RAL 5005)	050 White matt (RAL 9003)	064 Olive green (RAL 6003)
018 Flame red (RAL 3000)	026 Vermillion (RAL 2002)	053 Light green (RAL 6027)	065 Raspberry red (RAL 3027)
019 Sky blue (RAL 5015)	027 Heather violet (RAL 4003)	054 Ocean blue (RAL 5020)	066 Pastel green (RAL 6019)

SA-0210-MP Special colours on request (Enter RAL code Upcharge

How to measure, active degree and back angle



SH = Seat height **LLL** = Lower leg length **SD** = Seat depth **BH** = Backrest height **SA** = Seat angle

(All measurements are without seat or back rest upholsteries)

Seat height and lower leg length

CESA Seat height and Lower leg length depending on frame size							
Frame size	SH	Foot rest versions					
		3221-41SB	3221-43SB	3221-45SB	3221-46SB	3221-47SB	3221-40SB
		Footrest system one piece angle adjustable fixed	Footrest system one piece angle adjustable foldable	Footrest system two piece angle adjustable foldable	Footrest system one piece angle adjustable for short LLL	Footrest system one piece angle adjustable for short LLL	Footrest system one piece fold backwards
CESA size 1 (SA 5°)	36	27 - 33	27 - 33	27 - 33	15 - 22	19 - 25	27 - 34
	38	27 - 35	27 - 35	27 - 35	15 - 24	19 - 27	27 - 36
	40	27 - 37	27 - 37	27 - 37	15 - 26	19 - 29	27 - 38
	42	27 - 39	27 - 39	27 - 39	15 - 28	19 - 31	27 - 40
CESA size 2 (SA 5°)	37	30 - 34	30 - 34	30 - 34	19 - 24	23 - 27	30 - 35
	39	30 - 36	30 - 36	30 - 36	19 - 26	23 - 29	30 - 37
	41	30 - 38	30 - 38	30 - 38	19 - 28	23 - 31	30 - 39
	43	30 - 40	30 - 40	30 - 40	19 - 30	23 - 33	30 - 41
	45	30 - 42	30 - 42	30 - 42	19 - 32	23 - 35	30 - 43
CESA size 3 (SA 5°)	41	34 - 40	34 - 40	34 - 40	23 - 28	27 - 32	34 - 41
	43	34 - 42	34 - 42	34 - 42	23 - 30	27 - 34	34 - 43
	45	34 - 44	34 - 44	34 - 44	23 - 32	27 - 36	34 - 45
	47	34 - 44	34 - 44	34 - 44	23 - 34	27 - 38	34 - 45
	49	34 - 46	34 - 46	34 - 46	23 - 36	27 - 40	34 - 47

Seat width, seat depth and active degree

CESA Seat depth range depending on seat width and active degree (AD)																							
	AD	Seat width (cm)																					
		24		26		28		30		32		34		36		38		40		42		44	
CESA size 1 (SH 38, SA 5°)	1	24-30		24-32		26-34		28-36															
	2	25-30		25-32		26-34		28-36															
	3	26-30		26-32		26-34		28-36															
	4	27-30		27-32		27-34		28-36															
	5	28-30		28-32		28-34		28-36															
	6	29-30		29-32		29-34		29-36															
	7	30-30		30-32		30-34		30-36															
	8			31-32		31-34		31-36															
CESA size 2 (SH 41, SA 5°)	1			26-32		26-34		28-36		30-38		32-40		34-42									
	2			26-32		26-34		28-36		30-38		32-40		34-42									
	3			27-32		27-34		28-36		30-38		32-40		34-42									
	4			28-32		28-34		28-36		30-38		32-40		34-42									
	5			29-32		29-34		29-36		30-38		32-40		34-42									
	6			30-32		30-34		30-36		30-38		32-40		34-42									
	7			31-32		31-34		31-36		31-38		32-40		34-42									
	8			32		32-34		32-36		32-38		32-40		34-42									
CESA size 3 (SH 45, SA 5°)	1							32-36		32-38		32-40		34-42		36-44		38-46		40-48		42-48	
	2							32-36		32-38		32-40		34-42		36-44		38-46		40-48		42-48	
	3							32-36		32-38		32-40		34-42		36-44		38-46		40-48		42-48	
	4							33-36		33-38		33-40		34-42		36-44		38-46		40-48		42-48	
	5							34-36		34-38		34-40		34-42		36-44		38-46		40-48		42-48	
	6							35-36		35-38		35-40		35-42		36-44		38-46		40-48		42-48	
	7							36		36-38		36-40		36-42		36-44		38-46		40-48		42-48	
	8									38		37-40		37-42		37-44		38-46		40-48		42-48	

Upholstery colours

346 Black	321 Spacer fabric, black	(special colour)	348 Grey Soft	(special colour)
	347 Berry Soft	(special colour)	349 Turquoise Soft	(special colour)

All HOGGI upholstery are made from flame retardant Trevira CS. Exception spacer fabric (**321**).

Spoke guard print design (for detailed images see attachment or brochure "spoke guard")

Colour code

201	HOGGI blue	205	Train
202	Wonderland	206	Dragon
203	Mouse	207	Circus
204	Water world		



CONFIGURATIONS // ACCESSORIES // CESA



Seat panel



Back base, rigid



Back base, foldable and angle adjustable



Back panel in frame color



Nylon back cover



Pipe bracket footrest standard



Footrest system, one piece



Footrest system, foldable



Footrest system, two piece



Footrest system, one piece, for short LLL



One piece foot plate, fold backwards



Calf strap with upholstery



Tie down kit



Push-bar



Push bar „wide“



Push-handles



Push handles displaced to the back



Single push handles, foldable



Teleskopik tubes



Single push handles



Single push handles 90°



One hand push-handle delta form



HOGGI light sport brake



Fender




Fender with integrated brake system




Fender with arm rests



Fender with arm rests & integrated brake system



Therapy tray



Aluminium fender



Aluminium fender with arm rests

CONFIGURATIONS // ACCESSORIES // CESA



Aluminium fender with integr. brake system



Aluminium fender with brake and arm rests



Rear wheel HOGGI light



Rear wheel HOGGI light with drum brakes



Push rims aluminium anodised, standard



Push rims aluminium anodised, high



Push rims stainless steel



Push rims stainless, high



HOGGI Extra light



Spoke guard, transparent



Spoke guard, printed



Wheel with rubber tire



Skater wheel with LED




PU-wheel with plastic design rim



Anti tip



Tip assist



Bumper wheels



Push-rim cover



Ankle hugger



Lap belt




4-Point lap belt



Edge protector



Abduction block



Seat cushion




Back cushion



MDP pad holder swing away and detachable



MDP pad holder fixed



MDP trunk support



MDP trunk support curved




Universal headrest bracket

CONFIGURATIONS // ACCESSORIES // CESA



Headrest with occiput upholstery



Belt routing as a T-strap



Contoured seat cushion with cover



„HOGGI align“ Contoured back cushion



Aluminum hand rim with plastic anti-slip coating



Aluminum hand rim with plastic anti-slip coating



Closed ergonomic hand rims



HOGGI Grip Plus



Maxgrepp Original push rims



Maxgrepp Ergo Para push rims




Maxgrepp Ergo Para CP push rims



Preparation for electric drive „TWION“



Preparation for one-hand drive „ProActive“



Locking of the integrated brake



Carbon fender