



CLEO Active-wheelchair for children

Crash tested according to ISO 7176-19 and RESNA WC-4:2012: Section 19

Order form incl. Exportprices (€) valid until revoked, release date 04/2024

Basic configuration inclusive:

- Seat- and back panel
- Clothes cover
- Footrest hanger 90°
- **HOGGI** light rear wheels
- **HOGGI** light brake **CLEO**

Price:

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<input type="text"/>	<input type="text"/>
Invoice address	
Company	Customer number
<input type="text"/>	<input type="text"/>
Street	City, Zip code
<input type="text"/>	<input type="text"/>
E-Mail	Telephone
<input type="text"/>	<input type="text"/>
Fax	
<input type="text"/>	
Delivery address	
Company	
<input type="text"/>	
Street	City, Zip code
<input type="text"/>	<input type="text"/>



Basic configuration

- Powder coated aluminium frame with growth capacity
- **3 frame sizes:**
- Size 1: Seat width (SW) 24 - 30 cm in 2 cm steps
- Size 2: Seat width (SW) 24 - 34 cm in 2 cm steps
- Size 3: Seat width (SW) 26 - 34 cm in 2 cm steps
- Seat heights: 29 cm - 46 cm
- Seat angle: 0,5° - 11,5°
- Footrest hanger 90° inclusive ABS-footrest
- Backrest angle: -10° to +10°
- Seat- and back panel
- Clothes cover
- **HOGGI** light rear wheels inclusive push rims
- **HOGGI** light wheel locks
- Camber 6°, 9° or 12°
- Load capacity: 60 kg

Minimum information (Without this information your order can not be handled.)

Frame colour (Please see table „frame colours“)

Seat width in cm 24 26 28 30 32 34

Camber in degree 6° (3274-1292) 9° (3274-1293) 12° (3274-1294)

Seat height * (29 - 46 cm)

Seat angle * (0,5° - 11,5°)

* Please see attachment of "seat measurements"

Presettings

(Without this information your order can not be handled.)

Seat depth (SD) ¹ (22 - 46 cm)

Lower leg length (LLL) (13 - 36 cm)

Backrest angle ² (-10 to +10)

Active degree ³ (0 - 8 cm) Comfort degree (-5 to 3 cm)

¹ Please see table seat measurements

² Angle compared to the vertical line

³ Please see drawing "active degree"

At missing specification we will adjust to middle sized measurements

Basic variation

with/or: Please select only one version! SW=Seat width

	Order no.	Colour	Description	Qty.	Price €
	<input type="checkbox"/> 3274-00SW		CLEO size 1 (20")	1 pc.	
with	<input type="checkbox"/> 3274-1420		Base plates, aluminum size 1, black coated	01 1 set	incl.
or	<input type="checkbox"/> 3274-1410		Base plates, carbon * (Upcharge)	02 1 set	
or	<input type="checkbox"/> 3274-1430		Base plates alu "Comfort" long wheelbase **	03 1 set	incl.
	<input type="checkbox"/> 3275-00SW		CLEO size 2 (22")	1 pc.	
with	<input type="checkbox"/> 3275-1420		Base plates, aluminum size 2, black coated	01 1 set	incl.
or	<input type="checkbox"/> 3275-1410		Base plates, carbon * (Upcharge)	02 1 set	
or	<input type="checkbox"/> 3275-1430		Base plates alu "Comfort" long wheelbase **	03 1 set	incl.
	<input type="checkbox"/> 3276-00SW		CLEO size 3 (24")	1 pc.	
with	<input type="checkbox"/> 3276-1420		Base plates, aluminum size 3, black coated	01 1 set	incl.
or	<input type="checkbox"/> 3276-1410		Base plates, carbon * (Upcharge)	02 1 set	
or	<input type="checkbox"/> 3276-1430		Base plates alu "Comfort" long wheelbase **	03 1 set	incl.
* can not be combined with tie down kit					
** reduction of the active degree by 5 cm					
with	<input type="checkbox"/> 3274-25SW		Back base "rigid" with back panel, back height 20	04 1 set	incl.
or	<input type="checkbox"/> 3274-26SW		Back base "rigid" with back panel, back height 25	1 set	incl.
or	<input type="checkbox"/> 3274-27SW		Back base "rigid" with back panel, back height 30	1 set	incl.
or	<input type="checkbox"/> 3274-28SW		Back base "rigid" with back panel, back height 35	1 set	incl.
with	<input type="checkbox"/> 3276-26SW		Back base "reclinable", foldable with back panel, back height 25 (Upcharge)	05 1 set	
or	<input type="checkbox"/> 3276-27SW		Back base "reclinable", foldable with back panel, back height 30 (Upcharge)	1 set	
or	<input type="checkbox"/> 3276-28SW		Back base "reclinable", foldable with back panel, back height 35 (Upcharge)	1 set	
or	<input type="checkbox"/> 3276-41SW		Back base "reclinable", foldable with back panel, back height 40 (Upcharge)	1 set	
please specify angle adjustment range:					
□ 75°-95° □ 85°-105° □ 95° - 115°					
with	<input type="checkbox"/> 3274-7300		Tie down kit (minimum back rest height 30) *	06 1 set	
or	<input type="checkbox"/> 3274-7301		Back rest extension	07 1 set	
	<input type="checkbox"/> 3274-7370	346	Cushion for back rest extension (Upcharge)	08 1 pc	
* Use with back rest extension or head support if there is no support for shoulder and head mounted in the transport vehicle.					
with	<input type="checkbox"/> 3275-52SW		Push-bar	09 1 pc.	
or	<input type="checkbox"/> 3275-54SW		Push bar "wide"	10 1 pc.	
or	<input type="checkbox"/> 3275-53SW		Push-handles	11 1 pc	
or	<input type="checkbox"/> 3275-5502		One hand push-handle delta form (max. SW 36)	12 1 pc	

Brakes, fender

(Without this information your order can not be handled.)

	Order no.	Colour	Description	Qty.	Price €
with	<input type="checkbox"/> 3275-1750		HOGGI light sport brake	13 1 set	incl.
or	<input type="checkbox"/> 3221-1700		Standard brakes CLEO	14 1 set	incl.
with			Clothes cover	1 set	incl.
or	<input type="checkbox"/> 3275-160x		Fender	15 1 set	
or	<input type="checkbox"/> 3275-169x		Fender with integrated brake system	16 1 set	
with	<input type="checkbox"/> 3276-160x		Fender with height adjustable and detachable arm rests	17 1 set	
or	<input type="checkbox"/> 3276-169x		Fender with integrated brake system and with height adjustable and detachable arm rests	18 1 set	
with	<input type="checkbox"/> 7000-0051		PU-Armrest pads, 235 mm long	1 pair	incl.
or	<input type="checkbox"/> 7000-0052		PU-Armrest pads, 270 mm long	1 pair	incl.
	<input type="checkbox"/> 3240-8202		Therapy tray size 2, SW 24 and 26 cm *	19 1 pc	
	<input type="checkbox"/> 3240-8203		Therapy tray size 3, SW 28, 30 and 32 cm *	1 pc	
	<input type="checkbox"/> 3240-8204		Therapy tray size 4, SW 34 *	1 pc	
* only in conjunction with arm rests					

Aluminium fender

	Order no.	Colour	Description	Qty.
with	<input type="checkbox"/> 3277-1604		Wheel guards aluminum in frame color, 20-24"	20 1 set
or	<input type="checkbox"/> 3278-1604		Wheel guards aluminum in frame color, 20-24" with height adj. and detachable arm rests	21 1 set
or	<input type="checkbox"/> 3277-1694		Wheel guards aluminum in frame color 20-24" with integrated wheel lock levers	22 1 set
or	<input type="checkbox"/> 3278-1694		Wheel guards alu. frame color 20-24" with integrated brake system, height adjust. detachable arm rests	23 1 set

Rear wheels

with/or: Please select only one version!

	Order no.	Colour	Description	Qty.	Price €
with	<input type="checkbox"/> 1296-0xxx		Rear wheel HOGGI light	1 set	incl.
or	<input type="checkbox"/> 1297-1xxx		Rear wheel HOGGI light with drum brakes* (Upcharge)	1 set	
*in combination with push-bar or push-handles					
with	<input type="checkbox"/> 129x-x2xx		20" wheel diameter	1 set	incl.
or	<input type="checkbox"/> 129x-x3xx		22" wheel diameter	1 set	incl.
or	<input type="checkbox"/> 129x-x4xx		24" wheel diameter	1 set	incl.
with	<input type="checkbox"/> 129x-xx1x		Pneumatic tires 1"	24	1 set incl.
or	<input type="checkbox"/> 129x-xx2x		Flat free tires (20x1 3/8", 22x1" and 24x1")	25	1 set incl.
with	<input type="checkbox"/> 129x-xxx0		Without push-rims		
or	<input type="checkbox"/> 129x-xxx1		Push rims aluminium anodised, "standard" mounting	26	1 set incl.
or	<input type="checkbox"/> 129x-xxx2		Push rims aluminium anodised, "high" mounting	27	1 set incl.
or	<input type="checkbox"/> 129x-xxx3		Push rims stainless steel, "standard" mounting (Upcharge)	28	1 set
or	<input type="checkbox"/> 129x-xxx4		Push rims stainless steel, "high" mounting (Upcharge)	29	1 set
with		<input type="checkbox"/>	Push-rim mounting "narrow"		
or		<input type="checkbox"/>	Push-rim mounting "wide"		
or	<input type="checkbox"/> 1204-1155		Closed ergonomic hand rims 22", black anodized, for Hoggi-light wheels (Upcharge)	64	1 Set
or	<input type="checkbox"/> 1204-1156		Closed ergonomic hand rims 24", black anodized, for Hoggi-light wheels (Upcharge)		1 Set
with	<input type="checkbox"/> 1204-106X		"HOGGI Grip Plus" hand rim inserts black, for closed ergonomic hand rims (Upcharge)	65	1 Set
	<input type="checkbox"/> 1204-114X		Aluminum hand rim with plastic anti-slip coating, black, standard mounting * (Upcharge)	62	1 pair
	<input type="checkbox"/> 1204-114X		Aluminum hand rim with plastic anti-slip coating, black, high mounting * (Upcharge)		1 pair
* only in combination with HOGGI light rear wheels					
with	<input type="checkbox"/> 1296-4211		Pneumatic wheel 20x1" HOGGI Extra light * (Upcharge)	30	1 set
with	<input type="checkbox"/> 1296-4311		Pneumatic wheel 22x1" HOGGI Extra light * (Upcharge)		1 set
or	<input type="checkbox"/> 1296-4411		Pneumatic wheel 24x1" HOGGI Extra light * (Upcharge)		1 set
* incl. black alu push-rim					

Spoke guards // Old motifs // Sale

	Order no.	Colour	Description	Qty.	Price €
	<input type="checkbox"/> 1204-101x		Spoke guard, transparent *	31	1 set
	<input type="checkbox"/> 1204-101x		Spoke guard, transparent, for drum brakes wheels *		1 set
	<input type="checkbox"/> 1204-101x		Spoke guard, printed ** (Sale)	32	1 set
	<input type="checkbox"/> 1204-101x		Spoke guard, printed, for drum brakes ** (Sale)		1 set
* Motifs 2021 only possible in combination with 1204-2001-XXX					
** please enter code of your print design see attachment or brochure "spoke guard" // Only while stocks last!					
*** printed spoke guards are not available for 20" wheels					

Spoke guards // New motifs 2021 (with 10% dealer discount)

	Order no.	Colour	Description	Qty.	Price €
	<input type="checkbox"/> 1204-2001-000		Spoke guard Special interior foiling, single-coloured according to customer requirements ^{1 2}	1 set	
	<input type="checkbox"/> 1204-2001-2XX		Spoke guard Special interior foiling, printed ^{1 3}	1 set	
	<input type="checkbox"/> 1204-2001-SA		Spoke guard Special interior foiling, printed according to customer requirements ¹		Price on request
¹ only in combination with transparent spoke guard					
² please specify RAL tone					
³ Spoke guard motif selectable in "Order form spoke-guards" (available on the website)					

Front castors

with/or: Please select only one version!

	Order no.	Colour	Description	Qty.	Price €
with	<input type="checkbox"/> 1250-1002S		Wheel with rubber tire, black 100 x 32 mm	33	1 set incl.
or	<input type="checkbox"/> 1250-1003S		Wheel with rubber tire, black 125 x 32 mm		1 set incl.
or	<input type="checkbox"/> 1250-1005S		Wheel with rubber tire, black 140 x 40 mm		1 set incl.
or	<input type="checkbox"/> 1249-1005		Skater wheel with LED 100 x 24 mm (Upcharge)	34	1 set
or	<input type="checkbox"/> 1249-1006		Skater wheel with LED 125 x 24 mm (Upcharge)		1 set
or	<input type="checkbox"/> 1221-1003		PU-wheel with plastic design rim 100 x 34 mm (Upcharge)	35	1 set
or	<input type="checkbox"/> 1221-1004		PU-wheel with plastic design rim 125 x 34 mm (Upcharge)		1 set

Accessories

with/or: Please select only one version! SW=Seat width

	Order no.	Colour	Description	Qty.	Price €
	<input type="checkbox"/>	3275-2801	Footrest hanger 90°, lower leg length 13-26 cm	36 1 set	incl.
	<input type="checkbox"/>	3275-2802	Footrest hanger 90°, lower leg length 13-36 cm *	1 set	incl.
	<input type="checkbox"/>	3276-2900	Transfer support (not compatible with SA-0399)	37 1 pc.	
	<input type="checkbox"/>	3278-2801	Footrest hanger adj. 90-155°, LLL 19 - 25cm (Upcharge)	38 1 Set	
	<input type="checkbox"/>	3278-2802	Footrest hanger adj. 90-155°, LLL 19 - 35cm (Upcharge)	1 Set	
	<input type="checkbox"/>	SA-0430	Reinforced steel footplate-carrier (about 0,3 kg heavier)	1 pc.	incl.
			* applicable up from seat height 39 cm		
with	<input type="checkbox"/>	3275-2910	Footrest bracket angle adjustable, flip up type	39 1 set	incl.
or	<input type="checkbox"/>	3275-2960	Footrest bracket angle adjustable, with interlocking (Upcharge)	40 1 set	
with	<input type="checkbox"/>	3275-31SW	CLEO footrest ABS, flat version *	41 1 pc.	incl.
or	<input type="checkbox"/>	3275-30SW	CLEO footrest ABS, with heel stop * (Upcharge)	42 1 pc.	
with	<input type="checkbox"/>	3275-32SW	CLEO footrest aluminium, flat version in frame colour * (Upcharge)	43 1 pc.	
or	<input type="checkbox"/>	3275-33SW	CLEO footrest carbon, flat version * (Upcharge)	44 1 pc.	
			* The standard footrest width is 4 cm narrower than the seat width and fits almost all settings.		
			* Wider footrests are available on request. This can lead restrictions to changes in the settings after delivery.		
	<input type="checkbox"/>	3275-180x	Anti tip, swingaway, assembly left <input type="checkbox"/> =1 or right <input type="checkbox"/> =2	45 1 pc.	
	<input type="checkbox"/>	3275-180x	Tip assist, assembly left <input type="checkbox"/> =3 or right <input type="checkbox"/> =4	46 1 pc.	
	<input type="checkbox"/>	3275-1910	Bumper wheels	47 1 pc.	
	<input type="checkbox"/>	3231-7800	346 Abduction block size 1, swing away, detachable, (9 x 6 x 9 cm)	48 1 pc.	
	<input type="checkbox"/>	3231-7900	346 Abduction block size 2, swing away, detachable, (11 x 8x 10 cm)	1 pc.	
	<input type="checkbox"/>	32XX-71SW	346 Seat cushion; 3 cm thick	49 1 pc.	
	<input type="checkbox"/>	3275-7100	346 Back cushion, trough shaped, 2,5 cm thick	50 1 pc.	
	<input type="checkbox"/>	3276-7100	346 Back cushion for "reclinable" back, trough shaped, 2,5 cm thick	1 pc.	
	<input type="checkbox"/>	3201-3830	Universal headrest bracket (3245-7800 or third-party product) *	51 1 pc.	
	<input type="checkbox"/>	3245-7800	Headrest with occiput upholstery	52 1 pc.	
			* starting from backrest height 30 cm		
	<input type="checkbox"/>	3200-0401	Belt routing as a T-strap above the shoulders	53 1 pc.	
	<input type="checkbox"/>	1204-101x	Push-rim cover, silicone black fluted	1 set	
	<input type="checkbox"/>	3275-15SW	0XX Seat width enlarging set (enlarging range about 2 cm)	54 1 set	
	<input type="checkbox"/>	3276-15SW	0XX Seat width enlarging set (2 cm) for CLEO with "reclinable" back	55 1 set	
	<input type="checkbox"/>	3201-6400	Ankle hugger, size 1	56 1 pair	
	<input type="checkbox"/>	3203-6400	Ankle hugger, size 2	1 pair	
	<input type="checkbox"/>	3274-6100	Lap belt, black	57 1 pc.	
	<input type="checkbox"/>	3275-6702	4-Point lap belt, size S (SW 24 - 30)	58 1 pc.	
	<input type="checkbox"/>	3275-6703	4-Point lap belt, size M (SW 28 - 36)	1 pc.	

Custom-made

	Order no.	Colour	Description	Qty.	Price €
	<input type="checkbox"/>	3274-80SW	346 Contoured seat cushion with cover, black, SW 24 -32	59 1 pc.	
	<input type="checkbox"/>	3274-80SW	346 Contoured seat cushion with cover, black, SW 34	1 pc.	
	<input type="checkbox"/>	3275-7200	346 "HOGGI align" Contoured back cushion, black	60 1 pc.	
	<input type="checkbox"/>	SA-0286/2-KP	Abduction block 30 mm height extension kit *	1 set	
			* only in combination with abduction block and contoured seat cushion		
	<input type="checkbox"/>	SA-0210-MP	Special Frame colour, RAL code only (Upcharge)	1 pc.	
	<input type="checkbox"/>	6999-0014	060 Special Frame colour, Copper Gold (Upcharge)	1 pc.	
	<input type="checkbox"/>	6999-0015	CLEO component parts in various colours according to customer requirements *	Price per colour	
			* please note colours of all component parts below		
	<input type="checkbox"/>	6999-0017	Powder coated castor forks and both adapter for camber colour according to customer requirements *	1 set	
	<input type="checkbox"/>	SA-0279/4-KP	CLEO outdoor front end *	61 1 pc.	
			* additional costs for retrofitting		
			* only possible with two anti tippers (3275-1801 + 3275-1802).		
			* not compatible with 327X-1410 (base plates carbon)		
	<input type="checkbox"/>	6999-1026	Frame prepared for Outdoor front end	1 pc.	

Custom-made

Order no.	Colour	Description	Qty.	Price €
<input type="checkbox"/> SA-0390		Upholstery in special colour * * please enter colour code or upholstery code respectively		Price on request
<input type="checkbox"/> SA-0572-F-MP		Plastic anti-slip coating in custom colour: Pure white (RAL 9010), Silver grey (RAL 7001), Sky blue (RAL 5015), Traffic red (RAL 3020), Pure orange (RAL 2004)* (Upcharge)	63	1 pair
* only in combination with HOGGI light rear wheels				
<input type="checkbox"/> SA-0093-MP		Maxgrepp Original cell rubber coated push rims (Upcharge)	66	1 Set
<input type="checkbox"/> SA-0094-MP		Maxgrepp Ergo Para (closed push rims) (Upcharge)	67	1 Set
<input type="checkbox"/> SA-0095-MP		Maxgrepp Ergo Para CP (advanced closed push rims) (Upcharge)	68	1 Set
<input type="checkbox"/> SA-0384-MP		CLEO frame with parallel abduction (Upcharge)	69	1 Set
<input type="checkbox"/> SA-0394-MP		Coupling of two CLEOs one behind the other (Upcharge)	70	1 Set
<input type="checkbox"/> SA-0410-Min		Preparation for electric drive "TWION" from Alber		Price on request
<input type="checkbox"/> SA-0410/2-MP		Preparation for electric drive "e-fix" from Alber (Only together with armrests, for the joystick) (Upcharge)	71	1 pc.
<input type="checkbox"/> SA-0315/4-MP		Preparation for assembly of the one-hand drive of "ProActiv" (Upcharge)	72	1 pc.
<input type="checkbox"/> SA-0158-KP		Calf upholstery suitable to SW (Upcharge)		1 pc.
<input type="checkbox"/> SA-0427-MP		One-hand push bar with ratchet joint (Upcharge)	73	1 Set
<input type="checkbox"/> 3246-68SB		"CLEO Standard" trunk supports (Upcharge)	74	1 Set
<input type="checkbox"/> SA-0370/1-MP		Mounting adaption for CLEO angle and height adjustable trunk supports (Upcharge)		1 Set
<input checked="" type="checkbox"/> 3007-770X		CLEO angle and height adjustable trunk supports (Upcharge)	75	1 Set
	<input type="checkbox"/>	size 1 (8 x 10 x 2 cm)		
	<input type="checkbox"/>	size 2 (10 x 12 x 2 cm)		
	<input type="checkbox"/>	size 3 (12 x 13 x 2 cm)		
<input type="checkbox"/> SA-0370-MP		Mounting adaptation for Stealth "Nino" trunk supports (Upcharge)	76	Price on request
<input type="checkbox"/> SA-0505-MP		Locking of the integrated brake (Upcharge)	77	1 pc.
<input type="checkbox"/> SA-0594		Carbon fender (Upcharge)	78	Price on request

Custom-made Dynamic back base (not released for transport)


Order no.	Colour	Description	Qty.	Price €
<input type="checkbox"/> SA-0552/1-MP		Dynamic package in addition to back base angle adjustable from back height 30	79	1 Set
<input checked="" type="checkbox"/> SA-0431-MP		CLEO Back panel reinforced version	80	1 Set
				total
<input type="checkbox"/> SA-0552/2-MP		Dynamic package to CLEO Ti from back height 30	81	1 Set
<input checked="" type="checkbox"/> SA-0431-MP		CLEO Back panel reinforced version	81	1 Set
<input checked="" type="checkbox"/> SA-0555		CLEO Ti seat adapter reinforced steel version	82	1 pair
				total

Extra equipment

Order no.	Colour	Description	Qty.	Price €
<input type="checkbox"/> 3275-3100		CLEO trailer coupling for a "BOBBY-CAR" trailer		1 pc.
<input type="checkbox"/> SA-0366		Children name as film plot (23 x 95 mm) for sticking onto the back panel, Font: Calibri, film plot colour black or white		1 pc.
<input type="checkbox"/> 7000-0199		"Cosies" Thermofleece hand warmer for companion	83	1 set

Custom made 6999-0015 component parts in various colours

	Position	Colour	Component parts
<input type="checkbox"/>	A		Back panel
<input type="checkbox"/>	B		Back base
<input type="checkbox"/>	C		Sideplates, aluminium
<input type="checkbox"/>	D		Seat panel
<input type="checkbox"/>	E		Frame
<input type="checkbox"/>	F		Footrest hanger
<input type="checkbox"/>	G		Footrest, aluminium



Upholstery codes

346 Black	321 Spacer fabric, black	(special colour)	348 Grey Soft	(special colour)
	347 Berry Soft	(special colour)	349 Turquoise Soft	(special colour)


All HOGGI upholstery are made from flame retardant Trevira CS. Exception spacer fabric (**321**).

Frame colours

001 Black (RAL 9005)	020 Yellow green (RAL 6018)	028 Pearl violet (RAL 4011)	056 Turquoise blue (RAL 5018)
002 Silver	021 Grey aluminium (RAL 9007)	029 Pearl beige (RAL 1035)	057 Iron sparkling dark light
003 Pearl ruby red (RAL 3032)	022 White (RAL 9003)	032 Black dull matt (RAL 9005M)	058 Salmon orange (RAL 2012)
004 Pearl night blue (RAL 5026)	023 Pure orange (RAL 2004)	039 Signal violet (RAL 4008)	059 Pastel blue (RAL 5024)
005 Pearl gentian blue (RAL 5025)	024 Sparkling iron effect	040 Night blue (RAL 5022)	061 Light pink (RAL 3015)
017 Traffic yellow (RAL 1023)	025 Signal blue (RAL 5005)	050 White matt (RAL 9003)	064 Olive green (RAL 6003)
018 Flame red (RAL 3000)	026 Vermillion (RAL 2002)	053 Light green (RAL 6027)	065 Raspberry red (RAL 3027)
019 Sky blue (RAL 5015)	027 Heather violet (RAL 4003)	054 Ocean blue (RAL 5020)	066 Pastel green (RAL 6019)

SA-0210-MP Special colours on request (Enter RAL code **Upcharge**

How to measure



SH = Seat height **LLL** = Lower leg length **SD** = seat depth **BH** = Backrest height

Seat measurements

Seat width (SW)	24	26	28	30	32	34
Seat depth standard (SD)	22-30	24-32	26-34	28-36	30-38	32-40
Seat depth long (SD)	28-36	30-38	32-40	34-42	36-44	38-46

Frame size 1 20" rear wheels					Frame size 2 22" rear wheels					Frame size 3 24" rear wheels				
SH	Seat angle				SH	Seat angle				SH	Seat angle			
29,0	2,5	3,5			32,0	2,5	3,5			36,0	2,5	3,0		
29,5	0,5	3,5	4,5		32,5	0,5	3,5	4,0		36,5	0,5	3,0	4,0	
30,0	3,5	6,5	7,0		33,0	3,5	6,0	6,5		37,0	3,0	6,0	6,5	
31,5	3,5	8,0	10,5	11,0	34,5	3,5	7,5	10		38,5	3,0	7,0	9,5	10,0
32,5	0,5	3,5	8,5		35,5	0,5	3,5	7,5	11,5	39,5	0,5	3,0	7,5	11,0
33,5	2,5	3,5			36,5	2,5	3,5			40,5	2,5	3,0		
34,0	0,5	3,5	4,5		37,0	0,5	3,5	4,0		41,0	0,5	3,0	4,0	
35,0	3,5	6,5	7,0		38,0	3,5	6,0	6,5		42,0	3,0	6,0	6,5	
36,0	3,5	8,0	10,5	11,0	39,0	3,5	7,5	10,0		43,0	3,0	7,0	9,5	10,0
37,5	0,5	3,5	8,0		40,5	3,5	7,5	11,5		44,5	0,5	3,0	7,5	11,0
38,5	2,0	3,5	7,0		41,5	2,0	3,5	7,0	11,5	45,5	2,0	3,0	6,5	10,5
39,0	3,5	5,0	9,0		42,0	3,5	5,0	8,5		46,0	3,0	4,5	8,0	































Spoke guard print design



Colour code

201	HOGGI blue	205	Train
202	Wonderland	206	Dragon
203	Mouse	207	Circus
204	Water world		

CONFIGURATIONS // ACCESSORIES // CLEO

<p>01</p> 	<p>02</p> 	<p>03</p> 	<p>04</p> 	<p>05</p> 
Base plates, aluminium	Base plates, carbon	Base plates „Comfort“ long wheelbase	Back base with back panel, rigid	Back base with back panel, reclinable
<p>06</p> 	<p>07</p> 	<p>08</p> 	<p>09</p> 	<p>10</p> 
Tie down kit	Back rest extension	Cushion for back rest extension	Push-bar	Push bar „wide“
<p>11</p> 	<p>12</p> 	<p>13</p> 	<p>14</p> 	<p>15</p> 
Push-handles	One hand push-handle delta form	HOGGI light sport brake	Standard brakes	Fender
<p>16</p> 	<p>17</p> 	<p>18</p> 	<p>19</p> 	<p>20</p> 
Fender with integrated brake system	Fender with arm rests	Fender with arm rests and integr. brake system	Therapy tray	Aluminium fender
<p>21</p> 	<p>22</p> 	<p>23</p> 	<p>24</p> 	<p>24</p> 
Aluminium fender with arm rests	Aluminium fender with brake and arm rests	Aluminium fender with integr. brake system	Pneumatic tires	Pneumatic tires and drum brake
<p>25</p> 	<p>25</p> 	<p>26</p> 	<p>27</p> 	<p>28</p> 
Polyurethan tires	Polyurethan tires and drum brake	Push rims aluminium anodised, standard	Push rims aluminium anodised, high	Push rims stainless steel, standard

CONFIGURATIONS // ACCESSORIES // CLEO



29 Push rims stainless steel, high



30 HOGGI Extra light



31 Spoke guard, transparent



32 Spoke guard, printed



33 Wheel with rubber tire




34 Skater wheel with LED



35 PU-wheel with plastic design rim



36 Footrest hanger 90°



37 Transfer support




38 Adjustable footrest hanger 90-155°



39 Footrest bracket




40 Footrest bracket with interlocking



41 ABS, flat version



42 ABS with heel stop



43 Footrest: Aluminium in frame colour




44 Footrest: carbon



45 Anti tip




46 Tip assist



47 Bumper wheels



48 Abduction block



49 Seat cushion



50 Back cushion




51 Universal headrest bracket




52 Headrest with occiput upholstery



53 Belt routing as a T-strap



54 Seat width enlarging set



55 Seat width enlarging set for reclinable back



56 Ankle hugger



57 Lap belt



58 4-Point lap belt

CONFIGURATIONS // ACCESSORIES // CLEO



59

Contoured seat cushion with cover



60

„HOGGI align“ Contoured back cushion




61

Outdoor front end



62

Aluminum hand rim with plastic anti-slip coating



63

Aluminum hand rim with plastic anti-slip coating



64

Closed ergonomic hand rims



65

HOGGI Grip Plus



66

Maxgrepp Original push rims




67

Maxgrepp Ergo Para push rims




68

Maxgrepp Ergo Para CP push rims



69

CLEO frame with parallel abduction



70

Coupling of two CLEOs one behind the other



71

Electric drive „e-fix“




72

One-hand drive of „ProActiv“




73

One-hand push bar with ratchet join




74

Standard trunk supports



75

Angle and height adjustable trunk supports



76

Stealth „Nino“ trunk supports



77

Locking of the integrated brake




78

Carbon fender



79

Dynamic package CLEO



80

Back panel reinforced version



81

Dynamic package CLEO Ti



82

Back panel & seat adapter reinforced version



83

„Cosies“ Thermofleece hand warmer