

Partnering together to mobilise kids



#### CLEO Active-wheelchair for children

Crash tested according to ISO 7176-19 and RESNA WC-4:2012: Section 19 Order form incl. Exportprices (€) valid until revoked, release date 01/2025

#### Basic configuration inclusive:

- Seat- and back panel
- Clothes cover
- Footrest hanger 90°
- HOGGI light rear wheels
- HOGGI light brake CLEO

Price:

□ Offer	□ Order	
Order number / reference	Date	
Invoice address		
Company	Customer number	
Street	City, Zip code	
E-Mail	Telephone	
Fax		
Delivery address		
Company		
Street	City, Zip code	
Basic configuration		

- Powder coated aluminium frame with growth capacity
- 3 frame sizes:
- Size 1: Seat width (SW) 24 30 cm in 2 cm steps
- Size 2: Seat width (SW) 24 34 cm in 2 cm steps
- Size 3: Seat width (SW) 26 34 cm in 2 cm steps
- Seat heights: 29 cm 46 cm
- Seat angle: 0,5° 11,5°
- Footrest hanger 90° inclusive ABS-footrest

**Minimum information** (Without this information your order can not be handled.)

Frame colour		(Please see	table "fran	ne colours"	')				
Seat width in cm		□ 26	□ 28	□ 30	□ 32	□ 34	Seat height *		(29 - 46 cm)
Camber in degree	□ 6° (3	3274-1292)	□ 9° (3	274-1293)	□ 12° (3	3274-1294)	Seat angle *		(0,5° - 11,5°)
* Please see attachment of "seat measurements"									

# Presettings (Without this information your order can not be handled.) Seat depth (SD) 1 (22 - 46 cm) Lower leg length (LLL) (13 - 36 cm) Backrest angle 2 (-10 to +10) Comfort degree (-5 to 3 cm)

<sup>1</sup> Please see table seat measurements

<sup>2</sup> Angle compared to the vertical line

<sup>3</sup> Please see drawing "active degree"

At missing specification we will adjust to middle sized measurements

- Backrest angle: -10° to +10°
- Seat- and back panel
- Clothes cover
- HOGGI light rear wheels
- inclusive push rims
- HOGGI light wheel locks
- Camber 6°, 9° or 12°
- Load capacity: 60 kg

Basic	varia	tion		with/or: P	lease select only	y one version!	SW=Seat width
		Order no.	Colou	r Description		Qty.	Price €
		3274-00SW		CLEO size 1 (20")		1 pc.	
with		3274-1420		Base plates, aluminum size 1, black coated		01 1 set	incl.
or		3274-1410		Base plates, carbon *	(Upcharge)	02 1 set	
or		3274-1430		Base plates alu "Comfort" long wheelbase **		03 1 set	incl.
		3275-00SW		CLEO size 2 (22")		1 pc.	
with		3275-1420		Base plates, aluminum size 2, black coated		01 1 set	incl.
or		3275-1410		Base plates, carbon *	(Upcharge)	02 1 set	
or		3275-1430		Base plates alu "Comfort" long wheelbase **		03 1 set	incl.
		3276-00SW		CLEO size 3 (24")		1 pc.	
with		3276-1420		Base plates, aluminum size 3, black coated		01 1 set	incl.
or		3276-1410		Base plates, carbon *	(Upcharge)	02 1 set	
or		3276-1430		Base plates alu "Comfort" long wheelbase **		03 1 set	incl.
				* can not be combined with tie down kit			
				** reduction of the active degree by 5 cm			
with		3274-25SW		Back base "rigid" with back panel, back height 20		04 1 set	incl.
or		3274-26SW		Back base "rigid" with back panel, back height 25		1 set	incl.
or		3274-27SW		Back base "rigid" with back panel, back height 30		1 set	incl.
or		3274-28SW		Back base "rigid" with back panel, back height 35		1 set	incl.
				Back base "reclinable", foldable with back panel, back		05	
with		3276-26SW		height 25	(Upcharge)	05 1 set	
				Back base "reclinable", foldable with back panel, back			
or		3276-27SW		height 30	(Upcharge)	1 set	
				Back base "reclinable", foldable with back panel, back			
or		3276-28SW		height 35	(Upcharge)	1 set	
				Back base "reclinable", foldable with back panel, back			
or		3276-41SW		height 40	(Upcharge)	1 set	
				please specify angle adjustment range:			
				□ 75°-95° □ 85°-105° □ 95° - 115°		- <u> </u>	
with		3274-7300		Tie down kit (minimum back rest height 30) *		06 1 set	
or		3274-7301		Back rest extension		07 1 set	
		3274-7370	346	Cushion for back rest extention	(Upcharge)	08 1 pc	
				* Use with back rest extension or head support if there is			
				no support for shoulder and head mounted in the transport	t		
with		3275-52SW		vehicle.		09 1 pc.	
-		3275-523W		Push-bar		10 1 pc.	
or or		3275-543W		Push bar "wide"		11 1 pc	
				Push-handles			
or		3275-5502		One hand push-handle delta form (max. SW 36)		12 1 pc	

#### Brakes, fender

	Order no.	Colour Description		Qty.	Price #
with	3275-1750	HOGGI light sport brake	13	1 set	inc
or	3221-1700	Standard brakes CLEO	14	1 set	inc
with		Clothes cover		1 set	inc
or	3275-160x	Fender	15	1 set	
or	3275-169x	Fender with integrated brake system	16	1 set	
with	3276-160x	Fender with height adjustable and detachable arm rests	17	1 set	
or	3276-169x	Fender with integrated brake system and with height adjustable and detachable arm rests	18	1 set	
with	7000-0051	PU-Armrest pads, 235 mm long		1 pair	inc
or	7000-0052	PU-Armrest pads, 270 mm long		1 pair	inc
	3240-8202	Therapy tray size 2, SW 24 and 26 cm *	19	1 pc	
	3240-8203	Therapy tray size 3, SW 28, 30 and 32 cm *	_	1 pc	
	3240-8204	Therapy tray size 4, SW 34 *		1 pc	
		* only in conjunction with arm rests			

Alumi	nium	fender			
		Order no.	Colour Description		Qty.
with		3277-1604	Wheel guards aluminum in frame color, 20-24"	20	1 set
			Wheel guards aluminum in frame color, 20-24" with	24	
or		3278-1604	height adj. and detachable arm rests	21	1 set
			Wheel guards aluminum in frame color 20-24" with	22	
or		3277-1694	integrated wheel lock levers	22	1 set
			Wheel guards alu. frame color 20-24" with integrated	22	
r		3278-1694	brake system, height adjust. detachable arm rests	23	1 set

Rear	whee	ls		with/or: Pl	ease select only	y one version!
		Order no.	Colour Description		Qty.	Price €
with		1296-0xxx	Rear wheel HOGGI light		1 set	incl.
or		1297-1xxx	Rear wheel HOGGI light with drum brakes*	(Upcharge)	1 set	
			*in combination with push-bar or push-handles			
with		129x-x2xx	20" wheel diameter		1 set	incl.
or		129x-x3xx	22" wheel diameter		1 set	incl.
or		129x-x4xx	24" wheel diameter		1 set	incl.
with		129x-xx1x	Pneumatic tires 1"		24 1 set	incl.
or		129x-xx2x	Flat free tires (20x1 3/8", 22x1" and 24x1")		25 1 set	incl.
with		129x-xxx0	Without push-rims			
or		129x-xxx1	Push rims aluminium anodised, "standard" mounting		26 1 set	incl.
or		129x-xxx2	Push rims aluminium anodised, "high" mounting		27 1 set	incl.
or		129x-xxx3	Push rims stainless steel, "standard" mounting	(Upcharge)	28 1 set	
or		129x-xxx4	Push rims stainless steel, "high" mounting	(Upcharge)	29 1 set	
with			Push-rim mounting "narrow"			
or			Push-rim mounting "wide"			
			Closed ergonomic hand rims 22 ", black anodized, for		65	
or		1204-1155	Hoggi-light wheels	(Upcharge)	65 1 Set	
			Closed ergonomic hand rims 24 ", black anodized, for			
or		1204-1156	Hoggi-light wheels	(Upcharge)	1 Set	
			"HOGGI Grip Plus" hand rim inserts black, for closed		66 1 Set	
with		1204-106X	ergonomic hand rims	(Upcharge)	<sup>20</sup> 1 Set	
	_		Aluminum hand rim with plastic anti-slip coating, black,		63	
		1204-114X	standard mounting *	(Upcharge)	1 pair	
			Aluminum hand rim with plastic anti-slip coating, black,			
		1204-114X	high mounting *	(Upcharge)	1 pair	
			* only in combination with HOGGI light rear wheels			
with		1296-4211	Pneumatic wheel 20x1" HOGGI Extra light *	(Upcharge)	30 1 set	
with		1296-4311	Pneumatic wheel 22x1" HOGGI Extra light *	(Upcharge)	1 set	
or		1296-4411	Pneumatic wheel 24x1" HOGGI Extra light *	(Upcharge)	1 set	
			* incl. black alu push-rim			

Order no.	Colour Description			Qty.	Price	€
1204-101x	Spoke guard, transparent *		31	1 set		
1204-101x	Spoke guard, transparent, for drum brakes wheels *			1 set		
1204-101x	Spoke guard, printed **	(Sale)	32	1 set		
1204-101x	Spoke guard, printed, for drum brakes **	(Sale)		1 set		
	* Motifs 2021 only possible in combination with 1204-2001-XX	X				
	** please enter code of your print design see attachment or b	rochure				
	"spoke guard" // Only while stocks last!					
	*** printed spoke guards are not available for 20" wheels					

Spoke gua	rds // New mot	ifs 2021 (with 10% dealer discount)	
	Order no.	Colour Description	Qty. Price €
		Spoke guard Special interior foiling, single-coloured	
	1204-2001-00	0 according to customer requirements <sup>12</sup>	1 set
	1204-2001-2X		1 set
		Spoke guard Special interior foiling, printed according to	
	1204-2001-SA	Customer requirements <sup>1</sup>	Price on request
		<sup>1</sup> only in combination with transparent spoke guard	
		<sup>2</sup> please specify RAL tone	
		<sup>3</sup> Spoke guard motif selectable in "Order form spoke-	
		guards" (available on the website)	

Front castors with/or: Please select						
		Order no.	Colour Description		Qty.	Price €
with		1250-1002S	Wheel with rubber tire, black 100 x 32 mm	33	1 set	incl.
or		1250-1003S	Wheel with rubber tire, black 125 x 32 mm		1 set	incl.
or		1250-1005S	Wheel with rubber tire, black 140 x 40 mm		1 set	incl.
or		1249-1005	Skater wheel with LED 100 x 24 mm	(Upcharge) 34	1 set	
or		1249-1006	Skater wheel with LED 125 x 24 mm	(Upcharge)	1 set	
or		1221-1003	PU-wheel with plastic design rim 100 x 34 mm	(Upcharge) 35	1 set	
or		1221-1004	PU-wheel with plastic design rim 125 x 34 mm	(Upcharge)	1 set	

Acces			Colou	r Description	ase select only		Qty.	Price	
		3275-2801		Footrest hanger 90°, lower leg length 13-26 cm			1 set		incl
		3275-2802		Footrest hanger 90°, lower leg length 13-36 cm *			1 set		incl
		3276-2900		Transfer support (not compatible with SA-0399)			1 pc.		
		3278-2801			(Upcharge)		1 Set		
		3278-2802			(Upcharge)		1 Set		
		SA-0430		Reinforced steel footplate-carrier (about 0,3 kg heavier)	(openarge)		1 pc.		incl
		3A-0430		* applicable up from seat height 39 cm			i pc.	······	IIICI
ith		3275-2910		Footrest bracket angle adjustable, flip up type		39	1 set	i	incl
		3275-2960			(Upcharge)		1 set		
th		3275-31SW		CLEO footrest ABS, flat version *	(openaige)		1 pc.	i	incl
		3275-30SW			(Upcharge)				
th		3275-32SW			(Upcharge)				
		3275-33SW			(Upcharge)		1 pc.		
	_	0210 00011		* The standard footrest width is 4 cm narrower than the seat					
				fits almost all settings.	wiutii aliu				
				* Wider footrests are available on request. This can lead rest	rictions to				
				changes in the settings after delivery.					
		3275-180x		Anti tip, swingaway, assembly left $\Box$ =1 or right $\Box$ =2		45	1 pc.		
		3275-180x		Tip assist, assembly left $\Box$ =3 or right $\Box$ =4			1 pc.		
		3275-1910		Bumper wheels			1 pc.		
		3231-7800	346	Abduction block size 1, swing away, detachable, (9×6×9 cm	1)		1 pc.		
		3231-7900	346	Abduction block size 2, swing away, detachable, (1 x 8x 10 c			1 pc.		
		32XX-71SW	346	Seat cushion; 3 cm thick	411 <i>)</i>		1 pc.		
		3275-7100	346	Back cushion, trough shaped, 2,5 cm thick			1 pc.		
		3276-7100	346	Back cushion for "reclinable" back, trough shaped, 2,5 cm	thick		-		
		3201-3830	340	Universal headrest bracket (3245-7800 or third-party produ			1 pc. 1 pc.		
		3245-7800		Headrest with occiput upholstery			1 pc.		
		3243-7000		* starting from backrest height 30 cm		52	i pc.		
		3200-0401		Belt routing as a T-strap above the shoulders		52	1 pc.		
		1204-101x		Push-rim cover, silicone black fluted			1 set		
		3275-15SW	0XX				1 set		
		3276-15SW	0XX		hook		1 set		
		3201-6400	0//	<b>3 3 ( )</b>	DACK		1 pair		
		3201-0400		Ankle hugger, size 1			1 pair		
		3274-6100		Ankle hugger, size 2					
		3275-6702		Lap belt, black			1 pc. 1 pc.		
				4-Point lap belt, size S (SW 24 - 30)					
		3275-6703		4-Point lap belt, size M (SW 28 - 36)			1 pc.		
usto	om-ma		Colou	r Description			Otv	Price	:
				-			Qty.	FIICE	
		3274-80SW	346	Contoured seat cushion with cover, black, SW 24 -32			1 pc.		
		3274-80SW	346	Contoured seat cushion with cover, black, SW 34			1 pc.		
		3275-7200	346	"HOGGI align" Contoured back cushion, black		60	1 pc.		
		SA-0286/2-KP		Abduction block 30 mm height extension kit *			1 set		
				* only in combination with abduction block and contoured se	eat cushion				
		SA-0210-MP		Special Frame colour, RAL code only	(Upcharge)		1 pc.		
		6999-0014	060	• • •	(Upcharge)		1 pc.		
	_	0000 0015		CLEO component parts in various colours					
		6999-0015			Price per co	lour			
				* please note colours of all component parts below					
	_			Powder coated castor forks and both adapter for camber of	olour				
		6999-0017		according to customer requirements *			1 set		
		SA-0279/4-KP		CLEO outdoor front end *			1 pc.		

	according to customer requirements *	1 set
SA-0279/4-KP	CLEO outdoor front end *	61 1 pc.
	* additional costs for retrofitting	
	* only possible with two anti tippers (3275-1801 + 3275-1802).	
	* not compatible with 327X-1410 (base plates carbon)	
6999-1026	Frame prepared for Outdoor front end	1 pc.
	Mounting adapter for Outdoor front end aid for "phil &	<u></u>
SA-0619/2-KP	teds" sibling seat or universal car seat adapter	<sup>62</sup> 1 set

Custom-ma	ade							
	Order no.	Colou	r Description			Qty.	Price	€
	SA-0390		Upholstery in special colour *			Price	on reque	əst
			* please enter colour code or upholstery code respectively					
			Plastic anti-slip coating in custom colour: Pure white					
			(RAL 9010), Silver grey (RAL 7001), Sky blue (RAL		64			
	SA-0572-F-MF	)	5015), Traffic red (RAL 3020), Pure orange (RAL 2004)*	(Upcharge)		1 pair		
			* only in combination with HOGGI light rear wheels					
	SA-0093-MP		Maxgrepp Original cell rubber coated push rims	(Upcharge)	67	1 Set		
	SA-0094-MP		Maxgrepp Ergo Para (closed push rims)	(Upcharge)	68	1 Set		
	SA-0095-MP		Maxgrepp Ergo Para CP (advanced closed push rims)	(Upcharge)	69	1 Set		
	SA-0384-MP		CLEO frame with parallel abduction	(Upcharge)	70	1 Set		
	SA-0394-MP		Coupling of two CLEOs one behind the other	(Upcharge)	71	1 Set		
	SA-0410-Min		Preparation for electric drive "TWION" from Alber			Price	on reque	est
			Preparation for electric drive "e-fix" from Alber		70			
	SA-0410/2-MP	)	(Only together with armrests, for the joystick)	(Upcharge)	12	1 pc.		
			Preparation for assembly of the one-hand drive of		73			
	SA-0315/4-MP		"ProActiv"	(Upcharge)	13	1 pc.		
	SA-0158-KP		Calf upholstery suitable to SW	(Upcharge)		1 pc.		
	SA-0427-MP		One-hand push bar with ratchet joint	(Upcharge)	74	1 Set		
	3246-68SB		"CLEO Standard" trunk supports	(Upcharge)	75	1 Set		
			Mounting adaption for CLEO angle and height adjustable					
	SA-0370/1-MP		trunk supports	(Upcharge)		1 Set		
	🗷 3007-770X		CLEO angle and height adjustable trunk supports	(Upcharge)	76	1 Set		
			size 1 (8 x 10 x 2 cm)					
			size 2 (10 x 12 x 2 cm)					
			size 3 (12 x 13 x 2 cm)					
	SA-0370-MP		Mounting adaptation for Stealth "Nino" trunk supports	(Upcharge)	77	Price	on reque	əst
	SA-0505-MP		Locking of the integrated brake	(Upcharge)	78	1 pc.		
	SA-0594		Carbon fender	(Upcharge)	79	Price	on reque	est

## Custom-made Dynamic back base (not released for transport)

Order no. Co	lour Description	Qty.	Price	€
	Dynamic package in addition to back base angle	80		
SA-0552/1-MP	adjustable from back height 30	<sup>80</sup> 1 Set		
🗷 SA-0431-MP	CLEO Back panel reinforced version	81 1 Set		
		total		
SA-0552/2-MP	Dynamic package to CLEO Ti from back height 30	82 1 Set		
🗷 SA-0431-MP	CLEO Back panel reinforced version	83 1 Set		
 🗷 SA-0555	CLEO Ti seat adapter reinforced steel version	83 1 pair		
		total		

## Extra equipment

Order no.	Colour Description	Qty.	Price	€
3275-3100	CLEO trailer coupling for a "BOBBY-CAR" trailer	1 pc.		
	Children name as film plot (23 x 95 mm) for sticking onto the back			
SA-0366	panel, Font: Calibri, film plot colour black or white	1 pc.		
7000-0199	"Cosies" Thermofleece hand warmer for companion	84 1 set		

Cust	om mad	le 6999-0015	comp	onent parts in various co	lours		
	I	Position	Colou	r Component parts			
		А		Back panel			
						Α —	HOOGE
		В		Back base			
						в —	
		С		Sideplates, aluminium		C —	
						C	
		D		Seat panel		D —	V
						D	
		E		Frame			
						Е —	
		F		Footrest hanger		F —	
		G		Footrest, aluminium		G	
Link		odoo					
opno	olstery c	oues					
346	Black		321	Spacer fabric, black	(special colour)	348	Grey Soft (special colour)
			347	Berry Soft	(special colour)	349	Turquoise Soft (special colour)


All HOGGI upholstery are made from flame retardant Trevira CS. Exception spacer fabric (321).

## Frame colours

001 Black	(RAL 9005)	020 Yellow green	(RAL 6018)	028 Pearl violet	(RAL 4011)	056 Turqoise blue	(RAL 5018)
002 Silver		021 Greyaluminium	(RAL 9007)	029 Pearl beige	(RAL 1035)	057 Iron sparkling d	ark light
003 Pearl ruby red	(RAL 3032)	022 White	(RAL 9003)	032 Black dull matt	(RAL 9005M)	058 Salmon orange	(RAL 2012)
004 Pearl night blue	(RAL 5026)	023 Pure orange	(RAL 2004)	039 Signal violet	(RAL 4008)	059 Pastel blue	(RAL 5024)
005 Pearl gentiant blu	e (RAL 5025)	024 Sparkling iron ef	fect	040 Night blue	(RAL 5022)	061 Light pink	(RAL 3015)
017 Traffic yellow	(RAL 1023)	025 Signal blue	(RAL 5005)	050 White matt	(RAL 9003)	064 Olive green	(RAL 6003)
018 Flame red	(RAL 3000)	026 Vermillion	(RAL 2002)	053 Light green	(RAL 6027)	065 Raspberry red	(RAL 3027)
019 Skyblue	(RAL 5015)	027 Heather violet	(RAL 4003)	054 Ocean blue	(RAL 5020)	066 Pastel green	(RAL 6019)

SA-0210-MP Special colours on request (Enter RAL code Upcharge

## How to measure



SH = Seat height LLL = Lower leg length SD = seat depth BH = Backrest height

#### Seat measurements

Seat width (SW)	24	26	28	30	32	34
Seat depth standard (SD)	22-30	24-32	26-34	28-36	30-38	32-40
Seat depth long (SD)	28-36	30-38	32-40	34-42	36-44	38-46


	Frame					Frame size 2						e size		
	20" rea	ar whe	els			<u>22" rea</u>	ar whe	els			24" rea	ar whe	els	
SH		Seat	angle		SH		Seat	angle		SH		Seat	angle	
29,0	2,5	3,5			32,0	2,5	3,5			36,0	2,5	3,0		
29,5	0,5	3,5	4,5		32,5	0,5	3,5	4,0		36,5	0,5	3,0	4,0	
30,0	3,5	6,5	7,0		33,0	3,5	6,0	6,5		37,0	3,0	6,0	6,5	
31,5	3,5	8,0	10,5	11,0	34,5	3,5	7,5	10		38,5	3,0	7,0	9,5	10,0
32,5	0,5	3,5	8,5		35,5	0,5	3,5	7,5	11,5	39,5	0,5	3,0	7,5	11,0
33,5	2,5	3,5			36,5	2,5	3,5			40,5	2,5	3,0		
34,0	0,5	3,5	4,5		37,0	0,5	3,5	4,0		41,0	0,5	3,0	4,0	
35,0	3,5	6,5	7,0		38,0	3,5	6,0	6,5		42,0	3,0	6,0	6,5	
36,0	3,5	8,0	10,5	11,0	39,0	3,5	7,5	10,0		43,0	3,0	7,0	9,5	10,0
37,5	0,5	3,5	8,0		40,5	3,5	7,5	11,5		44,5	0,5	3,0	7,5	11,0
38,5	2,0	3,5	7,0		41,5	2,0	3,5	7,0	11,5	45,5	2,0	3,0	6,5	10,5
39,0	3,5	5,0	9,0		42,0	3,5	5,0	8,5		46,0	3,0	4,5	8,0	



Colour code

201	HOGGI blue	205	Train
202	Wonderland	206	Dragon
203	Mouse	207	Circus
204	Water world		

# CONFIGURATIONS // ACCESSORIES // CLEO



Polyurethan tires and drum brake

anodised, standard

Push rims aluminium anodised, high

Push rims stainless steel, standard

# CONFIGURATIONS // ACCESSORIES // CLEO



Seat width enlarging set

Seat width enlarging set for reclinable back

Ankle hugger

Lap belt

4-Point lap belt

# CONFIGURATIONS // ACCESSORIES // CLEO





"Cosies" Thermofleece hand warmer