



CLEO Ti Active-wheelchair for children

Crash tested according to ISO 7176-19 and RESNA WC-4:2012: Section 19

Basic configuration inclusive:

- Seat- and back panel
- Footrest hanger 90°
- Footrest ABS
- HOGGI light rear wheels
- HOGGI light brake CLEO

Price:

Order form incl. Exportprices (€) valid until revoked, release date 04/2024

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<hr/>	
Invoice address	
Company	Customer number
<hr/>	
Street	City, Zip Code
<hr/>	
E-Mail	Telephone
<hr/>	
Fax	
<hr/>	
Delivery address	
Company	
<hr/>	
Street	City, Zip Code
<hr/>	



Basic configuration

- Powder coated aluminium frame with growth capacity
- 3 frame sizes:
 - Size 1: Seat width (SW) 24 - 30 cm in 2 cm steps
 - Size 2: Seat width (SW) 24 - 34 cm in 2 cm steps
 - Size 3: Seat width (SW) 26 - 34 cm in 2 cm steps
- Convertible to other seat widths
- Seat heights front: 38 cm to 48 cm
- Footrest hanger 90° incl. ABS-footrest
- Backrest angle: -5° to +35°
- Seat- and back panel
- Clothes cover
- **HOGGI** light rear wheels inclusive push rims
- **HOGGI** light wheel locks
- Camber: 6°, 9° or 12°
- Load capacity: 60 kg

Minimum information (Without this information your order can not be handled.)

Frame colour (Please see table "frame colours")

Seat width (cm) 24 26 28 30 32 34

Camber 6° (3274-1292) 9° (3274-1293) 12° (3274-1294)

Seat height (cm) (38 - 48) see table „seat measurements“

Presettings

Seat depth (SD) ¹ (22 - 46 cm) **Lower leg length (LLL)** (13 - 36 cm)

Backrest angle ² (-5° to +35°)

Active degree ³ (0 - 8 cm)

¹ Please see table seat measurements ² Angle compared to the vertical line
At missing specification we will adjust to middle sized measurements

³ Please see drawing "active degree"

Basic variation

	Order no.	Colour	Description	Qty	Price €
	<input type="checkbox"/>	3277-00SW	CLEO Ti size 1 (20")	1 pc	
	<input checked="" type="checkbox"/>	SA-0399	CLEO Ti 45° tilt in space (-5° to +45°)	(Upcharge 01)	1 set
with	<input type="checkbox"/>	3277-1420	Base plates, aluminum size 1, black coated	02	1 set incl.
or	<input type="checkbox"/>	3277-1430	Base plates alu "Comfort" long wheelbase **	03	1 set incl.
	<input checked="" type="checkbox"/>	3275-1801	Anti tip, swivels away, assembly left	04	1 pc
	<input checked="" type="checkbox"/>	3275-1802	Anti tip, swivels away, assembly right		1 pc
	<input type="checkbox"/>	3278-00SW	CLEO Ti size 2 (22")	1 pc	
	<input checked="" type="checkbox"/>	SA-0399	CLEO Ti 45° tilt in space (-5° to +45°)	(Upcharge 04)	1 set
with	<input type="checkbox"/>	3278-1420	Base plates, aluminum size 2, black coated	02	1 set incl.
or	<input type="checkbox"/>	3278-1430	Base plates alu "Comfort" long wheelbase **	03	1 set incl.
	<input checked="" type="checkbox"/>	3275-1801	Anti tip, swivels away, assembly left	04	1 pc
	<input checked="" type="checkbox"/>	3275-1802	Anti tip, swivels away, assembly right		1 pc
	<input type="checkbox"/>	3279-00SW	CLEO Ti size 3 (24")	1 pc	
	<input checked="" type="checkbox"/>	SA-0399	CLEO Ti 45° tilt in space (-5° to +45°)	(Upcharge 01)	1 set
with	<input type="checkbox"/>	3279-1420	Base plates, aluminum size 3, black coated	02	1 set incl.
or	<input type="checkbox"/>	3279-1430	Base plates alu "Comfort" long wheelbase **	03	1 set incl.
			** reduction of the active degree by 5 cm		1 set
	<input checked="" type="checkbox"/>	3275-1801	Anti tip, swivels away, assembly left	04	1 pc
	<input checked="" type="checkbox"/>	3275-1802	Anti tip, swivels away, assembly right		1 pc
with	<input type="checkbox"/>	3279-26SW	Back base "angle adjustable" with back panel, back height 25	05	1 set incl.
or	<input type="checkbox"/>	3279-27SW	Back base "angle adjustable" with back panel, back height 30		1 set incl.
or	<input type="checkbox"/>	3279-28SW	Back base "angle adjustable" with back panel, back height 35		1 set incl.
or	<input type="checkbox"/>	3279-29SW	Back base "angle adjustable" with back panel, back height 40		1 set incl.
			please specify angle adjustment range:		
			<input type="checkbox"/> 75° - 95° <input type="checkbox"/> 85° - 105° <input type="checkbox"/> 95° - 115°		
with	<input type="checkbox"/>	3274-7300	Tie down kit (minimum back rest height 30) *	06	1 set
or	<input type="checkbox"/>	3274-7301	Back rest extension	07	1 set
	<input type="checkbox"/>	3274-7370	346 Cushion for back rest extension, black	(Upcharge 08)	1 pc
			* Use with back rest extension or head support if there is no support for shoulder and head mounted in the transport vehicle.		
with	<input type="checkbox"/>	3278-52SW	Push-bar	09	1 pc
or	<input type="checkbox"/>	3278-54SW	Push-bar "wide"	10	1 pc
or	<input type="checkbox"/>	3278-53SW	Push-handles	11	1 pc
or	<input type="checkbox"/>	3275-5502	One hand push-handle delta form (max. SW 36)	12	1 pc

Brakes, fender

	Order no.	Colour	Description	Qty	Price €
with	<input type="checkbox"/>	3275-1750	HOGGI light sport brake	13	1 set incl.
or	<input type="checkbox"/>	3221-1700	Standard brakes CLEO Ti	14	1 set incl.
with			Clothes cover		1 set incl.
or	<input type="checkbox"/>	3275-160x	Fender	15	1 set
or	<input type="checkbox"/>	3275-169x	Fender with integrated brake system	16	1 set
or	<input type="checkbox"/>	3276-160x	Fender with height adjustable and detachable arm rests	17	1 set
or	<input type="checkbox"/>	3276-169x	Fender with integrated brake system and with height adjustable and detachable arm rests	18	1 set
with	<input type="checkbox"/>	7000-0051	PU-Armrest pads, 235 mm long		1 pair incl.
or	<input type="checkbox"/>	7000-0052	PU-Armrest pads, 270 mm long		1 pair incl.
	<input type="checkbox"/>	3240-8202	Therapy tray size 2, SW 24 and 26 cm *	19	1 pc
	<input type="checkbox"/>	3240-8203	Therapy tray size 3, SW 28, 30 and 32 cm *		1 pc
	<input type="checkbox"/>	3240-8204	Therapy tray size 4, SW 34 *		1 pc
			* only in conjunction with arm rests		

Aluminum fender

	Order no.	Colour	Description	Qty	Price €
with	<input type="checkbox"/>	3277-1604	Wheel guards aluminum in frame color, 20-24"	20	1 set
or	<input type="checkbox"/>	3278-1604	Wheel guards aluminum in frame color, 20-24" with height adj. and detachable arm rests	21	1 set
or	<input type="checkbox"/>	3277-1694	Wheel guards aluminum in frame color 20-24" with integrated wheel lock levers	22	1 set
or	<input type="checkbox"/>	3278-1694	Wheel guards alu. frame color 20-24" with integrated brake system, height adjust. detachable arm rests	23	1 set

Rear wheels

	Order no.	Colour	Description	Qty	Price €
with	<input type="checkbox"/> 1296-0xxx		Rear wheel HOGGI light	1 set	incl.
or	<input type="checkbox"/> 1297-1xxx		Rear wheel HOGGI light with drum brakes* (Upcharge)	1 set	
			*in combination with push-bar or push-handles	1 set	
with	<input type="checkbox"/> 129x-xx1x		Pneumatic tires 1"	24 1 set	incl.
or	<input type="checkbox"/> 129x-xx2x		Flat free tires (20x1 3/8", 22x1" and 24x1")	25 1 set	incl.
with	<input type="checkbox"/> 129x-xxx0		Without push-rims	1 set	
or	<input type="checkbox"/> 129x-xxx1		Push rims aluminium anodised, "standard" mounting	26 1 set	incl.
or	<input type="checkbox"/> 129x-xxx2		Push rims aluminium anodised, "high" mounting	27 1 set	incl.
or	<input type="checkbox"/> 129x-xxx3		Push rims stainless steel, "standard" mounting	(Upcharge) 28 1 set	
or	<input type="checkbox"/> 129x-xxx4		Push rims stainless steel, "high" mounting	(Upcharge) 1 set	
with		<input type="checkbox"/>	Push-rim mounting "narrow"	1 set	
or		<input type="checkbox"/>	Push-rim mounting "wide"	1 set	
or	<input type="checkbox"/> 1204-1155		Closed ergonomic hand rims 22", black anodized, for Hoggi-light wheels	(Upcharge) 60 1 set	
or	<input type="checkbox"/> 1204-1156		Closed ergonomic hand rims 24", black anodized, for Hoggi-light wheels	(Upcharge) 1 set	
or	<input type="checkbox"/> 1204-106X		"HOGGI Grip Plus" hand rim inserts black, for closed ergonomic hand rims	(Upcharge) 61 1 set	
	<input type="checkbox"/> 1204-114X		Aluminum hand rim with plastic anti-slip coating, black, standard mounting *	(Upcharge) 58 1 pair	
	<input type="checkbox"/> 1204-114X		Aluminum hand rim with plastic anti-slip coating, black, high mounting *	(Upcharge) 1 pair	
with	<input type="checkbox"/> 1296-4211		Pneumatic wheel 20x1" HOGGI Extra light *	(Upcharge) 31 1 set	
with	<input type="checkbox"/> 1296-4311		Pneumatic wheel 22x1" HOGGI Extra light *	(Upcharge) 1 set	
or	<input type="checkbox"/> 1296-4411		Pneumatic wheel 24x1" HOGGI Extra light *	(Upcharge) 1 set	

* incl. black alu push-rim

Spoke guards // Old motifs // Sale

	Order no.	Colour	Description	Qty	Price €
	<input type="checkbox"/> 1204-101x		Spoke guard, transparent *	32 1 set	
	<input type="checkbox"/> 1204-101x		Spoke guard, transparent, for drum brakes wheels *	1 set	
	<input type="checkbox"/> 1204-101x		Spoke guard, printed **	(Sale) 33 1 set	
	<input type="checkbox"/> 1204-101x		Spoke guard, printed, for drum brakes **	(Sale) 1 set	

*** Motifs 2021 only possible in combination with 1204-2001-XXX**

**** please enter code of your print design see attachment or brochure "spoke guard" // Only while stocks last!**

***** printed spoke guards are not available for 20" wheels**

Spoke guards // New motifs 2021 (with 10% dealer discount)

	Order no.	Colour	Description	Qty	Price €
	<input type="checkbox"/> 1204-2001-000		Spoke guard Special interior foiling, single-coloured according to customer requirements ^{1 2}	1 set	
	<input type="checkbox"/> 1204-2001-2XX		Spoke guard Special interior foiling, printed ^{1 3}	1 set	
	<input type="checkbox"/> 1204-2001-SA		Spoke guard Special interior foiling, printed according to customer requirements ¹		Price on request

¹ **only in combination with transparent spoke guard**

² **please specify RAL tone**

³ **Spoke guard motif selectable in "Order form spoke-guards" (available on the website)**

Front castors

	Order no.	Colour	Description	Qty	Price €
with	<input type="checkbox"/> 1250-1002S		Wheel with rubber tire, black 100 x 32 mm	34 1 set	incl.
or	<input type="checkbox"/> 1250-1003S		Wheel with rubber tire, black 125 x 32 mm	1 set	incl.
or	<input type="checkbox"/> 1250-1005S		Wheel with rubber tire, black 140 x 40 mm	1 set	incl.
or	<input type="checkbox"/> 1249-1005		Skater wheel with LED 100 x 24 mm	(Upcharge) 35 1 set	
or	<input type="checkbox"/> 1249-1006		Skater wheel with LED 125 x 24 mm	(Upcharge) 1 set	
or	<input type="checkbox"/> 1221-1003		PU-wheel with plastic design rim 100 x 34 mm	(Upcharge) 36 1 set	
or	<input type="checkbox"/> 1221-1004		PU-wheel with plastic design rim 125 x 34 mm	(Upcharge) 1 set	

Accessories

	Order no.	Colour	Description	Qty	Price €
	<input type="checkbox"/> 3275-2801		Footrest hanger 90°, lower leg length 13-26 cm	37	1 set incl.
	<input type="checkbox"/> 3275-2802		Footrest hanger 90°, lower leg length 13-36 cm *		1 set incl.
	<input type="checkbox"/> 3278-2801		Footrest hanger adj. 90-155°, LLL 19 - 25cm	(Upcharge) 38	1 set
	<input type="checkbox"/> 3278-2802		Footrest hanger adj. 90-155°, LLL 19 - 35cm	(Upcharge)	1 set
	* applicable up from seat height 39 cm				
with	<input type="checkbox"/> 3275-2910		Footrest bracket angle adjustable, flip up type	39	1 set incl.
or	<input type="checkbox"/> 3275-2960		Footrest bracket angle adjustable, with interlocking	(Upcharge) 40	1 set
	<input type="checkbox"/> SA-0430		Reinforced steel footplate-carrier (about 0,3 kg heavier)		1 pc
with	<input type="checkbox"/> 3275-31SW		CLEO footrest ABS, flat version *	41	1 pc incl.
or	<input type="checkbox"/> 3275-30SW		CLEO footrest ABS, with heel stop *	(Upcharge) 42	1 pc
with	<input type="checkbox"/> 3275-32SW		CLEO footrest aluminium, flat version in frame colour *	(Upcharge) 43	1 pc
or	<input type="checkbox"/> 3275-33SW		CLEO footrest carbon, flat version *	(Upcharge) 44	1 pc
	* The standard footrest width is 4 cm narrower than the seat width and fits almost all settings. Wider footrests are available on request. This can lead restrictions to changes in the settings after delivery.				
	<input type="checkbox"/> 3275-1910		Bumper wheels	45	1 set
	<input type="checkbox"/> 3231-7800	346	Abduction block size 1, swing away, detachable, (9 x 6 x 9 cm)	46	1 pair
	<input type="checkbox"/> 3231-7900	346	Abduction block size 2, swing away, detachable, (11 x 8x 10 cm)		1 pair
	<input type="checkbox"/> 32XX-71SW	346	Seat cushion; 3 cm thick	47	1 pc
	<input type="checkbox"/> 3276-7100	346	Back cushion for "reclinable" back, trough shaped, 2,5 cm thick	48	1 pc
	<input type="checkbox"/> 3201-3830		Universal headrest bracket (3245-7800 or third-party product) *	49	1 pc
	<input type="checkbox"/> 3245-7800	346	Headrest with occiput upholstery	50	1 pc
	* starting from backrest height 30 cm				
	<input type="checkbox"/> 3200-0401		Belt routing as a T-strap above the shoulders	51	1 pc
	<input type="checkbox"/> 1204-101x		Push-rim cover, silicone black fluted		1 set
	<input type="checkbox"/> 3278-15SW		CLEO-Ti extension kit of 2 cm	52	1 set
	<input type="checkbox"/> 3201-6400		Ankle hugger, size 1	53	1 pair
	<input type="checkbox"/> 3203-6400		Ankle hugger, size 2		1 pair
	<input type="checkbox"/> 3274-6100		Lap belt, black	54	1 pc
	<input type="checkbox"/> 3275-6702		4-Point lap belt, size S (SW 24 - 30)	55	1 pc
	<input type="checkbox"/> 3275-6703		4-Point lap belt, size M (SW 28 - 36)		1 pc

Custom-made

	Order no.	Colour	Description	Qty	Price €
	<input type="checkbox"/> 3274-80SW	346	Contoured seat cushion with cover, black, SW 24 -32	56	1 pc
	<input type="checkbox"/> 3274-80SW	346	Contoured seat cushion with cover, black, SW 34		1 pc
	<input type="checkbox"/> 3275-7200	346	"HOGGI align" Contoured back cushion, black	57	1 pc
	<input type="checkbox"/> SA-0286/2-KP		Abduction block 30 mm height extension kit *		1 set
	* only in combination with abduction block and contoured seat cushion				
	<input type="checkbox"/> SA-0210-MP		Special Frame colour, RAL code only	(Upcharge)	1 pc
	<input type="checkbox"/> 6999-0014	060	Special Frame colour, Copper Gold	(Upcharge)	1 pc
	<input type="checkbox"/> 6999-0015		CLEO component parts in various colours according to customer requirements *		Per colour
	* please note colours of all component parts below				
	<input type="checkbox"/> 6999-0017		Powder coated castor forks and both adapter for camber colour according to customer requirements *		1 set
	<input type="checkbox"/> SA-0279/4-KP		CLEO outdoor front end *	58	1 pc
	* additional costs for retrofitting				
	<input type="checkbox"/> 6999-1026		Frame prepared for Outdoor front end		1 pc.
	<input type="checkbox"/> SA-0390		Upholstery in special colour *		
	<input type="checkbox"/> SA-0572-F-MP		Plastic anti-slip coating in custom colour: Pure white (RAL 9010), Silver grey (RAL 7001), Sky blue (RAL 5015), Traffic red (RAL 3020), Pure orange (RAL 2004)*	(Upcharge) 59	1 pair
	* only in combination with HOGGI light rear wheels				
	<input type="checkbox"/> SA-0093-MP		Maxgrepp Original cell rubber coated push rims	(Upcharge) 62	1 set
	<input type="checkbox"/> SA-0094-MP		Maxgrepp Ergo Para (closed push rims)	(Upcharge) 63	1 set
	<input type="checkbox"/> SA-0095-MP		Maxgrepp Ergo Para CP (advanced closed push rims)	(Upcharge) 64	1 set
	<input type="checkbox"/> SA-0384-MP		CLEO frame with parallel abduction	(Upcharge) 65	1 set
	<input type="checkbox"/> SA-0394-MP		Coupling of two CLEOs one behind the other	(Upcharge) 66	1 set
	<input type="checkbox"/> SA-0410-Min		Preparation for electric drive "TWION" from Alber	(Upcharge)	
	<input type="checkbox"/> SA-0410/3-MP		CLEO TI preparation for "e-fix"	(Upcharge) 67	1 pc
	<input type="checkbox"/> SA-0427-MP		CLEO One-hand push bar with ratchet joint	(Upcharge) 68	1 pc

Custom-made

Order no.	Colour	Description	Qty	Price €
<input type="checkbox"/> 3246-68SW		"CLEO Standard" trunk supports (Upcharge) 69	1 set	
<input type="checkbox"/> SA-0370/1-MP		Mounting adaption for CLEO angle and height adjustable trunk supports (Upcharge)	1 set	
<input checked="" type="checkbox"/> 3007-770X		CLEO angle and height adjustable trunk supports (Upcharge) 70	1 set	
	<input type="checkbox"/>	size 1 (8 x 10 x 2 cm)		
	<input type="checkbox"/>	size 2 (10 x 12 x 2 cm)		
	<input type="checkbox"/>	size 3 (12 x 13 x 2 cm)		
<input type="checkbox"/> SA-0370-MP		Mounting adaptation for Stealth "Nino" trunk supports (Upcharge) 71		
<input type="checkbox"/> SA-0505-MP		Locking of the integrated brake (Upcharge) 72	1 pc	
<input type="checkbox"/> SA-0594		Carbon fender (Upcharge) 73		

Custom-made Dynamic back base (not released for transport)


Order no.	Colour	Description	Qty.	Price €
<input type="checkbox"/> SA-0552/2-MP		Dynamic package to CLEO Ti from back height 30 74	1 set	
<input checked="" type="checkbox"/> SA-0431-MP		CLEO Back panel reinforced version 75	1 set	
<input checked="" type="checkbox"/> SA-0555		CLEO Ti seat adapter reinforced steel version 75	1 pair	
		total		

Extra equipment (with 10% dealer discount)

Order no.	Colour	Description	Qty.	Price €
<input type="checkbox"/> 3275-3100		CLEO trailer coupling for a "BOBBY-CAR" trailer	1 pc	
<input type="checkbox"/> SA-0366		Children name as film plot (23 x 95 mm) for sticking onto the back panel, Font: Calibri, film plot colour black or white	1 pc	
<input type="checkbox"/> 7000-0199		"Cosies" Thermofleece hand warmer for companion 76	1 set	

Custom made 6999-0015 component parts in various colours

Position	Colour	Component parts
<input type="checkbox"/> A		Back panel
<input type="checkbox"/> B		Back base
<input type="checkbox"/> C		Sideplates, aluminium
<input type="checkbox"/> D		Seat panel
<input type="checkbox"/> E		Frame
<input type="checkbox"/> F		Footrest hanger
<input type="checkbox"/> G		Footrest, aluminium



Upholstery codes

346 Black	321 Spacer fabric, black (special colour)	348 Grey Soft (special colour)
	347 Berry Soft (special colour)	349 Turquoise Soft (special colour)


All HOGGI upholstery are made from flame retardant Trevira CS. Exception spacer fabric **(321)**.

Frame colours

001 Black (RAL 9005)	020 Yellow green (RAL 6018)	028 Pearl violet (RAL 4011)	056 Turquoise blue (RAL 5018)
002 Silver	021 Grey aluminium (RAL 9007)	029 Pearl beige (RAL 1035)	057 Iron sparkling dark light
003 Pearl ruby red (RAL 3032)	022 White (RAL 9003)	032 Black dull matt (RAL 9005M)	058 Salmon orange (RAL 2012)
004 Pearl night blue (RAL 5026)	023 Pure orange (RAL 2004)	039 Signal violet (RAL 4008)	059 Pastel blue (RAL 5024)
005 Pearl gentian blue (RAL 5025)	024 Sparkling iron effect	040 Night blue (RAL 5022)	061 Light pink (RAL 3015)
017 Traffic yellow (RAL 1023)	025 Signal blue (RAL 5005)	050 White matt (RAL 9003)	064 Olive green (RAL 6003)
018 Flame red (RAL 3000)	026 Vermillion (RAL 2002)	053 Light green (RAL 6027)	065 Raspberry red (RAL 3027)
019 Sky blue (RAL 5015)	027 Heather violet (RAL 4003)	054 Ocean blue (RAL 5020)	066 Pastel green (RAL 6019)

SA-0210-MP Special colours on request (Enter RAL code Upcharge

Measuring and active degree



SH = Seat height LLL = Lower leg length SD = Seat depth BH = Backrest height
(All measures without seat and back cushions)

Seat width (SW)	24	26	28	30	32	34
Seat depth standard (SD)	22-30	24-32	26-34	28-36	30-38	32-40
Seat depth long (SD)	28-36	30-38	32-40	34-42	36-44	38-46

SH Frame size 1	38 cm	39 cm	40 cm	41 cm
SH Frame size 2	40,5 cm	41,5 cm	42,5 cm	43,5 cm
SH Frame size 3	45 cm	46 cm	47 cm	48 cm


CLEO 45° Seat tilt




Colour code

201	HOGGI blue	205	Train
202	Wonderland	206	Dragon
203	Mouse	207	Circus
204	Water world		

CONFIGURATIONS // ACCESSORIES // CLEO Ti



Tilt in space



Base plates, aluminium



Base plates „Comfort“ long wheelbase




Anti tip




Back base angle adjustable



Tie down kit



Back rest extension



Cushion for back rest extension



Push-bar



Push bar „wide“



Push-handles



One hand push-handle delta form



HOGGI light sport brake



Standard brakes



Fender




Fender with integrated brake system




Fender with arm rests



Fender with arm rests and integr. brake system



Therapy tray



Aluminium fender



Aluminium fender with arm rests



Aluminium fender with brake and arm rests



Aluminium fender with integr. brake system



Pneumatic tires



Pneumatic tires and drum brake



Polyurethan tires



Polyurethan tires and drum brake



Push rims aluminium anodised, standard



Push rims aluminium anodised, high




Push rims stainless steel, standard

CONFIGURATIONS // ACCESSORIES // CLEO Ti

<p>29</p> 	<p>30</p> 	<p>31</p> 	<p>32</p> 	<p>33</p> 
<p>Push rims stainless steel, high</p>	<p>HOGGI Extra light</p>	<p>Spoke guard, transparent</p>	<p>Spoke guard, printed</p>	<p>Wheel with rubber tire</p>
<p>34</p> 	<p>35</p> 	<p>36</p> 	<p>37</p> 	<p>38</p> 
<p>Skater wheel with LED</p>	<p>PU-wheel with plastic design rim</p>	<p>Footrest hanger 90°</p>	<p>Adjustable footrest hanger 90-155°</p>	<p>Footrest bracket</p>
<p>39</p> 	<p>40</p> 	<p>41</p> 	<p>42</p> 	<p>43</p> 
<p>Footrest bracket with interlocking</p>	<p>ABS, flat version</p>	<p>ABS with heel stop</p>	<p>Footrest: Aluminium in frame colour</p>	<p>Footrest: carbon</p>
<p>44</p> 	<p>45</p> 	<p>46</p> 	<p>47</p> 	<p>48</p> 
<p>Bumper wheels</p>	<p>Abduction block</p>	<p>Seat cushion</p>	<p>Back cushion</p>	<p>Universal headrest bracket</p>
<p>49</p> 	<p>50</p> 	<p>51</p> 	<p>52</p> 	<p>53</p> 
<p>Headrest with occiput upholstery</p>	<p>Belt routing as a T-strap</p>	<p>Seat width enlarging set</p>	<p>Ankle hugger</p>	<p>Lap belt</p>
<p>54</p> 	<p>55</p> 	<p>56</p> 	<p>57</p> 	<p>58</p> 
<p>4-Point lap belt</p>	<p>Contoured seat cushion with cover</p>	<p>„HOGGI align“ Contoured back cushion</p>	<p>Outdoor front end</p>	<p>Aluminum hand rim with plastic anti-slip coating</p>

CONFIGURATIONS // ACCESSORIES // CLEO Ti



59 Aluminum hand rim with plastic anti-slip coating



60 Closed ergonomic hand rims



61 HOGGI Grip Plus




62 Maxgrepp Original push rims




63 Maxgrepp Ergo Para push rims



64 Maxgrepp Ergo Para CP push rims



65 CLEO frame with parallel abduction




66 Coupling of two CLEOs one behind the other




67 Electric drive „e-fix“




68 One-hand push bar with ratchet joint



69 Standard trunk supports



70 Angle and height adjustable trunk supports



71 Stealth „Nino“ trunk supports



72 Locking of the integrated brake



73 Carbon fender



74 Dynamic package CLEO Ti



75 Back panel & seat adapter reinforced version



76 „Cosies“ Thermofleece hand warmer