

SHERPA

Release date: September 2025 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address

see invoice address

Company: _____

Street: _____

City, Zip code: _____



- Basic configuration
- Aluminium rigid frame
 - Push-bar or Push-handles
 - Knee lever wheel locks
 - System-Profile-Seat- & Back frame
 - HOGGI light rear wheels inclusive push rims

Frame colour

Black	Silver	Pearl ruby	Pearl night blue	Pearl gentian blue	Yellow green
Sparkling iron effect	Traffic yellow	Flame red	Sky blue	Pearl violet	Night blue
Grey aluminium	White	Pure orange	Signal blue	Pearl beige	White matt
Vermillion	Heather violet	Black dull matt	Signal violet	Light green	Ocean blue
Turquoise blue	Iron sparkling light	Salmon orange	Pastel blue	Light pink	Olive green
Raspberry red	Pastel green				


Special frame colour



Desired RAL colour:


Copper gold

Seat width and Back rest width



- Seat width: 36 cm
Back rest width: 34,8 cm
- Seat width: 40 cm
Back rest width: 38,8 cm
- Seat width: 44 cm
Back rest width: 42,8 cm
- Seat width: 48 cm
Back rest width: 46,8 cm

Seat height




Seat height	Front wheel diameter				
Rear wheel diameter	125		140		190
12"		48		48	48
16"	45	48	45	48	48
22"		45	45	50	50
24"		48		48	48

Camber (only for 22" & 24" wheels)



Seat-Profile length



Seat-Profile length


40 cm (approx. seat depth 29-39 cm)

45 cm (approx. seat depth 34-44 cm)


50 cm (approx. seat depth 39-49 cm)

55 cm (approx. seat depth 44-54 cm)


Universal-Back frame



Universal-Back frame
with physiological pivot point



Universal-Back frame for seat shells



Back guide with
length compensation

cross tube depth 130 mm

cross tube depth 55 mm

cross tube depth 130 mm

cross tube depth 55 mm



Notice:

Cannot be combined with armrests and a back panel.



Push-handles



Push-Bar, angle adjustable




Release of the backrest
angle adjustment




Notice:

Release via Bowden cable and
release lever on the push bar or
push handle.

System-Profile-Back frame



System-Profile-Back frame with physiological pivot point



System-Profile-Back frame for seat shells



Back-Profile length

40 cm (approx. BH RH 47-51 cm)

45 cm (approx. BH RH 52-56 cm)

50 cm (approx. BH 57-61 cm)



Back guide with length compensation



Back panel height adjustable, in frame colour



Back cushion 3 cm thick, black



Notice:

Cannot be combined with the back panel.



Notice:

Adjustment range from 42 cm to max. 10 cm above the length of the back profile possible.



Push-handles



Push-Bar, angle adjustable



Release of the backrest angle adjustment




Notice:

Release via Bowden cable and release lever on the push bar or push handle.


Seat panel



Seat panel



Seat height raiser,
adjustable to 2, 4 or 6 cm



Precision trapezoid adapter

size 1 (for SW 36) 175,00€

size 2 (for SW 40 & 48) 185,00€



Notice:

The minimum seat profile length must be observed.

SW 36 = 45 cm; SW 40 = 50 cm;

SW 44 + SW 48 = 55 cm

Seat cushion



Seat cushion 3 cm thick, black



Contoured seat cushion with cover



Notice:

Cannot be combined with hip supports.

Upholstery colours




Black




Spacer fabric, black



Berry Soft



Turquoise Soft



Grey Soft



Notice:

Black upholstery is the standard colour and is included in the basic configuration.

Head supports



Console with Headrest bracket
„universal“



Head support with occiput
upholstery incl. mounting bracket




Console for headrest bracket
(as for STEALTH headrest bracket)



Notice:

For head support with occiput
upholstery or for third-party
products.

Armrests



Armrests pivotably, angle adjustable
with 270 mm PU Pads



Pillar armrest with PU-pads,
removable, height adjustable



PU-armrest pads



Notice:

Cannot be combined with
the universal-back frame.

height: 22 - 31 cm

height: 27 - 36 cm

height: 38 - 48 cm



Notice:

Therapy tray only possible
in conjunction with special
construction SA-0586.

23,5 cm long

27 cm long

Brakes



Standard wheel locks



HOGGI light sport wheel locks




Integrated brake




Notice:

Only possible for 22" and 24" wheels.


Wheel guards



Plastic wheel guards



Plastic wheel guards with height adjustable and removable armrests



Wheel guards aluminum in frame colour



Wheel guards aluminum in frame colour with height adjustable and detachable arm rests



Vario wheel cover made of plastic



Vario wheel cover made of plastic with height adjustable and removable armrests



Vario wheel cover made of aluminium in frame colour



Carbon wheel cover with integrated brake



Notice:

Extension of the installation width to +6 or +7 cm possible.

Rear wheels

Select tires:

Pneumatic tires

Flat free tires

Select rear wheels:

Standard rear wheels

**Standard rear wheels
with drum brakes**

Select wheel diameter:



12" wheel diameter



16" wheel diameter



22" wheel diameter
with aluminium hand rim



24" wheel diameter
with aluminium hand rim



Push rim cover,
silicone black fluted

!
Notice:
Available for 22" and 24"
rear wheels.

Spoke guards

Spoke guard, transparent

Spoke guard Special interior foiling,
with motif

Spoke guard Special interior foiling,
**single-coloured according
to customer requirements:**

Spoke guard Special interior foiling,
printed according to customer requirements



Front castors



Wheel with rubber tire, black
125 x 32 mm with Aluminum
wheel fork



Wheel with rubber tire, black
140 x 40 mm with Aluminum
wheel fork



Wheel with PU tire
190 x 50 mm with Aluminum
wheel fork



Wheel with pneumatic tire
190 x 50 mm with Aluminum
wheel fork

Anti tip, Tip assist, Calf upholstery



Anti tipper, swing away

assembly left


assembly right



Anti tipper, swing away,
with tip assist

assembly left

assembly right



Calf upholstery

according to seat width

for single foot rests



Notice:

Cannot be combined with
Vario footrest system.

Foot rest hanger, Single foot rest hanger



Foot rest hanger 90°

LLL short 20 - 36 cm


LLL long 20 - 41 cm



Foot rest hanger angle adjustable
85° - 160°

LLL short 20 - 36 cm

LLL long 20 - 41 cm



Single foot rest hanger, separat
angle adjustable 85° - 160°


LLL short 20 - 36 cm

LLL long 20 - 41 cm

Single footrest, Knee angle retainers



Up and swing away single footrest
(with calf plates and heel strap)



Knee angle retainers standard
for 15 mm square tube

Foot rest bracket




Foot rest bracket angle adjustable,
flip up type



Foot rest bracket angle adjustable,
with interlocking


Foot rest




Flat plastic foot rest
SW 28 - 40 cm



Plastic foot rest with
heel stop, SW 28 - 40 cm



Alu Foot rest with heel stop,
SW 40 - 44 cm



Two single footplates
in boat shape

Vario Foot rest system

Vario Foot rest system 80°, anatomically angled abductable by 6 cm
Swivel frame in frame colour, can be swung away and removed



Plastic footboard standard, fold-up,
 LLL 33 - 44 cm

**Footb. surface W 125 x D 150 mm,
 total foot space 26 - 32 cm, (SW 36)**

**Footb. surface W 150 x D 150 mm,
 total foot space 31 - 37 cm, (SW 40)**

**Footb. surface W 150 x D 150 mm,
 total foot space 35 - 41 cm, (SW 44)**



Plastic footboard set angle-
 adjustable, fold-up, LLL 28 - 39 cm

**Footb. surface W 110 x D 160 mm,
 total foot space 23 - 29 cm, (SW 36)**

**Footb. surface W 120 x D 160 mm,
 total foot space 25 - 31 cm, (SW 40)**

**Footb. surface W 130 x D 160 mm,
 total foot space 27 - 33 cm, (SW 40)**

**Footb. surface W 145 x D 160 mm,
 total foot space 29 - 35 cm, (SW 44)**

**Footb. surface W 160 x D 160 mm,
 total foot space 31 - 37 cm, (SW 44)**

**Footb. surface W 160 x D 160 mm,
 total foot space 35 - 41 cm, (SW 48)**



Plastic footboard set angle-
 adjustable, fold-up, LLL 34 - 47 cm

**Footb. surface W 110 x D 160 mm,
 total foot space 23 - 29 cm, (SW 36)**

**Footb. surface W 120 x D 160 mm,
 total foot space 25 - 31 cm, (SW 40)**

**Footb. surface W 130 x D 160 mm,
 total foot space 27 - 33 cm, (SW 40)**

**Footb. surface W 145 x D 160 mm,
 total foot space 29 - 35 cm, (SW 44)**

**Footb. surface W 160 x D 160 mm,
 total foot space 31 - 37 cm, (SW 44)**

**Footb. surface W 160 x D 160 mm,
 total foot space 35 - 41 cm, (SW 48)**



Alu-footboard one piece with heel
 stop, angle-adjustable, fold-up,
 LLL 26 - 37 cm

**Footb. surface W 30-36 x D 18 cm,
 total foot space 30 - 36 cm, (SW 36)**

**Footb. surface W 34-38 x D 18 cm,
 total foot space 34 - 40 cm, (SW 40)**

**Footb. surface W 38-44 x D 18 cm,
 total foot space 38 - 44 cm, (SW 44)**

**Footb. surface W 42-48 x D 18 cm,
 total foot space 42 - 48 cm, (SW 48)**



Alu-footboard one piece with heel
 stop, angle-adjustable, fold-up,
 LLL 32 - 45 cm

**Footb. surface W 30-36 x D 18 cm,
 total foot space 30 - 36 cm, (SW 36)**

**Footb. surface W 34-38 x D 18 cm,
 total foot space 34 - 40 cm, (SW 40)**

**Footb. surface W 38-44 x D 18 cm,
 total foot space 38 - 44 cm, (SW 44)**

**Footb. surface W 42-48 x D 18 cm,
 total foot space 42 - 48 cm, (SW 48)**

Trunk supports



Trunk support, without upholstery

size 1 (8 x 10 x 2 cm)

size 2 (10 x 12 x 2 cm)

size 3 (12 x 13 x 2 cm)

Trunk support, with upholstery,
black

size 1 (8 x 10 x 2 cm)

size 2 (10 x 12 x 2 cm)

size 3 (12 x 13 x 2 cm)




Swing away trunk support

85 x 130 mm


100 x 165 mm

MDP trunk supports (can be ordered separately)

MDP pad holder swing away
and detachable

left

right

MDP pad holder fixed
left / right

MDP trunk support

8 x 10 x 2 cm

10 x 12 x 2 cm

12 x 13 x 2 cm



MDP trunk support curved

6 x 11,5 x 2,5 cm

8 x 14 x 2,5 cm

as a pair

as a pair

as a pair



Notice:

Only possible in combination with the universal back frame.

Thigh Supports



Thigh support
on the seat profile tubes

size 1 (290 x 140 mm)

size 2 (340 x 160 mm)

size 3 (390 x 180 mm)

**Part set for two customer
thigh support (without
pad)**



Thigh Support for pillar armrests

size 1 (290 x 140 mm)

size 2 (340 x 160 mm)

size 3 (390 x 180 mm)

**Part set for two customer
thigh support on pillar
armrests (without pad)**



Upholstery for Thigh Support, black




Notice:

Can be removed with the
pillar armrests.

Lap belt




Lap belt with Lap belt retainer



Lap belt retainer
(for third-party products)


Abduction block, Therapy tray



Abduction block,
swing away, detachable

LWH 11 x 8 x 10 cm

LWH 13 x 9 x 10 cm




Abduction block 60 mm height
extension kit



Notice:

Only in combination with abduction block and contoured seat cushion.



Therapy tray



Notice:

Only in conjunction with arm rests.

Tie down kit, Container for accessories




Tie down kit




Container for accessories


System Profile, Cap for system profile, Nut insert



System Profile 4N 30x30,
black anodized, length according
to customer specification



Cap for system profile, 30x30



Nut insert for system profile, M6

Outdoor front end, Cosies



Outdoor front end



„Cosies“ Thermofleece hand warmer
for companion

Remark / Customer request / Individual special custom-made request