

# SWINGBO® 2

Release date: September 2025 // valid until revoked

## Order

### Offer

Order number / reference: \_\_\_\_\_

Date: \_\_\_\_\_

## Invoice address

Customer number: \_\_\_\_\_

Company: \_\_\_\_\_

Street: \_\_\_\_\_

City, Zip code: \_\_\_\_\_

## Delivery address

see invoice address

Company: \_\_\_\_\_

Street: \_\_\_\_\_

City, Zip code: \_\_\_\_\_



## Basic configuration

- Aluminium rigid frame
- Knee lever wheel locks
- Push-bar or Push-handles
- HOGGI light rear wheels inclusive push rims

## Frame colour

|                       |                      |                 |                  |                    |              |
|-----------------------|----------------------|-----------------|------------------|--------------------|--------------|
|                       |                      |                 |                  |                    |              |
| Black                 | Silver               | Pearl ruby      | Pearl night blue | Pearl gentian blue | Yellow green |
|                       |                      |                 |                  |                    |              |
| Sparkling iron effect | Traffic yellow       | Flame red       | Sky blue         | Pearl violet       | Night blue   |
|                       |                      |                 |                  |                    |              |
| Grey aluminium        | White                | Pure orange     | Signal blue      | Pearl beige        | White matt   |
|                       |                      |                 |                  |                    |              |
| Vermillion            | Heather violet       | Black dull matt | Signal violet    | Light green        | Ocean blue   |
|                       |                      |                 |                  |                    |              |
| Turquoise blue        | Iron sparkling light | Salmon orange   | Pastel blue      | Light pink         | Olive green  |
|                       |                      |                 |                  |                    |              |
| Raspberry red         | Pastel green         |                 |                  |                    |              |

## Special frame colour



Copper gold

Desired RAL colour:

## Special paint finish (optional)




Powder coated castor forks  
in desired colour:  
(standard colour black)


**Notice:**

Please specify HOGGI frame colour.

## Seat width and seat depth



Seat width: 24 cm (actual mounting width: 23,5 cm)  
enter seat depth: (26 - 30 cm)

Seat width: 26 cm (actual mounting width: 25,5 cm)  
enter seat depth: (26 - 32 cm)

Seat width: 28 cm (actual mounting width: 27,5 cm)  
enter seat depth: (26 - 34 cm)


Seat width: 30 cm (actual mounting width: 29,5 cm)  
enter seat depth: (28 - 36 cm)

Seat width: 32 cm (actual mounting width: 31,5 cm)  
enter seat depth: (30 - 38 cm)

Seat width: 34 cm (actual mounting width: 33,5 cm)  
enter seat depth: (32 - 40 cm)

Seat width: 36 cm (actual mounting width: 35,5 cm)  
enter seat depth: (34 - 42 cm)

## Back height and Lower leg length





| Seat width                                      | 24 cm    | 26 cm    | 28 cm    | 30 cm    | 32 cm    | 34 cm    | 36 cm    |
|-------------------------------------------------|----------|----------|----------|----------|----------|----------|----------|
| <b>Backrest height</b><br>back trough<br>shaped | 36-44 cm | 37-45 cm | 38-48 cm | 39-49 cm | 40-50 cm | 42-54 cm | 43-55 cm |
| back panel<br>angle adjustable                  | 33-42 cm | 35-44 cm | 37-46 cm | 39-48 cm | 41-50 cm | 43-52 cm | 45-54 cm |

enter back height:

enter lower leg length: (16 - 37 cm)

## Active degree



enter active degree: (0 - 4 in 2,5 cm steps)

## Camber



## Seat height

| SWINGBO 2              | Front castors with<br>100 mm wheel diameter                                                                                                                              | Front castors with<br>125 mm wheel diameter                                                                                                                              | Front castors with<br>140 mm wheel diameter                                                                                                                              |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>20" rear wheels</b> | Camber: 3° 6° 9°<br>SH: 36 cm - x x<br>SH: 37 cm x - x<br>SH: 38 cm x x x<br>SH: 39 cm x x -<br>SH: 40 cm - x x<br>SH: 41 cm x - x<br>SH: 42 cm x x x<br>SH: 43 cm x - - | Camber: 3° 6° 9°<br>SH: 36 cm - x -<br>SH: 37 cm x - x<br>SH: 38 cm x x x<br>SH: 39 cm x x -<br>SH: 40 cm - x -<br>SH: 41 cm x x x<br>SH: 42 cm x - x<br>SH: 43 cm x - - | Camber: 3° 6° 9°<br>SH: 37 cm - x -<br>SH: 38 cm x x x<br>SH: 39 cm x - x<br>SH: 40 cm x - x<br>SH: 41 cm x x -<br>SH: 42 cm - x x<br>SH: 43 cm x - -                    |
|                        | Camber: 3° 6° 9°<br>SH: 38 cm - x x<br>SH: 39 cm - x x<br>SH: 40 cm x - x<br>SH: 41 cm x - x<br>SH: 42 cm x x -<br>SH: 43 cm x x x<br>SH: 44 cm - x -<br>SH: 45 cm x - - | Camber: 3° 6° 9°<br>SH: 38 cm - x x<br>SH: 39 cm - x x<br>SH: 40 cm x - x<br>SH: 41 cm x x x<br>SH: 42 cm x x -<br>SH: 43 cm - x x<br>SH: 44 cm - x -<br>SH: 45 cm x - - | Camber: 3° 6° 9°<br>SH: 38 cm - x x<br>SH: 39 cm x x -<br>SH: 40 cm - x x<br>SH: 41 cm - x x<br>SH: 42 cm x - x<br>SH: 43 cm x x x<br>SH: 44 cm x x -                    |
|                        | Camber: 3° 6° 9°<br>SH: 40 cm - - x<br>SH: 41 cm - x x<br>SH: 42 cm x x -<br>SH: 43 cm x x x<br>SH: 44 cm - x x<br>SH: 45 cm x x x<br>SH: 46 cm x - -                    | Camber: 3° 6° 9°<br>SH: 40 cm - - x<br>SH: 41 cm - x -<br>SH: 42 cm x x x<br>SH: 43 cm - x x<br>SH: 44 cm - - x<br>SH: 45 cm x x x<br>SH: 46 cm x - x<br>SH: 47 cm x x - | Camber: 3° 6° 9°<br>SH: 40 cm - - x<br>SH: 41 cm - x x<br>SH: 42 cm x - x<br>SH: 43 cm x x x<br>SH: 44 cm x x -<br>SH: 45 cm - x x<br>SH: 46 cm - x x<br>SH: 47 cm x x - |

x possible; - not possible

## Push handle / Push bar- and backrest variants


### Angle adjustable and foldable version:



Push bar, foldable and angle adjustable



Push bar „wide“, foldable and angle adjustable



Top view // reduction of seat width by 1,5 cm

Select the backrest angle adjustment range:

**75 - 115°**

**80 - 120°**

Extension push-bar

**5 cm extension**

**10 cm extension**

### Back, trough shaped:



With fastening material for „HOGGI align“ back cushions



Backrest and push bar folded down



Back angle adjustment

## Push handle / Push bar- and backrest variants

### Foldable variants:



Push handle, foldable



Push bar, foldable



Push bar „wide“, foldable

Select backrest angle:

**90°**                      **100°**

**95°**                      **105°**

Extension push-bar

**5 cm extension**

**10 cm extension**

### Back, trough shaped:



Backrest and push bar folded down



Reinforced locking claws

## Push handle / Push bar- and backrest variants

### Detachable variants:



Push handles, detachable



Push bar, detachable



Push bar „wide“, detachable

Extension push-bar

5 cm extension

10 cm extension

### Back panel angle- and height adjustable:




75 - 115°

80 - 120°



Back folded down



Top view // reduction of  
seat width by 1,5 cm

## Back cushion

Back cushion 2,5 cm

Through shaped back  
black

Back cushion, black  
SW 24 - 32 cm

Back cushion, black  
SW 34 - 36 cm



„HOGGI align“ Contoured  
back cushion, black

## Seat panel & seat cushion



Seat panel



Seat cushion 3 cm



Contoured seat cushion,  
with abduction

**Notice:**  
Cannot be  
combined  
with hip  
supports.

## Upholstery colours




Black




Spacer fabric, black



Berry Soft



Turquoise Soft



Grey Soft

## Headrests




Head support standard  
with upholstery



Headrest bracket universal for  
„Head support standard with  
upholstery“ or for thirdparty  
products




Head support with occiput  
upholstery incl. mounting  
bracket




Belt routing as a T-strap  
above the shoulders




## Armrests



Armrests

**height: 14 - 21 cm****height: 19 - 26 cm**

Armrests with tilt balancing

**height: 15 - 23 cm****height: 18 - 26 cm**

PU-armrest pads


**23,5 cm long****27 cm long****Notice:**

Depending on the configuration of the wheelchair, seat tilt may be restricted.

## Brakes



Knee lever wheel locks



Integrated brake

## Wheel guards



Plastic wheel guards



Carbon wheel cover

## Rear wheels



Standard rear wheel with aluminium hand rim



Standard rear wheel with drum brake and aluminium hand rim



Rear wheel HOGGI light



Rear wheel with drum brake

Select tires:

**Pneumatic tires**  
(1")

**Flat free tires**  
(20" x 1 3/8", 22" x 1" and 24" x 1")



**Push rims** (only selectable in combination with HOGGI light rear wheels)

Without push rims


Push rims aluminium

Push rims aluminium with **anti-slip coating**

Push rims stainless steel

**Anti-slip coating:**

Black



Pure white



Silver grey



Sky blue



Traffic red



Pure orange

**Push rim variant:**

standard



high

**Push rim mounting:**

wide

narrow

Push rim cover,  
silicone black fluted



**Push rims** (only selectable in combination with HOGGI light rear wheels)



Closed ergonomic hand rims



„HOGGI Grip Plus“ hand rim inserts  
black, for closed ergonomic  
hand rims



Maxgrepp Original



Maxgrepp Ergo Para



Maxgrepp Ergo Para CP

**Spoke guards**

Spoke guard, transparent

Spoke guard Special interior foiling,  
**with motif**

Spoke guard Special interior foiling,  
**single-coloured according  
to customer requirements:**

Spoke guard Special interior foiling,  
**printed according to customer requirements**




## Front castors



Wheel with rubber tire

**100 x 32 mm****125 x 32 mm****140 x 40 mm**

Skater wheel with LED

**100 x 24 mm****125 x 24 mm**

PU-wheel with plastic design rim

**100 x 34 mm****125 x 34 mm**


## Anti tip, Tip assist, Calf upholstery



Anti tip

**assembly left****assembly right**

Tip assist

**assembly left****assembly right**Calf upholstery  
suitable to seat width

## Foot rest hanger



Foot rest hanger 90°

**LLL short 16 - 27 cm****LLL long 16 - 37 cm**Foot rest hanger angle adjustable  
85° - 160°**LLL short 19 - 27 cm****LLL long 19 - 37 cm**

Single foot rest hanger 85° - 160°

**LLL short 19 - 27 cm****LLL long 19 - 37 cm**

## Foot rest bracket




Foot rest bracket angle adjustable, flip up type



Foot rest bracket angle adjustable, with interlocking

Reinforced **steel** footplate-carrier (about 0,3 kg heavier)

## Foot rest



Flat plastic foot rest



Plastic foot rest with heel stop



Foot rest aluminum in frame color



Two single footplates in boat shape

## Supports



Thoracic support, particular adjustable, size approx. 12 x 13 x 2 cm

**for back panel angle- and height adjustable**

**for back panel trough shaped**




Chest-/ thoracic support

**chest size 33 - 43 cm**

**chest size 39 - 49 cm**

**chest size 45 - 55 cm**



Straight hip supports

**22,5 x 10 cm, SW 28 - 34**

**25 x 11 cm, SW 28 - 34**

**25 x 11 cm, SW 36 - 44**



### Notice:

Only in conjunction with backrest angle and height adjustment. In conjunction with a back panel trough shaped, the mounting bracket angle must also be adjusted.



### Notice:

Hip pads reduce the maximum seat width by 3 cm.

**MDP trunk supports** (can be ordered separately)




MDP pad holder swing away  
and detachable

**left**

**right**



MDP pad holder fixed  
left / right




MDP trunk support

**8 x 10 x 2 cm**

**10 x 12 x 2 cm**

**12 x 13 x 2 cm**



MDP trunk support curved

**6 x 11,5 x 2,5 cm**

**8 x 14 x 2,5 cm**

**as a pair**

**as a pair**

**as a pair**

**Belt systems**



Lap belt



4-Point lap belt



Ankle hugger


**SW 22 - 30 cm**

**SW 28 - 34 cm**

**Ø 4 - 6,5 cm**

**Ø 6 - 8,5 cm**


## Abduction block, Therapy tray



Abduction block,  
swing away, detachable

**LWH 9 x 6 x 9 cm**

**LWH 11 x 8 x 10 cm**




Abduction block 30 mm height  
extension kit



**Notice:**

Only in combination with abduction block and contoured seat cushion.



Therapy tray

**SW 24 and 26 cm**

**SW 28 to 32 cm**

**SW 34 and 36 cm**



**Notice:**

Only in conjunction with arm rests.

## Tie down kit



Tie down kit



with reinforced locking claws




**Notice:**

Required for push bars/push handles that can be folded down when used in conjunction with tie down kit.

## Gripper clamps, Cosies



Four gripper clamps D25  
with fasteners



„Cosies“ Thermofleece hand warmer  
for companion



## Canopy, One hand drive



Canopy incl. adapter  
and rain protection



Reconstruction for ProAktiv One  
hand drive

**Notice:**

In combination with the  
trough-shaped back panel  
(with low back height),  
restrictions may occur.

## Outdoor front end, Sibling seat



Outdoor front end



Outdoor front end  
as device carrier variant



Mounting adapter for Outdoor front  
end aid for „phil & teds“ sibling seat  
or universal car seat adapter

## Remark / Customer request / Individual special custom-made request