

SWINGBO® 2 XL

Release date: September 2025 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address

see invoice address

Company: _____

Street: _____

City, Zip code: _____



Basic configuration

- Aluminium rigid frame
- Knee lever wheel locks
- Double anti tippers
- Push-bar or Push-handles
- Standard rear wheel with aluminium hand rim

Frame colour

Black	Silver	Pearl ruby	Pearl night blue	Pearl gentian blue	Yellow green
Sparkling iron effect	Traffic yellow	Flame red	Sky blue	Pearl violet	Night blue
Grey aluminium	White	Pure orange	Signal blue	Pearl beige	White matt
Vermillion	Heather violet	Black dull matt	Signal violet	Light green	Ocean blue
Turquoise blue	Iron sparkling light	Salmon orange	Pastel blue	Light pink	Olive green
Raspberry red	Pastel green				

enter lower leg length: (20 - 41 cm)

x possible; - not possible

Push handle / Push bar- and backrest variants


Angle adjustable and foldable version:



Push bar, foldable
and angle adjustable



Push bar „wide“, foldable
and angle adjustable



Top view // reduction of seat
width by 1,5 cm

Select the backrest angle
adjustment range:

75 - 115°

80 - 120°

Extension push-bar

5 cm extension

10 cm extension

Back, trough shaped:



With fastening material for
„HOGGI align“ back cushions



Backrest and push bar
folded down



Back angle adjustment

Push handle / Push bar- and backrest variants

Foldable variants:



Push handle, foldable



Push bar, foldable



Push bar „wide“, foldable

Select backrest angle:

90° 100°

95° 105°

Extension push-bar

5 cm extension

10 cm extension

Back, trough shaped:



Backrest and push bar
folded down



Reinforced locking claws

Push handle / Push bar- and backrest variants

Detachable variants:



Push handles, detachable



Push bar, detachable



Push bar „wide“, detachable


Extension push-bar

5 cm extension**10 cm extension**

Back panel angle- and height adjustable:

**75 - 115°****80 - 120°**

Back folded down

Top view // reduction of
seat width by 1,5 cm

Tilt mechanism support kit


Tilt mechanism support kit (recommended for load over 40 kg)

Back cushion


Back cushion 2,5 cm

**Through shaped back
SW 32 cm****Through shaped back
SW 36, 40, 44 cm****Back cushion, black
SW 32 cm****Back cushion, black
SW 36, 40, 44 cm**„HOGGI align“ Contoured
back cushion, black

Seat panel & seat cushion



Seat panel



Seat cushion 3 cm



Contoured seat cushion,
with abduction

Notice:
Cannot be
combined
with hip
supports.

Upholstery colours




Black




Spacer fabric, black



Berry Soft



Turquoise Soft



Grey Soft

Headrests



Head support standard
with upholstery



Headrest bracket universal for
„Head support standard with
upholstery“ or for thirdparty
products



Head support with occiput
upholstery incl. mounting
bracket



Belt routing as a T-strap
above the shoulders

Armrests



Armrests


height: 14 - 21 cm

height: 19 - 26 cm



Notice:

Depending on the configuration of the wheelchair, seat tilt may be restricted.




Armrests with tilt balancing

height: 15 - 23 cm

height: 18 - 26 cm

**height: 18 - 26 cm
reinforced version**



PU-armrest pads


23,5 cm long

27 cm long

Brakes



Knee lever wheel locks



Integrated brake

Wheel guards



Plastic wheel guards



Carbon wheel cover

Rear wheels



Standard rear wheel with aluminium hand rim



Standard rear wheel with drum brake and aluminium hand rim



Rear wheel HOGGI light



Rear wheel with drum brake

Select tires:

Pneumatic tires
(1")

Flat free tires
(20" x 1 3/8", 22" x 1" and 24" x 1")



Push rims (only selectable in combination with HOGGI light rear wheels)

Without push rims

Push rims aluminium

Push rims aluminium with **anti-slip coating**

Push rims stainless steel

Anti-slip coating:



Black



Pure white



Silver grey



Sky blue



Traffic red



Pure orange

Push rim variant:




standard



high

Push rim mounting:



wide

narrow




Push rim cover, silicone black fluted

Push rims (only selectable in combination with HOGGI light rear wheels)



Closed ergonomic hand rims



„HOGGI Grip Plus“ hand rim inserts
black, for closed ergonomic
hand rims



Maxgrepp Original



Maxgrepp Ergo Para



Maxgrepp Ergo Para CP

Spoke guards

Spoke guard, transparent

Spoke guard Special interior foiling,
with motif

Spoke guard Special interior foiling,
**single-coloured according
to customer requirements:**

Spoke guard Special interior foiling,
printed according to customer requirements



Front castors



Wheel with rubber tire

125 x 32 mm

140 x 40 mm



Wheel with PU tire
190 x 50 mm



Skater wheel with LED
125 x 24 mm



PU-wheel with
plastic design rim

125 x 34 mm


Anti tip, Tip assist, Calf upholstery



Double anti tippers



Double anti tippers with tip assist



Calf upholstery
suitable to seat width

assembly left

assembly right

Foot rest hanger, Footrests



Foot rest hanger 90°

LLL short 20 - 36 cm

LLL long 20 - 41 cm



Foot rest hanger angle
adjustable 85° - 160°

LLL short 20 - 36 cm


LLL long 20 - 41 cm



Single foot rest hanger
85° - 160°

LLL short 20 - 36 cm

LLL long 20 - 41 cm



Adjustment to the SHERPA
Vario Footrest system



Notice:

Can be used from a seat
width of 36 cm.

Foot rest bracket




Foot rest bracket angle adjustable,
flip up type




Foot rest bracket angle adjustable,
with interlocking

Foot rest




Flat plastic foot rest



Plastic foot rest with heel stop



Foot rest aluminum



Two single footplates
in boat shape

in frame color

in black

with heel plate



Notice:

Can be used from a seat
height of 45 cm.

Supports



Thoracic support, particular adjustable, size approx. 12 x 13 x 2 cm

for back panel angle- and height adjustable

for back panel trough shaped




Chest-/ thoracic support

chest size 33 - 43 cm

chest size 39 - 49 cm

chest size 45 - 55 cm



Straight hip supports



Notice:

Hip pads reduce the maximum seat width by 3 cm.



Notice:

Only in conjunction with backrest angle and height adjustment.
In conjunction with a back panel trough shaped, the mounting bracket angle must also be adjusted.

MDP trunk supports (can be ordered separately)




MDP pad holder swing away and detachable

left
right



MDP pad holder fixed left / right



MDP trunk support

8 x 10 x 2 cm
10 x 12 x 2 cm
12 x 13 x 2 cm



MDP trunk support curved

6 x 11,5 x 2,5 cm
8 x 14 x 2,5 cm

as a pair

as a pair

as a pair

Belt systems



Lap belt



4-Point lap belt



Ankle hugger


SW 28 - 34 cm

SW 34 - 42 cm

Ø 4 - 6,5 cm


Ø 6 - 8,5 cm

Abduction block, Therapy tray


Abduction block,
swing away, detachable

LWH 11 x 8 x 10 cm

LWH 13 x 9 x 10 cm

Abduction block 30 mm height
extension kit**Notice:**

Only in combination with abduction block and contoured seat cushion.



Therapy tray

SW 32 cm

SW 36 cm

SW 40 - 44 cm

**Notice:**

Only in conjunction with arm rests.

Tie down kit



Tie down kit



with reinforced locking claws

**Notice:**

Required for push bars/push handles that can be folded down when used in conjunction with tie down kit.

Gripper clamps, Cosies



Four gripper clamps D25
with fasteners



„Cosies“ Thermofleece hand warmer
for companion

Canopy, One hand drive



Canopy incl. adapter
and rain protection



Reconstruction for ProAktiv One
hand drive



Notice:

In combination with the
trough-shaped back panel
(with low back height),
restrictions may occur.

Outdoor front end, Sibling seat



Outdoor front end



Outdoor front end
as device carrier variant



Mounting adapter for Outdoor front
end aid for „phil & teds“ sibling seat
or universal car seat adapter

Remark / Customer request / Individual special custom-made request