

HOGGI

SPOKE GUARDS



 Designed and
made in Germany



As colourful as life!

www.hoggi.de

HOGGI SPOKE GUARDS

// NEW DESIGNS

©hoggi 2021

As colourful as life!

From now on we offer self-designed modern and colourful spoke guards.

The different designs match perfectly with many HOGGI wheelchair colours.

The two designs „Butterfly“ and „Paper plane“ can also be individualised with your own name and a desired background colour which gives the wheelchair a **personal touch!**

We are continuously working on new designs. Please visit our homepage to find our latest designs: www.hoggi.de



NORTH POLE



POLAR BEAR AVA



POLAR BEAR HENRY



PENGUIN



FLAMINA FLAMINGO



ASTRONAUT BLUE



ASTRONAUT BLACK



SHIP AHOY



ONE WORLD



FLOWER POWER TURQUOISE



FLOWER POWER PINK



WATERCOLOUR



MEXICO



SOCIAL MEDIA



TRANSPARENT



INTO THE JUNGLE



GRAFFITI



FAIRY DUST

// SO INDIVIDUAL

As unique as you!

The two designs „Butterfly“ and „Paper plane“ can be customised according to your desire.

The colour of the spoke guards can be perfectly matched to the frame colour of your wheelchair and your name clearly demonstrates:


This wheelchair belongs to you!



BUTTERFLY



PAPER PLANE



YOUR CLUB

// YOUR CLUB

Out and about with your favorite club.

We would be happy to create the logo of your favorite club on the spoke guards.

For this we need a license for one-time private use in advance, which can be requested from the clubs.

HOGGI SPOKE GUARDS

// ATTACHMENT AND PRINT



Fastening clips black



Drum brake wheel
with transparent
fastening clips

The HOGGI **spoke guard** is available for all HOGGI rear wheels in the sizes 22" and 24".

By using **fastening clips**, the spoke guard is attached precisely to the wheel.

Depending on the chosen design, the clips in **black** or **transparent** allow a subtle attachment to the wheel.

Thanks to the development and use of special production tools we are able to offer for every type of rear wheel a perfect solution with an optimal fit.

The HOGGI **spoke guard** therefore needs **no cable ties** for fastening. The assembly on the wheel with the help of the clips is simple and clean.

Printing and foiling:

The new designs shine in brilliant colours thanks to **high-resolution digital printing** on air duct film. Due to the foiling on the inside the print is protected and thus ensures a **long service life** of your spoke guard.



Spoke guard transparent on HOGGI light rear wheel 22"



Spoke guard transparent on HOGGI light rear wheel 24"



Spoke guard transparent on rear wheel with drum brake 22" and 24"



Print area & central cutout (see the mark)


// Specifications


Material	Polycarbonate // 1,5 mm // colourless
Outside diameter 24"	approx. 505 mm
Outside diameter 22"	approx. 455 mm
Central cutout	Ø 140 mm
Fixing	black and transparent fastening clips

HOGGI

Partnering together to mobilise kids.

HOGGI GmbH

 Eulerstraße 27
D-56235 Ransbach-Baumbach

 Telefon: (+49) 2623 / 92 499-0
Telefax: (+49) 2623 / 92 499-99

 info@hoggi.de

www.hoggi.de