

Order**Offer**

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address

see invoice address

Company: _____

Street: _____

City, Zip code: _____

**Basic configuration**

- Aluminium rigid frame
- Anti tipper left and right
- Push-bar or Push-handles
- Wheel locks integrated in wheel guards
- HOGGI light rear wheels inclusive push rims

Frame colour

Black	Silver	Pearl ruby	Pearl night blue	Pearl gentian blue	Yellow green
Sparkling iron effect	Traffic yellow	Flame red	Sky blue	Pearl violet	Night blue
Grey aluminium	White	Pure orange	Signal blue	Pearl beige	White matt
Vermillion	Heather violet	Black dull matt	Signal violet	Light green	Ocean blue
Turquoise blue	Iron sparkling light	Salmon orange	Pastel blue	Light pink	Olive green
Raspberry red	Pastel green				

Special frame colour



Copper gold

Desired RAL colour:

Special paint finish (optional)




Powder coated castor forks
in desired colour:
(standard colour black)


Notice:

Please specify HOGGI frame colour.

Seat width and seat depth



Seat width: 24 cm (actual mounting width: 23,5 cm)
enter seat depth: (24 - 30 cm)

Seat width: 26 cm (actual mounting width: 25,5 cm)
enter seat depth: (24 - 32 cm)

Seat width: 28 cm (actual mounting width: 27,5 cm)
enter seat depth: (26 - 34 cm)


Seat width: 30 cm (actual mounting width: 29,5 cm)
enter seat depth: (28 - 36 cm)

Seat width: 32 cm (actual mounting width: 31,5 cm)
enter seat depth: (30 - 38 cm)

Seat width: 34 cm (actual mounting width: 33,5 cm)
enter seat depth: (32 - 40 cm)

Seat width: 36 cm (actual mounting width: 35,5 cm)
enter seat depth: (34 - 42 cm)


Back height and Lower leg length



Seat width	24 cm	26 cm	28 cm	30 cm	32 cm	34 cm	36 cm
Backrest height back trough shaped	36-44 cm	37-45 cm	38-48 cm	39-49 cm	40-50 cm	42-54 cm	43-55 cm
back panel angle adjustable	33-42 cm	35-44 cm	37-46 cm	39-48 cm	41-50 cm	43-52 cm	45-54 cm

enter back height:

enter lower leg length: (16 - 37 cm)



Camber



Seat height

SWINGBO VTi	Front castors with 100 mm wheel diameter	Front castors with 125 mm wheel diameter	Front castors with 140 mm wheel diameter
20" rear wheels	Camber: 3° 6° 9° SH: 37 cm x x x SH: 38 cm x x x SH: 39 cm x x x SH: 40 cm x x x SH: 41 cm x x x SH: 42 cm x x x SH: 43 cm x - -	Camber: 3° 6° 9° SH: 37 cm x x x SH: 38 cm x x x SH: 39 cm x x x SH: 40 cm x x x SH: 41 cm x x x SH: 42 cm x x x SH: 43 cm x - -	Camber: 3° 6° 9° SH: 38 cm x x x SH: 39 cm x x - SH: 40 cm x x x SH: 41 cm x x - SH: 42 cm x x - SH: 43 cm - x -
	Camber: 3° 6° 9° SH: 39 cm - x x SH: 40 cm x x x SH: 41 cm x x x SH: 42 cm x x x SH: 43 cm x x x SH: 44 cm x x - SH: 45 cm x - -	Camber: 3° 6° 9° SH: 39 cm - x x SH: 40 cm x x x SH: 41 cm x x x SH: 42 cm x x x SH: 43 cm x x x SH: 44 cm x x - SH: 45 cm x - -	Camber: 3° 6° 9° SH: 39 cm - - x SH: 40 cm x x x SH: 41 cm x x x SH: 42 cm x x x SH: 43 cm x x x SH: 44 cm x x x SH: 45 cm x - -
	Camber: 3° 6° 9° SH: 42 cm x x x SH: 43 cm - x x SH: 44 cm x - x SH: 45 cm - x - SH: 46 cm x x x SH: 47 cm x x x	Camber: 3° 6° 9° SH: 42 cm x - x SH: 43 cm x x x SH: 44 cm x x x SH: 45 cm - x x SH: 46 cm x x x	Camber: 3° 6° 9° SH: 42 cm x x x SH: 43 cm x x x SH: 44 cm x x x SH: 45 cm x x x SH: 46 cm - x x SH: 47 cm x x x

x possible; - not possible

Push handle / Push bar- and backrest variants


Angle adjustable and foldable version:



Push bar, foldable
and angle adjustable



Push bar „wide“, foldable
and angle adjustable



Top view // reduction of seat
width by 1,5 cm

Select the backrest angle
adjustment range:

75 - 115°

80 - 120°

Extension push-bar

5 cm extension

10 cm extension


Back, trough shaped:



With fastening material for
„HOGGI align“ back cushions



Backrest and push bar
folded down



Back angle adjustment

Push handle / Push bar- and backrest variants

Foldable variants:



Push handle, foldable



Push bar, foldable



Push bar „wide“, foldable

Select backrest angle:

90° 100°

95° 105°

Extension push-bar

5 cm extension

10 cm extension

Back, trough shaped:



Backrest and push bar
folded down



Reinforced locking claws

Push handle / Push bar- and backrest variants

Detachable variants:



Push handles, detachable



Push bar, detachable



Push bar „wide“, detachable

Extension push-bar

5 cm extension

10 cm extension

Back panel angle- and height adjustable:




75 - 115°

80 - 120°



Back folded down



Top view // reduction of seat width by 1,5 cm

Back cushion

Back cushion 2,5 cm

**Through shaped back
SW 24 - 32 cm**

**Through shaped back
SW 34 - 36 cm**


**Back cushion, black
SW 24 - 32 cm**

**Back cushion, black
SW 34 - 36 cm**




„HOGGI align“ Contoured
back cushion, black

Seat panel & seat cushion



Seat panel



Seat cushion 3 cm



Contoured seat cushion,
with abduction

Notice:
Cannot be
combined
with hip
supports.

Upholstery colours




Black




Spacer fabric, black



Berry Soft



Turquoise Soft



Grey Soft



Notice:

Upholstery colour black is the standard colour and is included in the basic equipment.

Headrests



Head support standard
with upholstery



Headrest bracket universal for
„Head support standard with
upholstery“ or for thirdparty
products



Head support with occiput
upholstery incl. mounting
bracket




Belt routing as a T-strap
above the shoulders


Armrests



Armrests

height: 14 - 21 cm**height: 19 - 26 cm**

Armrests with tilt balancing

height: 15 - 23 cm**height: 18 - 26 cm**

PU-Armrest pads

23,5 cm long**27 cm long****Notice:**

Depending on the configuration of the wheelchair, seat tilt may be restricted.

Wheel guards, Brake system




Plastic wheel guards with integrated wheel lock levers



Wheel guards aluminum in frame colour with integrated wheel lock levers



Carbon wheel cover with integrated wheel lock levers



Locking of the integrated brake

Rear wheels



Rear wheel HOGGI light



Rear wheel with drum brake

Select tires:

Pneumatic tires
(1")**Flat free tires**
(20" x 1 3/8", 22" x 1" and 24" x 1")

Push rims

Without push rims

Push rims aluminium

Push rims aluminium
with **anti-slip coating**Push rims
stainless steel

Anti-slip coating:



Black



Pure white



Silver grey



Sky blue



Traffic red



Pure orange


Push rim variant:



standard



high




wide

narrow

Push rim cover,
silicone black fluted

Closed ergonomic hand rims

„HOGGI Grip Plus“ hand rim inserts
black, for closed ergonomic
hand rims

Maxgrepp Original



Maxgrepp Ergo Para



Maxgrepp Ergo Para CP

Spoke guards

Spoke guard, transparent

Spoke guard Special interior foiling,
with motif

Spoke guard Special interior foiling,
**single-coloured according
to customer requirements:**

Spoke guard Special interior foiling,
printed according to customer requirements



Front castors



Wheel with rubber tire

100 x 32 mm

125 x 32 mm


140 x 40 mm



Skater wheel with LED

100 x 24 mm

125 x 24 mm




PU-wheel with plastic design rim


100 x 34 mm

125 x 34 mm

Bumper wheels, Tip assist, Calf upholstery




Bumper wheels



Tip assist

assembly left

assembly right



Calf upholstery
suitable to seat width

Foot rest hanger



Foot rest hanger 90°

LLL short 16 - 27 cm

LLL long 16 - 37 cm

Foot rest hanger angle adjustable
85° - 160°

LLL short 19 - 27 cm

LLL long 19 - 37 cm



Single foot rest hanger 85° - 160°

LLL short 19 - 27 cm

LLL long 19 - 37 cm

Foot rest bracket

Foot rest bracket angle adjustable,
flip up typeFoot rest bracket angle adjustable,
with interlockingReinforced **steel** footplate-carrier
(about 0,3 kg heavier)

Foot rest



Flat plastic foot rest



Plastic foot rest with heel stop

Foot rest aluminum
in frame colorTwo single footplates
in boat shape

Supports



Thoracic support, particular adjustable, size approx. 12 x 13 x 2 cm

**for back panel
angle- and height
adjustable**

**for back panel
trough shaped**




Chest-/ thoracic support

chest size 33 - 43 cm

chest size 39 - 49 cm

chest size 45 - 55 cm



Straight hip supports

22,5 x 10 cm, SW 28 - 34

25 x 11 cm, SW 28 - 34

25 x 11 cm, SW 36 - 44



Notice:

Only in conjunction with backrest angle and height adjustment.

In conjunction with a back panel trough shaped, the mounting bracket angle must also be adjusted.



Notice:

Hip pads reduce the maximum seat width by 3 cm.

MDP trunk supports (can be ordered separately)




MDP pad holder swing away and detachable

left

right



MDP pad holder fixed left / right




MDP trunk support

8 x 10 x 2 cm

10 x 12 x 2 cm

12 x 13 x 2 cm



MDP trunk support curved

6 x 11,5 x 2,5 cm


8 x 14 x 2,5 cm

as a pair

as a pair

as a pair

Belt systems



Lap belt



4-Point lap belt



Ankle hugger


SW 22 - 30 cm

SW 28 - 34 cm

Ø 4 - 6,5 cm


Ø 6 - 8,5 cm

Abduction block, Therapy tray


Abduction block,
swing away, detachable

LWH 9 x 6 x 9 cm

LWH 11 x 8 x 10 cm

Abduction block 30 mm height
extension kit**Notice:**

Only in combination with abduction block and contoured seat cushion.



Therapy tray

SW 24 and 26 cm

SW 28 to 32 cm

SW 34 and 36 cm

**Notice:**

Only in conjunction with arm rests.

Tie down kit



Tie down kit



with reinforced locking claws

**Notice:**

Required for push bars/push handles that can be folded down when used in conjunction with tie down kit.

Gripper clamps, Cosies



Four gripper clamps D25 with fasteners



„Cosies“ Thermofleece hand warmer for companion



Notice:

Required for attaching a third-party back panel.

Canopy, One hand drive



Canopy incl. adapter and rain protection



Reconstruction for ProAktiv One hand drive



Notice:

In combination with the trough-shaped back panel (with low back height), restrictions may occur.

Outdoor front end, Sibling seat

Outdoor front end

Outdoor front end
as device carrier variantMounting adapter for Outdoor front
end aid for „phil & teds“ sibling seat
or universal car seat adapter**Remark / Customer request / Individual special custom-made request**