

SWINGBO® VTi XL

Release date: September 2025 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address

see invoice address

Company: _____

Street: _____

City, Zip code: _____



Basic configuration

- Aluminium rigid frame
- Anti Tipper left and right
- Push-handle or push-bar
- HOGGI light- & standard brake
- Rear wheels HOGGI light incl. push rims

Frame colour

Black	Silver	Pearl ruby	Pearl night blue	Pearl gentian blue	Yellow green
Sparkling iron effect	Traffic yellow	Flame red	Sky blue	Pearl violet	Night blue
Grey aluminium	White	Pure orange	Signal blue	Pearl beige	White matt
Vermillion	Heather violet	Black dull matt	Signal violet	Light green	Ocean blue
Turquoise blue	Iron sparkling light	Salmon orange	Pastel blue	Light pink	Olive green
Raspberry red	Pastel green				

Special frame colour



Copper gold

Desired RAL colour:


Special paint finish (optional)


Notice:

Please specify HOGGI frame colour.

Powder coated castor forks
in desired colour:
(standard colour black)

Seat width and seat depth




Seat width: 36 cm (actual mounting width: 35,5 cm)
enter seat depth: (34 - 46 cm)

Seat width: 40 cm (actual mounting width: 39,5 cm)
enter seat depth: (38 - 50 cm)

Seat width: 44 cm (actual mounting width: 43,5 cm)
enter seat depth: (38 - 50 cm)

Back height and Lower leg length





Seat width	36 cm	40 cm	44 cm
Backrest height			
back trough shaped	43-55 cm	45-57 cm	48-60 cm
back panel angle adjustable	45-54 cm	50-59 cm	54-63 cm

enter back height:

enter lower leg length:

(20 - 41 cm)

Active degree



enter active degree: (0 - 10 cm)

Camber



Seat height

SWINGBO VTi XL	Front castors with 100 mm wheel diameter	Front castors with 125 mm wheel diameter	Front castors with 140 mm wheel diameter
22" rear wheels Base frame high 	Camber: 3° 6° 9° SH: 44 cm x x x SH: 45 cm x x x SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x	Camber: 3° 6° 9° SH: 45 cm x x x SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x SH: 49 cm x - -	Camber: 3° 6° 9° SH: 45 cm - x - SH: 46 cm x x x SH: 47 cm x x - SH: 48 cm x x x SH: 49 cm x - -
24" rear wheels Base frame high 	Camber: 3° 6° 9° SH: 45 cm x x x SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x SH: 49 cm x x x SH: 50 cm x x x SH: 51 cm x x -	Camber: 3° 6° 9° SH: 45 cm x x x SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x SH: 49 cm x x x SH: 50 cm x x - SH: 51 cm x x -	Camber: 3° 6° 9° SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x SH: 49 cm x x x SH: 50 cm x x x SH: 51 cm x x -
22" rear wheels Base frame low 	Camber: 3° 6° 9° SH: 41 cm - - x SH: 42 cm - - x SH: 43 cm - x x SH: 44 cm x x x SH: 45 cm x x x SH: 46 cm x x x	Camber: 3° 6° 9° SH: 41 cm - - x SH: 42 cm - - x SH: 43 cm - x x SH: 44 cm x x x SH: 45 cm x x x SH: 46 cm x x x	Camber: 3° 6° 9° SH: 41 cm - - x SH: 42 cm - - x SH: 43 cm - x x SH: 44 cm x x x SH: 45 cm x x x SH: 46 cm x x x
24" rear wheels Base frame low 	Camber: 3° 6° 9° SH: 44 cm - - x SH: 45 cm x x x SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x	Camber: 3° 6° 9° SH: 44 cm - - x SH: 45 cm x x x SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x	Camber: 3° 6° 9° SH: 45 cm - - x SH: 46 cm x x x SH: 47 cm x x x SH: 48 cm x x x

x possible; - not possible

Push handle / Push bar- and backrest variants


Angle adjustable and foldable version:



Push bar, foldable
and angle adjustable



Push bar „wide“, foldable
and angle adjustable



Top view // reduction of seat
width by 1,5 cm

Select the backrest angle
adjustment range:

75 - 115°

80 - 120°

Extension push-bar

5 cm extension

10 cm extension


Back, trough shaped:



With fastening material for
„HOGGI align“ back cushions



Backrest and push bar
folded down



Back angle adjustment

Push handle / Push bar- and backrest variants

Foldable variants:



Push handle, foldable



Push bar, foldable



Push bar „wide“, foldable

Select backrest angle:

90° 100°

95° 105°

Extension push-bar

5 cm extension

10 cm extension

Back, trough shaped:



Backrest and push bar
folded down



Reinforced locking claws

Push handle / Push bar- and backrest variants

Detachable variants:



Push handles, detachable



Push bar, detachable



Push bar „wide“, detachable


Extension push-bar

5 cm extension**10 cm extension**

Back panel angle- and height adjustable:

**75 - 115°****80 - 120°**

Back folded down

Top view // reduction of
seat width by 1,5 cm

Tilt mechanism support kit

Tilt mechanism support kit (recommended for load over 40 kg)

Back cushion

Back cushion 2,5 cm

**Through shaped back
SW 36, 40, 44 cm****Back cushion, black
SW 36, 40, 44 cm**„HOGGI align“ Contoured
back cushion, black

Seat panel & seat cushion



Seat panel



Seat cushion 3 cm

Contoured seat cushion,
with abduction

Notice:
Cannot be combined
with hip supports.

Upholstery colours




Black




Spacer fabric, black



Berry Soft



Turquoise Soft




Grey Soft



Notice:
Upholstery colour black is the standard colour and is included in the basic equipment.

Headrests

Head support standard
with upholsteryHeadrest bracket universal for
„Head support standard with
upholstery“ or for thirdparty
productsHead support with occiput
upholstery incl. mounting
bracketBelt routing as a T-strap
above the shoulders

Armrests



Armrests

height: 14 - 21 cm**height: 19 - 26 cm**

Armrests with tilt balancing

height: 15 - 23 cm**height: 18 - 26 cm****height: 18 - 26 cm
reinforced version**

PU-Armrest pads

23,5 cm long**27 cm long****Notice:**

Depending on the configuration of the wheelchair, seat tilt may be restricted.


Brake variants



HOGGI light sport brake



Standard brakes



Integrated brake

Wheel guards, Locking of the integrated brake




Plastic wheel guards



Wheel guards aluminum in frame colour for integrated wheel lock levers



Carbon wheel cover



Locking of the integrated brake

Rear wheels



Rear wheel HOGGI light



Rear wheel with drum brake

Select tires:

Pneumatic tires
(7")

Flat free tires
(22" x 1" and 24" x 1")

Push rims

Without push rims

Push rims aluminium


Push rims aluminium
with **anti-slip coating**

Push rims
stainless steel

Anti-slip coating:



Black



Pure white



Silver grey



Sky blue



Traffic red



Pure orange

Push rim variant:



standard



high



wide

narrow

Push rim cover,
silicone black fluted

Push rim mounting:



Closed ergonomic hand rims

„HOGGI Grip Plus“ hand rim inserts
black, for closed ergonomic
hand rims

Push rims for HOGGI light rear wheels

Maxgrepp Original



Maxgrepp Ergo Para



Maxgrepp Ergo Para CP

Spoke guards

Spoke guard, transparent

Spoke guard Special interior foiling,
with motif

Spoke guard Special interior foiling,
**single-coloured according
to customer requirements:**

Spoke guard Special interior foiling,
printed according to customer requirements

**Front castors**

Wheel with rubber tire

100 x 32 mm

125 x 32 mm

140 x 40 mm



Skater wheel with LED

100 x 24 mm

125 x 24 mm




PU-wheel with plastic design rim


100 x 34 mm

125 x 34 mm


Bumper wheels, Tip assist, Calf upholstery



Bumper wheels



Tip assist



Calf upholstery
suitable to seat width

assembly left

assembly right

Foot rest hanger, Footrest system



Foot rest hanger 90°

LLL short 20 - 36 cm

LLL long 20 - 41 cm



Foot rest hanger angle
adjustable 85° - 160°

LLL short 20 - 36 cm


LLL long 20 - 41 cm



Single foot rest hanger
85° - 160°

LLL short 20 - 36 cm

LLL long 20 - 41 cm



Adjustment to the SHERPA
Vario Footrest system

Foot rest bracket




Foot rest bracket angle adjustable,
flip up type



Foot rest bracket angle adjustable,
with interlocking

Foot rest




Flat plastic foot rest



Plastic foot rest with heel stop



Foot rest aluminum



Two single footplates
in boat shape

in frame color

in black

with heel plate

!
Notice:
Suitable for seat heights of
45 cm and above.

Supports



Thoracic support, particular
adjustable, size approx. 12 x 13 x 2 cm

**for back panel
angle- and height
adjustable**

**for back panel
trough shaped**




Chest-/ thoracic support

chest size 33 - 43 cm

chest size 39 - 49 cm


chest size 45 - 55 cm

!
Notice:
Only in conjunction with
backrest angle and height
adjustment.
In conjunction with a back
panel trough shaped, the
mounting bracket angle
must also be adjusted.




Straight hip supports (250 x 100 mm)

!
Notice:
Hip pads reduce the maximum
seat width by 3 cm.

MDP trunk supports (can be ordered separately)MDP pad holder swing away
and detachableMDP pad holder fixed
left / right

MDP trunk support



MDP trunk support curved

left**right****8 x 10 x 2 cm****10 x 12 x 2 cm****12 x 13 x 2 cm****6 x 11,5 x 2,5 cm****8 x 14 x 2,5 cm****as a pair****as a pair****as a pair****Belt systems**

Lap belt




4-Point lap belt



Ankle hugger

SW 28 - 34 cm**SW 34 - 42 cm****Ø 4 - 6,5 cm****Ø 6 - 8,5 cm**


Abduction block, Therapy tray



Abduction block,
swing away, detachable

LWH 11 x 8 x 10 cm

LWH 13 x 9 x 10 cm




Abduction block 30 mm height
extension kit



Notice:

Only in combination with abduction block and contoured seat cushion.



Therapy tray



Notice:

Only in conjunction with arm rests.

Tie down kit



Tie down kit



with reinforced locking claws



Notice:

Required for push bars/push handles that can be folded down when used in conjunction with tie down kit.

Gripper clamps, Cosies



Four gripper clamps D25
with fasteners



„Cosies“ Thermofleece hand warmer
for companion

Canopy, One hand drive



Canopy incl. adapter
and rain protection



Reconstruction for ProAktiv One
hand drive



Notice:

In combination with the
trough-shaped back panel
(with low back height),
restrictions may occur.

Outdoor front end, Sibling seat



Outdoor front end



Outdoor front end
as device carrier variant



Mounting adapter for Outdoor front
end aid for „phil & teds“ sibling seat
or universal car seat adapter

Remark / Customer request / Individual special custom-made request