

www.hoggi.de
HOGGI®

Made in Germany
COBRA
by HOGGI 

COBRA



Indoor base for BINGO Evolution



- Gas pressure spring for height adjustment & 45° tilt

COBRA - Always on top

COBRA is developed especially as a indoor base with a 45° tilt mechanism and a continuous adjustable gas pressure spring for the height adjustment. It is made for different application areas:

- As **BINGO Evolution home under frame** to pick up the **BINGO Evolution** seat unit.
- As **COBRA seat shell frame** with a support for seat shells to pick up orthopaedic seat systems and seat shells (Seat width 34 cm or 39 cm).

COBRA can be equipped with 75 mm or 125 mm wheel diameter.



COBRA ist:

- stable
- space-saving
- 45° tilt
- continuous height adjustable and adjustable on load

COBRA

Technical data

Seat-/ lift height	
Caster diameter 75 mm	27 - 61 cm
Caster diameter 125 mm	33,5 - 67,5 cm
Total width	61 cm
Total length	71 cm
Caster diameter	75 or 125 mm
Tilt	0° to 45°
Total weight	
Caster diameter 75 mm	11,3 kg
Caster diameter 125 mm	13,3 kg
Load capacity	50 kg

HOGGI®

Partnerschaft für mobile Kinder

HOGGI GmbH • Eulerstraße 27 • D-56235 Ransbach-Baumbach • Internet: www.hoggi.de

E-mail: info@hoggi.de • Telefon: (+49) 2623 / 92 499-0 • Telefax: (+49) 2623 / 92 499-99