



BINGO Evolution rehab pushchair

Size 1 for children at the age of approx. 1 - 6 years

Crash tested according to ISO 7176-19 and ANSI/RESNA WC/Vol.1. - Section 19

Order form incl. Exportprices (€) valid until revoked, release date 04/2024

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<input type="text"/>	<input type="text"/>
Invoice address	
Company	Customer number
<input type="text"/>	<input type="text"/>
Street	City, Zip code
<input type="text"/>	<input type="text"/>
E-Mail	Telephone
<input type="text"/>	<input type="text"/>
Fax	
<input type="text"/>	
Delivery address	
Company	
<input type="text"/>	
Street	City, Zip code
<input type="text"/>	<input type="text"/>



Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable, rotatable

Seat unit

- IUS - Individual Upholstery System
- Seat depth 16 - 34 cm
- Seat width 18 - 30 cm
- Lower leg length 16 - 37 cm
- Back height 45 - 63 cm
- Recline 90° - 160°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 50 kg

Minimum information (Without this information your order can not be handled.)

Frame colour:	<input type="checkbox"/> Black dull matt 032	<input type="checkbox"/> Silver 002	
Upholstery single-coloured:	<input type="checkbox"/> Black 346	<input type="checkbox"/> Berry soft 347	<input type="checkbox"/> Grey soft 348
with accessories colour black	<input type="checkbox"/> Turquoise soft 349		
Upholstery two-coloured:	<input type="checkbox"/> Berry soft / Black 350	<input type="checkbox"/> Grey soft / Black 351	<input type="checkbox"/> Turquoise soft / Black 352
with accessories colour black			

Seat depth: (16 -34 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height: (45 - 63 cm)

Lower leg length: (16 - 37 cm)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colouir	Description	Qty.	Price
	<input type="checkbox"/> 3232-2000		BINGO Evolution seat unit, size 1	1 pc.	
with	<input type="checkbox"/> 3231-2900		Thigh supports, size 1, wide- and depth-adjustable * minimum seat width 20 cm	01 1 pair	incl.
or	<input type="checkbox"/> 3231-2600		Hip supports, size 1, wide- and depth-adjustable, removable	(upcharge) 02 1 pair	
	<input type="checkbox"/> 3232-5100		Upholstery, size 1, adaptable (IUS)	03 1 set	incl.
	<input type="checkbox"/> 3232-5190		Calf pad support plate and foam insert	04 1 set	
with	<input type="checkbox"/> 3232-2500		Footrest hanger 90°	05 1 pc.	incl.
or	<input type="checkbox"/> 3232-2402		Footrest hanger angle adjustable 85 - 170° * * minimum lower leg length 19 cm	(upcharge) 06 1 pc.	
or	<input type="checkbox"/> 3232-2504		Footrest hanger 90° (short LLL) with bracket and footrest flip up lock *	(upcharge) 1 pc.	
or	<input type="checkbox"/> 3232-2404		Footrest hanger angle adjustable 85-170° (short LLL) with bracket and footrest flip up lock * * lower leg length 15-36 cm	(upcharge) 1 pc.	
with	<input type="checkbox"/> 3231-249x		Footplate with heel stop, width: <input type="checkbox"/> 1=27 cm, <input type="checkbox"/> 2=29 cm, <input type="checkbox"/> 3=31 cm, <input type="checkbox"/> 4=33 cm	07 1 pc.	incl.
or	<input type="checkbox"/> 3231-2551		Footplate in boat form, width: approx. 30 cm	(upcharge) 07 1 pc.	
or	<input type="checkbox"/> 3231-2552	015	Alu foot plate in boat form, size 1	(upcharge) 1 pc.	
or	<input type="checkbox"/> 3231-2463	015	Alu-footplate "footprint" with heel stop, width 30cm	(upcharge) 08 1 pc.	
with	<input type="checkbox"/> 3240-2950		Footrest flip up lock	(upcharge) 09 1 pc.	
with	<input type="checkbox"/> SA-0430		Reinforced steel footplate-carrier (about 0,3 kg heavier)	10 1 pc.	incl.

Accessories - seat unit

	Orderno.	Colouir	Description	Qty.	Price
	<input type="checkbox"/> 3231-7600		Headrest pads low version size 1	11 1 set	
	<input type="checkbox"/> 3201-7600		Headrest pads, size 1	12 1 set	
	<input type="checkbox"/> 3201-7601		Headrest pads with spacer fabric for speakers, size 1	13 1 set	
	<input type="checkbox"/> 3231-6900		Anatomical neck-headrest with length compens, size 1	14 1 set	
	<input type="checkbox"/> 3231-7700		Trunk supports, size 1, stepless height adjustable (8 x 10 x 2 cm)	15 1 pair	
	<input type="checkbox"/> 3233-7700		Trunk supports, size 2, stepless height adjustable (10 x 12 x 2 cm)	1 pair	
	<input type="checkbox"/> 3201-6100		Lap belt	16 1 pc.	
	<input type="checkbox"/> 3231-6701		4-point lap belt, size XS (for seat width 18 - 24)	17 1 pc.	
	<input type="checkbox"/> 3231-6702		4-point lap belt, size S (for seat width 22 - 30)	1 pc.	
	<input type="checkbox"/> 3201-6200		Five-point harness	18 1 pc.	
	<input type="checkbox"/> 3201-6300		Foot straps	19 1 pair	
	<input type="checkbox"/> 3201-6400		Ankle huger, size 1	20 1 pair	
	<input type="checkbox"/> 3231-5300		Groin strap, size 1	21 1 pc.	
	<input type="checkbox"/> 3201-5800		Fixation jacket, size 0	22 1 pc.	
	<input type="checkbox"/> 3201-5400		Fixation jacket, size 1	1 pc.	
	<input type="checkbox"/> 3201-5500		Fixation jacket, size 2	1 pc.	
	<input type="checkbox"/> 3201-5600		Chest-shoulder harness, size 1	23 1 pc.	
	<input type="checkbox"/> 3201-5601		Chest-shoulder harness, size 1 with 4 buckels	24 1 pc.	
	<input type="checkbox"/> 3201-5700		Chest-shoulder harness, size 2	1 pc.	
	<input type="checkbox"/> 3201-5701		Chest-shoulder harness, size 2 with 4 buckels	1 pc.	
	<input type="checkbox"/> 3231-7800		Abduction block, size 1, swing away and removable	25 1 pc.	
	<input type="checkbox"/> 3231-7801		Abduction block, size 1, fixed	26 1 pc.	
	<input type="checkbox"/> 3231-7900		Abduction block, size 2, swing away and removable	1 pc.	
	<input type="checkbox"/> 3231-7901		Abduction block, size 2, fixed	1 pc.	
	<input type="checkbox"/> 3232-2203	032	Retainer for accessories, grab rail or therapy tray	27 1 pair	
	<input type="checkbox"/> 3231-8100		Grap rail with upholstery, size 1 *	28 1 pc.	
	<input type="checkbox"/> 3231-8200		Therapy tray, size 1, transparent ^{1) 2) 3)}	29 1 pc.	
	<input type="checkbox"/> 3231-8202		Therapy tray, size 1, transparent, height adjustable ^{1) 2)}	30 1 pc.	
	<input type="checkbox"/> 3231-8201		Therapy tray, size 1, grey ^{1) 2) 3)}	1 pc.	
	<input type="checkbox"/> 3231-8203		Therapy tray, size 1, grey, height adjustable ^{1) 2)}	1 pc.	
			¹⁾ only possible with retainer for accessories		
			²⁾ only possible with seat in driving direction		
			³⁾ not possible when using hip supports		
	<input type="checkbox"/> 3201-8300	306	Winter warmer, size 1, black	31 1 pc.	
	<input type="checkbox"/> 3201-8400	306	Winter warmer, size 2, black	1 pc.	
	<input type="checkbox"/> 7000-0018		Sheep skin insert, size 1	32 1 pc.	
	<input type="checkbox"/> 7000-0019		Sheep skin insert, size 2	1 pc.	
	<input type="checkbox"/> 3201-9200	306	Rain cape, size 1, black	33 1 pc.	
	<input type="checkbox"/> 3201-9300	306	Canopy incl. adapter and rain cover, black	34 1 pc.	

Accessories - seat unit

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3002-0010	032	BASIC home underframe size 1 with 75 mm twin castors	35	1 pc.
	<input type="checkbox"/> 3002-0011	032	BASIC home underframe size 1 with 125 mm castors	36	1 pc.
	<input type="checkbox"/> 3002-9500	022	Side covers in standard color white	37	
	<input type="checkbox"/> 3002-9500		Side covers matching the upholstery colour	38	
	<input type="checkbox"/> 3002-9500	032	Side covers in black		
	<input type="checkbox"/> 3002-5400		Push handle for BASIC home under frame		1 pc.
	<input type="checkbox"/> 3004-0000	032	COBRA indoor base, size 1, with 75 mm twin castors	39	1 pc.
	<input type="checkbox"/> 3004-0001	032	COBRA indoor base, size 1, with 125 mm castors	36	1 pc.
Presetting of the gas spring height adjustment: <input type="checkbox"/> up to 20kg, <input type="checkbox"/> up to 35kg, <input type="checkbox"/> up to 50kg					

Basic variations - mobility base swivelling front wheels

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3232-1000		BINGO Evolution mobility base, size 1, swivelling front wheels		1 pc.
	<input type="checkbox"/> 3238-1000		BINGO Evolution mobility base, size 1, swivelling front wheels, for respiratory equipment		1 pc.
with	<input type="checkbox"/> 1281-1027		Plastic front forks	40	1 set incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	41	1 set
or	<input type="checkbox"/> 1281-1028		Aluminum front forks	(upcharge) 42	1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks		1 set
with	<input type="checkbox"/> 3231-1301		Without suspension system		1 set incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable	(upcharge) 43	1 set
	<input type="checkbox"/> 1131-9017	022	Spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017	032	Spring set (black) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		Spring set (matching the upholstery color) for load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	44	1 pair incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	45	1 pair incl.
with	<input type="checkbox"/> 1219-0004S		Pneumatic wheels with 6-spoke rims	46	1 set incl.
or	<input type="checkbox"/> 1219-0030		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge) 47	1 set
or	<input type="checkbox"/> 1229-0004S		Polyurethan wheels with 6-spoke rims	48	1 set incl.
or	<input type="checkbox"/> 1229-0009		Polyurethan wheels with 7-spoke rims	49	1 set incl.
or	<input type="checkbox"/> 1229-0030		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge) 50	1 set
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use	51	1 set
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	52	1 pc. incl.
or	<input type="checkbox"/> 3231-9402		Accessory bag, size 1 (max. 8 kg) *	(upcharge) 53	1 pc.
* do not combine with tie down kit, use with seat in driving direction only					
with	<input type="checkbox"/> 3231-1501		Without tie down kit		1 set incl.
or	<input type="checkbox"/> 3231-7300		Tie down kit according to ISO 7176-19	(upcharge) 54	1 set

Basic variations - mobility base CROSS with Aluminium front forks

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3232-1030		BINGO Evolution CROSS, size 1 mobility base	1 pc.	
	<input type="checkbox"/> 3238-1030		BINGO Evo CROSS, size 1 mobility base for additional equipment	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	43	1 set
	<input type="checkbox"/> 1131-9017	022	Spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017	032	Spring set (black) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		Spring set (matching the upholstery color) for load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	44	1 pair incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	45	1 pair incl.
with	<input type="checkbox"/> 1219-0050		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	47	1 set incl.
or	<input type="checkbox"/> 1229-0050		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	50	1 set incl.
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	52	1 pc. incl.
or	<input type="checkbox"/> 3231-9402		Accessory bag, size 1 (max. 8 kg) * (upcharge)	53	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3231-1501		Without tie down kit		1 set incl.
or	<input type="checkbox"/> 3231-7300		Tie down kit according to ISO 7176-19 (upcharge)	54	1 set
with	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks		1 set
with	<input type="checkbox"/> 3234-7100		Anti tip left and right *		1 set
			* not possible in combination with frame for respiratory equipment		

! BINGO Evolution Cross: Custom-made only on request !

Basic variations - mobility base rigid front wheels

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3231-1000		BINGO Evolution mobility base, size 1, rigid front wheels	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	43	1 set
	<input type="checkbox"/> 1131-9017	022	Spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017	032	Spring set (black) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		Spring set (matching the upholstery color) for load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	44	1 pair incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	45	1 pair incl.
with	<input type="checkbox"/> 1219-0025S		Pneumatic wheels with 6-spoke rims	46	1 set incl.
or	<input type="checkbox"/> 1219-0029		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	47	1 set
or	<input type="checkbox"/> 1229-0021S		Polyurethan wheels with 6-spoke rims	48	1 set incl.
or	<input type="checkbox"/> 1229-0029		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	50	1 set
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use	51	1 set
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	52	1 pc. incl.
or	<input type="checkbox"/> 3231-9402		Accessory bag, size 1 (max. 8 kg) * (upcharge)	53	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3231-1501		Without tie down kit		1 set incl.
or	<input type="checkbox"/> 3231-7300		Tie down kit according to ISO 7176-19 (upcharge)	54	1 set

Basic variations - mobility base TRIKE (custom-made product)

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/>	3236-1000	BINGO Evolution mobility base TRIKE, size 1	1 pc.	
with	<input type="checkbox"/>	3231-1301	Without suspension system	1 set	incl.
or	<input type="checkbox"/>	3231-7400	Frame suspension system, fourfold adjustable (upcharge)	43	1 set
	<input type="checkbox"/>	1131-9017 022	Spring set (white) for load up to 15 kg		
	<input type="checkbox"/>	1131-9017 032	Spring set (black) for load up to 15 kg		
	<input type="checkbox"/>	1131-9017	Spring set (matching the upholstery color) for load up to 15 kg		
	<input type="checkbox"/>	1131-9018 022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/>	1131-9018 032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/>	1131-9018	Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/>	1131-9019 022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/>	1131-9019 032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/>	1131-9019	Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/>	3231-9500	Side covers matching the upholstery color (Standard)	44	1 pair incl.
or	<input type="checkbox"/>	3231-9500	Side covers matching the frame color	45	1 pair incl.
with	<input type="checkbox"/>	1219-0011S	Pneumatic wheels with 6-spoke rims	46	1 set incl.
or	<input type="checkbox"/>	1219-0031	Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	47	1 set
or	<input type="checkbox"/>	1229-0011S	Polyurethan wheels with 6-spoke rims	48	1 set incl.
or	<input type="checkbox"/>	1229-0031	Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	50	1 set
	<input type="checkbox"/>	3231-7201	Attendant brake, two hand use	51	1 set
with	<input type="checkbox"/>	3232-9401	Soft accessory bag, size 1 (max. 3 kg)	52	1 pc. incl.
or	<input type="checkbox"/>	3231-9402	Accessory bag, size 1 (max. 8 kg) * (upcharge)	53	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input checked="" type="checkbox"/>	3231-1501	Without tie down kit		1 set incl.

Accessories - mobility base

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/>	3202-7100	Swivel lock for Plastic front forks	41	1 set
	<input type="checkbox"/>	3232-7100	Swivel lock for Aluminum front forks		1 set
	<input type="checkbox"/>	3231-7400E	Frame suspension system, fourfold adjustable, retrofit (upcharge)	43	1 set
	<input type="checkbox"/>	1131-9017 022	Spring set (white) for load up to 15 kg		
	<input type="checkbox"/>	1131-9017 032	Spring set (black) for load up to 15 kg		
	<input type="checkbox"/>	1131-9017	Spring set (matching the upholstery color) for load up to 15 kg		
	<input type="checkbox"/>	1131-9018 022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/>	1131-9018 032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/>	1131-9018	Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/>	1131-9019 022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/>	1131-9019 032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/>	1131-9019	Spring set (matching the upholstery color) for load above to 25 kg		
	<input type="checkbox"/>	3231-9402E	Accessory bag, size 1 (max. 8 kg), retrofit	53	1 pc.
	<input type="checkbox"/>	3231-7300E	Tie down kit according to ISO 7176-19	54	1 set
	<input type="checkbox"/>	3231-9500E	Side covers matching the upholstery colour or the frame color * * please enter colour code		1 pair
	<input type="checkbox"/>	3231-2552E 015	Alu foot plate in boat form, size 1		1 pc.

Custom-made

Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/> 3202-4100		Tray for additional equipment, size 1 (mobility base 3238-10X0 only)	55 1 pc.	
<input type="checkbox"/> 3238-4400		BINGO Evolution Cross tray for additional equipment size 1 (max. 25kg)	1 pc.	
<input type="checkbox"/> 3231-4200		Oxygen bottle holder (d=100 mm) for bottles up to 75 cm *	56 1 pc.	
		* assembly on the right side only		
<input type="checkbox"/> SA-0225-MP		Adjustable Angle brackets for contoured oxygen units (upcharge)	57 1 set	
<input type="checkbox"/> 3231-4201		Oxygen bottle holder (d=100 mm) for short bottles up to 40 cm *		
		* assembly on the right side only		
<input type="checkbox"/> SA-0031-KP		Infusions bottle holder	58 1 pc.	
<input type="checkbox"/> SA-0016/3-MP		Reinforced rear axle for mobility base max. load 70 kg (upcharge)	59 1 pc.	
<input type="checkbox"/> SA-0298-MP		Seat adapter 10 cm heightened * (upcharge)	1 set	
		* with seat in driving direction only		
<input type="checkbox"/> SA-0321-MP		Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base * (upcharge)	60 1 set	
		* only possible with 2 mobility bases of one size, not for base 3230-1000		
		* not possible in combination with aluminium front forks & mobility base Cross		
<input type="checkbox"/> SA-0441-MP		Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution mobility base 3232-1000 or 3238-1000 rear * (upcharge)	1 set	
		* not possible in combination with aluminium front forks & mobility base Cross		

Extra equipment

Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/> 7000-0010		Buggy-Board	61 1 pc.	
<input type="checkbox"/> 7000-0012	306	Sun shade, black	62 1 pc.	
<input type="checkbox"/> 7000-0006		Headrest speakers, for headrest pads 3201-7601	63 1 set	
<input type="checkbox"/> 7000-0202	306	Changing backpack, black	64 1 pc.	
<input type="checkbox"/> 7000-0199		"Cosies" Thermofleece hand warmer for companion	65 1 pair	

Pre-assembly of pads, upholstery parts & belts

Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/> 6999-1012		BINGO Evolution Pre-assembly service (without settings)	1 pc.	


CONFIGURATIONS // ACCESSORIES // BINGO Evolution




01 Thigh supports



02 Hip supports



03 IUS - Upholstery



04 Calf pad support plate and foam insert




05 Footrest hanger 90°



06 Footrest hanger angle adjustable 85 - 170°



07 Footplate with heel stop or boat form



08 Alu-footplate „foot-print“ with heel stop



09 Footrest flip up lock



10 Reinforced steel footplate-carrier



11 Headrest pads low version



12 Headrest pads



13 Headrest pads with spacer fabric for speakers



14 Anatomical neck-headrest



15 Trunk supports



16 Lap belt



17 4-point lap belt



18 Five-point harness



19 Foot straps



20 Ankle hugger



21 Groin strap



22 Fixation jacket



23 Chest-shoulder harness




24 Chest-shoulder harness with 4 buckles



25 Abduction block, swing away and removable




26 Abduction block, fixed




27 Retainer for accessories



28 Grap rail with upholstery



29 Therapy tray



30 Therapy tray, height adjustable

CONFIGURATIONS // ACCESSORIES // BINGO Evolution



31 Winter warmer



32 Sheep skin insert



33 Rain cape



34 Canopy



35 BASIC home underframe, 75 mm twin castors



36 125 mm castors



37 Side covers in standard color white



38 Side covers matching the upholstery colour



39 COBRA indoor base, 75 mm twin castors



40 Plastic front forks



41 Swivel lock



42 Aluminum front forks



43 Frame suspension system



44 Side covers matching the upholstery color



45 Side covers matching the frame color



46 Pneumatic wheel, front wheel



47 Pneumatic wheel, rear wheel



48 Pneumatic tyres front wheel, Drum brake wheel



49 Pneumatic tyres rear wheel, Drum brake wheel



50 Polyurethan wheel, front wheel



51 Polyurethan wheel, rear wheel



52 Polyurethan wheel front wheel, black rim



53 Polyurethan wheel rear wheel, black rim



54 Polyurethan tyres front wheel, Drum brake wheel



55 Polyurethan tyres rear wheel, Drum brake wheel



56 Attendant brake, two hand use




57 Soft accessory bag



58 Accessory bag



59 Tie down kit according to ISO 7176-19



60 Tray for additional equipment

CONFIGURATIONS // ACCESSORIES // BINGO Evolution




56 Oxygen bottle holder



57 Adjustable Angle brackets



58 Infusions bottle holder



59 Reinforced rear axle



60 Connecting parts for mobility bases



61 Buggy-Board



62 Sun shade, grey



63 Headrest speakers



64 Changing backpack



65 „Cosies“ Thermofleece hand warmer