



## FLUX posterior walker

Size 1 - 4

Order form incl. Exportprices (€) valid until revoked, release date 04/2024

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<input type="text"/>	<input type="text"/>
<b>Invoice address</b>	
Company	Customer number
<input type="text"/>	<input type="text"/>
Street	City, Zip code
<input type="text"/>	<input type="text"/>
E-Mail	Telephone
<input type="text"/>	<input type="text"/>
Fax	
<input type="text"/>	
<b>Delivery address</b>	
Company	
<input type="text"/>	
Street	City, Zip code
<input type="text"/>	<input type="text"/>



### Basic configuration

#### Durable, compact frame in a light aluminium construction

- Grab rail height- and depth adjustable
- Front swivel wheels
- Polyurethan tires
- Engageable and disengageable free-wheeling system on rear wheels

#### Stock- Frame colours: (available immediately)

- Size 1: **003** Pearl ruby red (RAL 3032)
- Size 2: **004** Pearl night blue (RAL 5026)
- Size 3: **005** Pearl gentian blue (RAL 5025)
- Size 4: **021** Grey aluminium (RAL 9007)

**Upholstery:** Black **306**

### Minimum information ( Without this information your order can not be handled.)

Frame colour  ( Please see table „frame colours“)

### Technical data: Measures (cm) and weight (kg)

		<b>FLUX size 1</b>	<b>FLUX size 2</b>	<b>FLUX size 3</b>	<b>FLUX size 4</b>
<b>Grab rail standard</b>	Grip height	48 - 61	57 - 81	69 - 87	82 - 102
	Grip width	41	46	51	51
<b>Grab rail with universal grips</b>	Grip height	43 - 68	51 - 88	62 - 94	75 - 102
	Grip width	29 - 41	34 - 46	39 - 51	39 - 51
<b>Frame</b>	Total width	64	69	74	78
	Load capacity	60 kg	75 kg	75 kg	100 kg
<b>Wheel diameter</b>	Front / rear	15 / 20	15 / 20	15 / 20	20 / 30

**Basic variations**

with/or: Please choose one version, only!

	Orderno.	Colouir	Description	Qty.	Price
	<input type="checkbox"/> 3301-0000		<b>FLUX, size 1</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3301-1500		Grab rail standard	<b>01</b>	1 pc. incl.
or	<input type="checkbox"/> 3301-5100		Grab rail with universal handle	<b>(upcharge) 02</b>	1 pc.
with	<input type="checkbox"/> 3301-6700		Low friction reverse roll lock	<b>03</b>	1 set
	<input type="checkbox"/> 3302-0000		<b>FLUX, size 2</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3302-1500		Grab rail standard	<b>01</b>	1 pc. incl.
or	<input type="checkbox"/> 3302-5100		Grab rail with universal handle	<b>(upcharge) 02</b>	1 pc.
with	<input type="checkbox"/> 3301-6700		Low friction reverse roll lock	<b>03</b>	1 set
	<input type="checkbox"/> 3303-0000		<b>FLUX, size 3</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3303-1500		Grab rail standard	<b>01</b>	1 pc. incl.
or	<input type="checkbox"/> 3303-5100		Grab rail with universal handle	<b>(upcharge) 02</b>	1 pc.
	<input type="checkbox"/> 3313-0000		<b>FLUX, size 4</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3313-1600		Grab rail standard	<b>01</b>	1 pc. incl.
or	<input type="checkbox"/> 3313-5200		Grab rail with universal handle	<b>(upcharge) 02</b>	1 pc.
with	<input type="checkbox"/> 1229-0008S		Polyurethan tires with 6-spoke rims	<b>04</b>	1 set incl.
or	<input type="checkbox"/> 1229-0009		Polyurethan wheels with 7-spoke rims	<b>05</b>	1 set incl.

**Accessories**

	Orderno.	Colouir	Description	Qty.	Price
	<input type="checkbox"/> 3301-5200		Lower arm support with handles, 14 x 9,4 cm	<b>06</b>	1 set
	<input type="checkbox"/> 3301-5300		Lower arm support with handles, 3-axial adjustable, 14 x 9,4 cm		1 set
	<input type="checkbox"/> 3302-5200		Lower arm support with handles, 22 x 10,4 cm	<b>07</b>	1 set
	<input type="checkbox"/> 3302-5300		Lower arm support with handles, 3-axial adjustable, 22 x 10,4 cm		1 set
	<input type="checkbox"/> 3301-6600		Hand disk	<b>08</b>	1 set
	<input type="checkbox"/> 3301-5400		Bumper wheels	<b>09</b>	1 set
	<input type="checkbox"/> 3301-5500		Drag brake, fits to size 1 - 3	<b>10</b>	1 set
	<input type="checkbox"/> 3313-5500		Drag brake, fits to size 4*		1 set
	<input type="checkbox"/> 3301-5600		Brake with wheel lock, fits to size 1 - 3	<b>11</b>	1 set
	<input type="checkbox"/> 3313-5600		Brake with wheel lock, fits to size 4*		1 set
			<b>* FLUX 4 "Drag brake" and "Brake with wheel lock" can not be mounted at once</b>		
	<input type="checkbox"/> 3301-5700		Anti-tips, fits to size 1 - 3	<b>12</b>	1 set
	<input type="checkbox"/> 3301-5701		Anti-tips with castors fits to Flux size 1-3	<b>13</b>	1 set
	<input type="checkbox"/> 3313-5700		Anti-tips, fits to size 4	<b>14</b>	1 set
	<input type="checkbox"/> 3313-5701		Anti-tips with castors fits to Flux size 4	<b>15</b>	1 set
	<input type="checkbox"/> 3301-5800		Seat, fits to size 1 + 2	<b>16</b>	1 pc.
	<input type="checkbox"/> 3313-5800		Seat, fits to size 3 + 4		1 pc.
	<input type="checkbox"/> 3301-5801		Seat folds down, suitable for size 1	<b>17</b>	1 pc.
	<input type="checkbox"/> 3302-5801		Seat folds down, suitable for size 2		1 pc.
	<input type="checkbox"/> 3301-5900		Hip-/bottom pads, multi adjustable	<b>18</b>	1 pc.
	<input type="checkbox"/> 3301-5970		Lap belt for Hip-/bottom pads	<b>19</b>	1 pc.
	<input type="checkbox"/> 3301-6000		Dynamic hip control kit for orthoses *	<b>20</b>	1 set
	<input type="checkbox"/> 3301-6001		Dynamic hip control, suspension soft (yellow) *	<b>21</b>	1 set
	<input type="checkbox"/> 3301-6002		Dynamic hip control, suspension hard (red) *		1 set
			<b>* Hip belt is not included</b>		
	<input type="checkbox"/> 3301-6300		Hip belt for dynamic hip control, hip size 53 - 83 cm	<b>22</b>	1 pc.
	<input type="checkbox"/> 3302-6300		Hip belt for dynamic hip control, hip size 65 - 115 cm		1 pc.
	<input type="checkbox"/> 3301-6100		Carrier, fits size 1 and size 2 *	<b>23</b>	1 pc.
	<input type="checkbox"/> 3303-6100		Carrier, fits size 3 and size 4 *		1 pc.
			<b>* combination with anti-tips only</b>		
	<input type="checkbox"/> 3301-6200		Basket suitable for carrier 3301-6100 *	<b>24</b>	1 pc.
	<input type="checkbox"/> 3303-6200		Basket suitable for carrier 3303-6100 *		1 pc.
			<b>* combination with carrier only</b>		
	<input type="checkbox"/> 3202-7100		Swivel lock	<b>25</b>	1 set

## Accessories

	Orderno.	Colouir	Description	Qty.	Price
<input type="checkbox"/>	3302-6400		Sling seat, hip size 57 - 77 cm	26	1 pc.
<input type="checkbox"/>	3303-6400		Sling seat, hip size 73 - 103 cm		1 pc.
<input type="checkbox"/>	3301-6500		Sling seat rail, size 1	27	1 pc.
<input type="checkbox"/>	3302-6500		Sling seat rail, size 2		1 pc.
<input type="checkbox"/>	3303-6500		Sling seat rail, size 3, fits to size 3 and size 4		1 pc.
<input type="checkbox"/>	3301-6800		Trunk support size 1, multi adjustable	28	1 pc.
<input type="checkbox"/>	3302-6800		Trunk support size 2, multi adjustable		1 pc.
<input type="checkbox"/>	3303-6800		Trunk support size 3, multi adjustable		1 pc.

## Extra equipment

	Order no.	Colouir	Description	Qty.	Price
<input type="checkbox"/>	SA-0210-MP		Special Frame colour, RAL code only	(upcharge)	1 pc.
<input type="checkbox"/>	6999-0014	060	Special Frame colour, Copper Gold	(upcharge)	1 pc.
<input type="checkbox"/>	SA-0548/1-KP		FLUX attendant-push bar size 1	(upcharge) 29	1 pc.
<input type="checkbox"/>	SA-0548/2-KP		FLUX attendant-push bar size 2	(upcharge)	1 pc.
<input type="checkbox"/>	SA-0548/3-KP		FLUX attendant-push bar size 3+4	(upcharge)	1 pc.


## Frame colours

<b>001</b> Black (RAL 9005)	<b>020</b> Yellow green (RAL 6018)	<b>028</b> Pearl violet (RAL 4011)	<b>056</b> Turquoise blue (RAL 5018)
<b>002</b> Silver	<b>021</b> Grey aluminium (RAL 9007)	<b>029</b> Pearl beige (RAL 1035)	<b>057</b> Iron sparkling dark light
<b>003</b> Pearl ruby red (RAL 3032)	<b>022</b> White (RAL 9003)	<b>032</b> Black dull matt (RAL 9005M)	<b>058</b> Salmon orange (RAL 2012)
<b>004</b> Pearl night blue (RAL 5026)	<b>023</b> Pure orange (RAL 2004)	<b>039</b> Signal violet (RAL 4008)	<b>059</b> Pastel blue (RAL 5024)
<b>005</b> Pearl gentiant blue (RAL 5025)	<b>024</b> Sparkling iron effect	<b>040</b> Night blue (RAL 5022)	<b>061</b> Light pink (RAL 3015)
<b>017</b> Traffic yellow (RAL 1023)	<b>025</b> Signal blue (RAL 5005)	<b>050</b> White matt (RAL 9003)	<b>064</b> Olive green (RAL 6003)
<b>018</b> Flame red (RAL 3000)	<b>026</b> Vermillion (RAL 2002)	<b>053</b> Light green (RAL 6027)	<b>065</b> Raspberry red (RAL 3027)
<b>019</b> Sky blue (RAL 5015)	<b>027</b> Heather violet (RAL 4003)	<b>054</b> Ocean blue (RAL 5020)	<b>066</b> Pastel green (RAL 6019)

\* The following frame colors are in stock and available immediately:

Pearl ruby red for size 1, Pearl night blue for size 2, Pearl gentiant blue for size 3, Grey aluminium for size 4

# CONFIGURATIONS // ACCESSORIES // FLUX



01 Grab rail, standard



02 Grab rail with universal handle



03 Reverse roll lock



04 Polyurethan front wheel with 6-spoke rims



Polyurethan rear wheel with 6-spoke rims



05 Polyurethan front wheel with 7-spoke rims



06 Polyurethan rear wheel with 7-spoke rims



07 Lower arm support



08 Lower arm support



09 Hand disk



10 Bumper wheels



11 Drag brake



12 Brake with wheel lock



13 Anti-tips, fits to size 1 - 3



14 Anti-tips with castors, fits to size 1 - 3



15 Anti-tips, fits to size 4



16 Anti-tips with castors fits to size 4



17 Seat



18 Seat folds down



19 Hip-/bottom pads



20 Lap belt for Hip-/bottom pads



21 Dynamic hip control



22 Dynamic hip control, suspension soft



23 Hip belt for dynamic hip control



24 Carrier



25 Basket



26 Swivel lock



27 Sling seat



28 Sling seat rail



Trunk support

## CONFIGURATIONS // ACCESSORIES // FLUX

29



Attendant-push bar