

# FLUX

Release date: April 2026 // valid until revoked

### Order

### Offer

Order number / reference: \_\_\_\_\_

Date: \_\_\_\_\_



### Invoice address

Customer number: \_\_\_\_\_

Company: \_\_\_\_\_

Street: \_\_\_\_\_

City, Zip code: \_\_\_\_\_

### Delivery address

see invoice address

Company: \_\_\_\_\_

Street: \_\_\_\_\_

City, Zip code: \_\_\_\_\_

### Basic configuration

- Grab rail height- and depth adjustable
- Front swivel wheels
- Polyurethan tires
- Engageable and disengageable free-wheeling system on rear wheels

## Frame colour



\*The following frame colors are in stock and available immediately:

Perl ruby for size 1, Perl night blue for size 2, Perl gentian blue for size 3, Grey aluminium for size 4

## Special frame colour



Copper gold



Pastel violet, matt



Pastel blue, matt



Pastel green, matt

Desired RAL colour:

## Technical data

		size 1	size 2	size 3	size 4
<b>Grab rail standard</b>	Grip height	48 - 61	57 - 81	69 - 87	82 - 102
	Grip width	41	46	51	51
<b>Grab rail with universal grips</b>	Grip height	43 - 68	51 - 88	62 - 94	75 - 102
	Grip width	29 - 41	34 - 46	39 - 51	39 - 51
<b>Frame</b>	Total width	64	69	74	78
	Load capacity	60 kg	75 kg	75 kg	100 kg
<b>Wheel diameter</b>	Front / rear	15 / 20	15 / 20	15 / 20	20 / 30



FLUX size 1



FLUX size 2




FLUX size 3



FLUX size 4


## Grab rail



Grab rail standard



Grab rail with universal handle



FLUX attendant-push bar

Lower arm support




Lower arm support with handles

14 x 9,4 cm  
22 x 10,4 cm



Lower arm support with handles, 3-axial adjustable


14 x 9,4 cm  
22 x 10,4 cm



Rotatable lower arm support with angle and depth adjustable handles, Ø 30 mm



Rotatable lower arm support with angle and depth adjustable handles, Ø 18 mm



Rotatable lower arm support with depth adjustable 40 mm ball handle

**Notice:**  
Cannot be combined with the driving and parking brake.

Hand disk, Bumper wheels, Drag brake, Brake with wheel lock



Hand disk



Bumper wheels



Drag brake



Brake with wheel lock

**Notice:**  
Drag brake and brake with wheel lock cannot be combined on FLUX size 4.

Reverse roll lock



Reverse roll lock



Low friction reverse roll lock

**!** Notice:  
Only possible for FLUX sizes 1 & 2.

Tires (Only selectable for FLUX size 4); Swivel lock



Polyurethan tires with 6-spoke rims



Polyurethan wheels with 7-spoke rims



Swivel lock

Anti-tips



Anti-tips, fits to size 1 - 3



Anti-tips with castors, fits to size 1 - 3



Anti-tips, fits to size 4



Anti-tips with castors, fits to size 4

Seat, Hip-/bottom pads



Seat

seat height size 1 37 cm  
 seat height size 2 43 cm  
 seat height size 3 58,5 cm  
 seat height size 4 57,5 cm



Seat folds down

seat height size 1 37 cm  
 seat height size 2 43 cm




Hip-/bottom pads, multi adjustable



Lap belt for Hip-/bottom pads

**!** Notice:  
 Only possible for FLUX sizes 1 & 2.

Dynamic hip control kit, Hip belt



Dynamic hip control kit

for orthoses  
 suspension soft  
 suspension hard



Hip belt for dynamic hip control

53 - 83 cm  
 65 x 115 cm



Trunk support, multi adjustable

54 - 78 cm  
 64 - 92 cm  
 74 - 105 cm

**!** Notice:  
 The hip belt is not included in the scope of delivery and must be selected separately.

Carrier, Basket, Sling seat rail, Sling seat



Carrier



Basket



Sling seat rail



Sling seat

**!** **Notice:**  
Only possible in combination with anti-tip.

**!** **Notice:**  
Only possible in combination with carrier.

**hip size:**  
**57 - 77 cm**

**hip size:**  
**73 - 103 cm**

**!** **Notice:**  
Only possible in combination with sling seat rail.

Remark / Customer request / Individual special custom-made request

All prices are subject to change and non-binding. Errors, misprints, and price changes are reserved.  
Binding prices are only those stated in our written offer or order confirmation.