



SUPRA 2.0 Abdu Active-wheelchair for kids & adolescents

Crash tested according to ISO 7176-19

Basic configuration:

- Foot rest standard
- Adjustable backrest with adjustable back cover
- Rear wheels **HOGGI** light
- **HOGGI** light sports brake
- Seat cover

Price:

Order form incl. Exportprices (€) Release date 04/2024, valid until revoked

☐ Offer	☐ Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



Basic configuration

- Powder coated aluminium rigid frame active wheelchair with 6 cm abduction growing seat depth
- 3 size variations:
 - Size 1: SW 24-34cm in 2cm-steps, SD 23-39 cm
 - Size 2: SW 28-38cm in 2cm-steps, SD 24-44 cm
 - Size 3: SW 28-38cm in 2cm-steps, SD 29-49 cm
- Rear wheels with hollow chamber rims and low-resistance tires incl. handrims, removable via Quick Release Axis
- Wheel camber optional 3°, 6° or 9°
- Backrest angle adjustable from 70° to 115°
- Adjustable back cover or back panel
- Seat height front: 36 cm to 51 cm
- Seat angle: 0° to 10°
- Convertible to other seat widths
- Foot rest: Standard
- Maximum load: 60 kg size 1, 75 kg size 2, 100 kg size 3

Minimum information (Without this information your order can not be handled.)

Frame colour									
Seat width	in cm	☐ 24	☐ 26	☐ 28	☐ 30	☐ 32	☐ 34	☐ 36	☐ 38
Camber		☐ 3° (3282-1220)		☐ 6° (3282-1230)		☐ 9° (3282-1240)			

Presettings

Seat angle	0° - 10°	Seat height *	(36 - 51 cm)
Seat depth (SD) **	(23 - 49 cm)	Lower Leg Length (LLL) ***	(16 - 46 cm)
Backrest angle	☐ 70° ☐ 75° ☐ 80° ☐ 85° ☐ 90° ☐ 95° ☐ 100° ☐ 105° ☐ 110° ☐ 115°		
Backrest height	(20 - 35 cm) partly only with telescopic tubes/single push handle	AD	(0 - 10 cm)

* see table "Seat height"

** Possible SD frame size 1: 19 - 34 cm; size 2: 20 - 39 cm; size 3: 25 - 44 cm, see table "Seat depth"

*** Possible LLL frame size see table „Lower leg length“

At missing specifications, we will adjust to:

SD = SW, seat angle = 4°, seat height = mid dimension, LLL = seat depth, back height = minimum, back angle = 90°, active degree = 3
rear wheel: **HOGGI** light with aluminium design rim, push rim mounting "standard", front wheel: wheel with rubber tire

Basic variation

with/or: Please select only one version! SW=Seat width

	Orderno.	Colour	Description	Qty.	Price	€
	☐ 3284-00SW		SUPRA 2.0 Abdu size 1, SD 23 - 39 cm	1 pc.		
with	☐ 3284-22SW	320	Seat cover size 1	01	1 set	incl.
or	☐ 3281-23SW	052	Seat panel, size 1	02	1 set	incl.
	☐ 3285-00SW		SUPRA 2.0 Abdu size 2, SD 24 - 44 cm	1 pc.		
with	☐ 3285-22SW	320	Seat cover size 2	01	1 set	incl.
or	☐ 3282-23SW	052	Seat panel, size 2	02	1 set	incl.
	☐ 3286-00SW		SUPRA 2.0 Abdu size 3, SD 29 - 49 cm	1 pc.		
with	☐ 3286-22SW	320	Seat cover size 3	01	1 set	incl.
or	☐ 3283-23SW	052	Seat panel, size 3	02	1 set	incl.
or	☐ 3281-2299		Without seat cover		1 set	
with	☐ 3281-35SW	000	Back rest frame, back height 20 cm	03	1 set	incl.
with	☐ 3281-36SW	000	Back rest frame, back height 25 cm		1 set	incl.
with	☐ 3281-37SW	000	Back rest frame, back height 30 cm		1 set	incl.
with	☐ 3281-38SW	000	Back rest frame, back height 35 cm		1 set	incl.
with	☐ 3281-24SW	320	Adaptable back cover back height 20 - 25 cm	04	1 set	incl.
or	☐ 3281-25SW	320	Adaptable back cover back height 25 - 35 cm		1 set	incl.
or	☐ 3281-26SW	320	Adaptable back cover back height 30 - 40 cm		1 set	incl.
or	☐ 3281-27SW	320	Adaptable back cover back height 35 - 40 cm		1 set	incl.
			* Height adjustment only with 3251-2175 or 3261-5102			
or	☐ 3282-52SB	000	Back panel in frame colour, back height 25 (max. SW 36)	05	1 set	incl.
or	☐ 3282-53SB	000	Back panel in frame colour, back height 30		1 set	incl.
or	☐ 3282-54SB	000	Back panel in frame colour, back height 35		1 set	incl.
or	☐ 3282-55SB	000	Back panel in frame colour, back height 40 (with back rest frame BH 35)		1 set	incl.
or	☐ 3281-2499		Without adaptable back cover or back panel		1 set	
with	☐ 3251-2175		Telescopic tubes for back rest height adjustment (from BH 25)	06	1 set	
or	☐ 3261-5102		Single push handle for back rest height adjustment (from BH 25)	07	1 set	
			* single push handles are plugged into the back rest frame tubes			
or	☐ 3281-5201		Standard single push handles 90°, size 1	08	1 set	
or	☐ 3281-5202		Standard single push handles 90°, size 2		1 set	
or	☐ 3281-55SW		Push bar SUPRA 2.0	09	1 pc.	
or	☐ 3281-53SW		Push-handles SUPRA 2.0	10	1 set	
or	☐ 3281-5200		Push handles displaced to the back, height adjustable	11	1 set	
or	☐ 3281-5203		Single push handles, foldable	12	1 pc.	
or	☐ 3275-5502		One hand push-handle delta form (max. SW 36)	13	1 pc.	
with	☐ 3281-3401		Back rest fixed	14	1 set	incl.
or	☐ 3281-34SW		Back rest frame with additional relax position (15°) Retractable and lockable in transport position	15	1 set	
with	☐ 3281-7300		Tie down kit according to ISO 7176-19 *	16	1 set	
with	☐ 3274-7301		Back rest extension (only with back panel)	17	1 set	
with	☐ 3274-7370	346	Cushion for back rest extension, black check	18	1 pc.	
			* If there is no shoulder and head support mounted in the handicapped transport vehicle, SUPRA 2.0 needs a back panel and back extension.			
with	☐ 3284-47SW		Tube frame foot rest	19	1 pc.	incl.
with	☐ 3284-41SW		Footrest one piece aluminium footplate, angle adjustable	20	1 set	
with	☐ 3284-40SW		Footrest one piece aluminium footplate, fold backwards	21	1 set	
or	☐ 3284-46SW		Footrest one piece aluminium footplate, angle adjustable, for short LLL reduced by 11 cm	22	1 set	
or	☐ 3284-48SW		Footrest one piece aluminium footplate, angle adjustable, for short LLL reduced by 7 cm	23	1 set	
or	☐ 3284-43SW		Footrest one piece aluminium footplate, angle adjustable, foldable	24	1 set	
or	☐ 3284-45SW		Single foot plate footrest, angle adjustable, foldable	24	1 pair	

Brakes, fender

with/or: Please select only one version!

	Orderno.	Colour Description	Qty.	Price	€
with ☐	3281-1750	HOGGI light sport brake	25 1 set		incl.
or ☐	3285-1700	Standard brakes SUPRA 2.0 Abdu	1 set		incl.
with ☐	3283-1610	Standard side guards suitable for frame size 1	1 pair		incl.
or ☐	3283-1630	Standard side guards suitable for frame size 2+3 (SH 38-45)	26 1 pair		incl.
or ☐	3283-1640	Standard side guards suitable for frame size 2+3 (SH 44-52)	27 1 pair		incl.
or ☐	3283-1650	Standard side guards with armrests suitable for frame size 1 (upcharge)	1 pair		
or ☐	3283-1670	Standard side guards with armrests suitable for frame size 2+3 (SH 38-45) (upcharge)	28 1 pair		
or ☐	3283-1680	Standard side guards with armrests suitable for frame size 2+3 (SH 44-52) (upcharge)	29 1 pair		
or ☐	3282-160x	Side guards plane in frame colour (0=20", 2=22", 4=24") (upcharge)	1 pair		
or ☐	3283-160x	Side guards plane in fr. colour with armrests (0=20", 2=22", 4=24") (upcharge)	30 1 pair		
with ☐	7000-0051	PU-Armrest pads, 235 mm long	1 pair		incl.
or ☐	7000-0052	PU-Armrest pads, 270 mm long	1 pair		incl.
or ☐	3282-165y-000	Wheel guards aluminum in frame colours, 20-24" mounting: ☐ fixed ☐ detachable (upcharge)	31 1 set		
or ☐	3282-166y-000	Wheel guards aluminum in frame colour with integrated wheel lock levers, mounting: fixed (upcharge)	32 1 set		
or ☐	3282-167y-000	Wheel guards aluminum in frame colour with height adjustable, detachable arm rests, 20-24" mounting: ☐ fixed ☐ detachable (upcharge)	33 1 set		
or ☐	3282-168y-000	Wheel guards aluminum in frame colour with integrated wheel lock levers, height adjustable, detachable arm rests, 20-24", mounting: fixed (upcharge)	34 1 set		
with ☐	7000-0051	PU-Armrest pads, 235 mm long	1 pair		incl.
or ☐	7000-0052	PU-Armrest pads, 270 mm long	1 pair		incl.
☐	3240-8202	Therapy tray size 2, SW 24 and 26 cm *	35 1 pc.		
☐	3240-8203	Therapy tray size 3, SW 28, 30 and 32 cm *	1 pc.		
☐	3240-8204	Therapy tray size 4, SW 34, 36 and 38 cm *	1 pc.		
☐	3240-8205	Therapy tray size 5, SW 40, 42 and 44 cm *	1 pc.		
* in combination with arm rests only					

Rear wheels

with/or: Please select only one version!

	Orderno.	Colour Description	Qty.	Price	€
with ☐	1293-0xxx	Rear wheel HOGGI light	1 set		incl.
or ☐	1293-2xxx	Rear wheel HOGGI light with drum brakes* (upcharge)	1 set		
*in combination with push-bar or push-handles only					
with ☐	1293-x2xx	20" wheel diameter	1 set		incl.
or ☐	1293-x3xx	22" wheel diameter	1 set		incl.
or ☐	1293-x4xx	24" wheel diameter	1 set		incl.
with ☐	1293-xx1x	Pneumatic tires 1"	36 1 set		incl.
or ☐	1293-xx2x	Flat free tires (20x1 3/8", 22x1" and 24x1")	37 1 set		incl.
with ☐	1293-xxx0	Without push-rims			
or ☐	1293-xxx1	Push rims aluminium anodised, "standard" mounting	38 1 set		incl.
or ☐	1293-xxx2	Push rims aluminium anodised, "high" mounting	39 1 set		incl.
or ☐	1293-xxx3	Push rims stainless steel, "standard" mounting (upcharge)	40 1 set		
or ☐	1293-xxx4	Push rims stainless steel, "high" mounting (upcharge)	41 1 set		
with ☐		☐ Push-rim mounting "narrow"			
or ☐		☐ Push-rim mounting "wide"			
or ☐	1204-1155	Closed ergonomic hand rims 22 ", black anodized, for Hoggi-light wheels (Upcharge)	69 1 Set		
or ☐	1204-1156	Closed ergonomic hand rims 24 ", black anodized, for Hoggi-light wheels (Upcharge)	1 Set		
with ☐	1204-106X	"HOGGI Grip Plus" hand rim inserts black, for closed ergonomic hand rims (Upcharge)	70 1 Set		

Rear wheels				with/or: Please select only one version!		
	Orderno.	Colour	Description		Qty.	Price €
☐	1204-114X		Aluminum hand rim with plastic anti-slip coating, black, standard mounting *	(Upcharge)	71 1 pair	
☐	1204-114X		Aluminum hand rim with plastic anti-slip coating, black, high mounting *	(Upcharge)	1 pair	
* only in combination with HOGGI light rear wheels						
or	☐	1296-4211	Pneumatic wheel 20x1" HOGGI Extra light *	(upcharge)	42 1 set	
or	☐	1296-4311	Pneumatic wheel 22x1" HOGGI Extra light *	(upcharge)	1 set	
or	☐	1296-4411	Pneumatic wheel 24x1" HOGGI Extra light *	(upcharge)	1 set	
* incl. black alu push-rim						

Spoke guards // Old motifs // Sale

	Order no.	Colour	Description		Qty.	Price €
☐	1204-101x		Spoke guard, transparent *		43 1 set	
☐	1204-101x		Spoke guard, transparent, for drum brakes wheels *		1 set	
☐	1204-101x		Spoke guard, printed **	(Sale)	44 1 set	
☐	1204-101x		Spoke guard, printed, for drum brakes **	(Sale)	1 set	
* Motifs 2021 only possible in combination with 1204-2001-XXX						
** please enter code of your print design see attachment or brochure "spoke guard" // Only while stocks last!						
*** printed spoke guards are not available for 20" wheels						

Spoke guards // New motifs 2021 (with 10% dealer discount)

	Order no.	Colour	Description		Qty.	Price €
☐	1204-2001-000		Spoke guard Special interior foiling, single-coloured according to customer requirements ^{1 2}		1 set	
☐	1204-2001-2XX		Spoke guard Special interior foiling, printed ^{1 3}		1 set	
☐	1204-2001-SA		Spoke guard Special interior foiling, printed according to customer requirements ¹		Price on request	
¹ only in combination with transparent spoke guard						
² please specify RAL tone						
³ Spoke guard motif selectable in "Order form spoke-guards" (available on the website)						

Front castors				with/or: Please select only one version!		
	Orderno.	Colour	Description		Qty.	Price €
with	☐	1250-1002	Wheel with rubber tire, black 100 x 32 mm		45 1 set	incl.
or	☐	1250-1003	Wheel with rubber tire, black 125 x 32 mm		1 set	incl.
or	☐	1250-1005	Wheel with rubber tire, black 140 x 40 mm		1 set	incl.
or	☐	1249-1005	Skater wheel with LED 100 x 24 mm	(upcharge)	46 1 set	
or	☐	1249-1006	Skater wheel with LED 125 x 24 mm	(upcharge)	1 set	
or	☐	1221-1003	PU-wheel with plastic design rim 100 x 34 mm	(upcharge)	47 1 set	
or	☐	1221-1004	PU-wheel with plastic design rim 125 x 34 mm	(upcharge)	1 set	

Accessories

SW=Seat width

	Orderno.	Colour	Description		Qty.	Price €
☐	3281-880x		Calf strap SUPRA 2.0 with upholstery, black		48 1 set	
☐	3281-1901		SUPRA 2.0 Anti tip 20", swing away, left		49 1 pc.	
☐	3281-1902		SUPRA 2.0 Anti tip 20", swing away, right		1 pc.	
☐	3282-1901		SUPRA 2.0 Anti tip 22 + 24", swing away, left		1 pc.	
☐	3282-1902		SUPRA 2.0 Anti tip 22 + 24", swing away, right		1 pc.	
☐	3281-190x		Tip assist, left ☐ =3 or right ☐ =4		50 1 pc.	
☐	3251-71SW	346	Seat cushion; 3 cm thick, size 1		51 1 pc.	
☐	3252-71SW	346	Seat cushion; 3 cm thick, size 2		1 pc.	
☐	3253-71SW	346	Seat cushion; 3 cm thick, size 3		1 pc.	
☐	3276-7100	346	Back cushion for back panle, 2,5 cm thick (SB und RH)		52 1 pc.	

Accessories				SW=Seat width		
	Orderno.	Colour	Description	Qty.	Price	€
with ☐	2240-1505		MDP pad holder swing away and detachable for attachment to round tube D25 mm, left side	53	1 pc.	
with ☐	2240-1506		MDP pad holder swing away and detachable for attachment to round tube D25 mm, right side		1 pc.	
or ☐	2240-1405		MDP pad holder fixed for attachment to round tube D25 mm, left/right side	54	1 pc.	
with ☐	2240-1461	346	MDP trunk support size 1 (8 x 10 x 2 cm), cover black	55	1 pc.	
or ☐	2240-1462	346	MDP trunk support size 2 (10 x 12 x 2 cm), cover black		1 pc.	
or ☐	2240-1463	346	MDP trunk support size 3 (12 x 13 x 2 cm), cover black		1 pc.	
or ☐	2240-1471		MDP trunk support curved size 1 (6 x 11,5 x 2,5 cm), PU w/o cover	56	1 pc.	
or ☐	2240-1472		MDP trunk support curved size 2 (8 x 14 x 2,5 cm), PU w/o cover		1 pc.	
☐	3201-3830		Universal headrest bracket (for 3245-7800 or third-party products) *	57	1 pc.	
☐	3245-7800	306	Headrest with occiput upholstery	58	1 pc.	
* starting from backrest height 30 cm						
☐	3200-0401		Belt routing as a T-strap above the shoulders	59	1 pc.	
☐	3251-77SW		Nylon back cover to cover the back upholstery belts	60	1 pc.	
☐	3251-7600		Edge protector for back rest tubes and telescopic tubes	61	1 set	
☐	1204-101x		Push-rim cover, silicone black fluted	62	1 set	
☐	3281-6100		Lap belt	63	1 pc.	
☐	3201-6400		Ankle hugger, size 1	64	1 pair	
☐	3203-6400		Ankle hugger, size 2		1 pair	
☐	3220-6702		4-Point lap belt, size S (SW 24 - 30)	65	1 pc.	
☐	3220-6703		4-Point lap belt, size M (SW 28 - 36)		1 pc.	
☐	3220-6704		4-Point lap belt, size L (SW 34 - 42)		1 pc.	
☐	3231-7800	346	Abduction block size 1, swing away, detachable, (9 x 6 x 9 cm) *	66	1 set	
☐	3231-7900	346	Abduction block size 2 swing away, detachable, (11 x 8 x 10 cm) *		1 set	
☐	3233-7800	346	Abduction block size 3 swing away, detachable, (12 x 9 x 10 cm) *		1 set	
* only with seat panel 328X-23SW						

Custom-made				SW=Seat width		
	Orderno.	Colour	Description	Qty.	Price	€
☐	3281-80SW	346	Contoured seat cushion with cover, black, SW 24 - 32	67	1 pc.	
☐	3281-80SW	346	Contoured seat cushion with cover, black, SW 34 - 44		1 pc.	
☐	3275-7200	346	"HOGGI align" Contoured back cushion, black	68	1 pc.	
☐	SA-0286/2-KP		Abduction block 30 mm height extension kit *		1 set	
* only in combination with abduction block and contoured seat cushion						
☐	6999-0014	060	Special Frame colour, Copper Gold (Upcharge)		1 pc.	
☐	6999-0021/1		Castor forks powder coated according to customer requirements * (no glitter) (Upcharge)		1 set	
☐	6999-0021/2		Back retainer plates powder coated, color according to customer specification (Upcharge)		1 set	
☐	6999-0021/3		Alu-wheel guards powder coated, color according to customer specification (Upcharge)		1 set	
☐	SA-0390		Upholstery in special colour *		Price on request	
* please enter colour code or upholstery code respectively						
☐	SA-0572-F-MP		Plastic anti-slip coating in custom colour: Pure white (RAL 9010), Silver grey (RAL 7001), Sky blue (RAL 5015), Traffic red (RAL 3020), Pure orange (RAL 2004)* (Upcharge)	72	1 pair	
* only in combination with HOGGI light rear wheels						
☐	SA-0093-MP		Maxgrepp Original cell rubber coated push rims (Upcharge)	73	1 set	
☐	SA-0094-MP		Maxgrepp Ergo Para (closed push rims) (Upcharge)	74	1 set	
☐	SA-0095-MP		Maxgrepp Ergo Para CP (advanced closed push rims) (Upcharge)	75	1 set	
☐	SA-0210		Special Frame colour, RAL code only (Upcharge)		1 pc.	
☐	SA-0505-MP		Locking of the integrated brake	76	1 pc.	
☐	SA-0594		Carbon wheel cover (Upcharge)	77	Price on request	

Upholstery colours

346	Black	321	Spacer fabric, black (special colour)	348	Grey Soft (special colour)
		347	Berry Soft (special colour)	349	Turquoise Soft (special colour)


All HOGGI upholstery are made from flame retardant Trevira CS. Exception spacer fabric (321).

Restriction of the footrest systems

SUPRA 2.0 Footrest systems	Restriction footrest systems*			
	Minimum seat angle			Front castors
	1°	4°	6°	125 140
Tube frame footrest (standard)	x			x
Footrest one piece aluminium footplate, angle adjustable	x			x
Footrest one piece aluminium footplate, fold backwards	x			x
Footrest one piece aluminium footplate, angle adjustable, for short LLL reduced by 11 cm		x		x
Footrest one piece aluminium footplate, angle adjustable, for short LLL reduced by 7 cm		x		x
Footrest one piece aluminium footplate, angle adjustable, foldable		x		x
Single foot plate footrest, angle adjustable, foldable	x			x
Single foot plate footrest, angle adjustable, foldable			x	x

* concerns SUPRA 2.0 (SW 30-44) and SUPRA 2.0 Abdu

How to measure and active degree



SH = Seat height LLL = Lower leg length SD = Seat depth BH = Back height SW = Seat width AD = Active degree
(All measurements are without seat or back upholstery)

Frame colours

001 Black (RAL 9005)	020 Yellow green (RAL 6018)	028 Pearl violet (RAL 4011)	056 Turquoise blue (RAL 5018)
002 Silver	021 Grey aluminium (RAL 9007)	029 Pearl beige (RAL 1035)	057 Iron sparkling dark light
003 Pearl ruby red (RAL 3032)	022 White (RAL 9003)	032 Black dull matt (RAL 9005M)	058 Salmon orange (RAL 2012)
004 Pearl night blue (RAL 5026)	023 Pure orange (RAL 2004)	039 Signal violet (RAL 4008)	059 Pastel blue (RAL 5024)
005 Pearl gentian blue (RAL 5025)	024 Sparkling iron effect	040 Night blue (RAL 5022)	061 Light pink (RAL 3015)
017 Traffic yellow (RAL 1023)	025 Signal blue (RAL 5005)	050 White matt (RAL 9003)	064 Olive green (RAL 6003)
018 Flame red (RAL 3000)	026 Vermillion (RAL 2002)	053 Light green (RAL 6027)	065 Raspberry red (RAL 3027)
019 Sky blue (RAL 5015)	027 Heather violet (RAL 4003)	054 Ocean blue (RAL 5020)	066 Pastel green (RAL 6019)

SA-0210-MP Special colours on request (Enter RAL code , Upcharge

Colour code	
201	HOGGI blue
202	Wonderland
203	Mouse
204	Water world
205	Train
206	Dragon
207	Circus



CONFIGURATIONS // ACCESSORIES // SUPRA 2.0 Abdu



Seat cover




Seat panel



Back rest frame



Adaptable back cover



Back panel



Telescopic tubes



Single push handle



Single push handles 90°



Push bar



Push-handles




Push handles
displaced to the back



Single push handles,
foldable



One hand push-handle
delta form




Back rest fixed




Back rest, retractable



Tie down kit



Back rest extension



Cushion for
back rest extension



Tube frame footrest



Footrest one piece,
angle adjustable



Footrest one piece,
fold backwards



Footrest, for short LLL



Footrest,
foldable




Single foot plate footrest




HOGGI light sport brake



Side guards plane
with armrests



Standard side guards



Standard side guards




Standard side guards
with armrests



Standard side guards
with armrests

CONFIGURATIONS // ACCESSORIES // SUPRA 2.0 Abdu



31

Wheel guards aluminum



32

Wheel guards with
wheel lock levers




33

Wheel guards, height
adjustable, detachable



34

Wheel guards, arm rests,
wheel lock levers



35

Therapy tray



36

Pneumatic tires



Pneumatic tires
with drum brakes



37

Flat free tires



Flat free tires
with drum brakes



38

Push rims, standard



39

Push rims, high



40

Push rims stainless steel,
standard



41

Push rims stainless steel,
high



42

HOGGI Extra light



43

Spoke guard, transparent



44

Spoke guard, printed



45

Wheel with rubber tire



46

Skater wheel



47

PU-wheel with plastic
design rim



48

Calf strap



49

Anti tip



50

Tip assist



51

Seat cushion



52

Back cushion



53

MDP pad holder swing
away and detachable



54

MDP pad holder fixed



55

MDP trunk support



56

MDP trunk support
curved



57


Universal headrest
bracket



58

Headrest with occiput
upholstery

CONFIGURATIONS // ACCESSORIES // SUPRA 2.0 Abdu



Belt routing as a T-strap



Nylon back cover



Edge protector



Push-rim cover



Lap belt



Ankle hugger



4-Point lap belt



Abduction block



Contoured seat cushion with cover



„HOGGI align“ Contoured back cushion



Closed ergonomic hand rims



HOGGI Grip Plus



Aluminum hand rim with plastic anti-slip coating



Aluminum hand rim with plastic anti-slip coating




Maxgrepp Original push rims



Maxgrepp Ergo Para push rims



Maxgrepp Ergo Para CP push rims



Locking of the integrated brake



Carbon fender