

www.hoggi.de

Made in Germany  
by HOGGI 

---

# SHARKY

Multifunctional bath- and toilet wheelchair.



Designed and hand built in Germany.

# SHARKY

Product description.

## SHARKY - Multifunctional bath- and toilet wheelchair.

In the development of **SHARKY**, the main focus has been on multifunctional applications and practical handling.

Unique in this product segment is its easy foldability to a very compact folding size, with a height of just 28 cm, while the seat is attached!

The tilt in space mechanism with locking bolts guarantees a safe 45° tilt in space function - perfect for showering.

The extremely slim design of the removable seat makes it very easy to use in the bathtub near the tub floor. The effect: Low water consumption and quickly reached water level.

The easy and comfortable height adjustment works with two stainless steel gas pressure springs.



**SHARKY:** Backrest recline



**SHARKY:** 45° tilt

# SHARKY

Technical data.

Specifications	Size 1	Size 2
Seat width	24, 27, 30 cm	29, 32, 35 cm
Seat depth	27 - 36 cm	31 - 40 cm
Backrest height	45 cm	55 cm
Maximum traversable toilet height	53 cm	53 cm
Seat angle	-10° to 45°	-10° to 45°
Backrest angle	72° - 120°	72° - 120°
Knee angle	78° - 138°	78° - 138°
Fußwinkel	80° - 100°	80° - 100°
Armrest height	14,5 - 20,5 cm	14,5 - 20,5 cm
Wheel diameter	75, 100 mm	75, 100 mm
Weight limit	60 kg	60 kg



SHARKY: Side view



SHARKY: Very compact folding size

## SHARKY overview:

- 45° tilt for high safety
- Simple and practical use
- Very compact folding size (even with attached seat)
- Extremely slim, removable seat offers water-saving use in the bathtub
- Easy and comfortable height adjustment with two gas pressure springs

---

# SHARKY

Multifunctional bath- and toilet wheelchair.

## **HOGGI**

Partnerschaft für mobile Kinder.

---

### **HOGGI GmbH**

 Eulerstraße 27  
D-56235 Ransbach-Baumbach

 Telefon: (+49) 2623 / 92 499-0  
Telefax: (+49) 2623 / 92 499-99

 [info@hoggi.de](mailto:info@hoggi.de)

[www.hoggi.de](http://www.hoggi.de)