



## SHARKY Shower and toilet wheelchair for kids

Order form incl. Exportprices Release date 04/2024, valid until revoked

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<b>Invoice address</b>	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
<b>Delivery address</b>	
Company	
Street	City, Zip code



### Basic configuration

#### Mobility base

- Compact folding
- Height adjustable, toilet height for Ø75mm rollers 43-53 cm
- Tilt in space -10° to 45°

#### Seat unit

- removable from the base (quick release)
- PU padding on the back, seat and hip supports
- sealable hygiene opening
- Back is angle adjustable, foldable and removable

	size 1	size 2
- Seat width	24, 27, 30 cm	29, 32, 35 cm
- Seat depth	27-36 cm	31-40 cm
- Back height	45 cm	55 cm
- Back angle	72-120°	
- Maximum load	60 kg	

Frame colour: Signal white 022 (RAL 9003) PU-Pads colour: Sky blue 019 (similar RAL 5015)

### Presettings

Seat depth

## Basic variation

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> 3401-2000		<b>SHARKY Seat unit size 1</b>	<b>01</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3401-2161	019	PU-hip support pads for seat width 24 cm	<b>02</b>	1 pair	incl.
or	<input type="checkbox"/> 3401-2162	019	PU-hip support pads for seat width 27 cm	<b>02</b>	1 pair	incl.
or	<input type="checkbox"/> 3401-2163	019	PU-hip support pads for seat width 30 cm	<b>02</b>	1 pair	incl.
with	<input type="checkbox"/> 3401-2602	019	Back insert, long, BH 45 cm	<b>03</b>	1 pc.	incl.
or	<input type="checkbox"/> 3401-2603	019	Back insert, short, BH 35 cm	<b>04</b>	1 pc.	incl.
	<input type="checkbox"/> 3401-1000		<b>SHARKY mobility base size 1</b>	<b>05</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3401-1150		Double castors 75 mm	<b>06</b>	1 set	incl.
or	<input type="checkbox"/> 3401-1160		Double castors 100 mm	<b>(upcharge) 07</b>	1 set	

## Basic variation

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> 3402-2000		<b>SHARKY Seat unit size 2</b>	<b>01</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3401-2161	019	PU-hip support pads for seat width 29 cm	<b>02</b>	1 pair	incl.
or	<input type="checkbox"/> 3401-2162	019	PU-hip support pads for seat width 32 cm	<b>02</b>	1 pair	incl.
or	<input type="checkbox"/> 3401-2163	019	PU-hip support pads for seat width 35 cm	<b>02</b>	1 pair	incl.
with	<input type="checkbox"/> 3402-2602	019	Back insert, long, BH 55 cm	<b>03</b>	1 pc.	incl.
or	<input type="checkbox"/> 3402-2603	019	Back insert, short, BH 45 cm	<b>04</b>	1 pc.	incl.
	<input type="checkbox"/> 3402-1000		<b>SHARKY mobility base size 2</b>	<b>05</b>	<b>1 pc.</b>	
with	<input type="checkbox"/> 3401-1150		Double castors 75 mm	<b>06</b>	1 set	incl.
or	<input type="checkbox"/> 3401-1160		Double castors 100 mm	<b>(upcharge) 07</b>	1 set	

## Accessories

	Orderno.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> 3401-2360		Push handle unit	<b>08</b>	1 pc.	
	<input type="checkbox"/> 3401-2500	022	Footrest system size 1, angle adjustable, removable (LLL 22 - 40 cm)	<b>09</b>	1 set	
	<input type="checkbox"/> 3402-2500	022	Footrest system size 2, angle adjustable, removable (LLL 22 - 40 cm)	<b>09</b>	1 set	
	<input type="checkbox"/> 3401-2380	019	Back rest extension with PU-pad (not in combination with headrest)	<b>10</b>	1 Set	
	<input type="checkbox"/> 3401-7600	019	Headrest with PU-pad, multi-directionally adjustable, incl. retainer	<b>11</b>	1 set	
	<input type="checkbox"/> 3401-7700	019	Lateral PU pads (11 x 6 cm), width and height adjustable, swing away	<b>12</b>	1 pair	
	<input type="checkbox"/> 3402-7700	019	Lateral PU pads (14 x 8 cm), width and height adjustable, swing away		1 pair	
	<input type="checkbox"/> 3401-5601	019	Chest-shoulder harness neopren, size 1	<b>13</b>	1 pc.	
	<input type="checkbox"/> 3401-5602	019	Chest-shoulder harness neopren, size 2		1 pc.	
	<input type="checkbox"/> 3401-6200		Lap belt	<b>14</b>	1 pc.	
	<input type="checkbox"/> 3401-6301	019	Lap belt with neopren pads size 1 (4 x 13 cm)	<b>15</b>	1 pc.	
	<input type="checkbox"/> 3401-6302	019	Lap belt with neopren pads size 2 (4 x 17 cm)		1 pc.	
	<input type="checkbox"/> 3401-6401		A2J Pelvis-retraction belt 3D size XS (L 46 - 68 cm)		1 pc.	
	<input type="checkbox"/> 3401-6402		A2J Pelvic-retraction belt 3D size S (L 48 - 71 cm)	<b>16</b>	1 pc.	
	<input type="checkbox"/> 3401-6403		A2J Pelvic-retraction belt 3D size M (L 56 - 76 cm)		1 pc.	
	<input type="checkbox"/> 3401-6100	019	Armrests with PU pads, removable, height adjustable (14.5-20.5 cm)	<b>17</b>	1 pair	
	<input type="checkbox"/> 3401-8100		Therapy table size 1, depth adjustable, removable *	<b>18</b>	1 pc.	
	<input type="checkbox"/> 3402-8100		Therapy table size 2, depth adjustable, removable *	<b>18</b>	1 pc.	
			<b>* only in combination with 3401-6100-019 armrests with PU pads</b>			
	<input type="checkbox"/> 3401-7800	019	Abduction block size 1 (LWH 100x85x80 cm), swing away and removable	<b>19</b>	1 pc.	
	<input type="checkbox"/> 3402-7800	019	Abduction block size 2 (LWH 140x100x80 cm), swing away and removable	<b>19</b>	1 pc.	
	<input type="checkbox"/> 3401-6500	019	Calf support neopren size 1	<b>20</b>	1 pc.	
	<input type="checkbox"/> 3402-6500	019	Calf support neopren size 2	<b>20</b>	1 pc.	
	<input type="checkbox"/> 3401-6601	019	Ankle huggers neoprene size 1	<b>21</b>	1 pair	
	<input type="checkbox"/> 3401-6602	019	Ankle huggers neoprene size 2		1 pair	
	<input type="checkbox"/> 3401-5100		Toilet bucket, incl. retainer	<b>22</b>	1 pc.	
	<input type="checkbox"/> 3401-5110	019	Splash guard PU, pluggable	<b>23</b>	1 pc.	
	<input type="checkbox"/> 3401-2266	019	PU seat infill pad for hygiene opening	<b>24</b>	1 pc.	
	<input type="checkbox"/> 3401-3000		Bath subframe with seat angle adjustment, size 1	<b>25</b>	1 pc.	
	<input type="checkbox"/> 3402-3000		Bath subframe with seat angle adjustment, size 2	<b>25</b>	1 pc.	
	<input type="checkbox"/> 3401-3100		Toilet adapter (to mount on the toilet) size 1	<b>26</b>	1 pc.	
	<input type="checkbox"/> 3402-3100		Toilet adapter (to mount on the toilet) size 2	<b>26</b>	1 pc.	
	<input type="checkbox"/> 3401-2161E	019	PU-hip support pads for seat width 24/29 cm (replacement)	<b>02</b>	1 pair	
	<input type="checkbox"/> 3401-2162E	019	PU-hip support pads for seat width 27/32 cm (replacement)	<b>02</b>	1 pair	
	<input type="checkbox"/> 3401-2163E	019	PU-hip support pads for seat width 30/35 cm (replacement)	<b>02</b>	1 pair	

# CONFIGURATIONS // ACCESSORIES // SHARKY



01 Seat unit



02 PU-hip support pads




03 Back insert, long



04 Back insert, short



05 Mobility base



06 Double castors 75 mm




07 Double castors 100 mm



08 Push handle unit



09 Footrest system, angle adjustable



10 Back rest extension with PU-pad (without headrest)



11 Headrest



12 Lateral PU pads




13 Chest-shoulder harness



14 Lap belt



15 Lap belt with neopren pads



16 A2J Pelvic-retraction belt



17 Armrests



18 Therapy table



19 Abduction block



20 Calf support



21 Ankle huggers



22 Toilet bucket




23 Splash guard



24 PU seat infill pad for hygiene opening



25 Bath subframe



26 Toilet adapter