



SHARKY Shower and toilet wheelchair for kids

Order form incl. Exportprices Release date 08/2021, valid until 31.12.2021

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



Basic configuration

Mobility base

- Compact folding
- Height adjustable, toilet height 43-52 cm
- Tilt in space -10° to 45°

Seat unit

- removable from the base (quick release)
- PU padding on the back, seat and hip supports
- sealable hygiene opening
- Back is angle adjustable, foldable and removable

	size 1	size 2 in progress
- Seat width	24, 27, 30 cm	
- Seat depth	27-36 cm	
- Back height	45 cm	
- Back angle	72-120°	
- Maximum load	60 kg	

Frame colour: Signal white 022 (RAL 9003) PU-Pads colour: Sky blue 019 (similar RAL 5015)

Presettings

Seat depth



Basic variation

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> 3401-2000		SHARKY Seat unit size 1	01	1 pc.	
with	<input type="checkbox"/> 3401-2161	019	PU-hip support pads for seat width 24 cm	02	1 pair	incl.
or	<input type="checkbox"/> 3401-2162	019	PU-hip support pads for seat width 27 cm	02	1 pair	incl.
or	<input type="checkbox"/> 3401-2163	019	PU-hip support pads for seat width 30 cm	02	1 pair	incl.
	<input type="checkbox"/> 3401-1000		SHARKY mobility base size 1	03	1 pc.	
with	<input type="checkbox"/> 3401-1150		Double castors 75 mm	04	1 set	incl.
or	<input type="checkbox"/> 3401-1160		Double castors 100 mm	(upcharge) 05	1 set	

Accessories

	Orderno.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> 3401-2360		Push handle unit	06	1 pc.	
	<input type="checkbox"/> 3401-2400		Footrest system size 1, angle adjustable, removable (LLL 22 - 40 cm)	07	1 set	
	<input type="checkbox"/> 3401-7600	019	Headrest with PU-pad, multi-directionally adjustable, incl. retainer	08	1 set	
	<input type="checkbox"/> 3401-7700	019	Lateral PU pads (11 x 6 cm), width and height adjustable, swing away	09	1 pair	
	<input type="checkbox"/> 3402-7700	019	Lateral PU pads (14 x 8 cm), width and height adjustable, swing away		1 pair	
	<input type="checkbox"/> 3401-5601	019	Chest-shoulder harness neopren, size 1	10	1 pc.	
	<input type="checkbox"/> 3401-5602	019	Chest-shoulder harness neopren, size 2		1 pc.	
	<input type="checkbox"/> 3401-6200		Lap belt	11	1 pc.	
	<input type="checkbox"/> 3401-6301	019	Lap belt with neopren pads size 1 (4 x 13 cm)	12	1 pc.	
	<input type="checkbox"/> 3401-6302	019	Lap belt with neopren pads size 1 (4 x 17 cm)		1 pc.	
	<input type="checkbox"/> 3401-6402		A2J Pelvic-retraction belt 3D size "S" (SW 24 - 30)	13	1 pc.	
	<input type="checkbox"/> 3401-6403		A2J Pelvic-retraction belt 3D size "M" (SW 24 - 30)		1 pc.	
	<input type="checkbox"/> 3401-6100	019	Armrests with PU pads, removable, height adjustable (14.5-20.5 cm)	14	1 pair	
	<input type="checkbox"/> 3401-8100		Therapy table, depth adjustable, removable *	15	1 pc.	
			* available from 2nd quarter			
	<input type="checkbox"/> 3401-7800	019	Abduction block size 1 (LWH 100x85x80 cm), swing away and removable	16	1 pc.	
	<input type="checkbox"/> 3401-6500	019	Calf support neopren	17	1 pc.	
	<input type="checkbox"/> 3401-6601	019	Ankle huggers neoprene size 1	18	1 pair	
	<input type="checkbox"/> 3401-6602	019	Ankle huggers neoprene size 2		1 pair	
	<input type="checkbox"/> 3401-5100		Toilet bucket, incl. retainer	19	1 pc.	
	<input type="checkbox"/> 3401-5110	019	Splash guard PU, pluggable	20	1 pc.	
	<input type="checkbox"/> 3401-2266	019	PU seat infill pad for hygiene opening *	21	1 pc.	
			* available from 2nd quarter			
	<input type="checkbox"/> 3401-3000		Bath subframe with seat angle adjustment	22	1 pc.	
	<input type="checkbox"/> 3401-3100		Toilet adapter (to mount on the toilet)		1 pc.	
	<input type="checkbox"/> 3401-2161E	019	PU-hip support pads for seat width 24 cm (replacement)	02	1 pair	
	<input type="checkbox"/> 3401-2162E	019	PU-hip support pads for seat width 27 cm (replacement)	02	1 pair	
	<input type="checkbox"/> 3401-2163E	019	PU-hip support pads for seat width 30 cm (replacement)	02	1 pair	

CONFIGURATIONS // ACCESSORIES // SHARKY




01 Seat unit



02 PU-hip support pads



03 Mobility base



04 Double castors 75 mm



05 Double castors 100 mm



06 Push handle unit



07 Footrest system, angle adjustable



08 Headrest



09 Lateral PU pads



10 Chest-shoulder harness



11 Lap belt



12 Lap belt with neopren pads



13 A2J Pelvic-retraction belt



14 Armrests



15 Therapy table



16 Abduction block



17 Calf support



18 Ankle huggers



19 Toilet bucket



20 Splash guard



21 PU seat infill pad for hygiene opening



22 Bath subframe