



TIMO Therapy chair

HMV-Nr.: 26.11.05.004X TIMO Size 2

Order form incl. Exportprices (€), release date 01/2025, valid until revoked

☐ Offer	☐ Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



Basic configurations

Baseframe

- Tilt in space -10° to 45°
- Height adjustment from 26cm to 67cm
- 4 pieces transport rollers 75mm, fixable
- Frame parts selectable in two colors
- Maximum load 70kg

Seat unit

- Back base angle adjustable
- Armrests height, angle, width and depth adjustable, removable
- Operation for height adjustment and tilting on the push handle
- Aluminum footboard in frame color, foldable
- Seat cushion and cushion back base selectable in 4 colours
- Maximum load 60 kg

Adjustment ranges

- Seat width 23-35
- Seat depth 24-43
- Seat tilt (-)10° - 45°
- Back height short 40-50
- Back height medium 46-56
- Back height long 52-62
- Back tilt 80°-125°
- Lower leg length 18-45
- Armrest width 29, 32, 37, 40

Minimum information (Without this information your order can not be handled.)

- | | | |
|------------------------------------|---|---|
| Frame colour: | ☐ Signal white 022 (RAL 9003) | ☐ Black dull matt 032 |
| Upholstery single-coloured: | ☐ Black 346 | ☐ Berry soft 347 |
| with accessories colour black | ☐ Turquoise soft 349 | ☐ Grey soft 348 |

Presettings (At missing specifications we will adjust to medium sizes.)

- | | | | |
|---------------------|--|--------------------------|--|
| Seat depth: | <input type="text" value="24"/> (24 - 43 cm) | Lower leg length: | <input type="text" value="18"/> (18 - 45 cm) (20-47 with contoured seat cushion) |
| Back height: | <input type="text" value="40"/> (40 - 62 cm) | Armrest width: | <input type="text" value="29"/> (clear width 29, 32, 37, 40) |

Basic variations

with/or: Please select only one version!

	Article no.	Colour	Description	Picture	Qty.	Price
	☐ 4112-0000		TIMO therapy chair size Gr. 2		1 pc.	
with	☐ 7000-0051		Armrest pads 235mm	01	1 pair	incl.
or	☐ 7000-0052		Armrest pads 270mm		1 pair	incl.
	☐ 4111-6140		Clamping lever for height adjustment armrests	02	1 set	
	☐ 4111-6150		Star handle for height adjustment armrests		1 set	
with	☐ 4112-2400		Footrest standard (fixed to the seat frame) size 2	03	1 set	incl.
or	☐ 4112-2500		Removable footrest size 2	(upcharge) 04	1 set	
with	☐ 3245-4302		Foot rest bracket angle adjustable, flip up type	05	1 set	incl.
or	☐ 3245-4303		Foot rest bracket angle adjustable, with interlocking	(upcharge) 06	1 set	
	☐ 4112-2520		Footrest padding size 2	(upcharge) 07	1 pc.	
with	☐ 4111-5100		Push handle standard (fixed to the back base)	07	1 set	incl.
or	☐ 4111-5101		Removable push handle	(upcharge) 08	1 set	
with	☐ 4112-2380		Cushion for back base TIMO size 2	09	1 set	incl.
or	☐ 4112-2200		Cushion for back base size 2 with adjustable hip pads	(upcharge) 10	1 set	
with	☐ 4112-7100		Seat cushion standard TIMO size 2	11	1 set	incl.
or	☐ 4112-8001		Contoured seat cushion TIMO size 2, SD 24-28	(upcharge) 12	1 pc.	
or	☐ 4112-8002		Contoured seat cushion TIMO size 2, SD 29-33	(upcharge) 13	1 pc.	
or	☐ 4112-8003		Contoured seat cushion TIMO size 2, SD 34-38	(upcharge) 14	1 pc.	
or	☐ 4112-8004		Contoured seat cushion TIMO size 2, SD 39-43	(upcharge) 15	1 pc.	
	☐ 4112-7110		Incontinence cover for seat cushion standard size 2	(upcharge) 16	1 pc.	
	☐ SA-0313/2-MP		Additional A-Tex layer on the foam cushion as incontinence protection for contoured seat cushion	(upcharge) 17	1 pc.	
with	☐ 4111-1150		Double castors 75 mm	13	1 set	incl.
or	☐ 4111-1160		Double castors 100 mm	(upcharge) 14	1 set	

Accessories

	Article no.	Colour	Description		Qty.	Price
	☐ 4112-2601		Attachment back short (back height 41-51) size 2	15	1 set	
and	☐ 4112-2660		Cushion for back top short TIMO size 2	15	1 pc.	incl.
	☐ 4112-2602		Attachment back short (back height 47-57)	16	1 set	
and	☐ 4112-2670		Cushion for back top medium TIMO size 2	16	1 pc.	incl.
	☐ 4112-2603		Attachment back short (back height 53-63)	17	1 set	
and	☐ 4112-2680		Cushion for back top long TIMO size 2	17	1 pc.	incl.
	☐ 4111-2699		Back displacement of the attachable back by 1.5cm		1 set	incl.
	☐ 4112-2901		Upper arm guide incl. cushion size 2	18	1 pair	
	☐ 4111-2820		Universal headrest bracket asym. (for 3245-7800 or third-party products)	19	1 set	
	☐ 3245-7800	306	Head rest with mounting bracket	20	1 pc.	
	☐ 4111-6802		Thigh supports, wide- and depth-adjustable, abductible, incl. cushion size 2	21	1 pair	
	☐ 4111-6850		Clamping lever for adjustment thigh supports	22	1 set	
	☐ 4111-6860		Star handle for adjustment thigh supports		1 set	
	☐ 4112-6801		Thigh supports, wide- and depth-adjustable, abductible, incl. cushion size 3	21	1 pair	
	☐ 4111-6850		Clamping lever for adjustment thigh supports	22	1 set	
	☐ 4111-6860		Star handle for adjustment thigh supports		1 set	
	☐ 4111-5702		Fixation jacket size 2	23	1 pc.	
	☐ 4111-5703		Fixation jacket size 3		1 pc.	
	☐ 4111-5602		Chest-shoulder harness size 2	24	1 pc.	
	☐ 4111-5603		Chest-shoulder harness size 3		1 pc.	
	☐ 4111-7702	346	MDP trunk support, swing away, size 2	25	1 set	
	☐ 4111-7703	346	MDP trunk support, swing away, size 3		1 set	
	☐ 4111-7802	346	MDP trunk support, standard size 2	26	1 set	
	☐ 4111-7803	346	MDP trunk support, standard size 3		1 set	
	☐ 3201-6100		Hip belt	27	1 pc.	
	☐ 4111-6702		4-point lap belt size S (for seat width 22 - 30)	28	1 pc.	
	☐ 4111-6703		4-point lap belt size M (for seat width 28 - 36 cm)		1 pc.	
	☐ 3233-5300		Groin strap size 2	29	1 pc.	
	☐ 4111-6710E		Belt guide on the TIMO seat profiles with folding buckles 25 mm (retrofit)	30	1 pc.	
	☐ 4111-5610E		Shoulder strap guide with folding buckles 25 mm (retrofit)	31	1 pc.	

Accessories

	Article no.	Colour	Description		Qty.	Price
☐	3231-7900	346	Abduction block size 2 swing away and detachable		32 1 pc.	
☐	3231-7901	346	Abduction block size 2 fixed		33 1 pc.	
☐	3233-7800	346	Abduction block size 3 swing away and detachable		1 pc.	
☐	3233-7801	346	Abduction block size 3 fixed		1 pc.	
☐	3201-6300		Foot straps		34 1 pair	
☐	3203-6400		Ankle huger size 2		35 1 pair	
☐	3290-4538	015	Heel plate alu, black		36 1 pc.	
☐	4111-8200		Therapy table, transparent		37 1 pc.	
☐	4111-8201		Therapy table, anthracite		38 1 pc.	
☐	4111-8300		Wooden Therapy table		1 pc.	

Custom-made**(separate discount rate)**

	Article no.	Colour	Description		Qty.	Price
☐	SA-0210		Upcharge for custom colour	(upcharge)	1 pc.	

CONFIGURATIONS // ACCESSORIES // TIMO



Armrest pads



Clamping lever for height adjustment armrests



Footrest standard



Removable footrest



Foot rest bracket, flip up type



Foot rest bracket, with interlocking



Push handle standard



Removable push handle



Cushion for back base



Cushion for back base with adjustable hip pads



Seat cushion standard



Contoured seat cushion



Double castors 75 mm



Double castors 100 mm



Attachment back short



Back top medium



Attachment back long



Upper arm guide



Universal headrest bracket asym.



Head rest with mounting bracket



Thigh supports



Clamping lever for adjustment thigh supports



Fixation jacket



Chest-shoulder harness



MDP trunk support, swing away



MDP trunk support, standard



Hip belt



4-point lap belt



Groin strap



Belt guide with folding buckles

CONFIGURATIONS // ACCESSORIES // TIMO



Shoulder strap guide
with folding buckles



Abduction block swing
away and detachable



Abduction block fixed



Foot straps



Ankle hugger



Heel plate aluminium



Therapy table,
transparent



Therapy table,
anthracite