

Designed and  
made in Germany

# CRUISER



So small, so light, so agile.

[www.hoggi.de](http://www.hoggi.de)

# NEW: ANTERIOR WALKER

## // ACTIVE ON THE MOVE

The **HOGGI** product family (category: walking aids) grows with **CRUISER**.

**CRUISER** offers even the smallest beginners **safety** and **support** for their first walking attempts.

The compact **CRUISER** design also impresses with its **agility** and **low weight**.

The **newly designed brake** can be operated easily and with minimal effort, even by small children.



Thanks to the „**scissor folding mechanism**“, the **CRUISER** can be folded quickly and easily to a **very compact folding size**.

The included nylon seat cushion folds with it.

On the way, it serves as a place to sit to take a deep breath and recharge your batteries, so that the remaining meters can be mastered with fun and energy.

Accessories such as a push bar, underarm supports, a nylon storage bag and a horizontal grab rail enable a wide range for everyday use.

Thanks to larger wheels (option: outdoor wheel set), **CRUISER** can also be safely used outdoors.

**CRUISER** is available in all well-known **HOGGI** wheelchair colours.



Scan for detailed product information.

