




## CRUISER Anterior walking aid

Size 0 - 2

Order form incl. Exportprices (\$) valid until revoked, release date 01/2025

|                                |                                |
|--------------------------------|--------------------------------|
| <input type="checkbox"/> Offer | <input type="checkbox"/> Order |
| Order number / reference       | Date                           |
| <input type="text"/>           | <input type="text"/>           |
| <b>Invoice address</b>         |                                |
| Company                        | Customer number                |
| <input type="text"/>           | <input type="text"/>           |
| Street                         | City, Zip code                 |
| <input type="text"/>           | <input type="text"/>           |
| E-Mail                         | Telephone                      |
| <input type="text"/>           | <input type="text"/>           |
| Fax                            |                                |
| <input type="text"/>           |                                |
| <b>Delivery address</b>        |                                |
| Company                        |                                |
| <input type="text"/>           |                                |
| Street                         | City, Zip code                 |
| <input type="text"/>           | <input type="text"/>           |



### Basic configuration

**Compact foldable frame in very stable, lightweight aluminum construction**

- Handles individually adjustable in height and width
- Polyurethane tires
- Swivelling front wheels
- Folding seat cover

**Upholstery:**

All sizes: Black

**320**

### Minimum information (Without this information your order can not be handled.)

**Frame colour:**  Enter color code, color code and color selection see page 2

**Basic variation**

with/or: Please select only one version! SW=Seat width

|      | Order no.                          | Colour | Description                                | Qty.         | Price | €     |
|------|------------------------------------|--------|--|--------------|-------|-------|
|      | <input type="checkbox"/> 3305-0000 |        | <b>CRUISER Gr. 0</b>                       | <b>1 pc.</b> |       |       |
| with | <input type="checkbox"/> 3305-2401 |        | Wheel set 150/200                          | 1 set        |       | incl. |
| or   | <input type="checkbox"/> 3305-2402 |        | Wheel set 150/200 with swivel lock         | (Upcharge)   | 1 set |       |
| or   | <input type="checkbox"/> 3305-1901 |        | Outdoor wheel set 200/250                  | (Upcharge)   | 1 set |       |
| or   | <input type="checkbox"/> 3305-1902 |        | Outdoor wheel set 200/250 with swivel lock | (Upcharge)   | 1 set |       |
| with | <input type="checkbox"/> 3305-1370 | 320    | Seat cover size 0                          | 1 Stk        |       | incl. |
| or   | <input type="checkbox"/> 3305-1330 | 320    | Seat cover size 0 with back strap          | (Upcharge)   | 1 Stk |       |
|      | <input type="checkbox"/> 3306-0000 |        | <b>CRUISER size 1</b>                      | <b>1 pc.</b> |       |       |
| with | <input type="checkbox"/> 3305-2401 |        | Wheel set 150/200                          | 1 set        |       | incl. |
| or   | <input type="checkbox"/> 3305-2402 |        | Wheel set 150/200 with swivel lock         | (Upcharge)   | 1 set |       |
| or   | <input type="checkbox"/> 3305-1901 |        | Outdoor wheel set 200/250                  | (Upcharge)   | 1 set |       |
| or   | <input type="checkbox"/> 3305-1902 |        | Outdoor wheel set 200/250 with swivel lock | (Upcharge)   | 1 set |       |
| with | <input type="checkbox"/> 3306-1370 | 320    | Seat cover size 1                          | 1 Stk        |       | incl. |
| or   | <input type="checkbox"/> 3306-1330 | 320    | Seat cover size 1 with back strap          | (Upcharge)   | 1 Stk |       |
|      | <input type="checkbox"/> 3307-0000 |        | <b>CRUISER size 2</b>                      | <b>1 pc.</b> |       |       |
| with | <input type="checkbox"/> 3305-2401 |        | Wheel set 150/200                          | 1 set        |       | incl. |
| or   | <input type="checkbox"/> 3305-2402 |        | Wheel set 150/200 with swivel lock         | (Upcharge)   | 1 set |       |
| or   | <input type="checkbox"/> 3305-1901 |        | Outdoor wheel set 200/250                  | (Upcharge)   | 1 set |       |
| or   | <input type="checkbox"/> 3305-1902 |        | Outdoor wheel set 200/250 with swivel lock | (Upcharge)   | 1 set |       |
| with | <input type="checkbox"/> 3307-1370 | 320    | Seat cover size 2                          | 1 pc.        |       | incl. |
| or   | <input type="checkbox"/> 3307-1330 | 320    | Seat cover size 2 with back strap          | (Upcharge)   | 1 pc. |       |

**Accessories**

|  | Order no.                          | Colour | Description   | Qty.  | Price | € |
|--|------------------------------------|--------|---|-------|-------|---|
|  | <input type="checkbox"/> 3305-2560 |        | Lower arm support rotatable with angle and depth-adjustable 22 mm handle      | 1 set |       |   |
|  | <input type="checkbox"/> 3305-2570 |        | Lower arm support rotatable with angle and depth-adjustable 15 mm handle      | 1 set |       |   |
|  | <input type="checkbox"/> 3305-2580 |        | Lower arm support rotatable with central depth-adjustable 40 mm ball handle   | 1 set |       |   |
|  | <input type="checkbox"/> 3305-2100 |        | Drag brake  | 1 set |       |   |
|  | <input type="checkbox"/> 3305-1500 |        | Driving and parking brake, suitable for size 0 + 1 *                          | 1 set |       |   |
|  | <input type="checkbox"/> 3305-1501 |        | Driving and parking brake, suitable for size 0 + 1, Outdoor wheel set *       | 1 set |       |   |
|  | <input type="checkbox"/> 3307-1500 |        | Driving and parking brake, suitable for size 2 *                              | 1 set |       |   |
|  | <input type="checkbox"/> 3307-1501 |        | Driving and parking brake, suitable for size 2, Outdoor wheel set *           | 1 set |       |   |
|  |                                    |        | <b>* not compatible with the lower arm support with ball handle 3305-2580</b> |       |       |   |
|  | <input type="checkbox"/> 3305-2700 |        | Backstops for both rear wheels (size 0-2)                                     | 1 set |       |   |
|  | <input type="checkbox"/> 3305-2701 |        | Backstops for both outdoor rear wheels (size 0-2)                             | 1 set |       |   |
|  | <input type="checkbox"/> 3305-2200 | 320    | Storage bag size 0  | 1 pc. |       |   |
|  | <input type="checkbox"/> 3306-2200 | 320    | Storage bag size 1  | 1 pc. |       |   |
|  | <input type="checkbox"/> 3307-2200 | 320    | Storage bag size 2  | 1 pc. |       |   |
|  | <input type="checkbox"/> 3305-2600 |        | Push bar for accompanying person size 0-2, retainer in frame color            | 1 pc. |       |   |
|  | <input type="checkbox"/> 3305-2620 |        | Push bar for accompanying person size 0-2 Outdoor, retainer in frame color    | 1 pc. |       |   |

**Custom-made**

|  | Order no.                           | Colour | Description                           | Qty.       | Price | € |
|--|-------------------------------------|--------|---------------------------------------|------------|-------|---|
|  | <input type="checkbox"/> SA-0210-MP |        | Special Frame colour, RAL code only * | (Upcharge) | 1 pc. |   |
|  |                                     |        | <b>* please specify RAL color</b>     |            |       |   |
|  | <input type="checkbox"/> 6999-0014  | 060    | Special Frame colour, Copper Gold     | (Upcharge) | 1 pc. |   |

**Frame colours**

|  |                                      |  |                                      |
|--|--------------------------------------|--|--------------------------------------|
| <b>001</b> Black (RAL 9005)              | <b>020</b> Yellow green (RAL 6018)   | <b>028</b> Pearl violet (RAL 4011)     | <b>056</b> Turquoise blue (RAL 5018) |
| <b>002</b> Silver                        | <b>021</b> Grey aluminium (RAL 9007) | <b>029</b> Pearl beige (RAL 1035)      | <b>057</b> Iron sparkling dark light |
| <b>003</b> Pearl ruby red (RAL 3032)     | <b>022</b> White (RAL 9003)          | <b>032</b> Black dull matt (RAL 9005M) | <b>058</b> Salmon orange (RAL 2012)  |
| <b>004</b> Pearl night blue (RAL 5026)   | <b>023</b> Pure orange (RAL 2004)    | <b>039</b> Signal violet (RAL 4008)    | <b>059</b> Pastel blue (RAL 5024)    |
| <b>005</b> Pearl gentian blue (RAL 5025) | <b>024</b> Sparkling iron effect     | <b>040</b> Night blue (RAL 5022)       | <b>061</b> Light pink (RAL 3015)     |
| <b>017</b> Traffic yellow (RAL 1023)     | <b>025</b> Signal blue (RAL 5005)    | <b>050</b> White matt (RAL 9003)       | <b>064</b> Olive green (RAL 6003)    |
| <b>018</b> Flame red (RAL 3000)          | <b>026</b> Vermillion (RAL 2002)     | <b>053</b> Light green (RAL 6027)      | <b>065</b> Raspberry red (RAL 3027)  |
| <b>019</b> Sky blue (RAL 5015)           | <b>027</b> Heather violet (RAL 4003) | <b>054</b> Ocean blue (RAL 5020)       | <b>066</b> Pastel green (RAL 6019)   |

**Technical data** (dimensions in cm and weights in kg)

|   |                             | Size 0   |      | Size 1   |      | Size 2   |      |
|---|-----------------------------|----------|------|----------|------|----------|------|
|   |                             | Min.     | Max. | Min.     | Max. | Min.     | Max. |
| A | Handle height [cm]          | 32       | 59   | 32       | 59   | 54       | 83   |
| a | Height forearm support [cm] |          |      |          |      |          |      |
| B | Handle distance [cm]        | 18       | 32   | 22       | 36   | 36       | 50   |
| C | Seat height [cm]            | 25       |      | 33       |      | 45       |      |
| D | Total width [cm]            | 47       |      | 51       |      | 65       |      |
| E | Total length [cm]           | 55       |      | 55       |      | 58       |      |
| F | Folding size (HxLxD) [cm]   | 42x55x29 |      | 50x55x27 |      | 68x58x29 |      |
|   | max. load [kg]              | 30       |      | 50       |      | 65       |      |
|   | Weight [kg]                 | 5,4      |      | 5,6      |      | 5,9      |      |