



CRUISER Anterior walking aid

Size 0 - 2

Order form incl. Exportprices (\$) valid until revoked, release date 05/2024

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<input type="text"/>	<input type="text"/>
Invoice address	
Company	Customer number
<input type="text"/>	<input type="text"/>
Street	City, Zip code
<input type="text"/>	<input type="text"/>
E-Mail	Telephone
<input type="text"/>	<input type="text"/>
Fax	
<input type="text"/>	
Delivery address	
Company	
<input type="text"/>	
Street	City, Zip code
<input type="text"/>	<input type="text"/>



Basic configuration

Compact foldable frame in very stable, lightweight aluminum construction

- Handles individually adjustable in height and width
- Polyurethane tires
- Swivelling front wheels
- Folding seat cover

Upholstery:

All sizes: Black

320

Minimum information (Without this information your order can not be handled.)

Frame colour: Enter color code, color code and color selection see page 2

Basic variation

with/or: Please select only one version! SW=Seat width

	Order no.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> 3305-0000		CRUISER Gr. 0	1 pc.		
with	<input type="checkbox"/> 3305-2401		Wheel set 150/200	1 set		incl.
or	<input type="checkbox"/> 3305-2402		Wheel set 150/200 with swivel lock	(Upcharge)	1 set	
or	<input type="checkbox"/> 3305-1901		Outdoor wheel set 200/250	(Upcharge)	1 set	
or	<input type="checkbox"/> 3305-1902		Outdoor wheel set 200/250 with swivel lock	(Upcharge)	1 set	
with	<input type="checkbox"/> 3305-1370	320	Seat cover size 0	1 Stk		incl.
or	<input type="checkbox"/> 3305-1330	320	Seat cover size 0 with back strap	(Upcharge)	1 Stk	
	<input type="checkbox"/> 3306-0000		CRUISER size 1	1 pc.		
with	<input type="checkbox"/> 3305-2401		Wheel set 150/200	1 set		incl.
or	<input type="checkbox"/> 3305-2402		Wheel set 150/200 with swivel lock	(Upcharge)	1 set	
or	<input type="checkbox"/> 3305-1901		Outdoor wheel set 200/250	(Upcharge)	1 set	
or	<input type="checkbox"/> 3305-1902		Outdoor wheel set 200/250 with swivel lock	(Upcharge)	1 set	
with	<input type="checkbox"/> 3306-1370	320	Seat cover size 1	1 Stk		incl.
or	<input type="checkbox"/> 3306-1330	320	Seat cover size 1 with back strap	(Upcharge)	1 Stk	
	<input type="checkbox"/> 3307-0000		CRUISER size 2	1 pc.		
with	<input type="checkbox"/> 3305-2401		Wheel set 150/200	1 set		incl.
or	<input type="checkbox"/> 3305-2402		Wheel set 150/200 with swivel lock	(Upcharge)	1 set	
or	<input type="checkbox"/> 3305-1901		Outdoor wheel set 200/250	(Upcharge)	1 set	
or	<input type="checkbox"/> 3305-1902		Outdoor wheel set 200/250 with swivel lock	(Upcharge)	1 set	
with	<input type="checkbox"/> 3307-1370	320	Seat cover size 2	1 pc.		incl.
or	<input type="checkbox"/> 3307-1330	320	Seat cover size 2 with back strap	(Upcharge)	1 pc.	

Accessories

	Order no.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> 3305-2100		Drag brake	1 set		
	<input type="checkbox"/> 3305-1500		Driving and parking brake, suitable for size 0 + 1	1 set		
	<input type="checkbox"/> 3305-1501		Driving and parking brake, suitable for size 0 + 1, Outdoor wheel set	1 set		
	<input type="checkbox"/> 3307-1500		Driving and parking brake, suitable for size 2	1 set		
	<input type="checkbox"/> 3307-1501		Driving and parking brake, suitable for size 2, Outdoor wheel set	1 set		
	<input type="checkbox"/> 3305-2200	320	Storage bag size 0	1 pc.		
	<input type="checkbox"/> 3306-2200	320	Storage bag size 1	1 pc.		
	<input type="checkbox"/> 3307-2200	320	Storage bag size 2	1 pc.		

Custom-made

	Order no.	Colour	Description	Qty.	Price	€
	<input type="checkbox"/> SA-0210-MP		Special Frame colour, RAL code only *	(Upcharge)	1 pc.	
			* please specify RAL color			
	<input type="checkbox"/> 6999-0014	060	Special Frame colour, Copper Gold	(Upcharge)	1 pc.	

Frame colours

001 Black (RAL 9005)	020 Yellow green (RAL 6018)	028 Pearl violet (RAL 4011)	056 Turquoise blue (RAL 5018)
002 Silver	021 Grey aluminium (RAL 9007)	029 Pearl beige (RAL 1035)	057 Iron sparkling dark light
003 Pearl ruby red (RAL 3032)	022 White (RAL 9003)	032 Black dull matt (RAL 9005M)	058 Salmon orange (RAL 2012)
004 Pearl night blue (RAL 5026)	023 Pure orange (RAL 2004)	039 Signal violet (RAL 4008)	059 Pastel blue (RAL 5024)
005 Pearl gentian blue (RAL 5025)	024 Sparkling iron effect	040 Night blue (RAL 5022)	061 Light pink (RAL 3015)
017 Traffic yellow (RAL 1023)	025 Signal blue (RAL 5005)	050 White matt (RAL 9003)	064 Olive green (RAL 6003)
018 Flame red (RAL 3000)	026 Vermillion (RAL 2002)	053 Light green (RAL 6027)	065 Raspberry red (RAL 3027)
019 Sky blue (RAL 5015)	027 Heather violet (RAL 4003)	054 Ocean blue (RAL 5020)	066 Pastel green (RAL 6019)

Technical data (dimensions in cm and weights in kg)

		Size 0		Size 1		Size 2	
		Min.	Max.	Min.	Max.	Min.	Max.
A	Handle height [cm]	32	59	32	59	54	83
a	Height forearm support [cm]						
B	Handle distance [cm]	18	32	22	36	36	50
C	Seat height [cm]	25		33		45	
D	Total width [cm]	47		51		65	
E	Total length [cm]	55		55		58	
F	Folding size (HxLxD) [cm]	42x55x29		50x55x27		68x58x29	
	max. load [kg]	30		50		65	
	Weight [kg]	5,4		5,6		5,9	