

CRUISER

Release date: February 2026 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address see invoice address

Company: _____

Street: _____

City, Zip code: _____



Basic configuration

- Handles individually adjustable in height & width
- Polyurethane tires
- Swivelling front wheels
- Folding seat cover

Frame colour



Special frame colour



Copper gold



Pastel violet, matt



Pastel blue, matt



Pastel green, matt

Desired RAL colour:

Technical data

	size 0		size 1		size 2	
	min.	max.	min.	max.	min.	max.
Handle height (cm)	32	59	40	69	44	83
Height forearm support (cm)						
Handle distance (cm)	18	32	22	36	36	50
Seat height (cm)	25		33		45	
Total width (cm)	47		51		65	
Total length (cm)	55		55		58	
Folding size (HxLxD) (cm)	42 x 55 x 29		50 x 55 x 27		68 x 58 x 29	
max. load (kg)	30		50		65	
Weight (kg)	5,4		5,6		5,9	



CRUISER size 0




CRUISER size 1




CRUISER size 2

Seat cover



Seat cover



Seat cover with back strap

Tires



Wheel set 150/200



Wheel set 150/200 with swivel lock




Outdoor wheel set 200/250



Outdoor wheel set 200/250 with swivel lock

Brakes



Drag brake




Driving and parking brake




Driving and parking brake, Outdoor wheel set

!
Notice:
Not compatible with the lower arm support with ball handle.

Backstops




Backstops for both rear wheels




Backstops for both outdoor rear wheels


Forearm rests



Rotatable forearm rests with angle and depth adjustable handles, tube/handle diameter 22/29 mm




Rotatable forearm rests with angle and depth adjustable handles, tube/handle diameter 15/18 mm



Rotatable forearm rests with central depth-adjustable 40 mm ball handle

! **Notice:**
Cannot be combined with the driving and parking brake.


Push bar, Storage bag



Push bar for accompanying person, retainer in frame color




Push bar for accompanying person Outdoor, retainer in frame color




Storage bag

Handles, Handle bar



Handles, tube/handle diameter 22/29 mm



Handles, tube/handle diameter 15/18 mm



Handle bar width adjustable and removable

Remark / Customer request / Individual special custom-made request

All prices are subject to change and non-binding. Errors, misprints, and price changes are reserved.
Binding prices are only those stated in our written offer or order confirmation.