



BINGO SWITCH Mini rehab pushchair

Size Mini for children at the age of approx. 1 - 3 years

Order form incl. Exportprices (€) valid until revoked, release date 01/2025

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



*The illustration shows the BINGO SWITCH in size 1

Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable, rotatable
- Maximum load 60 kg

Seat unit

- IUS – Individual Upholstery System
- Seat depth 12 - 27 cm
- Seat width 13 - 25 cm
- Lower leg length 16 - 36 cm
- Back height 37 - 51 cm
- Should height 25,5 – 39,5 cm
- Recline 80° - 180°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 25 kg

Minimum information (Without this information your order can not be handled.)

Frame colour:	<input type="checkbox"/> Black dull matt 032
Upholstery single-coloured:	<input type="checkbox"/> Black 346 <input type="checkbox"/> Berry soft 347 <input type="checkbox"/> Grey soft 348
with accessories colour black	<input type="checkbox"/> Turquoise soft 349
Upholstery two-coloured:	<input type="checkbox"/> Berry soft / Black 350 <input type="checkbox"/> Grey soft / Black 351 <input type="checkbox"/> Turquoise soft / Black 352
with accessories colour black	

Seat depth: (12 -27 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height: (37 - 51 cm) **Lower leg length:** (16 - 36 cm)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3230-2000	032	BINGO SWITCH Mini seat unit incl. foot rest hanger angle adjustable	1 pc.	
with	<input type="checkbox"/> 3230-2900		Thigh supports wide- and depth-adjustable * * minimum seat width 16 cm	01 1 pair	incl.
or	<input type="checkbox"/> 3230-2600		Hip supports, wide- and depth-adjustable, removable	(upcharge) 02 1 pair	
with	<input type="checkbox"/> 3230-5200		Upholstery, size 1, adaptable (IUS)	03 1 set	incl.
	<input type="checkbox"/> 3230-5190		Calf pad support plate and foam insert	04 1 set	
with	<input type="checkbox"/> 3231-2810		Foot rest bracket angle adjustable, flip up type	05 1 pc.	incl.
or	<input type="checkbox"/> 3231-2820		Foot rest bracket angle adjustable, flip up type, with flip up lock	(upcharge) 1 pc.	
or	<input type="checkbox"/> 3231-2830		Reinforced foot rest bracket "Steel SA-0430" angle adjustable, flip up type	06 1 pc.	incl.
or	<input type="checkbox"/> 3231-2840		Reinforced foot rest bracket "Steel SA-0430" angle adjustable, flip up type, with flip up lock * Extra weight ca. 0,3 kg	(upcharge) 1 pc.	
or	3231-2850		Foot rest bracket angle adjustable, with flip up lock, for short lower leg lengths (LLL reduction by 4 cm)	(upcharge) 1 pc.	
with	<input type="checkbox"/> 3231-249x		Footplate with heel stop width: <input type="checkbox"/> 1=27 cm, <input type="checkbox"/> 2=29 cm, <input type="checkbox"/> 3=31 cm	07 1 pc.	incl.
or	<input type="checkbox"/> 3231-248x		Evo foot rest with heel stop for short LLL width: <input type="checkbox"/> 1=27 cm <input type="checkbox"/> 2=29 cm <input type="checkbox"/> 3=31 cm	07 1 pc.	incl.
or	<input type="checkbox"/> 3231-2551		Foot plate in boat form, width ca. 30 cm	(upcharge) 07 1 pc.	
or	<input type="checkbox"/> 3231-2553		Bingo-Evo foot plate in boat form size 1, for short LLL	(upcharge) 07 1 pc.	
or	<input type="checkbox"/> 3231-2552	015	Alu foot plate in boat form, size 1 (not for short LLL)	(upcharge) 1 pc.	
or	<input type="checkbox"/> 3231-2463	015	Alu-footplate "footprint" with heel stop, width ca. 30 cm	(upcharge) 08 1 pc.	

Accessories - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3230-7601		Headrest pads low version size 1	09 1 set	
	<input type="checkbox"/> 3230-7602		Headrest pads, size 1	10 1 set	
	<input type="checkbox"/> 3230-6900		Anatomical neck-headrest with length compens, size 1	11 1 set	
	<input type="checkbox"/> 3230-7700		Trunk supports, size 1, stepless height adjustable (8 x 10 x 2 cm)	12 1 pair	
	<input type="checkbox"/> 3230-6100		Lap belt Mini	13 1 pc.	
	<input type="checkbox"/> 3231-6701		4-point lap belt, size XS (for seat width 18 - 24cm)	14 1 pc.	
	<input type="checkbox"/> 3230-6200		Five-point harness Mini	15 1 pc.	
	<input type="checkbox"/> 3201-6300		Foot straps	16 1 pair	
	<input type="checkbox"/> 3201-6400		Ankle huger, size 1	17 1 pair	
	<input type="checkbox"/> 3230-5300		Groin strap, Mini	18 1 pc.	
	<input type="checkbox"/> 3201-5800		Fixation jacket, size 0	19 1 pc.	
	<input type="checkbox"/> 3201-5400		Fixation jacket, size 1	1 pc.	
	<input type="checkbox"/> 3230-5600		Chest-shoulder harness, size 0	20 1 pc.	
	<input type="checkbox"/> 3201-5600		Chest-shoulder harness, size 1	1 pc.	
	<input type="checkbox"/> 3201-5601		Chest-shoulder harness, size 1 with 4 buckels	21 1 pc.	
	<input type="checkbox"/> 3201-5700		Chest-shoulder harness, size 2	1 pc.	
	<input type="checkbox"/> 3201-5701		Chest-shoulder harness, size 2 with 4 buckels	1 pc.	
	<input type="checkbox"/> 3231-7800		Abduction block, size 1, swing away and removable	22 1 pc.	
	<input type="checkbox"/> 3231-7801		Abduction block, size 1, fixed	23 1 pc.	
	<input type="checkbox"/> 3230-2205	032	Retainer for accessories, grab rail or therapy tray	24 1 pair	
	<input type="checkbox"/> 3230-8102	306	Grab rail with imitation leather, pivotable, size Mini, black ¹⁾	25 1 pc.	
	<input type="checkbox"/> 3230-8400		Therapy tray, size 1, transparent ^{1) 2) 3)}	26 1 pc.	
	<input type="checkbox"/> 3230-8402		Therapy tray, size 1, transparent, height adjustable ^{1) 2)}	27 1 pc.	
	<input type="checkbox"/> 3230-8401		Therapy tray, size 1, grey ^{1) 2) 3)}	28 1 pc.	
	<input type="checkbox"/> 3230-8403		Therapy tray, size 1, grey, height adjustable ^{1) 2)}	1 pc.	
			¹⁾ only possible with retainer for accessories		
			²⁾ only possible with seat in driving direction		
			³⁾ not possible when using hip supports		

Accessories - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	3230-8300	306	Winter warmer Mini, black (up to SD 20)	29	1 pc.
<input type="checkbox"/>	3201-8300	306	Winter warmer, size 1, black (from SD 18)		1 pc.
<input type="checkbox"/>	7000-0018		Sheep skin insert, size 1	30	1 pc.
<input type="checkbox"/>	3201-9200	306	Rain cape, size 1, black	31	1 pc.
<input type="checkbox"/>	3230-9300	306	Canopy incl. adapter and rain cover, black	32	1 pc.
<input type="checkbox"/>	3002-0012	032	BASIC home underframe size Mini with 75 mm twin castors	33	1 pc.
<input type="checkbox"/>	3002-0013	032	BASIC home underframe size Mini with 125 mm castors	34	1 pc.
<input type="checkbox"/>	3002-9500	022	Side covers in standard color white	35	
<input type="checkbox"/>	3002-9500		Side covers matching the upholstery colour	36	
<input type="checkbox"/>	3002-9500	032	Side covers in black		
<input type="checkbox"/>	3002-5400		Push handle for BASIC home under frame		1 pc.
<input type="checkbox"/>	3004-0002	032	COBRA indoor base, size Mini , with 75 mm twin castors	37	1 pc.
<input type="checkbox"/>	3004-0003	032	COBRA indoor base, size Mini , with 125 mm castors	34	1 pc.
Presetting of the gas spring height adjustment: <input type="checkbox"/> up to 20kg, <input type="checkbox"/> up to 35kg, <input type="checkbox"/> up to 50kg					

Basic variations - mobility base with steering wheel attachment

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	3030-1000	032	BINGO SWITCH Mini mobility base with steering wheel attachment	1 pc.	
<input type="checkbox"/>	3031-1000	032	BINGO SWITCH mini mobility base for additional devices with steering wheel attachment	1 pc.	
with	<input type="checkbox"/>	1281-1027	Plastic front forks	38	1 set incl.
	<input type="checkbox"/>	3202-7100	Swivel lock for Plastic front forks	39	1 set
or	<input type="checkbox"/>	1281-1028	Aluminum front forks	(upcharge) 40	1 set
	<input type="checkbox"/>	3232-7100	Swivel lock for Aluminum front forks		1 set
with	<input type="checkbox"/>	3031-1301	032 Without suspension system		1 set incl.
or	<input type="checkbox"/>	3031-7400	Frame suspension system, fourfold adjustable	(upcharge) 41	1 set
	<input type="checkbox"/>	1131-9017	022 spring set (white) for load up to 15 kg		
	<input type="checkbox"/>	1131-9017	spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/>	1131-9018	022 spring set (white) for load up to 25 kg		
	<input type="checkbox"/>	1131-9018	spring set (in upholstery color) foar load up to 25 kg		
with	<input type="checkbox"/>	3231-9500	Side covers matching the upholstery color (Standard)	42	1 pair incl.
or	<input type="checkbox"/>	3231-9500	032 Side covers matching the frame color	43	1 pair incl.
with	<input type="checkbox"/>	1219-0004S	Pneumatic wheels with 6-spoke rims	44	1 set incl.
or	<input type="checkbox"/>	1219-0030	Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge) 45	1 set
or	<input type="checkbox"/>	1229-0004S	Polyurethan wheels with 6-spoke rims	46	1 set incl.
or	<input type="checkbox"/>	1229-0009	Polyurethan wheels with 7-spoke rims	47	1 set incl.
or	<input type="checkbox"/>	1229-0030	Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge) 48	1 set
	<input type="checkbox"/>	3231-7201	Attendant brake, two hand use	49	1 set
with	<input type="checkbox"/>	3232-9401	Soft accessory bag, size 1 (max. 3 kg)	50	1 pc. incl.
or	<input type="checkbox"/>	3030-9402	Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge) 51	1 pc.
* do not combine with tie down kit, use with seat in driving direction only					
or	<input type="checkbox"/>	3202-4100	Tray for additional equipment size 1 (only with 3033-10xx) (Custom-made)	(upcharge) 52	1 pc.
with	<input type="checkbox"/>	3230-1501	032 Without tie down kit		1 set incl.
or	<input type="checkbox"/>	3230-7300	032 Tie down kit according to ISO 7176-19 (from 18 kg body weight) *		1 set
* only possible with chest-shoulder harness or fixation jacket as positioning aid					
with	<input type="checkbox"/>	3231-7100	Anti tip left and right *		1 set
* not possible in conjunction with frames for additional devices					

Basic variations - mobility base with center steering wheel attachment					with/or: Please select only one version!	
	Orderno.	Colour	Description	Qty.	Price	
	<input type="checkbox"/> 3030-1020	032	BINGO SWITCH Mini mobility base with center steering wheel attachment			
	<input type="checkbox"/> 3031-1020	032	BINGO SWITCH Mini mobility base for additional devices with center steering wheel attachment			
with	<input type="checkbox"/> SA-0149/2-KP		Swivellock for center steering wheel			
with	<input type="checkbox"/> 3031-1301	032	Without suspension system			incl.
or	<input type="checkbox"/> 3031-7400		Frame suspension system, fourfold adjustable	(upcharge)	41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg			
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg			
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg			
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg			
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)		42	1 pair incl.
or	<input type="checkbox"/> 3231-9500	032	Side covers matching the frame color		43	1 pair incl.
with	<input type="checkbox"/> 1219-0056		Pneumatic wheels with 6-spoke rims		44	1 set incl.
or	<input type="checkbox"/> 1219-0057		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge)	45	1 set
or	<input type="checkbox"/> 1229-0056		Polyurethan wheels with 6-spoke rims		46	1 set incl.
or	<input type="checkbox"/> 1229-0057		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge)	48	1 set
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use		49	1 set
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)		50	1 pc. incl.
or	<input type="checkbox"/> 3030-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge)	51	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only			
or	<input type="checkbox"/> 3202-4100		Tray for additional equipment size 1 (only with 3033-10xx) (Custom-made)	(upcharge)	52	1 pc.
with	<input checked="" type="checkbox"/> 3230-1501	032	Without tie down kit			1 set incl.
with	<input type="checkbox"/> 3231-7100		Anti tip left and right *			1 set
			* not possible in conjunction with frames for additional devices			

Basic variations - mobility base CROSS with Aluminium front forks					with/or: Please select only one version!	
	Orderno.	Colour	Description	Qty.	Price	
	<input type="checkbox"/> 3030-1030	032	BINGO SWITCH CROSS Mini mobility base with steering wheel attachment			
	<input type="checkbox"/> 3031-1030	032	BINGO SWITCH CROSS Mini mobility base for additional devices with steering wheel attachment			
with	<input type="checkbox"/> 3031-1301		Without suspension system			1 set incl.
or	<input type="checkbox"/> 3031-7400		Frame suspension system, fourfold adjustable	(upcharge)	41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg			
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg			
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg			
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg			
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)		42	1 pair incl.
or	<input type="checkbox"/> 3231-9500	032	Side covers matching the frame color		43	1 pair incl.
with	<input type="checkbox"/> 1219-0050		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function			1 set incl.
or	<input type="checkbox"/> 1229-0050		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function			1 set incl.
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)		50	1 pc. incl.
or	<input type="checkbox"/> 3030-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge)	51	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only			
or	<input type="checkbox"/> 3238-4400		Tray for additional equipment size 1 (max. 25kg) only with 3031-1030/40 (Custom-made)	(Aufpreis)	52	1 pc.
with	<input type="checkbox"/> 3230-1501	032	Without tie down kit			1 set incl.
or	<input type="checkbox"/> 3230-7300	032	Tie down kit according to ISO 7176-19 (from 18 kg body weight) *			1 set
			* only possible with chest-shoulder harness or fixation jacket as positioning aid			
with	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks			1 set
with	<input type="checkbox"/> 3234-7100		Anti tip left and right *			1 set
			* not possible in conjunction with frames for additional devices			

Basic variations - mobility base CROSS with Aluminium front forks				with/or: Please select only one version!	
	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3030-1040	032	BINGO SWITCH CROSS Mini mobility base with center steering wheel attachment	1 pc.	
	<input type="checkbox"/> 3031-1040	032	BINGO SWITCH CROSS Mini mobility base for additional devices with center steering wheel attachmentsteering	1 pc.	
with	<input type="checkbox"/> 3031-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3031-7400	032	Frame suspension system, fourfold adjustable (upcharge)	41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	42	1 pair incl.
or	<input type="checkbox"/> 3231-9500	032	Side covers matching the frame color	43	1 pair incl.
with	<input type="checkbox"/> 1219-0058		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function		1 set incl.
or	<input type="checkbox"/> 1229-0058		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function		1 set incl.
with	<input type="checkbox"/> 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	50	1 pc. incl.
or	<input type="checkbox"/> 3030-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge)	51	1 pc.
			* do not combine with tie down kit, use with seat in driving direction only		
or	<input type="checkbox"/> 3238-4400		Tray for additional equipment size 1 (max. 25kg) only with 3031-1030/40 (Custom-made) (upcharge)	52	1 pc.
with	<input checked="" type="checkbox"/> 3230-1501	032	Without tie down kit		1 set incl.
with	<input type="checkbox"/> SA-0149/2-KP		Swivellock for center steering wheel (Custom-made)		1 set
with	<input type="checkbox"/> 3234-7100		Anti tip left and right *		1 set
			* not possible in conjunction with frames for additional devices		

! BINGO Evolution Cross: Custom-made only on request !

Accessories - mobility base				with/or: Please select only one version!	
	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	39	1 set
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks		1 set
	<input type="checkbox"/> SA-0149/2-KP		Swivellock for center steering wheel		1 set
	<input type="checkbox"/> 3031-7400E	032	Frame suspension system, fourfold adjustable	41	1 set
	<input type="checkbox"/> 1131-9017	022	spring set (white) for load up to 15 kg		
	<input type="checkbox"/> 1131-9017		spring set (in upholstery color) foar load up to 15 kg		
	<input type="checkbox"/> 1131-9018	022	spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		spring set (in upholstery color) foar load up to 25 kg		
	<input type="checkbox"/> 3230-9402E		Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge)	51	1 pc.
			* only possible with chest-shoulder harness or fixation jacket as positioning aid		
	<input type="checkbox"/> 3231-9500E		Side covers matching the upholstery or frame colour		1 pair

Custom-made with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	3202-4100E		Tray for additional equipment size 1 (only with 3031-10xx)	52	1 pc.
<input type="checkbox"/>	3238-4400E		CROSS tray for additional equipment size 1 (max. 25kg) only with 3031-1030/40	52	1 pc.
<input type="checkbox"/>	3030-4200		Oxygen bottle holder (d=100 mm) long, right hand attachment especially for bottles up to 75 cm *	53	1 pc.
<input type="checkbox"/>	3030-4200_CROSS		Oxygen bottle holder (d=100 mm) long, right hand attachment especially for bottles up to 75 cm *	53	1 pc.
<input type="checkbox"/>	3030-4201		Oxygen bottle holder short version (d=100 mm), right hand attachment especially for 2 L bottles *		1 pc.
<input type="checkbox"/>	3030-4201_CROSS		Oxygen bottle holder short version (d=100 mm), right hand attachment especially for 2 L bottles *		1 pc.
<input type="checkbox"/>	SA-0225-MP		Two adjustable angle holders for contoured oxygen bottle with velcro fasteners (upcharge)	54	1 set
<input type="checkbox"/>	SA-0031-KP		Infusions bottle holder	55	1 pc.
<input type="checkbox"/>	SA-0605-MP		Footrest carrier for short LLL		1 pc.
<input type="checkbox"/>	SA-0016/3-MP		Reinforced rear axle (10 kg more load capacity) (upcharge)	56	1 pc.
<input type="checkbox"/>	SA-0298-MP		Seat adapter 10 cm heightened (upcharge)		1 set
<input type="checkbox"/>	SA-0321-MP		Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base * (upcharge)	57	1 set
* only possible with 2 mobility bases of one size					
* mobility base BINGO Evolution Mini only possible at the front					
* not possible in combination with aluminium front forks & mobility base Cross					
<input type="checkbox"/>	SA-0441-MP		Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution mobility base 3238-1011 rear * (upcharge)		1 set
* not possible in combination with aluminium front forks & mobility base Cross					
<input type="checkbox"/>	SA-0619/1-KP		Mounting adapter for Outdoor front end aid for "phil & teds" sibling seat or universal car seat adapter	58	1 set

Change attachments with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	3032-1300E	032	Steering wheel change attachment BINGO SWITCH		1 pair
with	<input type="checkbox"/>	1281-1027	Plastik front forks	38	1 set incl.
	<input type="checkbox"/>	3202-7100	Swivel lock for Plastic front forks	39	1 set
or	<input type="checkbox"/>	1281-1028	Aluminium front forks (upcharge)	40	1 set
	<input type="checkbox"/>	3232-7100	Swivel lock for Aluminium front forks		1 set
with	<input type="checkbox"/>	1219-0055	Steering wheels with pneumatic tyres on 3-spoke rims	44	1 set incl.
or	<input type="checkbox"/>	1229-0055	Steering wheels with PU tyres on 3-spoke rims	46	1 set incl.
or	<input type="checkbox"/>	1229-0054	Steering wheels with PU tyres on 7-spoke rims	47	1 set incl.
<input type="checkbox"/>	3032-1501E	032	Center steering wheel interchangeable attachment BINGO SWITCH (Mini/1), pneumatic tires		1 set
<input type="checkbox"/>	3032-1502E	032	Center steering wheel interchangeable attachment BINGO SWITCH (Mini/1), PU tires		1 set
<input type="checkbox"/>	SA-0149/2-KP		Swivellock for center steering wheel		1 set

Extra equipment with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	7000-0010		Buggy-Board	59	1 pc.
<input type="checkbox"/>	7000-0012	306	Sun shade, black	60	1 pc.
<input type="checkbox"/>	7000-0006		Headrest speakers, for headrest pads 3201-7601		1 set
<input type="checkbox"/>	7000-0202	306	Changing backpack, black	61	1 pc.
<input type="checkbox"/>	7000-0199		"Cosies" Thermofleece hand warmer for companion	62	1 set

Pre-assembly of pads, upholstery parts & belts

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	6999-1012		BINGO Evolution Pre-assembly service (without settings)		1 pc.

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH Mini




Thigh supports



Hip supports



IUS - Upholstery



Calf pad support plate and foam insert




Footrest flip up lock



Reinforced steel footplate-carrier



Footplate with heel stop or boat form



Alu-footplate „foot-print“with heel stop



Headrest pads low version



Headrest pads high version



Anatomical neck-headrest



Trunk supports



Lap belt



4-point lap belt



Five-point harness



Foot straps



Ankle huger



Groin strap



Fixation jacket



Chest-shoulder harness



Chest-shoulder harness with 4 buckles



Abduction block, swing away and removable




Abduction block, fixed




Retainer for accessories



Grap rail



Therapy tray



Therapy tray, height adjustable



Therapy tray, grey



Winter warmer



Sheep skin insert

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH Mini



31 Rain cape



32 Canopy



33 BASIC home underframe, 75 mm twin castors



34 125 mm castors



35 Side covers in standard color white



36 Side covers matching the upholstery colour



37 COBRA indoor base, 75 mm twin castors



38 Plastic front forks



39 Swivel lock



40 Aluminum front forks



41 Frame suspension system



42 Side covers matching the upholstery color



43 Side covers matching the frame color



44 Pneumatic wheel, front wheel



45 Pneumatic wheel, rear wheel



46 Pneumatic tyres front wheel, Drum brake wheel



47 Pneumatic tyres rear wheel, Drum brake wheel



48 Polyurethan wheel, front wheel



49 Polyurethan wheel, rear wheel



50 Polyurethan wheel front wheel, black rim



51 Polyurethan wheel rear wheel, black rim



52 Polyurethan tyres front wheel, Drum brake wheel



53 Polyurethan tyres rear wheel, Drum brake wheel




54 Attendant brake, two hand use



55 Soft accessory bag



56 Accessory bag



57 Tray for additional equipment



58 Oxygen bottle holder




59 Adjustable Angle brackets



60 Infusions bottle holder

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH Mini



Reinforced rear axle



Connecting parts for mobility bases



Mounting adapter for Outdoor front end



Buggy-Board



Sun shade, grey



Changing backpack



„Cosies“ Thermofleece hand warmer