

HOGGI®

Partnering together to mobilise kids



BINGO SWITCH Mini rehab pushchair

Size Mini for children at the age of approx. 1 - 3 years

Order form incl. Exportprices (€) valid until revoked, release date 03/2024

☐ Offer	☐ Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



*The illustration shows the BINGO SWITCH in size 1

Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable, rotatable
- Maximum load 60 kg

Seat unit

- IUS – Individual Upholstery System
- Seat depth 12 - 27 cm
- Seat width 13 - 25 cm
- Lower leg length 16 - 36 cm
- Back height 37 - 51 cm
- Should height 25,5 – 39,5 cm
- Recline 80° - 180°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 25 kg

Minimum information (Without this information your order can not be handled.)

Frame colour:	☐ Black dull matt 032		
Upholstery single-coloured:	☐ Black 346	☐ Berry soft 347	☐ Grey soft 348
with accessories colour black	☐ Turquoise soft 349		
Upholstery two-coloured:	☐ Berry soft / Black 350	☐ Grey soft / Black 351	☐ Turquoise soft / Black 352
with accessories colour black			

Seat depth: (12 -27 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height: (37 - 51 cm)

Lower leg length: (16 - 36 cm)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
☐	3230-2000K	032	BINGO Evolution Mini seat unit			
			incl. foot rest hanger angle adjustable		1 pc.	
with	☐ 3230-2900		Thigh supports wide- and depth-adjustable *		01 1 pair	incl.
			* minimum seat width 16 cm			
or	☐ 3230-2600		Hip supports, wide- and depth-adjustable, removable	(upcharge)	02 1 pair	
with	☐ 3230-5200		Upholstery, size 1, adaptable (IUS)		03 1 set	incl.
	☐ 3230-5190		Calf pad support plate and foam insert		04 1 set	
with	☐ 3231-2810		Foot rest bracket angle adjustable, flip up type		05 1 pc.	incl.
or	☐ 3231-2820		Foot rest bracket angle adjustable, flip up type, with flip up lock	(upcharge)	1 pc.	
or	☐ 3231-2830		Reinforced foot rest bracket "Steel SA-0430" angle adjustable, flip up type		06 1 pc.	incl.
or	☐ 3231-2840		Reinforced foot rest bracket "Steel SA-0430" angle adjustable, flip up type, with flip up lock	(upcharge)	1 pc.	
			* Extra weight ca. 0,3 kg			
or	3231-2850		Foot rest bracket angle adjustable, with flip up lock, for short lower leg lengths (LLL reduction by 4 cm)	(upcharge)	1 pc.	
with	☐ 3231-249x		Footplate with heel stop		07 1 pc.	incl.
			width: ☐ 1=27 cm, ☐ 2=29 cm, ☐ 3=31 cm			
or	☐ 3231-248x		Evo foot rest with heel stop for short LLL		07 1 pc.	incl.
			width: ☐ 1=27 cm ☐ 2=29 cm ☐ 3=31 cm			
or	☐ 3231-2551		Foot plate in boat form, width ca. 30 cm	(upcharge)	07 1 pc.	
or	☐ 3231-2553		Bingo-Evo foot plate in boat form size 1, for short LLL	(upcharge)	07 1 pc.	
or	☐ 3231-2552	015	Alu foot plate in boat form, size 1 (not for short LLL)	(upcharge)	1 pc.	
or	☐ 3231-2463	015	Alu-footplate "footprint" with heel stop, width ca. 30 cm	(upcharge)	08 1 pc.	

Accessories - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
☐	3230-7601		Headrest pads low version size 1		09 1 set	
☐	3230-7602		Headrest pads, size 1		10 1 set	
☐	3230-6900		Anatomical neck-headrest with length compens, size 1		11 1 set	
☐	3230-7700		Trunk supports, size 1, stepless height adjustable (8 x 10 x 2 cm)		12 1 pair	
☐	3230-6100		Lap belt Mini		13 1 pc.	
☐	3231-6701		4-point lap belt, size XS (for seat width 18 - 24cm)		14 1 pc.	
☐	3230-6200		Five-point harness Mini		15 1 pc.	
☐	3201-6300		Foot straps		16 1 pair	
☐	3201-6400		Ankle hugger, size 1		17 1 pair	
☐	3230-5300		Groin strap, Mini		18 1 pc.	
☐	3201-5800		Fixation jacket, size 0		19 1 pc.	
☐	3201-5400		Fixation jacket, size 1		1 pc.	
☐	3230-5600		Chest-shoulder harness, size 0		20 1 pc.	
☐	3201-5600		Chest-shoulder harness, size 1		1 pc.	
☐	3201-5601		Chest-shoulder harness, size 1 with 4 buckels		21 1 pc.	
☐	3201-5700		Chest-shoulder harness, size 2		1 pc.	
☐	3201-5701		Chest-shoulder harness, size 2 with 4 buckels		1 pc.	
☐	3231-7800		Abduction block, size 1, swing away and removable		22 1 pc.	
☐	3231-7801		Abduction block, size 1, fixed		23 1 pc.	
☐	3230-2201	032	Retainer for accessories, grab rail or therapy tray		24 1 pair	
☐	3230-8100		Grab rail with upholstery, size 1 ¹⁾		25 1 pc.	
☐	3230-8200		Therapy tray, size 1, transparent ^{1) 2) 3)}		26 1 pc.	
☐	3230-8202		Therapy tray, size 1, transparent, height adjustable ^{1) 2)}		27 1 pc.	
☐	3230-8201		Therapy tray, size 1, grey ^{1) 2) 3)}		1 pc.	
☐	3230-8203		Therapy tray, size 1, grey, height adjustable ^{1) 2)}		1 pc.	
			¹⁾ only possible with retainer for accessories			
			²⁾ only possible with seat in driving direction			
			³⁾ not possible when using hip supports			

Accessories - seat unit					with/or: Please select only one version!		
	Orderno.	Colour	Description		Qty.	Price	
	☐ 3230-8300	306	Winter warmer Mini, black (up to SD 20)		28	1 pc.	
	☐ 3201-8300	306	Winter warmer, size 1, black (from SD 18)			1 pc.	
	☐ 7000-0018		Sheep skin insert, size 1		29	1 pc.	
	☐ 3201-9200	306	Rain cape, size 1, black		30	1 pc.	
	☐ 3230-9300	306	Canopy incl. adapter and rain cover, black		31	1 pc.	
	☐ 3002-0012	032	BASIC home underframe size Mini with 75 mm twin castors		32	1 pc.	
	☐ 3002-0013	032	BASIC home underframe size Mini with 125 mm castors		33	1 pc.	
	☐ 3002-9500	022	Side covers in standard color white		34		
	☐ 3002-9500		Side covers matching the upholstery colour		35		
	☐ 3002-9500	032	Side covers in black				
	☐ 3002-5400		Push handle for BASIC home under frame			1 pc.	
	☐ 3004-0002	032	COBRA indoor base, size Mini , with 75 mm twin castors		36	1 pc.	
	☐ 3004-0003	032	COBRA indoor base, size Mini , with 125 mm castors		33	1 pc.	
	Presetting of the gas spring height adjustment: ☐ up to 20kg, ☐ up to 35kg, ☐ up to 50kg						
Basic variations - mobility base with steering wheel attachment					with/or: Please select only one version!		
	Orderno.	Colour	Description		Qty.	Price	
	☐ 3030-1000	032	BINGO SWITCH Mini mobility base with steering wheel attachment		1 pc.		
	☐ 3031-1000	032	BINGO SWITCH mini mobility base for additional devices with steering wheel attachment		1 pc.		
with	☐ 1281-1027		Plastic front forks		37	1 set	incl.
	☐ 3202-7100		Swivel lock for Plastic front forks		38	1 set	
or	☐ 1281-1028		Aluminum front forks	(upcharge)	39	1 set	
	☐ 3232-7100		Swivel lock for Aluminum front forks			1 set	
with	☐ 3031-1301	032	Without suspension system			1 set	incl.
or	☐ 3031-7400		Frame suspension system, fourfold adjustable	(upcharge)	40	1 set	
	☐ 1131-9017	022	spring set (white) for load up to 15 kg				
	☐ 1131-9017		spring set (in upholserly color) foar load up to 15 kg				
	☐ 1131-9018	022	spring set (white) for load up to 25 kg				
	☐ 1131-9018		spring set (in upholserly color) foar load up to 25 kg				
with	☐ 3231-9500		Side covers matching the upholstery color (Standard)		41	1 pair	incl.
or	☐ 3231-9500	032	Side covers matching the frame color		42	1 pair	incl.
with	☐ 1219-0004S		Pneumatic wheels with 6-spoke rims		43	1 set	incl.
or	☐ 1219-0030		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge)	44	1 set	
or	☐ 1229-0004S		Polyurethan wheels with 6-spoke rims		45	1 set	incl.
or	☐ 1229-0009		Polyurethan wheels with 7-spoke rims		46	1 set	incl.
or	☐ 1229-0030		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge)	47	1 set	
	☐ 3231-7201		Attendant brake, two hand use		48	1 set	
with	☐ 3232-9401		Soft accessory bag, size 1 (max. 3 kg)		49	1 pc.	incl.
or	☐ 3030-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge)	50	1 pc.	
	* do not combine with tie down kit, use with seat in driving direction only						
or	☐ 3202-4100		Tray for additional equipment size 1 (only with 3033-10xx) (Custom-made)	(upcharge)	51	1 pc.	
with	☐ 3230-1501	032	Without tie down kit			1 set	incl.
	* only possible with chest-shoulder harness or fixation jacket as positioning aid						
with	☐ 3231-7100		Anti tip left and right *			1 set	
	* not possible in conjunction with frames for additional devices						

Basic variations - mobility base with center steering wheel attachment					with/or: Please select only one version!		
	Orderno.	Colour	Description		Qty.	Price	
	☐ 3030-1020	032	BINGO SWITCH Mini mobility base with center steering wheel attachment		1 pc.		
	☐ 3031-1020	032	BINGO SWITCH Mini mobility base for additional devices with center steering wheel attachment		1 pc.		
with	☐ SA-0149/2-KP		Swivellock for center steering wheel		1 set		
with	☐ 3031-1301	032	Without suspension system		1 set	incl.	
or	☐ 3031-7400		Frame suspension system, fourfold adjustable	(upcharge)	40	1 set	
	☐ 1131-9017	022	spring set (white) for load up to 15 kg				
	☐ 1131-9017		spring set (in upholstery color) foar load up to 15 kg				
	☐ 1131-9018	022	spring set (white) for load up to 25 kg				
	☐ 1131-9018		spring set (in upholstery color) foar load up to 25 kg				
with	☐ 3231-9500		Side covers matching the upholstery color (Standard)		41	1 pair	
or	☐ 3231-9500	032	Side covers matching the frame color		42	1 pair	
with	☐ 1219-0056		Pneumatic wheels with 6-spoke rims		43	1 set	
or	☐ 1219-0057		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge)	44	1 set	
or	☐ 1229-0056		Polyurethan wheels with 6-spoke rims		45	1 set	
or	☐ 1229-0057		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge)	47	1 set	
	☐ 3231-7201		Attendant brake, two hand use		48	1 set	
with	☐ 3232-9401		Soft accessory bag, size 1 (max. 3 kg)		49	1 pc.	
or	☐ 3030-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge)	50	1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only				
or	☐ 3202-4100		Tray for additional equipment size 1 (only with 3033-10xx) (Custom-made)	(upcharge)	51	1 pc.	
with	☒ 3230-1501	032	Without tie down kit		1 set	incl.	
with	☐ 3231-7100		Anti tip left and right *		1 set		
			* not possible in conjunction with frames for additional devices				
Basic variations - mobility base CROSS with Aluminium front forks					with/or: Please select only one version!		
	Orderno.	Colour	Description		Qty.	Price	
	☐ 3030-1030	032	BINGO SWITCH CROSS Mini mobility base with steering wheel attachment		1 pc.		
	☐ 3031-1030	032	BINGO SWITCH CROSS Mini mobility base for additional devices with steering wheel attachment		1 pc.		
with	☐ 3031-1301		Without suspension system		1 set	incl.	
or	☐ 3031-7400		Frame suspension system, fourfold adjustable	(upcharge)	40	1 set	
	☐ 1131-9017	022	spring set (white) for load up to 15 kg				
	☐ 1131-9017		spring set (in upholstery color) foar load up to 15 kg				
	☐ 1131-9018	022	spring set (white) for load up to 25 kg				
	☐ 1131-9018		spring set (in upholstery color) foar load up to 25 kg				
with	☐ 3231-9500		Side covers matching the upholstery color (Standard)		41	1 pair	
or	☐ 3231-9500	032	Side covers matching the frame color		42	1 pair	
with	☐ 1219-0050		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function		1 set	incl.	
or	☐ 1229-0050		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function		1 set	incl.	
with	☐ 3232-9401		Soft accessory bag, size 1 (max. 3 kg)		49	1 pc.	
or	☐ 3030-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) *	(upcharge)	50	1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only				
or	☐ 3238-4400		Tray for additional equipment size 1 (max. 25kg) only with 3031-1030/40 (Custom-made)	(Aufpreis)	51	1 pc.	
with	☐ 3230-1501	032	Without tie down kit		1 set	incl.	
			* only possible with chest-shoulder harness or fixation jacket as positioning aid				
with	☐ 3232-7100		Swivel lock for Aluminum front forks		1 set		
with	☐ 3234-7100		Anti tip left and right *		1 set		
			* not possible in conjunction with frames for additional devices				

Basic variations - mobility base CROSS with Aluminium front forks					with/or: Please select only one version!	
	Orderno.	Colour	Description		Qty.	Price
	☐ 3030-1040	032	BINGO SWITCH CROSS Mini mobility base with center steering wheel attachment		1 pc.	
	☐ 3031-1040	032	BINGO SWITCH CROSS Mini mobility base for additional devices with center steering wheel attachmentsteering		1 pc.	
with	☐ 3031-1301		Without suspension system		1 set	incl.
or	☐ 3031-7400	032	Frame suspension system, fourfold adjustable	(upcharge) 40	1 set	
	☐ 1131-9017	022	spring set (white) for load up to 15 kg			
	☐ 1131-9017		spring set (in upholstery color) foar load up to 15 kg			
	☐ 1131-9018	022	spring set (white) for load up to 25 kg			
	☐ 1131-9018		spring set (in upholstery color) foar load up to 25 kg			
with	☐ 3231-9500		Side covers matching the upholstery color (Standard)	41	1 pair	incl.
or	☐ 3231-9500	032	Side covers matching the frame color	42	1 pair	incl.
with	☐ 1219-0058		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function		1 set	incl.
or	☐ 1229-0058		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function		1 set	incl.
with	☐ 3232-9401		Soft accessory bag, size 1 (max. 3 kg)	49	1 pc.	incl.
or	☐ 3030-9402		Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge)	50	1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only			
or	☐ 3238-4400		Tray for additional equipment size 1 (max. 25kg) only with 3031-1030/40 (Custom-made)	(upcharge) 51	1 pc.	
with	☒ 3230-1501	032	Without tie down kit		1 set	incl.
with	☐ SA-0149/2-KP		Swivellock for center steering wheel (Custom-made)		1 set	
with	☐ 3234-7100		Anti tip left and right *		1 set	
			* not possible in conjunction with frames for additional devices			
! BINGO Evolution Cross: Custom-made only on request !						

Accessories - mobility base					with/or: Please select only one version!	
	Orderno.	Colour	Description		Qty.	Price
	☐ 3202-7100		Swivel lock for Plastic front forks	38	1 set	
	☐ 3232-7100		Swivel lock for Aluminum front forks		1 set	
	☐ SA-0149/2-KP		Swivellock for center steering wheel		1 set	
	☐ 3031-7400E	032	Frame suspension system, fourfold adjustable	40	1 set	
	☐ 1131-9017	022	spring set (white) for load up to 15 kg			
	☐ 1131-9017		spring set (in upholstery color) foar load up to 15 kg			
	☐ 1131-9018	022	spring set (white) for load up to 25 kg			
	☐ 1131-9018		spring set (in upholstery color) foar load up to 25 kg			
	☐ 3230-9402E		Accessory bag, for "mini" mobility bases (max. 8 kg) * (upcharge)	50	1 pc.	
			* only possible with chest-shoulder harness or fixation jacket as positioning aid			
	☐ 3231-9500E		Side covers matching the upholstery or frame colour		1 pair	

Custom-made				with/or: Please select only one version!	
	Orderno.	Colour	Description	Qty.	Price
☐	3202-4100E		Tray for additional equipment size 1 (only with 3031-10xx)	51	1 pc.
☐	3238-4400E		CROSS tray for additional equipment size 1 (max. 25kg) only with 3031-1030/40	51	1 pc.
☐	3030-4200		Oxygen bottle holder (d=100 mm) long, right hand attachment especially for bottles up to 75 cm *	52	1 pc.
☐	3030-4200_CROSS		Oxygen bottle holder (d=100 mm) long, right hand attachment especially for bottles up to 75 cm *	52	1 pc.
☐	3030-4201		Oxygen bottle holder short version (d=100 mm), right hand attachment especially for 2 L bottles *		1 pc.
☐	3030-4201_CROSS		Oxygen bottle holder short version (d=100 mm), right hand attachment especially for 2 L bottles *		1 pc.
☐	SA-0225-MP		Two adjustable angle holders for contoured oxygen bottle with velcro fasteners (upcharge)	53	1 set
☐	SA-0031-KP		Infusions bottle holder	54	1 pc.
☐	SA-0605-MP		Footrest carrier for short LLL		1 pc.
☐	SA-0016/3-MP		Reinforced rear axle (10 kg more load capacity) (upcharge)	55	1 pc.
☐	SA-0298-MP		Seat adapter 10 cm heightened * (upcharge)		1 set
☐	SA-0321-MP		Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base * (upcharge)	56	1 set
* only possible with 2 mobility bases of one size					
* mobility base BINGO Evolution Mini only possible at the front					
* not possible in combination with aluminium front forks & mobility base Cross					
☐	SA-0441-MP		Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution mobility base 3238-1011 rear * (upcharge)		1 set
* not possible in combination with aluminium front forks & mobility base Cross					

Change attachments				with/or: Please select only one version!	
	Orderno.	Colour	Description	Qty.	Price
☐	3032-1300E	032	Steering wheel change attachment BINGO SWITCH		1 pair
with ☐	1281-1027		Plastik front forks	37	1 set incl.
☐	3202-7100		Swivel lock for Plastic front forks	38	1 set
or ☐	1281-1028		Aluminium front forks (upcharge)	39	1 set
☐	3232-7100		Swivel lock for Aluminium front forks		1 set
with ☐	1219-0055		Steering wheels with pneumatic tyres on 3-spoke rims	43	1 set incl.
or ☐	1229-0055		Steering wheels with PU tyres on 3-spoke rims	45	1 set incl.
or ☐	1229-0054		Steering wheels with PU tyres on 7-spoke rims	46	1 set incl.
☐	3032-1501E	032	Center steering wheel interchangeable attachment BINGO SWITCH (Mini/1), pneumatic tires		1 set
☐	3032-1502E	032	Center steering wheel interchangeable attachment BINGO SWITCH (Mini/1), PU tires		1 set
☐	SA-0149/2-KP		Swivellock for center steering wheel		1 set

Extra equipment				with/or: Please select only one version!	
	Orderno.	Colour	Description	Qty.	Price
☐	7000-0010		Buggy-Board	57	1 pc.
☐	7000-0012	306	Sun shade, black	58	1 pc.
☐	7000-0006		Headrest speakers, for headrest pads 3201-7601		1 set
☐	7000-0202	306	Changing backpack, black	59	1 pc.
☐	7000-0199		"Cosies" Thermofleece hand warmer for companion	60	1 set

Pre-assembly of pads, upholstery parts & belts					
	Orderno.	Colour	Description	Qty.	Price
☐	6999-1012		BINGO Evolution Pre-assembly service (without settings)		1 pc.


CONFIGURATIONS // ACCESSORIES // BINGO SWITCH Mini




Thigh supports



Hip supports



IUS - Upholstery



Calf pad support plate and foam insert




Footrest flip up lock



Reinforced steel footplate-carrier



Footplate with heel stop or boat form



Alu-footplate „foot-print“ with heel stop



Headrest pads low version



Headrest pads high version



Anatomical neck-headrest



Trunk supports



Lap belt



4-point lap belt



Five-point harness



Foot straps



Ankle huger



Groin strap



Fixation jacket



Chest-shoulder harness




Chest-shoulder harness with 4 buckles



Abduction block, swing away and removable




Abduction block, fixed




Retainer for accessories



Grap rail with upholstery



Therapy tray



Therapy tray, height adjustable



Winter warmer



Sheep skin insert



Rain cape

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH Mini



Canopy



BASIC home underframe,
75 mm twin castors



125 mm castors



Side covers in standard
color white



Side covers matching the
upholstery colour



COBRA indoor base,
75 mm twin castors



Plastic front forks



Swivel lock



Aluminum front forks



Frame suspension
system



Side covers matching the
upholstery color



Side covers matching the
frame color



Pneumatic wheel,
front wheel



Pneumatic wheel,
rear wheel



Pneumatic tyres front
wheel, Drum brake wheel



Pneumatic tyres rear
wheel, Drum brake wheel



Polyurethan wheel,
front wheel



Polyurethan wheel,
rear wheel



Polyurethan wheel front
wheel, black rim



Polyurethan wheel rear
wheel, black rim



Polyurethan tyres front
wheel, Drum brake wheel



Polyurethan tyres rear
wheel, Drum brake wheel




Attendant brake,
two hand use



Soft accessory bag



Accessory bag



Tray for additional
equipment




Oxygen bottle holder



Adjustable Angle
brackets



Infusions bottle holder



Reinforced rear axle

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH Mini

56



Connecting parts for
mobility bases

57



Buggy-Board

58



Sun shade, grey

59



Changing backpack

60



„Cosies“ Thermofleece
hand warmer