



BINGO SWITCH rehab pushchair

Size 2 XL for children at the age of approx. 8 - 12 years

Order form incl. Exportprices (€) valid until revoked, release date 03/2024

☐ Offer	☐ Order
Order number / reference	Date
Invoice address	
Company	Customer number
Street	City, Zip code
E-Mail	Telephone
Fax	
Delivery address	
Company	
Street	City, Zip code



*The illustration shows the BINGO SWITCH in size 1

Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable
- Load capacity 60 kg

Seat unit

- IUS - Individual Upholstery System
- Seat depth 30 - 45 cm
- Seat width 23 - 35 cm
- Lower leg length 16 - 40 cm
- Back height 62 - 83 cm *
- Recline 90° - 135°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 50 kg

* From a back height of 75 cm the long push handle is recommended

Minimum information (Without this information your order can not be handled.)

Frame colour:	☐ Black dull matt 032	☐ Silver 002	
Upholstery single-coloured:	☐ Black 346	☐ Berry soft 347	☐ Grey soft 348
with accessories colour black	☐ Turquoise soft 349		
Upholstery two-coloured:	☐ Berry soft / Black 350	☐ Grey soft / Black 351	☐ Turquoise soft / Black 352
with accessories colour black			

Seat depth: (30 - 45 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height: (62 - 83 cm)

Lower leg length: (16 - 40 cm)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
	☐ 3235-2000		BINGO Evolution seat unit, size 2 XL		1 pc.	
with	☐ 3235-2900		Thigh supports, size 2 XL, wide- and depth-adjustable	01	1 pair	incl.
			* minimum seat width 25 cm			
or	☐ 3235-2600		Hip supports, size 2 XL, wide- and depth-adjustable, removable	(upcharge) 02	1 pair	
with	☐ 3235-5200		Upholstery, size 2 XL, adaptable (IUS)	03	1 set	incl.
	☐ 3234-5190		Calf pad support plate and foam insert	04	1 set	
with	☐ 3234-2500K		Footrest hanger 90°	05	1 pc.	incl.
or	☐ 3234-2402K		Footrest hanger angle adjustable 85 - 170° *	(upcharge) 06	1 pc.	
			* minimum lower leg length 19 cm			
with	☐ 3231-2810		Foot rest bracket angle adjustable, flip up type	07	1 pc.	incl.
or	☐ 3231-2820		Foot rest bracket angle adjustable, flip up type, with flip up lock	(upcharge)	1 pc.	
or	☐ 3231-2830		Reinforced foot rest bracket "Steel SA-0430"angle adjustable, flip up type	08	1 pc.	incl.
or	☐ 3231-2840		Reinforced foot rest bracket "Steel SA-0430"angle adjustable, flip up type, with flip up lock	(upcharge)	1 pc.	
			* Extra weight ca. 0,3 kg			
or	3231-2850		Foot rest bracket angle adjustable, with flip up lock, for short lower leg lengths (LLL reduction by 4 cm)	(upcharge)	1 pc.	
with	☐ 3231-249x		Evo foot rest with heel stop	09	1 pc.	incl.
			width ☐ 3=31 ☐ 4=33 ☐ 5=35 (cm)			
or	☐ 3231-248x		Evo foot rest with heel stop for short LLL	09	1 pc.	incl.
			width ☐ 3=31 ☐ 4=33 ☐ 5=35 (cm)			
or	☐ 3233-2551		Foot plate in boat form, width ca. 35 cm	(upcharge) 09	1 pc.	
or	☐ 3233-2553		Bingo-Evo foot plate in boat form size 2, for short LLL	(upcharge) 09	1 pc.	
or	☐ 3233-2552	015	Alu foot plate in boat form, size 2 (not for short LLL)	(upcharge)	1 pc.	
or	☐ 3231-2465	015	Alu-footplate "footprint" with heel stop, width ca. 34 cm	(upcharge) 10	1 pc.	

Accessories - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
	☐ 3233-7600		Headrest pads low version, size 2	11	1 set	
	☐ 3203-7600		Headrest pads, size 2	12	1 set	
	☐ 3203-7601		Headrest pads with spacer fabric for speakers, size 2	13	1 set	
	☐ 3233-6900		Anatomical neck-headrest with length compens, size 2	14	1 set	
	☐ 3233-7700		Trunk supports, size 2, stepless height adjustable (10 x 12 x 2 cm)	15	1 pair	
	☐ 3234-7700		Trunk supports, size 3, stepless height adjustable (12 x 13 x 2 cm)		1 pair	
	☐ 3201-6100		Lap belt	16	1 pc.	
	☐ 3231-6702		4-point lap belt, size S (for seat width 22 - 30)	17	1 pc.	
	☐ 3231-6703		4-point lap belt, size M (for seat width 28 - 36)		1 pc.	
	☐ 3201-6200		Five-point harness	18	1 pc.	
	☐ 3201-6300		Foot straps	19	1 pair	
	☐ 3203-6400		Ankle hugger, size 2	20	1 pair	
	☐ 3233-5300		Groin strap, size 2	21	1 pc.	
	☐ 3201-5500		Fixation jacket, size 2	22	1 pc.	
	☐ 3203-5400		Fixation jacket, size 3		1 pc.	
	☐ 3201-5700		Chest-shoulder harness, size 2	23	1 pc.	
	☐ 3201-5701		Chest-shoulder harness, size 2 with 4 buckels	24	1 pc.	
	☐ 3203-5600		Chest-shoulder harness, size 3		1 pc.	
	☐ 3203-5601		Chest-shoulder harness, size 3 with 4 buckels		1 pc.	
	☐ 3231-7900		Abduction block, size 2, swing away and removable	25	1 pc.	
	☐ 3231-7901		Abduction block, size 2, fixed	26	1 pc.	
	☐ 3233-7800		Abduction block, size 3, swing away and removable		1 pc.	
	☐ 3233-7801		Abduction block, size 3, fixed		1 pc.	
	☐ 3232-2203	032	Retainer for accessories, grab rail or therapy tray	27	1 pair	
	☐ 3233-8100		Grab rail with upholstery, size 2 ¹⁾	28	1 pc.	

Accessories - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
☐	3233-8200		Therapy tray, size 2, transparent ^{1) 2) 3)}	29	1 pc.	
☐	3233-8202		Therapy tray, size 2, transparent, height adjustable ^{1) 2)}	30	1 pc.	
☐	3233-8201		Therapy tray, size 2, grey ^{1) 2) 3)}		1 pc.	
☐	3233-8203		Therapy tray, size 2, grey, height adjustable ^{1) 2)}		1 pc.	
			¹⁾ only possible with retainer for accessories			
			²⁾ only possible with seat in driving direction			
			³⁾ not possible when using hip supports			
☐	3201-8400	306	Winter warmer, size 2, black	31	1 pc.	
☐	3203-8300	306	Winter warmer, size 3, black		1 pc.	
☐	7000-0019		Sheep skin insert, size 2	32	1 pc.	
☐	7000-0020		Sheep skin insert, size 3		1 pc.	
☐	3203-9200	306	Rain cape, size 2, black	33	1 pc.	
☐	3235-9300	306	Canopy incl. adapter and rain cover, black	34	1 pc.	
☐	3003-0010	032	BASIC home underframe size 2 with 75 mm twin castors	35	1 pc.	
☐	3003-0011	032	BASIC home underframe size 2 with 125 mm castors	36	1 pc.	
☐	3002-9500	022	Side covers in standard color white	37		
☐	3002-9500		Side covers matching the upholstery colour	38		
☐	3002-9500	032	Side covers in black			
☐	3002-5400		Push handle for BASIC home under frame		1 pc.	
☐	3005-0000	032	COBRA indoor base size 2, with 75 mm twin castors	39	1 pc.	
☐	3005-0001	032	COBRA indoor base size 2, with 125 mm castors	36	1 pc.	
	Presetting of the gas spring height adjustment: ☐ up to 20kg, ☐ up to 35kg, ☐ up to 50kg					

Basic variations - mobility base with steering wheel attachment

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
☐	3036-1000		BINGO SWITCH size 2XL mobility base with steering wheel attachment		1 pc.	
☐	1014,75		BINGO SWITCH size 2XL mobility base for additional devices with steering wheel attachment		1 pc.	
with	☐	1281-1027	Plastic front forks	40	1 set	incl.
	☐	3202-7100	Swivel lock for Plastic front forks	41	1 set	
or	☐	1281-1028	Aluminum front forks	(upcharge) 42	1 set	
	☐	3232-7100	Swivel lock for Aluminum front forks		1 set	
	☐	SA-0042/1-MP	Push handle extended by 7 cm		1 pc.	
with	☐	3031-1301	Without suspension system		1 set	incl.
or	☐	3031-7400	Frame suspension system, fourfold adjustable	(upcharge) 43	1 set	
	☐	1131-9018	022 Spring set (white) for load up to 25 kg			
	☐	1131-9018	032 Spring set (black) for load up to 25 kg			
	☐	1131-9018	Spring set (matching the upholstery color) for load up to 25 kg			
	☐	1131-9019	022 Spring set (white) for load above 25 kg			
	☐	1131-9019	032 Spring set (black) for load above 25 kg			
	☐	1131-9019	Spring set (matching the upholstery color) for load above to 25 kg			
with	☐	3231-9500	Side covers matching the upholstery color (Standard)	44	1 pair	incl.
or	☐	3231-9500	Side covers matching the frame color	45	1 pair	incl.
with	☐	1219-0004S	Pneumatic wheels with 6-spoke rims	46	1 set	incl.
	☐		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	(upcharge) 47	1 set	
or	☐	1219-0030	Polyurethan wheels with 6-spoke rims	48	1 set	incl.
or	☐	1229-0004S	Polyurethan wheels with 7-spoke rims	49	1 set	incl.
	☐		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	(upcharge) 50	1 set	
or	☐	1229-0030	Attendant brake 2XL, two hand use	51	1 set	
with	☐	3234-9401	Soft accessory bag size 2 (max. 3 kg)	52	1 pc.	incl.
or	☐	3034-9402	Accessory bag size 2 (max. 8 kg) *	(upcharge) 53	1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only			
	☐	3204-4100	Tray for additional equipment size 2 (only with 3037-1000) (Custom-made)	(upcharge) 55	1 pc.	
or	☐	3231-1501	Without tie down kit		1 set	incl.
with	☐	3231-7100	Anti tip left and right *		1 set	
			* not possible in conjunction with frames for additional devices			

Basic variations - mobility base with center steering wheel attachment

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	☐ 3036-1020	032	BINGO SWITCH size 2XL mobility base with center steering wheel attachment	1 pc.	
	☐ 3037-1020	032	BINGO SWITCH size 2XL mobility base for additional devices with center steering wheel attachment	1 pc.	
	☐ SA-0149/2-KP		Swivellock for center steering wheel	1 set	
	☐ SA-0042/1-MP		Upper part push handle extended by 7 cm	1 pc.	
with	☐ 3031-1301	032	Without suspension system	1 set	incl.
or	☐ 3031-7400		Frame suspension system, fourfold adjustable (upcharge)	43 1 set	
	☐ 1131-9018	022	Spring set (white) for load up to 25 kg		
	☐ 1131-9018	032	Spring set (black) for load up to 25 kg		
	☐ 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	☐ 1131-9019	022	Spring set (white) for load above 25 kg		
	☐ 1131-9019	032	Spring set (black) for load above 25 kg		
	☐ 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	☐ 3231-9500		Side covers matching the upholstery color (Standard)	44 1 pair	incl.
or	☐ 3231-9500	032	Side covers matching the frame color	45 1 pair	incl.
with	☐ 1219-0056		Pneumatic wheels with 6-spoke rims	46 1 set	incl.
or	☐ 1219-0057		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	47 1 set	
or	☐ 1229-0056		Polyurethan wheels with 6-spoke rims	48 1 set	incl.
or	☐ 1229-0057		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	50 1 set	
	☐ 3235-7201		Attendant brake BINGO Evo 2XL, two hand use	51 1 set	
with	☐ 3234-9401		Soft accessory bag, size 2 (max. 3 kg)	52 1 pc.	incl.
or	☐ 3034-9402		Accessory bag, size 2 (max. 8 kg) * (upcharge)	53 1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only		
or	☐ 3204-4100		Tray for additional equipment size 2 (only with 3037-1000) (Custom-made) (upcharge)	55 1 pc.	
with	☒ 3231-1501	032	Without tie down kit	1 set	incl.
with	☐ 3231-7100		Anti tip left and right *	1 set	
			* not possible in conjunction with frames for additional devices		

Basic variations - mobility base CROSS with Aluminium front forks

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
	☐ 3036-1030		BINGO SWITCH CROSS size 2XL mobility base with steering wheel attachment		1 pc.	
	☐ 3037-1030		BINGO SWITCH CROSS size 2XL mobility base for additional devices with steering wheel attachment		1 pc.	
	☐ SA-0042/1-MP		Upper part push handle extended by 7 cm		1 pc.	
with	☐ 3031-1301		Without suspension system		1 set	incl.
or	☐ 3031-7400		Frame suspension system, fourfold adjustable	(upcharge)	43 1 set	
	☐ 1131-9018	022	Spring set (white) for load up to 25 kg			
	☐ 1131-9018	032	Spring set (black) for load up to 25 kg			
	☐ 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg			
	☐ 1131-9019	022	Spring set (white) for load above 25 kg			
	☐ 1131-9019	032	Spring set (black) for load above 25 kg			
	☐ 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg			
with	☐ 3231-9500		Side covers matching the upholstery color (Standard)		44 1 pair	incl.
or	☐ 3231-9500	032	Side covers matching the frame color		45 1 pair	incl.
with	☐ 1219-0050		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function		1 set	incl.
or	☐ 1229-0050		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function		1 set	incl.
with	☐ 3234-9401		Soft accessory bag, size 2 (max. 3 kg)		52 1 pc.	incl.
or	☐ 3034-9402		Accessory bag, size 2 (max. 8 kg) *	(upcharge)	53 1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only			
	☐ 3239-4400		BINGO Evolution Cross Tray for additional equipment size 2 (max. 25 kg) (only with 3037-1030/40) (Custom-made)		55 1 pc.	
or				(upcharge)		
with	☐ 3231-1501		Without tie down kit		1 set	incl.
with	☐ 3232-7100		Swivel lock for Aluminum front forks		1 set	
with	☐ 3234-7100		Anti tip left and right *		1 set	
			* not possible in conjunction with frames for additional devices			

Basic variations - mobility base CROSS with Aluminium front forks

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
☐	3036-1040		BINGO SWITCH CROSS size 2XL mobility base with center steering wheel attachment	1 pc.	
☐	3037-1040		BINGO SWITCH CROSS size 2XL mobility base for additional devices with center steering wheel attachmentsteering	1 pc.	
☐	SA-0042/1-MP		Upper part push handle extended by 7 cm	1 pc.	
with ☐	3031-1301		Without suspension system	1 set	incl.
or ☐	3031-7400		Frame suspension system, fourfold adjustable (upcharge) 43	1 set	
☐	1131-9018	022	Spring set (white) for load up to 25 kg		
☐	1131-9018	032	Spring set (black) for load up to 25 kg		
☐	1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
☐	1131-9019	022	Spring set (white) for load above 25 kg		
☐	1131-9019	032	Spring set (black) for load above 25 kg		
☐	1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with ☐	3231-9500		Side covers matching the upholstery color (Standard)	44	1 pair incl.
or ☐	3231-9500	032	Side covers matching the frame color	45	1 pair incl.
with ☐	1219-0058		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function		1 set incl.
or ☐	1229-0058		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function		1 set incl.
with ☐	3234-9401		Soft accessory bag, size 2 (max. 3 kg)	52	1 pc. incl.
or ☐	3034-9402		Accessory bag, size 2 (max. 8 kg) * (upcharge) 53	1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only		
☐	3239-4400		BINGO Evolution Cross Tray for additional equipment size 2 (max. 25 kg) (only with 3037-1030/40) (Custom-made) (upcharge) 55	1 pc.	
or ☒	3231-1501		Without tie down kit	1 set	incl.
with ☐	SA-0149/2-KP		Swivellock for center steering wheel (Custom-made)	1 set	
with ☐	3234-7100		Anti tip left and right *	1 set	
			* not possible in conjunction with frames for additional devices		

! BINGO Evolution Cross: Sonderbauten nur auf Anfrage !

Accessories - mobility base

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
☐	3202-7100		Swivel lock for Plastic front forks	41	1 set
☐	3232-7100		Swivel lock for Aluminum front forks		1 set
☐	SA-0149/2-KP		Swivellock for center steering wheel		1 set
☐	3031-7400E		Frame suspension system, fourfold adjustable	43	1 set
☐	1131-9018	022	Spring set (white) for load up to 25 kg		
☐	1131-9018	032	Spring set (black) for load up to 25 kg		
☐	1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
☐	1131-9019	022	Spring set (white) for load above 25 kg		
☐	1131-9019	032	Spring set (black) for load above 25 kg		
☐	1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
☐	3034-9402E		Accessory bag, size 2 (max. 8 kg) *	53	1 pc.
☐	3231-9500E		Side covers matching the upholstery or frame colour		1 pair
☐	3233-2552E	015	Aluminium foot plate in boat form, size 2		1 pc.

Custom-made

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
☐	3204-4100E		Tray for additional equipment size 2 (only with 3037-10xx)	55	1 pc.	
☐	3239-4400E		BINGO Evolution Cross tray for additional equipment size 2 (max. 25 kg) only with 3037-1030/40		1 pc.	
☐	3030-4200		Oxygen bottle holder (d=100 mm) long, right hand attachment especially for bottles up to 75 cm *	56	1 pc.	
☐	3030-4200_CROSS		Oxygen bottle holder (d=100 mm) long, right hand attachment especially for bottles up to 75 cm *	56	1 pc.	
☐	3030-4201		Oxygen bottle holder short version (d=100 mm), right hand attachment especially for 2 L bottles *		1 pc.	
☐	3030-4201_CROSS		Oxygen bottle holder short version (d=100 mm), right hand attachment especially for 2 L bottles *		1 pc.	
☐	SA-0225-MP		Two adjustable angle holders for contoured oxygen bottle with velcro fasteners	(upcharge) 57	1 set	
☐	SA-0031-KP		Infusions bottle holder	58	1 pc.	
☐	SA-0016/4-MP		Reinforced rear axle (10 kg more load capacity)	(upcharge) 59	1 pc.	
☐	SA-0298-MP		Seat adapter 10 cm heightened *	(upcharge)	1 set	
			* only possible with seat facing the direction of travel			
☐	SA-0355-MP		Mobility base and seat 5 cm wider *	(upcharge)	1 pc.	
☐	SA-0355/2-MP		Tray frame 5 cm wider *	(upcharge)	1 pc.	
			* only possible in frame colour black (032)			

Change attachments

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
☐	3032-1300E		Steering wheel change attachment BINGO SWITCH		1 pair	
with ☐	1281-1027		Plastik front forks	40	1 set	incl.
☐	3202-7100		Swivel lock for Plastic front forks	41	1 set	
or ☐	1281-1028		Aluminium front forks	(upcharge) 42	1 set	
☐	3232-7100		Swivel lock for Aluminium front forks		1 set	
with ☐	1219-0055		Steering wheels with pneumatic tyres on 3-spoke rims	46	1 set	incl.
or ☐	1229-0055		Steering wheels with PU tyres on 3-spoke rims	48	1 set	incl.
or ☐	1229-0054		Steering wheels with PU tyres on 7-spoke rims	49	1 set	incl.
☐	3034-1501E		Center steering wheel interchangeable attachment BINGO SWITCH (2 / 2XL), pneumatic tires		1 set	
☐	3034-1502E		Center steering wheel interchangeable attachment BINGO SWITCH (2 / 2XL), PU tires		1 set	
☐	SA-0149/2-KP		Swivellock for center steering wheel		1 set	

Extra equipment

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
☐	7000-0010		Buggy-Board *	60	1 pc.	
			* not possible with tie down kit			
☐	7000-0012	306	Sun shade, black	61	1 pc.	
☐	7000-0006		Headrest speakers, for headrest pads 3203-7601	62	1 set	
☐	7000-0202	306	Changing backpack, black	63	1 pc.	
☐	7000-0199		"Cosies" Thermofleece hand warmer for companion	64	1 pair	

Pre-assembly of pads, upholstery parts & belts

	Orderno.	Colour	Description		Qty.	Price
☐	6999-1012		BINGO Evolution Pre-assembly service (without settings)		1 pc.	

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH




Thigh supports



Hip supports



IPS - Individuelles Polstersystem



Calf pad support plate



Footrest hanger 90°



Footrest hanger angle adjustable 85 - 170°




Footrest flip up lock



Reinforced steel foot rest bracket



Footplate with heel stop or boat form



Alu-footplate „foot-print“ with heel stop



Headrest pads low version



Headrest pads



Headrest pads with spacer fabric for speakers



Anatomical neck-headrest



Trunk supports



Lap belt



4-point lap belt



Five-point harness



Foot straps



Ankle hugger



Groin strap



Fixation jacket



Chest-shoulder harness



Chest-shoulder harness with 4 buckles



Abduction block, swing away and removable




Abduction block, fixed




Retainer for accessories



Grap rail with upholstery



Therapy tray



Therapy tray, height adjustable

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH



Winter warmer



Sheep skin insert



Rain cape



Canopy



BASIC home underframe,
75 mm twin castors



125 mm castors



Side covers in standard
color white



Side covers matching the
upholstery colour



COBRA indoor base,
75 mm twin castors



Plastic front forks



Swivel lock



Aluminum front forks



Frame suspension
system



Side covers matching the
upholstery color



Side covers matching
the frame color



Pneumatic wheel,
front wheel



Pneumatic wheel,
rear wheel



Pneumatic tyres front
wheel, Drum brake wheel



Pneumatic tyres rear
wheel, Drum brake wheel



Polyurethan wheel,
front wheel



Polyurethan wheel,
rear wheel



Polyurethan wheel front
wheel, black rim



Polyurethan wheel rear
wheel, black rim



Polyurethan tyres front
wheel, Drum brake wheel



Polyurethan tyres rear
wheel, Drum brake wheel



Attendant brake,
two hand use




Soft accessory bag



Accessory bag



Tie down kit according to
ISO 7176-19



Tray for additional
equipment

CONFIGURATIONS // ACCESSORIES // BINGO SWITCH




Oxygen bottle holder



Adjustable Angle brackets



Infusions bottle holder



Reinforced rear axle



Buggy-Board



Sun shade, grey



Headrest speakers



Changing backpack



„Cosies“ Thermofleece hand warmer