

BINGO SWITCH

Release date: February 2026 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____



Delivery address

see invoice address

Company: _____

Street: _____

City, Zip code: _____

Basic configuration

- Individual Upholstery System (IUS)
- Thigh supports
- Footrest hanger 90°
- Angle adjustable Footplate
- Side covers

BINGO SWITCH Trendline



Trendline package

- Upholstery in reed green
- Canopy in reed green
- Frame colour black
- Hand loop



Frame colour black

Upholstery colour reed green

Frame colour



Black dull matt



Silver


Upholstery colour



Black



Berry soft



Berry soft // Black



Grey soft



Grey soft // Black



Turquoise soft



Turquoise soft // Black



Presettings // size 1

Seat width, size 1

enter seat width: (18 - 30 cm)

Back height, size 1

enter back height: (45 - 63 cm)

Seat depth, size 1

enter seat depth: (16 - 34 cm)

Lower leg length, size 1

enter lower leg length: (16 - 37 cm)



Notice:

If dimensions are missing, we adjust to medium dimensions.

Presettings // size 2

Seat width, size 2

enter seat width: (23 - 35 cm)

Back height, size 2

enter back height: (54 - 75 cm)*

Seat depth, size 2

enter seat depth: (22 - 40 cm)

Lower leg length, size 2

enter lower leg length: (16 - 40 cm)



Notice:

If dimensions are missing, we adjust to medium dimensions.

* From a back height of 65 cm, the long push handle is recommended.

Presettings // size 2 XL ! Considered as custom-made with corresponding accessories !

Seat width, size 2 XL

enter seat width: (23 - 35 cm)

Back height, size 2 XL

enter back height: (62 - 83 cm)**

Seat depth, size 2 XL

enter seat depth: (30 - 45 cm)

Lower leg length, size 2 XL

enter lower leg length: (16 - 40 cm)




Notice:

If dimensions are missing, we adjust to medium dimensions.

** From a back height of 75 cm, the long push handle is recommended.

Seat unit



BINGO SWITCH seat unit

su size 1

su size 2

su size 2 XL

Headrest pads



Headrest pads,
low version



Headrest pads



Headrest pads with spacer
fabric, for speakers



Anatomical neck-headrest
with length compens

Supports



Trunk supports

8 x 10 x 2 cm

10 x 12 x 2 cm

12 x 13 x 2 cm



Thigh supports, wide- and depth-
adjustable, removable



Hip supports, wide- and depth-
adjustable, removable



Notice:

Size 1: Minimum seat width
is 20 cm.

Size 2 / Size 2 XL: Minimum
seat width is 25 cm.

Supports



Standard side supports, removable, foldable with stop clips




PositionPlus side supports, height and angle adjustable, width adjustment by 10 cm


Calf pad support plate, Calf upholstery



Calf pad support plate and foam insert



Calf upholstery suitable to seat width



Calf upholstery for single foot rest hanger



Notice:

If calf upholsteries are fitted to the footrest hanger, the minimum LLL increases by 6 cm.

Footrest hanger



Footrest hanger 90°



Footrest hanger angle adjustable 85 - 175°



Single footrest hangers adjustable 90-160°

Footplates



Footplate with heel stop

1=27 cm
2=29 cm
3=31 cm
4=33 cm
5=35 cm

for short
lower leg lengths

1=27 cm
2=29 cm
3=31 cm
4=33 cm
5=35 cm



Alu footplate with heel stop

1=27 cm
2=29 cm
3=31 cm
4=33 cm
5=35 cm



Footplate in boat form

su width approx. 30 cm
size 1

su width approx. 35 cm
size 2
size 2 XL

for short
lower leg lengths

**Notice:**


Widths 1=27 cm and 2=29 cm are only possible for the size 1 seat unit. Widths 4=33 cm and 5=35 cm are only possible for size 2 and size 2 XL seating units.



Alu foot plate in boat form

su width approx. 30 cm
size 1

su width approx. 35 cm
size 2
size 2 XL

Alu-footplate „fingerprint“
with heel stop

su width approx. 30 cm
size 1

su width approx. 34 cm
size 2
size 2 XL



Single footplates in boat form

Footrest flip up lock, Steel footplate-carrier



Foot rest bracket angle adjustable,
flip up type

with flip up lock

**with flip up lock,
for short lower leg lengths
(reduction by 4 cm)**



Reinforced steel footplate-carrier
(about 0,3 kg heavier) angle
adjustable, flip up type

with flip up lock



Footboard holder for single footrests

Belting



Lap belt



4-point lap belt



Five-point harness

size XS
sw 18 - 24 cm

size S
sw 22 - 30 cm

size M
sw 28 - 36 cm

Belting



Groin strap



Foot straps



Ankle hugger

su	size 1
size 1	ø 4 - 6,5 cm
su	size 2
size 2	ø 6 - 8,5 cm
size 2 XL	



Fixation jacket

- size 0**
- size 1**
- size 2**
- size 3**



Chest-shoulder harness

- size 1**
- size 2**
- size 3**



Chest-shoulder harness
with 4 buckles

- size 1**
- size 2**
- size 3**

Abduction block



Abduction block,
swing away and removable

size 1

size 2

size 3



Abduction block, fixed

size 1

size 2

size 3



Notice:

The abduction block in size 3 is only for the seat unit in size 2 and size 2 XL possible.

Retainer for accessories, Grab rail



Retainer for accessories,
grab rail or therapy tray



Grab rail with imitation leather,
pivotable, black



Notice:

Only possible with retainer for accessories.


Therapy tray




Therapy tray, transparent



Therapy tray, transparent,
height adjustable



Therapy tray, anthracite



Therapy tray, anthracite,
height adjustable



Notice:

Only possible with retainer for accessories / with seat in the direction of travel.
Hip supports only possible in combination with the height adjustable therapy tray.

Winter warmer, Rain cape, Canopy



Winter warmer

size 1

size 2

size 3



Sheep skin insert

size 1

size 2

size 3



Rain cape, black

Canopy incl. adapter
and rain cover**Notice:**

If the **Trendline** package is selected, the canopy is supplied in reed green.

Mobility bases




Steering wheel attachment

Steering wheel attachment,
for additional devices


Center steering wheel attachment

Center steering wheel attachment,
for additional devicesCross with steering wheel
attachmentCross with steering wheel
attachment, for additional
devices

Mobility bases



Cross with center steering wheel attachment



Cross with center steering wheel attachment, for additional devices



Equipment carrier attachment



Equipment carrier attachment, for additional devices



Cross with equipment carrier attachment



Cross with equipment carrier attachment, for additional devices

Steering wheel change attachment



Steering wheel change attachment



Device carrier interchangeable attachment



Mounting adapter for device carrier attachment for „phil&teds“ sibling seat or universal car seat adapter



Plastik front forks



Swivel lock for Plastic front forks



Aluminium front forks



Swivel lock for Aluminium front forks



Steering wheels on 3-spoke rims



Steering wheels with PU tyres on 7-spoke rims

!
Notice:
Wheel removable via quick release.

Pneumatic wheels

Polyurethan wheels

Center steering wheel interchangeable attachment



Center steering wheel interchangeable attachment

Pneumatic wheels

Polyurethan wheels



Swivellock for center steering wheel

Side covers



Side covers matching the upholstery color



Side covers matching the frame color

Frame suspension system

without suspension system



Frame suspension system, fourfold adjustable



Please select the desired spring set:

Spring set white

Spring set black

Spring set matching the upholstery color

for load up to 15 kg



for load up to 25 kg

for load above to 25 kg



Notice:

Only available for the mobility base in **size 1**.

Front forks, Swivel lock



Plastic front forks



Aluminum front forks

**Notice:**


Cross mobility base only possible with aluminum front forks. These are included with the Cross mobility base.



Swivel lock for Plastic front forks



Swivel lock for Aluminum front forks



Swivellock for center steering wheel

Tires (not selectable for the Cross mobility base)



Standard wheels

Drum brake wheels with
tyres and with locking functionPolyurethan wheels with
7-spoke rims

Pneumatic wheels

Polyurethan wheels

**Notice:**

Wheel removable via quick release.

Pneumatic wheels

Polyurethan wheels

**Notice:**

Only possible on the mobility base with steering wheel attachment.

Wheel can be removed via quick-release fastener.

Tires (only selectable for the Cross mobility base)

Drum brake wheels with tyres on black 5-spoke rims, with locking function

Pneumatic wheels

Polyurethan wheels

**Notice:**

Wheel removable via quick release.

Brake variant

Attendant brake, two hand use

**Notice:**

Not for drum brake wheels.

Anti tip, Push handle, Tie down kit

Anti tip left and right



Push handle extended by 7 cm



Tie down kit according to ISO 7176-19

**Notice:**

Not possible with the mobility base for additional equipment.

**Notice:**

Not possible with the mobility base with center steering wheel attachment or equipment carrier attachment .

Accessory bags



Soft accessory bag (max. 3 kg)



Accessory bag (max. 8 kg)

**Notice:**

Not possible with tie down kit & only possible with seat facing the direction of travel.

Tray for additional equipment, Oxygen bottle holder



Tray for additional equipment



Oxygen bottle holder (d=100 mm)



Oxygen bottle holder for contoured oxygen devices

for bottles up to 75 cm

for 2 liter bottles
(short version)

for bottles up to 75 cm

for 2 liter bottles
(short version)

Attachment left

Attachment right



Attachment left

Attachment right

Infusions pole, Reinforced rear axle, Seat adapter, Mobility base and seat wider (only possible for the 2 XL mobility base)



Infusion pole with 2-armed holder, articulated joint

Reinforced rear axle
(10 kg more load capacity)Seat adapter 10 cm
heightenedMobility base and seat
5 cm wider**Notice:**

Not possible for the Cross mobility base.

**Notice:**

Only possible in black frame color.
Additional costs apply for all widened accessories.

Connecting parts



Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base



Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution swivelling or swivelling mobility base for additional equipment rear



Connecting parts for coupling a BINGO Evolution mobility base and a HOGGI wheelchair

**Notice:**

Only possible with 2 mobility bases of one size.
Not possible in combination with aluminium front forks & mobility base Cross.

**Notice:**

Not possible in combination with aluminium front forks & mobility base Cross.

SWINGBO VTi**CLEO**

Sun shade, Changing backpack, Headrest speakers



Sun shade, black



Changing backpack, black



Headrest speakers, for headrest pads

Cosies, Buggy-Board



„Cosies“ Thermofleece hand warmer for companion



Buggy-Board

BASIC home underframe



BASIC home underframe

**with 75 mm
twin castors**

**with 125 mm
twin castors**



Side covers

in standard color white

**matching the upholstery
colour**

in black



Push handle for BASIC
home under frame

COBRA indoor base



COBRA indoor base

**with 75 mm
twin castors**

**with 125 mm
twin castors**

Presetting of the gas spring height adjustment:

up to 20 kg

up to 35 kg

up to 50 kg

Remark / Customer request / Individual special custom-made request

All prices are subject to change and non-binding. Errors, misprints, and price changes are reserved.
Binding prices are only those stated in our written offer or order confirmation.