

www.hoggi.de

Made in Germany
by HOGGI 

BINGO EVOLUTION

Innovative and modern rehab pushchair.



Designed and hand built in Germany.

BINGO EVOLUTION CROSS

The rehab pushchair for all challenges in a proven concept.

BINGO Evolution Cross: The innovative and modern rehab push chair

HOGGI extends its range of rehab buggies for kids and adolescents with the addition of **BINGO Evolution Cross** - the new solid rehab buggy meeting the **challenging needs of daily applications**. Furthermore the version **BINGO Evolution Cross** convinces with its extreme stable back frame and large rear wheels for all outdoor situations.

The drum brake (basic configuration of **BINGO Evolution Cross**) ensures reliable control even on inclined terrain. Thanks to the compact construction of the mobility base, **BINGO Evolution Cross** provides extreme agility - an advantage in crowded pedestrian zones as well as in narrows of the terrain.

The adaptability and the ability to grow with the child, correspond to those of **BINGO Evolution**.



	Size 1	Size 2	Size 2 XL
Seat width	18 - 30 cm	23 - 35 cm	23 - 35 cm
Seat depth	16 - 34 cm	22 - 40 cm	30 - 45 cm
Tilt in space seat unit			
Seat unit front facing	45°	45°	45°
Seat unit rear facing	45°	45°	-
Back height	45 - 63 cm	54 - 75 cm	62 - 83
Recline			
Seat unit front facing	90° - 160°	90° - 160°	90° - 135°
Seat unit rear facing	90° - 160°	90° - 160°	-
Lower leg length	16 - 37 cm	16 - 40 cm	16 - 40 cm
Height of push bar	75 - 110 cm	75 - 110 cm	86 - 118 cm
Overall width	63 cm	68 cm	68 cm
Folding size (L x W x H)	78 x 62 x 38 cm	78 x 67 x 38 cm	78 x 67 x 38 cm
Weight limit			
Seat unit	50 kg	50 kg	50 kg
Mobility base	60 kg	60 kg	60 kg