



BINGO Evolution rehab pushchair

Size 2 for children at the age of approx. 4 - 10 years

Crash tested according to ISO 7176-19 and ANSI/RESNA WC/Vol.1. - Section 19

Order form incl. Exportprices (€) valid until revoked, release date 10/2023

<input type="checkbox"/> Offer	<input type="checkbox"/> Order
Order number / reference	Date
<input type="text"/>	<input type="text"/>
Invoice address	
Company	Customer number
<input type="text"/>	<input type="text"/>
Street	City, Zip code
<input type="text"/>	<input type="text"/>
E-Mail	Telephone
<input type="text"/>	<input type="text"/>
Fax	
<input type="text"/>	
Delivery address	
Company	
<input type="text"/>	
Street	City, Zip code
<input type="text"/>	<input type="text"/>



Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable, rotatable

Seat unit

- IUS - Individual Upholstery System
- Seat depth 22 - 40 cm
- Seat width 23 - 35 cm
- Lower leg length 16 - 40 cm
- Back height 54 - 75 cm *
- Recline 90° - 160°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 50 kg

*From a back height of 65 cm the long push handle is recommended

Minimum information (Without this information your order can not be handled.)

Frame colour:	<input type="checkbox"/> Black dull matt 032	<input type="checkbox"/> Silver 002	
Upholstery single-coloured:	<input type="checkbox"/> Black 346	<input type="checkbox"/> Berry soft 347	<input type="checkbox"/> Grey soft 348
with accessories colour black	<input type="checkbox"/> Turquoise soft 349		
Upholstery two-coloured:	<input type="checkbox"/> Berry soft / Black 350	<input type="checkbox"/> Grey soft / Black 351	<input type="checkbox"/> Turquoise soft / Black 352
with accessories colour black			

Seat depth: (22 - 40 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height: (54 - 75 cm)

Lower leg length: (16 - 40 cm)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
	<input type="checkbox"/> 3234-2000		BINGO Evolution seat unit, size 2		1 pc.	
with	<input type="checkbox"/> 3233-2900		Thigh supports, size 2, wide- and depth-adjustable		01 1 pair	incl.
			* minimum seat width 25 cm			
or	<input type="checkbox"/> 3233-2600		Hip supports, size 2, wide- and depth-adjustable, removable	(upcharge)	02 1 pair	
with	<input type="checkbox"/> 3234-5100		Upholstery, size 2, adaptable (IUS)		03 1 set	incl.
	<input type="checkbox"/> 3234-5190		Calf pad support plate and foam insert		04 1 set	
with	<input type="checkbox"/> 3234-2500		Footrest hanger 90°		05 1 pc.	incl.
or	<input type="checkbox"/> 3234-2402		Footrest hanger angle adjustable 85 - 170° *	(upcharge)	06 1 pc.	
			* minimum lower leg length 19 cm			
or	<input type="checkbox"/> 3232-2504		Footrest hanger 90° (short LLL) with bracket and footrest flip up lock *	(upcharge)	1 pc.	
or	<input type="checkbox"/> 3232-2404		Footrest hanger angle adjustable 85-170° (short LLL) with bracket and footrest flip up lock *	(upcharge)	1 pc.	
			* lower leg length 15-36 cm			
with	<input type="checkbox"/> 3231-249x		Footplate with heel stop, width: <input type="checkbox"/> 3=31, <input type="checkbox"/> 4=33, <input type="checkbox"/> 5=35 (cm)		07 1 pc.	incl.
or	<input type="checkbox"/> 3233-2551		Footplate in boat form, width: approx. 35 cm	(upcharge)	07 1 pc.	
or	<input type="checkbox"/> 3233-2552	015	Alu foot plate in boat form, size 2	(upcharge)	1 pc.	
	<input type="checkbox"/> 3240-2950		Footrest flip up lock	(upcharge)	08 1 pc.	
with	<input type="checkbox"/> SA-0430		Reinforced steel footplate-carrier (about 0,3 kg heavier)		09 1 pc.	incl.

Accessories - seat unit

	Orderno.	Colour	Description		Qty.	Price
	<input type="checkbox"/> 3233-7600		Headrest pads low version, size 2		10 1 set	
	<input type="checkbox"/> 3203-7600		Headrest pads, size 2		11 1 set	
	<input type="checkbox"/> 3203-7601		Headrest pads with spacer fabric for speakers, size 2		12 1 set	
	<input type="checkbox"/> 3233-6900		Anatomical neck-headrest with length compens, size 2		13 1 set	
	<input type="checkbox"/> 3233-7700		Trunk supports, size 2, stepless height adjustable (10 x 12 x 2 cm)		14 1 pair	
	<input type="checkbox"/> 3234-7700		Trunk supports, size 3, stepless height adjustable (12 x 13 x 2 cm)		1 pair	
	<input type="checkbox"/> 3201-6100		Lap belt		15 1 pc.	
	<input type="checkbox"/> 3231-6702		4-point lap belt, size S (for seat width 22 - 30)		16 1 pc.	
	<input type="checkbox"/> 3231-6703		4-point lap belt, size M (for seat width 28 - 36)		1 pc.	
	<input type="checkbox"/> 3201-6200		Five-point harness		17 1 pc.	
	<input type="checkbox"/> 3201-6300		Foot straps		18 1 pair	
	<input type="checkbox"/> 3203-6400		Ankle hugger, size 2		19 1 pair	
	<input type="checkbox"/> 3233-5300		Groin strap, size 2		20 1 pc.	
	<input type="checkbox"/> 3201-5500		Fixation jacket, size 2		21 1 pc.	
	<input type="checkbox"/> 3203-5400		Fixation jacket, size 3		1 pc.	
	<input type="checkbox"/> 3201-5700		Chest-shoulder harness, size 2		22 1 pc.	
	<input type="checkbox"/> 3201-5701		Chest-shoulder harness, size 2 with 4 buckels		23 1 pc.	
	<input type="checkbox"/> 3203-5600		Chest-shoulder harness, size 3		1 pc.	
	<input type="checkbox"/> 3203-5601		Chest-shoulder harness, size 3 with 4 buckels		1 pc.	
	<input type="checkbox"/> 3231-7900		Abduction block size 2, swing away and removable		24 1 pc.	
	<input type="checkbox"/> 3231-7901		Abduction block size 2, fixed		25 1 pc.	
	<input type="checkbox"/> 3233-7800		Abduction block size 3, swing away and removable		1 pc.	
	<input type="checkbox"/> 3233-7801		Abduction block size 3, fixed		1 pc.	
	<input type="checkbox"/> 3232-2203	032	Retainer for accessories, grab rail or therapy tray		26 1 pair	
	<input type="checkbox"/> 3233-8100		Grab rail with upholstery, size 2 ¹⁾		27 1 pc.	
	<input type="checkbox"/> 3233-8200		Therapy tray, size 2, transparent ^{1) 2) 3)}		28 1 pc.	
	<input type="checkbox"/> 3233-8202		Therapy tray, size 2, transparent, height adjustable ^{1) 2)}		29 1 pc.	
	<input type="checkbox"/> 3233-8201		Therapy tray, size 2, grey ^{1) 2) 3)}		1 pc.	
	<input type="checkbox"/> 3233-8203		Therapy tray, size 2, grey, height adjustable ^{1) 2)}		1 pc.	
			¹⁾ only possible with retainer for accessories			
			²⁾ only possible with seat in driving direction			
			³⁾ not possible when using hip supports			
	<input type="checkbox"/> 3201-8400	306	Winter warmer, size 2, black		30 1 pc.	
	<input type="checkbox"/> 3203-8300	306	Winter warmer, size 3, black		1 pc.	
	<input type="checkbox"/> 7000-0019		Sheep skin insert, size 2		31 1 pc.	
	<input type="checkbox"/> 7000-0020		Sheep skin insert, size 3		1 pc.	
	<input type="checkbox"/> 3203-9200	306	Rain cape, size 2, black		32 1 pc.	
	<input type="checkbox"/> 3201-9300	306	Canopy incl. adapter and rain cover, black		33 1 pc.	

Accessories - seat unit

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	3003-0010	032	BASIC home underframe size 2 with 75 mm twin castors	34	1 pc.
<input type="checkbox"/>	3003-0011	032	BASIC home underframe size 2 with 125 mm castors	35	1 pc.
<input type="checkbox"/>	3002-9500	022	Side covers in standard color white	36	
<input type="checkbox"/>	3002-9500		Side covers matching the upholstery colour	37	
<input type="checkbox"/>	3002-9500	032	Side covers in black		
<input type="checkbox"/>	3002-5400		Push handle for BASIC home under frame		1 pc.
<input type="checkbox"/>	3005-0000	032	COBRA indoor base, size 2, with 75 mm twin castors	38	1 pc.
<input type="checkbox"/>	3005-0001	032	COBRA indoor base, size 2, with 125 mm castors	35	1 pc.
Presetting of the gas spring height adjustment: <input type="checkbox"/> up to 20kg, <input type="checkbox"/> up to 35kg, <input type="checkbox"/> up to 50kg					

Basic variations - mobility base swivelling front wheels

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	3234-1000		BINGO Evolution mobility base, size 2, swivelling front wheels	1	pc.
<input type="checkbox"/>	3239-1000		BINGO Evolution mobility base, size 2, swivelling front wheels, for respiratory equipment	1	pc.
with <input type="checkbox"/>	1281-1027		Plastic front forks	39	1 set incl.
<input type="checkbox"/>	3202-7100		Swivel lock for Plastic front forks	40	1 set
or <input type="checkbox"/>	1281-1028		Aluminum front forks (upcharge)	41	1 set
<input type="checkbox"/>	3232-7100		Swivel lock for Aluminum front forks		1 set
<input type="checkbox"/>	SA-0042/1-MP		Push handle extended by 7 cm		1 pc.
with <input type="checkbox"/>	3231-1301		Without suspension system		1 set incl.
or <input type="checkbox"/>	3231-7400		Frame suspension system, fourfold adjustable (upcharge)	42	1 set
<input type="checkbox"/>	1131-9018	022	Spring set (white) for load up to 25 kg		
<input type="checkbox"/>	1131-9018	032	Spring set (black) for load up to 25 kg		
<input type="checkbox"/>	1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
<input type="checkbox"/>	1131-9019	022	Spring set (white) for load above 25 kg		
<input type="checkbox"/>	1131-9019	032	Spring set (black) for load above 25 kg		
<input type="checkbox"/>	1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with <input type="checkbox"/>	3231-9500		Side covers matching the upholstery color (Standard)	43	1 pair incl.
or <input type="checkbox"/>	3231-9500		Side covers matching the frame color	44	1 pair incl.
with <input type="checkbox"/>	1219-0004S		Pneumatic wheels with 6-spoke rims	45	1 set incl.
or <input type="checkbox"/>	1219-0030		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	46	1 set
or <input type="checkbox"/>	1229-0004S		Polyurethan wheels with 6-spoke rims	47	1 set incl.
or <input type="checkbox"/>	1229-0009		Polyurethan wheels with 7-spoke rims	48	1 set incl.
or <input type="checkbox"/>	1229-0030		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	49	1 set
<input type="checkbox"/>	3231-7201		Attendant brake, two hand use	50	1 set
<input type="checkbox"/>	3235-7201		Attendant brake 2, two hand use (for extended push handles)		1 set
with <input type="checkbox"/>	3234-9401		Soft accessory bag, size 2 (max. 3 kg)	51	1 pc. incl.
or <input type="checkbox"/>	3233-9402		Accessory bag, size 2 (max. 8 kg) * (upcharge)	52	1 pc.
* do not combine with tie down kit, use with seat in driving direction only					
with <input type="checkbox"/>	3231-1501		Without tie down kit		1 set incl.
or <input type="checkbox"/>	3231-7300		Tie down kit according to ISO 7176-19 (upcharge)	53	1 set

Basic variations - mobility base CROSS with Aluminum front forks

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3234-1030		BINGO Evolution CROSS, size 2 mobility base	1 pc.	
	<input type="checkbox"/> 3239-1030		BINGO Evo CROSS, size 2 mobility base for additional equipment	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	42 1 set	
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	43 1 pair	incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	44 1 pair	incl.
or	<input type="checkbox"/> 1219-0050		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	46 1 set	incl.
or	<input type="checkbox"/> 1229-0050		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	49 1 set	incl.
with	<input type="checkbox"/> 3234-9401		Soft accessory bag, size 2 (max. 3 kg)	51 1 pc.	incl.
or	<input type="checkbox"/> 3233-9402		Accessory bag, size 2 (max. 8 kg) * (upcharge)	52 1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3231-1501		Without tie down kit	1 set	incl.
or	<input type="checkbox"/> 3231-7300		Tie down kit according to ISO 7176-19 (upcharge)	53 1 set	
with	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set	
with	<input type="checkbox"/> 3234-7100		Anti tip left and right *	1 set	
			* not possible in combination with frame for respiratory equipment		

! BINGO Evolution Cross: Custom-made only on request !
Basic variations - mobility base rigid front wheels

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3233-1000		BINGO Evolution mobility base size 2, rigid front wheels	1 pc.	
	<input type="checkbox"/> SA-0042/1-MP		Push handle extended by 7 cm	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	42 1 set	
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	43 1 pair	incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	44 1 pair	incl.
with	<input type="checkbox"/> 1219-0025S		Pneumatic wheels with 6-spoke rims	45 1 set	incl.
or	<input type="checkbox"/> 1219-0029		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	46 1 set	
or	<input type="checkbox"/> 1229-0021S		Polyurethan wheels with 6-spoke rims	47 1 set	incl.
or	<input type="checkbox"/> 1229-0029		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	49 1 set	
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use	50 1 set	
	<input type="checkbox"/> 3235-7201		Attendant brake 2, two hand use (for extended push handles)	1 set	
with	<input type="checkbox"/> 3234-9401		Soft accessory bag, size 2 (max. 3 kg)	51 1 pc.	incl.
or	<input type="checkbox"/> 3233-9402		Accessory bag, size 2 (max. 8 kg) * (upcharge)	52 1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3231-1501		Without tie down kit	1 set	incl.
or	<input type="checkbox"/> 3231-7300		Tie down kit according to ISO 7176-19 (upcharge)	53 1 set	

Basic variations - mobility base TRIKE (custom-made product)

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3237-1000		BINGO Evolution mobility base TRIKE, size 2	1 pc.	
	<input type="checkbox"/> SA-0042/1-MP		Push handle extended by 7 cm	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	42 1 set	
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	43 1 pair	incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	44 1 pair	incl.
with	<input type="checkbox"/> 1219-0011S		Pneumatic wheels with 6-spoke rims	45 1 set	incl.
or	<input type="checkbox"/> 1219-0031		Drum brake wheels, pneumatic tyres on black 5-spoke rims (upcharge)	46 1 set	
or	<input type="checkbox"/> 1229-0011S		Polyurethan wheels with 6-spoke rims	47 1 set	incl.
or	<input type="checkbox"/> 1229-0031		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	49 1 set	
	<input type="checkbox"/> 3231-7201		Attendant brake, two hand use	50 1 set	
	<input type="checkbox"/> 3235-7201		Attendant brake 2, two hand use (for extended push handles)	1 set	
with	<input type="checkbox"/> 3234-9401		Soft accessory bag, size 2 (max. 3 kg)	51 1 pc.	incl.
or	<input type="checkbox"/> 3233-9402		Accessory bag, size 2 (max. 8 kg) * (upcharge)	52 1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only		
	<input checked="" type="checkbox"/> 3231-1501		Without tie down kit	1 set	incl.

Accessories - mobility base

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	40 1 set	
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set	
	<input type="checkbox"/> 3231-7400E		Frame suspension system, fourfold adjustable, retrofit	42 1 set	
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
	<input type="checkbox"/> 3233-9402E		Accessory bag, size 2 (max. 8 kg), retrofit	52 1 pc.	
	<input type="checkbox"/> 3231-7300E		Tie down kit according to ISO 7176-19	53 1 set	
	<input type="checkbox"/> 3231-9500E		Side covers matching the upholstery colour or the frame color *	1 pair	
			* please enter colour code		
	<input type="checkbox"/> 3233-2552E	015	Alu foot plate in boat form, size 2	1 pc.	

Custom-made

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	3204-4100		Tray for additional equipment, size 2 (mobility base 3239-10X0 only)	54 1 pc.	
<input type="checkbox"/>	3239-4400		BINGO Evolution Cross tray for additional equipment size 2 (max. 25kg)	1 pc.	
<input type="checkbox"/>	3231-4200		Oxygen bottle holder (d=100 mm) for bottles up to 75 cm *	55 1 pc.	
			* assembly on the right side only		
<input type="checkbox"/>	SA-0225-MP		Adjustable Angle brackets for contoured oxygen units (upcharge)	56 1 set	
<input type="checkbox"/>	3231-4201		Oxygen bottle holder (d=100 mm) for short bottles up to 40 cm *	1 pc.	
			* assembly on the right side only		
<input type="checkbox"/>	SA-0031-KP		Infusions bottle holder	57 1 pc.	
<input type="checkbox"/>	SA-0016/4-MP		Reinforced rear axle for mobility base max. load 70 kg (upcharge)	58 1 pc.	
<input type="checkbox"/>	SA-0298-MP		Seat adapter 10 cm heightened * (upcharge)	1 set	
			* with seat in driving direction only		
<input type="checkbox"/>	SA-0321-MP		Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base * (upcharge)	59 1 set	
			* only possible with 2 mobility bases of one size, not for base 3230-1000		
			* not possible in combination with aluminium front forks & mobility base Cross		
<input type="checkbox"/>	SA-0441-MP		Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution mobility base 3232-1000 or 3238-1000 rear * (upcharge)	1 set	
			* not possible in combination with aluminium front forks & mobility base Cross		

Extra equipment

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	7000-0010		Buggy-Board	60 1 pc.	
<input type="checkbox"/>	7000-0012	306	Sun shade, black	61 1 pc.	
<input type="checkbox"/>	7000-0006		Headrest speakers, for headrest pads 3201-7601	62 1 set	
<input type="checkbox"/>	7000-0202	306	Changing backpack, black	63 1 pc.	
<input type="checkbox"/>	7000-0199		"Cosies" Thermofleece hand warmer for companion	64 1 pair	

Pre-assembly of pads, upholstery parts & belts

	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	6999-1012		BINGO Evolution Pre-assembly service (without settings)	1 pc.	

CONFIGURATIONS // ACCESSORIES // BINGO Evolution



Thigh supports



Hip supports



IUS - Upholstery



Calf pad support plate and foam insert



Footrest hanger 90°



Footrest hanger angle adjustable 85 - 170°



Footplate with heel stop or boat form



Footrest flip up lock



Reinforced steel footplate-carrier



Headrest pads low version



Headrest pads



Headrest pads with spacer fabric for speakers



Anatomical neck-headrest



Trunk supports



Lap belt



4-point lap belt



Five-point harness



Foot straps



Ankle hugger



Groin strap



Fixation jacket



Chest-shoulder harness




Chest-shoulder harness with 4 buckles



Abduction block, swing away and removable




Abduction block, fixed



Retainer for accessories



Grap rail with upholstery



Therapy tray



Therapy tray, height adjustable



Winter warmer

CONFIGURATIONS // ACCESSORIES // BINGO Evolution




CONFIGURATIONS // ACCESSORIES // BINGO Evolution



Adjustable Angle
brackets



Infusions bottle holder



Reinforced rear axle



Connecting parts for
mobility bases



Buggy-Board



Sun shade, grey



Headrest speakers



Changing backpack



„Cosies“ Thermofleece
hand warmer