

HOGGI®

Partnering together to mobilise kids



BINGO Evolution rehab pushchair - custom-made

Size 2 XL for children at the age of approx. 8 - 12 years

Crash tested according to ISO 7176-19 and ANSI/RESNA WC/Vol.1. - Section 19
Order form incl. Exportprices (€) valid until revoked, release date 10/2023

☐ Offer

☐ Order

Order number / reference

Date

Invoice address

Company

Customer number

Street

City, Zip code

E-Mail

Telephone

Fax

Delivery address

Company

Street

City, Zip code



Basic configuration

Mobility base

- Compact folding
- 45° tilt in space
- Seat quickly detachable

Seat unit

- IUS - Individual Upholstery System
- Seat depth 30 - 45 cm
- Seat width 23 - 35 cm
- Lower leg length 16 - 40 cm
- Back height 62 - 83 cm *
- Recline 90° - 135°
- Ergonomic lap belt positioning
- Angle adjustable foot plate
- Maximum load 50 kg

* From a back height of 75 cm the long push handle is recommended

Minimum information (Without this information your order can not be handled.)

Frame colour:

☐ Black dull matt 032

☐ Silver 002

Upholstery single-coloured:

☐ Black 346

☐ Berry soft 347

☐ Grey soft 348

with accessories colour black

☐ Turquoise soft 349

Upholstery two-coloured:

☐ Berry soft / Black 350

☐ Grey soft / Black 351

☐ Turquoise soft / Black 352

with accessories colour black

Seat depth:

(30 - 45 cm)

Presettings (At missing specifications we will adjust to medium sizes. Side covers matching the upholstery color)

Back height:

(62 - 83 cm)

Lower leg length:

(16 - 40 cm)

Basic variations - seat unit

with/or: Please select only one version!

	Orderno.	Colour	Description		Qty.	Price
	<input type="checkbox"/> 3235-2000		BINGO Evolution seat unit, size 2 XL		1 pc.	
with	<input type="checkbox"/> 3235-2900		Thigh supports, size 2 XL, wide- and depth-adjustable		01 1 pair	incl.
			* minimum seat width 25 cm			
or	<input type="checkbox"/> 3235-2600		Hip supports, size 2 XL, wide- and depth-adjustable, removable	(upcharge)	02 1 pair	
	<input type="checkbox"/> 3235-5200		Upholstery, size 2 XL, adaptable (IUS)		03 1 set	incl.
	<input type="checkbox"/> 3234-5190		Calf pad support plate and foam insert		04 set	
with	<input type="checkbox"/> 3234-2500		Footrest hanger 90°		05 1 pc.	incl.
or	<input type="checkbox"/> 3234-2402		Footrest hanger angle adjustable 85 - 170° *	(upcharge)	06 1 pc.	
			* minimum lower leg length 19 cm			
			Footrest hanger 90° (short LLL) with bracket and footrest flip up lock *	(upcharge)	1 pc.	
or	<input type="checkbox"/> 3232-2504		Footrest hanger angle adjustable 85-170° (short LLL)	(upcharge)	1 pc.	
or	<input type="checkbox"/> 3232-2404		with bracket and footrest flip up lock *	(upcharge)	1 pc.	
			* lower leg length 15-36 cm			
with	<input type="checkbox"/> 3231-249x		Footplate with heel stop, width: <input type="checkbox"/> 3=31, <input type="checkbox"/> 4=33, <input type="checkbox"/> 5=35 (cm)		07 1 pc.	incl.
or	<input type="checkbox"/> 3233-2551		Footplate in boat form, width: approx. 35 cm	(upcharge)	07 1 pc.	
or	<input type="checkbox"/> 3233-2552	015	Alu foot plate in boat form, size 2	(upcharge)	1 pc.	
	<input type="checkbox"/> 3240-2950		Footrest flip up lock	(upcharge)	08 1 pc.	
			Reinforced steel foot rest bracket, (Additional weight approx. 0.3 kg)		09 1 pc.	incl.
with	<input type="checkbox"/> SA-0430					

Accessories - seat unit

	Orderno.	Colour	Description		Qty.	Price
	<input type="checkbox"/> 3233-7600		Headrest pads low version, size 2		10 1 set	
	<input type="checkbox"/> 3203-7600		Headrest pads, size 2		11 1 set	
	<input type="checkbox"/> 3203-7601		Headrest pads with spacer fabric for speakers, size 2		12 1 set	
	<input type="checkbox"/> 3233-6900		Anatomical neck-headrest with length compens, size 2		13 1 set	
	<input type="checkbox"/> 3233-7700		Trunk supports, size 2, stepless height adjustable (10 x 12 x 2 cm)		14 1 pair	
	<input type="checkbox"/> 3234-7700		Trunk supports, size 3, stepless height adjustable (12 x 13 x 2 cm)		1 pair	
	<input type="checkbox"/> 3201-6100		Lap belt		15 1 pc.	
	<input type="checkbox"/> 3231-6702		4-point lap belt, size S (for seat width 22 - 30)		16 1 pc.	
	<input type="checkbox"/> 3231-6703		4-point lap belt, size M (for seat width 28 - 36)		1 pc.	
	<input type="checkbox"/> 3201-6200		Five-point harness		17 1 pc.	
	<input type="checkbox"/> 3201-6300		Foot straps		18 1 pair	
	<input type="checkbox"/> 3203-6400		Ankle hugger, size 2		19 1 pair	
	<input type="checkbox"/> 3233-5300		Groin strap, size 2		20 1 pc.	
	<input type="checkbox"/> 3201-5500		Fixation jacket, size 2		21 1 pc.	
	<input type="checkbox"/> 3203-5400		Fixation jacket, size 3		1 pc.	
	<input type="checkbox"/> 3201-5700		Chest-shoulder harness, size 2		22 1 pc.	
	<input type="checkbox"/> 3201-5701		Chest-shoulder harness, size 2 with 4 buckels		23 1 pc.	
	<input type="checkbox"/> 3203-5600		Chest-shoulder harness, size 3		1 pc.	
	<input type="checkbox"/> 3203-5601		Chest-shoulder harness, size 3 with 4 buckels		1 pc.	
	<input type="checkbox"/> 3231-7900		Abduction block, size 2, swing away and removable		24 1 pc.	
	<input type="checkbox"/> 3231-7901		Abduction block, size 2, fixed		25 1 pc.	
	<input type="checkbox"/> 3233-7800		Abduction block, size 3, swing away and removable		1 pc.	
	<input type="checkbox"/> 3233-7801		Abduction block, size 3, fixed		1 pc.	
	<input type="checkbox"/> 3232-2203	032	Retainer for accessories, grab rail or therapy tray		26 1 pair	
	<input type="checkbox"/> 3233-8100		Grab rail with upholstery, size 2 ¹⁾		27 1 pc.	
	<input type="checkbox"/> 3233-8200		Therapy tray, size 2, transparent ^{1) 2) 3)}		28 1 pc.	
	<input type="checkbox"/> 3233-8202		Therapy tray, size 2, transparent, height adjustable ^{1) 2)}		29 1 pc.	
	<input type="checkbox"/> 3233-8201		Therapy tray, size 2, grey ^{1) 2) 3)}		1 pc.	
	<input type="checkbox"/> 3233-8203		Therapy tray, size 2, grey, height adjustable ^{1) 2)}		1 pc.	
			¹⁾ only possible with retainer for accessories			
			²⁾ only possible with seat in driving direction			
			³⁾ not possible when using hip supports			
	<input type="checkbox"/> 3201-8400	306	Winter warmer, size 2, black		30 1 pc.	
	<input type="checkbox"/> 3203-8300	306	Winter warmer, size 3, black		1 pc.	
	<input type="checkbox"/> 7000-0019		Sheep skin insert, size 2		31 1 pc.	
	<input type="checkbox"/> 7000-0020		Sheep skin insert, size 3		1 pc.	
	<input type="checkbox"/> 3203-9200	306	Rain cape, size 2, black		32 1 pc.	
	<input type="checkbox"/> 3235-9300	306	Canopy incl. adapter and rain cover, black		33 1 pc.	

Accessories - seat unit

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3003-0010	032	BASIC home underframe size 2 with 75 mm twin castors	34 1 pc.	
	<input type="checkbox"/> 3003-0011	032	BASIC home underframe size 2 with 125 mm castors	35 1 pc.	
	<input type="checkbox"/> 3002-9500	022	Side covers in standard color white	36	
	<input type="checkbox"/> 3002-9500		Side covers matching the upholstery colour	37	
	<input type="checkbox"/> 3002-9500	032	Side covers in black		
	<input type="checkbox"/> 3002-5400		Push handle for BASIC home under frame		1 pc.
	<input type="checkbox"/> 3005-0000	032	COBRA indoor base size 2, with 75 mm twin castors	38 1 pc.	
	<input type="checkbox"/> 3005-0001	032	COBRA indoor base size 2, with 125 mm castors	35 1 pc.	
Presetting of the gas spring height adjustment: <input type="checkbox"/> up to 20kg, <input type="checkbox"/> up to 35kg, <input type="checkbox"/> up to 50kg					

Basic variations - mobility base swivelling front wheels

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3235-1000		BINGO Evolution mobility base, size 2-XL, swivelling front wheels	1 pc.	
	<input type="checkbox"/> 3239-1006		BINGO Evolution mobility base, size 2-XL, swivelling front wheels, for respiratory equipment	1 pc.	
with	<input type="checkbox"/> 1281-1027		Plastic front forks	39 1 set	incl.
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	40 1 set	
or	<input type="checkbox"/> 1281-1028		Aluminum front forks (upcharge)	41 1 set	
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set	
	<input type="checkbox"/> SA-0042/1-MP		Push handle extended by 7 cm	1 pc.	
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	42 1 set	
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	43 1 pair	incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	44 1 pair	incl.
with	<input type="checkbox"/> 1219-0004S		Pneumatic wheels with 6-spoke rims	45 1 set	incl.
or	<input type="checkbox"/> 1219-0030		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function (upcharge)	46 1 set	
or	<input type="checkbox"/> 1229-0004S		Polyurethan wheels with 6-spoke rims	47 1 set	incl.
or	<input type="checkbox"/> 1229-0009		Polyurethan wheels with 7-spoke rims	48 1 set	incl.
or	<input type="checkbox"/> 1229-0030		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function (upcharge)	49 1 set	
	<input type="checkbox"/> 3235-7201		Attendant brake 2XL, two hand use	50 1 set	
with	<input type="checkbox"/> 3234-9401		Soft accessory bag size 2 (max. 3 kg)	51 1 pc.	incl.
or	<input type="checkbox"/> 3233-9402		Accessory bag size 2 (max. 8 kg) * (upcharge)	52 1 pc.	
* do not combine with tie down kit, use with seat in driving direction only					
with	<input type="checkbox"/> 3231-1501		Without tie down kit	1 set	incl.
or	<input type="checkbox"/> 3231-7300		Tie down kit according to ISO 7176-19 (upcharge)	53 1 set	

Basic variations - mobility base CROSS with Aluminum front forks

with/or: Please select only one version!

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3235-1030		BINGO Evolution CROSS Gr.2 XL mobility base	1 pc.	
	<input type="checkbox"/> 3239-1008		BINGO Evo CROSS, size 2 XL mobility base for additional equipment		
with	<input type="checkbox"/> 3231-1301		Without suspension system	1 set	incl.
or	<input type="checkbox"/> 3231-7400		Frame suspension system, fourfold adjustable (upcharge)	42 1 set	
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
with	<input type="checkbox"/> 3231-9500		Side covers matching the upholstery color (Standard)	43 1 pair	incl.
or	<input type="checkbox"/> 3231-9500		Side covers matching the frame color	44 1 pair	incl.
or	<input type="checkbox"/> 1219-0050		Drum brake wheels with pneumatic tyres on black 5-spoke rims, with locking function	46 1 set	incl.
or	<input type="checkbox"/> 1229-0050		Drum brake wheels with polyurethan tyres on black 5-spoke rims, with locking function	49 1 set	incl.
with	<input type="checkbox"/> 3234-9401		Soft accessory bag size 2 (max. 3 kg)	51 1 pc.	incl.
or	<input type="checkbox"/> 3233-9402		Accessory bag size 2 (max. 8 kg) * (upcharge)	52 1 pc.	
			* do not combine with tie down kit, use with seat in driving direction only		
with	<input type="checkbox"/> 3231-1501		Without tie down kit	1 set	incl.
or	<input type="checkbox"/> 3231-7300		Tie down kit according to ISO 7176-19 (upcharge)	53 1 set	
with	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set	
with	<input type="checkbox"/> 3234-7100		Anti tip left and right *	1 set	
			* not possible in combination with frame for respiratory equipment		

! BINGO Evolution Cross: Custom-made only on request !
Accessories - mobility base

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3202-7100		Swivel lock for Plastic front forks	40 1 set	
	<input type="checkbox"/> 3232-7100		Swivel lock for Aluminum front forks	1 set	
	<input type="checkbox"/> 3231-7400E		Frame suspension system, fourfold adjustable, retrofit	42 1 set	
	<input type="checkbox"/> 1131-9018	022	Spring set (white) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018	032	Spring set (black) for load up to 25 kg		
	<input type="checkbox"/> 1131-9018		Spring set (matching the upholstery color) for load up to 25 kg		
	<input type="checkbox"/> 1131-9019	022	Spring set (white) for load above 25 kg		
	<input type="checkbox"/> 1131-9019	032	Spring set (black) for load above 25 kg		
	<input type="checkbox"/> 1131-9019		Spring set (matching the upholstery color) for load above to 25 kg		
	<input type="checkbox"/> 3233-9402E		Accessory bag size 2 (max. 8 kg), retrofit	52 1 pc.	
	<input type="checkbox"/> 3231-7300E		Tie down kit according to ISO 7176-19	53 1 set	
	<input type="checkbox"/> 3231-9500E		Side covers matching the upholstery colour or the frame color *	1 pair	
			* please enter colour code		
	<input type="checkbox"/> 3233-2552E	015	Alu foot plate in boat form, size 2	1 pc.	

Custom-made

	Orderno.	Colour	Description	Qty.	Price
	<input type="checkbox"/> 3204-4100		Tray for additional equipment, size 2 (mobility base 3239-100X only)	54 1 pc.	
	<input type="checkbox"/> 3239-4400		BINGO Evolution Cross tray for additional equipment size 2 (max. 25kg)	1 pc.	
	<input type="checkbox"/> 3231-4200		Oxygen bottle holder (d=100 mm) for bottles up to 75 cm *	55 1 pc.	
			* assembly on the right side only		
	<input type="checkbox"/> SA-0225-MP		Adjustable Angle brackets for contoured oxygen units (upcharge)	56 1 set	
	<input type="checkbox"/> 3231-4201		Oxygen bottle holder (d=100 mm) for short bottles up to 40 cm *	1 pc.	
			* assembly on the right side only		
	<input type="checkbox"/> SA-0031-KP		Infusions bottle holder	57 1 pc.	
	<input type="checkbox"/> SA-0016/4-MP		Reinforced rear axle for mobility base max. load 70 kg (upcharge)	58 1 pc.	
	<input type="checkbox"/> SA-0298-MP		Seat adapter 10 cm heightened * (upcharge)	1 set	
			* with seat in driving direction only		
	<input type="checkbox"/> SA-0355-MP		Mobility base and seat 5 cm wider * (upcharge)	1 pc.	
	<input type="checkbox"/> SA-0355/2-MP		Tray frame 5 cm wider * (upcharge)	1 pc.	
			* only possible in frame colour black (032)		

Extra equipment					
	Orderno.	Colour	Description	Qty.	Price
<input type="checkbox"/>	7000-0010		Buggy-Board	59 1 pc.	
<input type="checkbox"/>	7000-0012	306	Sun shade, black	60 1 pc.	
<input type="checkbox"/>	7000-0006		Headrest speakers, for headrest pads 3203-7601	61 1 set	
<input type="checkbox"/>	7000-0202	306	Changing backpack, black	62 1 pc.	
<input type="checkbox"/>	7000-0199		"Cosies" Thermofleece hand warmer for companion	63 1 pair	

Pre-assembly of pads, upholstery parts & belts					
	Orderno.	Colour	Description	Qty.	
<input type="checkbox"/>	6999-1012		BINGO Evolution Pre-assembly service (without settings)	1 pc.	

CONFIGURATIONS // ACCESSORIES // BINGO Evolution




Thigh supports



Hip supports



IPS - Individuelles Polstersystem



Calf pad support plate



Footrest hanger 90°



Footrest hanger angle adjustable 85 - 170°



Footplate with heel stop or boat form



Footrest flip up lock



Reinforced steel foot rest bracket



Headrest pads low version



Headrest pads



Headrest pads with spacer fabric for speakers



Anatomical neck-headrest



Trunk supports



Lap belt



4-point lap belt



Five-point harness



Foot straps



Ankle hugger



Groin strap



Fixation jacket



Chest-shoulder harness




Chest-shoulder harness with 4 buckles



Abduction block, swing away and removable




Abduction block, fixed




Retainer for accessories



Grap rail with upholstery



Therapy tray



Therapy tray, height adjustable



Winter warmer

CONFIGURATIONS // ACCESSORIES // BINGO Evolution

31		32		33		34		35	
	Sheep skin insert		Rain cape		Canopy		BASIC home underframe, 75 mm twin castors		125 mm castors
36		37		38		39		40	
	Side covers in standard color white		Side covers matching the upholstery colour		COBRA indoor base, 75 mm twin castors		Plastic front forks		Swivel lock
41		42		43		44		45	
	Aluminum front forks		Frame suspension system		Side covers matching the upholstery color		Side covers matching the frame color		Pneumatic wheels front wheel
45		46				47			
	Pneumatic wheel, rear wheel		Pneumatic tyres front wheel, Drum brake wheel		Pneumatic tyres rear wheel, Drum brake wheel		Polyurethan wheel, front wheel		Polyurethan wheel, rear wheel
48				49				50	
	Polyurethan wheel front wheel, black rim		Polyurethan wheel rear wheel, black rim		Polyurethan tyres front wheel, Drum brake wheel		Polyurethan tyres rear wheel, Drum brake wheel		Attendant brake, two hand use
51		52		53		54		55	
	Soft accessory bag		Accessory bag		Tie down kit according to ISO 7176-19		Tray for additional equipment		Oxygen bottle holder


CONFIGURATIONS // ACCESSORIES // BINGO Evolution



Adjustable Angle brackets



Infusions bottle holder



Reinforced rear axle



Buggy-Board



Sun shade, grey



Headrest speakers



Changing backpack



„Cosies“ Thermofleece hand warmer