

BINGO Evolution

Release date: February 2026 // valid until revoked

Order

Offer

Order number / reference: _____

Date: _____

Invoice address

Customer number: _____

Company: _____

Street: _____

City, Zip code: _____

Delivery address see invoice address

Company: _____

Street: _____

City, Zip code: _____



- Basic configuration**
- Individual Upholstery System (IUS)
 - Thigh supports
 - Footrest hanger 90°
 - Angle adjustable Footplate
 - Side covers

BINGO Evolution Trendline



- Trendline package**
- Upholstery in reed green
 - Canopy in reed green
 - Frame colour black
 - Hand loop
 - Whitewall rims (optional)



Frame colour black

Upholstery colour reed green

Frame colour



Black dull matt



Silver

Upholstery colour



Black



Berry soft



Berry soft // Black



Grey soft



Grey soft // Black



Turquoise soft



Turquoise soft // Black



Presettings // size 1

Seat width, size 1

enter seat width: (18 - 30 cm)

Back height, size 1

enter back height: (45 - 63 cm)

Seat depth, size 1

enter seat depth: (16 - 34 cm)

Lower leg length, size 1

enter lower leg length: (16 - 37 cm)



Notice:

If dimensions are missing, we adjust to medium dimensions.

Presettings // size 2

Seat width, size 2

enter seat width: (23 - 35 cm)

Back height, size 2

enter back height: (54 - 75 cm)*

Seat depth, size 2

enter seat depth: (22 - 40 cm)

Lower leg length, size 2

enter lower leg length: (16 - 40 cm)



Notice:

If dimensions are missing, we adjust to medium dimensions.

* From a back height of 65 cm, the long push handle is recommended.

Presettings // size 2 XL ! Considered as custom-made with corresponding accessories !

Seat width, size 2 XL

enter seat width: (23 - 35 cm)

Back height, size 2 XL

enter back height: (62 - 83 cm)**

Seat depth, size 2 XL

enter seat depth: (30 - 45 cm)

Lower leg length, size 2 XL

enter lower leg length: (16 - 40 cm)




Notice:

If dimensions are missing, we adjust to medium dimensions.

** From a back height of 75 cm, the long push handle is recommended.

Seat unit



BINGO Evolution seat unit

su size 1

su size 2

su size 2 XL

Headrest pads



Headrest pads,
low version



Headrest pads



Headrest pads with spacer
fabric, for speakers



Anatomical neck-headrest
with length compens

Supports



Trunk supports

8 x 10 x 2 cm

10 x 12 x 2 cm

12 x 13 x 2 cm



Thigh supports, wide- and depth-
adjustable, removable



Notice:

Size 1: Minimum seat width
is 20 cm.

Size 2 / Size 2 XL: Minimum
seat width is 25 cm.



Hip supports, wide- and depth-
adjustable, removable

Supports



Standard side supports, removable, foldable with stop clips




PositionPlus side supports, height and angle adjustable, width adjustment by 10 cm


Calf pad support plate, Calf upholstery



Calf pad support plate and foam insert



Calf upholstery suitable to seat width



Calf upholstery for single foot rest hanger

! **Notice:**
If calf upholsteries are fitted to the footrest hanger, the minimum LLL increases by 6 cm.

Footrest hanger



Footrest hanger 90°



Footrest hanger angle adjustable 85 - 175°



Single footrest hangers adjustable 90-160°

Footplates



Footplate with heel stop

- 1=27 cm
- 2=29 cm
- 3=31 cm
- 4=33 cm
- 5=35 cm

for short lower leg length

- 1=27 cm
- 2=29 cm
- 3=31 cm
- 4=33 cm
- 5=35 cm



Alu footplate with heel stop

- 1=27 cm
- 2=29 cm
- 3=31 cm
- 4=33 cm
- 5=35 cm



Footplate in boat form

su width approx. 30 cm
size 1

su width approx. 35 cm
size 2
size 2 XL

for short lower leg length



Notice:


Widths 1=27 cm and 2=29 cm are only possible for the size 1 seat unit. Widths 4=33 cm and 5=35 cm are only possible for size 2 and size 2 XL seating units.



Alu footplate in boat form

su width approx. 30 cm
size 1

su width approx. 35 cm
size 2
size 2 XL



Alu footplate „footprint“ with heel stop

su width approx. 30 cm
size 1

su width approx. 34 cm
size 2
size 2 XL



Single footplates in boat form

Footrest flip up lock, Steel footplate-carrier



Footrest bracket angle adjustable, flip up type

with flip up lock

with flip up lock, for short LLL (reduction by 4 cm)



Reinforced footrest bracket steel angle adjustable, flip up type (extra weight approx. 0,3 kg)

with flip up lock



Footboard holder for single footrests

Belting



Lap belt



4-point lap belt



Five-point harness

- size XS**
sw 18 - 24 cm
- size S**
sw 22 - 30 cm
- size M**
sw 28 - 36 cm

Belting



Groin strap



Foot straps



Ankle huger

su **size 1**
size 1 ø 4 - 6,5 cm

su **size 2**
size 2 ø 6 - 8,5 cm
size 2 XL



Fixation jacket

size 0
size 1
size 2
size 3



Chest-shoulder harness

size 1
size 2
size 3



Chest-shoulder harness
 with 4 buckles

size 1
size 2
size 3

Abduction block



Abduction block, swing away and removable

size 1

size 2

size 3



Abduction block, fixed

size 1

size 2

size 3

Notice:
The abduction block in size 3 is only for the seat unit in size 2 and size 2 XL possible.

Retainer for accessories, Grab rail



Retainer for accessories, grab rail or therapy tray



Grab rail with imitation leather, pivotable, black

Notice:
Only possible with retainer for accessories.


Therapy tray




Therapy tray, transparent



Therapy tray, transparent, height adjustable



Therapy tray, anthracite



Therapy tray, anthracite, height adjustable

Notice:
Only possible with retainer for accessories / with seat in the direction of travel.
Hip supports only possible in combination with the height adjustable therapy tray.

Winter warmer, Rain cape, Canopy



Winter warmer

- size 1
- size 2
- size 3



Sheep skin insert

- size 1
- size 2
- size 3




Rain cape, black



Canopy incl. adapter and rain cover

! Notice:
If the **Trendline** package is selected, the canopy is supplied in reed green.

Mobility bases



Rigid




Swivelling



Swivelling for additional equipment



Cross



Cross for additional equipment

Side covers



Side covers matching the upholstery color



Side covers matching the frame color

Frame suspension system

Without suspension system



Frame suspension system, fourfold adjustable



Please select the desired spring set:

Spring set white

Spring set black

Spring set matching the upholstery color

for load up to 15 kg



for load up to 25 kg

for load above to 25 kg



Notice:

Only available for the mobility base in **size 1**.

Front forks, Swivel lock (selectable for the Swivelling and the Cross mobility base)



Plastic front forks



Swivel lock for Plastic front forks




Aluminum front forks



Swivel lock for Aluminum front forks

! **Notice:**
Cross mobility base only possible with aluminum front forks. These are included with the Cross mobility base.


Tires (selectable for the swivelling, the rigid mobility base)




Standard wheels



Drum brake wheels with tyres and with locking function



Polyurethane wheels with 7-spoke rims



Polyurethane wheels with whitewall rims

incl.

Pneumatic wheels

Polyurethane wheels

Pneumatic wheels

Polyurethane wheels

! **Notice:**
Wheel removable via quick release.

! **Notice:**
Only possible with the swivelling mobility base.
Wheel can be removed via quick-release fastener.

! **Notice:**
Only possible with the **Trendline** version of the swivelling mobility base.
Wheel can be removed via quick-release fastener.

Tires (only selectable for the Cross mobility base)



Drum brake wheels with tyres on black 5-spoke rims, with locking function

Pneumatic wheels

Polyurethan wheels

!
Notice:
Wheel removable via quick release.

Brake variant



Attendant brake, two hand use

!
Notice:
Not for drum brake wheels.

Anti tip, Push handle, Tie down kit



Anti tip left and right



Push handle extended by 7 cm



Tie down kit according to ISO 7176-19

!
Notice:
Can only be selected for the Cross mobility base. Not possible with the mobility base for additional equipment

Accessory bags



Soft accessory bag (max. 3 kg)



Accessory bag (max. 8 kg)

!
Notice:
Not possible with tie down kit & only possible with seat facing the direction of travel.

Tray for additional equipment, Oxygen bottle holder



Tray for additional equipment



Oxygen bottle holder (d=100 mm)



Oxygen bottle holder for contoured oxygen devices

for bottles up to 75 cm

for 2 liter bottles
(short version)

for bottles up to 75 cm

for 2 liter bottles
(short version)

Attachment left

Attachment right

Attachment left


Attachment right

!
Notice:
Only possible with the swivelling & CROSS mobility base to accommodate the tray for additional equipment.

Infusions bottle holder, Reinforced rear axle, Seat adapter



Infusion pole with 2-armed holder, articulated joint



Reinforced rear axle
(10 kg more load capacity)



Seat adapter 10 cm heightened

Mobility base and seat wider (only possible for the 2 XL mobility base)



Notice:
Only possible in black frame color.
Additional costs apply for all widened accessories.

Mobility base and seat 5 cm wider


Connecting parts



Connecting parts for 2 BINGO Evolution mobility bases incl. swivel lock on rear mobility base



Connecting parts for Gesslein baby stroller Model F4 Air+ / S4 Air+ ahead with BINGO Evolution swivelling or swivelling mobility base for additional equipment rear



Connecting parts for coupling a BINGO Evolution mobility base and a HOGGI wheelchair



Notice:
Only possible with 2 mobility bases of one size.
Not possible in combination with aluminium front forks & mobility base Cross.



Notice:
Not possible in combination with aluminium front forks & mobility base Cross.

SWINGBO VTi

CLEO

Sun shade, Changing backpack, Headrest speakers



Sun shade, black



Changing backpack, black



Headrest speakers, for headrest pads

Cosies, Buggy-Board



„Cosies“ Thermofleece hand warmer
for companion



Buggy-Board

BASIC home underframe



BASIC home underframe

**with 75 mm
twin castors**

**with 125 mm
twin castors**



Side covers

in standard color white

**matching the upholstery
colour**

in black



Push handle for BASIC
home under frame

COBRA indoor base



COBRA indoor base

**with 75 mm
twin castors**

**with 125 mm
twin castors**

Presetting of the gas spring height adjustment:

up to 20 kg

up to 35 kg

up to 50 kg

Remark / Customer request / Individual special custom-made request

All prices are subject to change and non-binding. Errors, misprints, and price changes are reserved.
Binding prices are only those stated in our written offer or order confirmation.