

Active drives and add on drives

	BINGO Evolution	ZIP	DURO	SWINGBO 2	SWINGBO 2 XL	SUPRA 2.0	SWINGBO VTi	SWINGBO VTi XL	CESA	CLEO	CLEO Ti	CLEO Teen	SHERPA
Ulrich Alber													
TWION (only for 24") e-motion (22" or 24") Power assist add-on drive	–	–	–	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
e-fix, Electrical add-on drive	–	–	–	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Scalamobil, Stair climber	–	–	–	✗	✗	✗	✓	✓	–	–	–	–	✓
Viamobil, Push- and brake aid	–	–	–	✗	V25	✗	✓	✓	–	–	–	–	V25
Alber via Go	–	–	–	–	–	–	✓	✓	–	–	–	–	–
e-Pilot	–	–	–	–	–	–	–	–	✓	–	–	–	–
SMOOV	–	–	–	–	–	✓	–	–	✓	–	–	–	–
AAT													
MAX-E, Electriv drive	–	–	–	✓	✓	–	✓	✓	–	–	–	–	✓
V-MAX, Push- and brake aid	–	–	–	✓	✓	–	✓	✓	–	–	–	–	✓
V-MAX mini, Push aid	✓	–	✓	–	–	✓	–	–	–	–	–	–	–
MOVILINO, Push- and brake aid	–	–	–	✗	✗	–	–	–	–	–	–	–	–
S-MAX, Stair climber	✓	✓	–	✗	✓	–	–	–	–	–	–	–	✓
SOLO, Wheel hub drive	–	–	–	–	✓	–	✓	✓	–	✓	–	✓	✓
SERVO, Power-assisted auxiliary drive	–	–	–	–	✓	✓	✓	✓	–	✓	–	✓	✓
SANO													
LIFTKAR PT uni 130, Stair climber	✗	–	–	✓	✓	–	✓	✓	–	✓	✓	✓	✓
SUNRISE MEDICAL													
Empulse R20, Push aid	✓	–	✓	✓	✓	–	✓	✓	–	✓	✓	✓	✓
Empulse, Wheel drive	✗	✗	✗	✗	–	✓	–	✓	✗	–	–	✓	✗
AMF-Bruns													
Wheelchair restraint system (AMF)	✓	–	–	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DECON													
Easy Wheel	–	–	–	✓	–	–	✓	–	–	✓	✓	–	✓
E-Drive	–	–	–	✓	–	–	✓	–	–	✓	✓	–	✓
E-Drive 16"	–	–	–	✗	–	–	✗	–	–	✗	✗	–	✓
E-Move	–	–	–	✓	–	–	✓	–	–	✓	✓	–	✓

This table is intended as a general guide. Whether a combination is possible depends on the individual wheelchair configuration and requires consultation with our development department.